



Yes, you can.®



Invacare® **Alegio™**

User Manual (EN)

Bruksanvisning (SV)

Gebruiksaanwijzing (NL)

Manuel d'utilisation (FR)

Manual del usuario (ES)


Manuale d'uso (IT)

Manual de Utilização (PT)



Quality Declaration

Congratulations with your new bed Invacare®Alegio™ from Invacare®.

Your new bed is -marked in accordance with directive 93/42/EEC concerning medical devices. Invacare®Alegio™ is developed and constructed with consideration for the user and others handling or assisting at the bed.

Furthermore, the bed is developed in accordance with the European Standard EN 1970/A1.

Invacare®Alegio™ is supervised and quality controlled throughout the entire production process, and the finished bed is inspected by our finished goods control.

Identification label are located on the bed confirming that the finished goods control has approved the bed.

Please read the entire user's manual before using the bed.

Invacare®France Operation SAS is certified according to ISO 9001 and ISO 13485.

DK

KVALITETSDEKLARATION

Tillykke med Deres nye seng Invacare®Alegio™ fra Invacare® France Operation SAS.

Deres nye seng er CE - mærket og lever op til alle krav i henhold til direktiv 93/42/EOF om medicinske anordninger.

Sengen er udviklet og konstrueret under størst mulig hensyntagen til brugeren samt alle andre, der enten håndterer sengen eller hjælper til ved sengen.

Invacare®Alegio™ er udviklet under hensyntagen til de sikkerhedsmæssige krav i den Europæiske Standard EN 1970/A1. Hver seng har gennem hele produktionsforløbet været overvåget og kontrolleret, og den færdige seng er blevet inspiceret af vores færdigvarekontrol.

Typiske påsættelse sengen som dokumentation for at færdigvarekontrollen har godkendt sengen. Før ibrugtagning af Deres seng skal De gennemlæse brugermanualen grundigt.

Invacare® France Operation SAS er certificeret i henhold til ISO 9001 og ISO 13485.

DE

QUALITÄTSDOKUMENTATION

Herlichen Glückwunsch zu Ihrem neuen Pflegebett Invacare®Alegio™ von Invacare®France Operation SAS.

Ihr neuer Pflegebett ist gemäß der Richtlinie 93/42/EWG für Medizinprodukte CE-gekennzeichnet.

Der Invacare®Alegio™ wurde unter Berücksichtigung der Bedürfnisse von Benutzer und Pflegepersonal entwickelt und in Übereinstimmung mit der Europäischen Norm EN 1970/A1 konzipiert und hergestellt.

Während des gesamten Herstellungsprozesses unterliegt der Invacare®Alegio™ einer ständigen Qualitätskontrolle und wird im Anschluss nochmals durch unsere Endkontrolle geprüft. Das Typenschild sind am Bett angebracht, um die Abnahme durch unsere Qualitätskontrolle zu bestätigen.

Bitte lesen die gesamte Bedienungsanleitung, bevor Sie den Bett in Gebrauch nehmen.

Invacare® France Operation SAS ist zertifiziert nach ISO 9001 und ISO 13485.

NL

KWALITEITSGARANTIE

Gefeliciteerd met de aanschaf van uw nieuwe bed Alegio™ van Invacare® France Operation SAS.

Invacare®Alegio™ is CE - gecertificeerd en goedgekeurd conform richtlijn 93/42/EEC betreffende medische hulpmiddelen.

Invacare®Alegio™ is ontwikkeld en geconstrueerd met inachtneming van de behoeften van de gebruiker en verzorgers.

Invacare®Alegio™ is ontwikkeld in overeenstemming met de Europese Standaard EN 1970/A1.

Gedurende het gehele productieproces is het Alegio™ bed gecontroleerd op kwaliteitsaspecten en het complete bed is na productie zorgvuldig gecontroleerd. Het productlabel zijn op het bed bevestigd om aan te tonen dat de bedden zijn gecontroleerd en goed zijn bevonden door de afdeling productiecontrole. Wij verzoeken u vriendelijk de gehele gebruikershandleiding te lezen voordat u het bed in gebruik neemt.

Invacare® France Operation SAS is ISO 9001 en ISO 13485 gecertificeerd.

IT

DICHIARAZIONE DI QUALITÀ

Complimenti per aver scelto il letto Invacare®Alegio™ prodotto da Invacare® France Operation SAS.

Il vostro nuovo letto è marcato CE - ai sensi della Direttiva 93/42/EEC relativa ai dispositivi medici.

Il letto è stato progettato e costruito con un occhio di riguardo per gli utilizzatori e per i loro assistenti. Invacare®Alegio™ è stato realizzato nel rispetto della Normativa Europea EN 1970/A1.

Il letto è stato oggetto di accurate verifiche qualitative durante l'intero processo produttivo; una volta completato è stato controllato dal nostro servizio prodotti finiti.

Un'apposita targhetta sono state apposte sul letto a conferma dell'avvenuta verifica e accettazione del prodotto da parte del servizio qualità. Prima di utilizzare il letto vi invitiamo a leggere integralmente il manuale d'uso.

Invacare® France Operation SAS è un'azienda certificata ai sensi della Normativa ISO 9001 e ISO 13485.

ES

DECLARACIÓN DE CALIDAD

Enhorabuena por su nueva cama Invacare®Alegio™ de Invacare® France Operation SAS.

Su nueva cama cuenta con el marcado CE - de acuerdo a la directiva 93/42/EEC que hace referencia a los aparatos médicos.

Invacare®Alegio™ ha sido diseñada y fabricada teniendo en cuenta a los usuarios y las personas que los asisten.

La cama Invacare®Alegio™ ha sido diseñada de acuerdo a la normativa Europea EN 1970/A1.

El proceso de producción del modelo Invacare®Alegio™ ha sido supervisado en su totalidad y su calidad inspeccionada, por nuestro control de producto acabado.

La cama cuenta con una placa y marca CE que confirma que la misma ha sido inspeccionada por el control de producto acabado. Por favor lea el manual antes de utilizar la cama.

Invacare® France Operation SAS es una empresa certificada ISO 9001 e ISO 13485.

FR

DÉCLARATION DE QUALITÉ

Félicitations! Vous avez choisi votre nouveau lit Invacare®Alegio™.

Votre nouveau lit est marqué CE conformément à la directive 93/42/EEC concernant les dispositifs médicaux.

Invacare®Alegio™ a été développé et construit en considérant systématiquement les besoins de l'utilisateur et des tierces personnes lors de la manipulation du lit ou de son utilisation.

Invacare®Alegio™ a été développé conformément au Standard Européen NF EN 1970/A1.

Invacare®Alegio™ a été supervisé et contrôlé tout au long du processus de fabrication et le lit achevé a été inspecté par le contrôle des produits finis. La plaque d'identification sur le lit attestant que le contrôle des produits finis a approuvé le lit. Le lit répond aux exigences de l'analyse de risques de la norme NF EN 14971.

Nous vous remercions de lire le Manuel de l'Utilisateur dans son intégralité avant d'utiliser le lit.

Invacare® est certifiée ISO 9001 et ISO 13485.

PT

DECLARAÇÃO DE QUALIDADE

Parabéns pela sua nova cama Invacare®Alegio™ da Invacare® France Operation SAS.

A cama tem a marca CE - em conformidade com a directiva 93/42/EEC referente a aparelhos médicos.

A cama foi concebida e desenhada, tendo em consideração o seu utilizador e o seu(s) assistente(s), que o ajudarão a manipular a cama. A cama foi concebida em conformidade com o Standard Europeu EN 1970/A1.

Durante todo o processo de fabrico e produção, a cama Invacare®Alegio™ foi supervisionada e a sua qualidade controlada, sendo o produto final inspeccionado e testado pelo nosso controlo de qualidade. As etiquetas são colocadas na cama, após a aprovação final da cama, confirmando e garantindo a conformidade com o nosso controlo de qualidade. Por favor, leia atentamente este manual de utilizador antes de utilizar a sua cama.

A Invacare® France Operation SAS está certificada em conformidade com ISO 9001 e ISO 13485.

Invacare® Alegio™

USER MANUAL (EN)

EN

BRUKSANVISNING (SV)

SV

GEBRUIKSAANWIJZING (NL)

NL

MANUEL D'UTILISATION (FR)

FR

MANUAL DEL USUARIO (ES)

ES

MANUALE D'USO (IT)

IT

MANUAL DE UTILIZAÇÃO (PT)

PT

Table of contents

User's part	8
1. In general	8
2. Operating the Invacare®Alegio™	9
Technical part	11
3. Safety guidelines	11
4. Dismounting/mounting the Invacare®Alegio™	12
5. Removal of cables from the control unit	15
6. Mounting/dismounting the accessories	15
7. Emergency lowering of the backrest and/or thigh section .	16
8. Operating the accessories	17
9. Accessories	18
10. Cleaning	19
11. Maintenance and check-ups	19
12. Maintenance chart	20
13. Disposal	20
14. Trouble-shooting the electrical system	21
15. Technical specifications	22
16. Electrical data	23
17. Weights	23
18. Electro Magnetic Compliance (EMC)	24

User's part

Congratulations with your choice of the *Invacare*® **Alegio**™ nursing bed.

- The *Invacare*® **Alegio**™ is designed for patients over the age of 12 years especially for home care.
- The *Invacare*® **Alegio**™ combines excellent stability and ergonomic design with easy disassembly and operation.
- Max. weight: 170 kg.
Max. patient weight: 135 kg (provided that the weight of the mattress and the accessories do not exceed 35 kg).
Important! The max. load of the bed must not be exceeded.

In order to optimise the patient's comfort, *Invacare*® France Operation SAS recommends using a 16 cm mattress.

1. In general

Please read the entire User's part carefully before using the bed.

All indications of right and left are based on a patient lying on the back in the bed.



Check whether the bed shows any signs of damage. If the bed is damaged, see terms of delivery.



Please contact your *Invacare*® supplier if, contrary to our expectations, a problem should arise in connection with the delivered product.



Invacare® France Operation SAS accepts no liability for any use, change or assembly of the product other than as stated in this User's Manual.



The bed must not be used by patients under 12 years of age, or by patients with body size equivalent to an average 12-year-old or smaller.



Electromagnetic interference between the bed and other electrical products can occur.
To reduce or eliminate such electromagnetic interference, increase the distance between the bed and the products or switch them off.
This medical bed can be used together with medical electrical equipment connected to the heart (intracardially) or the blood vessels (intravascularly), provided that following points are respected:
- The medical bed should be equipped with means for potential equalization connection marked out by a symbol shown in the back of this manual.
- Medical electrical equipment should not be fixed on the bed's metallic accessories such as side rails, lifting pole, drip rod, bed ends ect.
In addition, the medical electrical equipment power supply cord should be kept clear of the accessories or any other moving part of the bed.



If the bed is used by restless (spasms) or confused persons, an ACP box can be used for blocking the hand control functions.



Ensure the bed is adjusted to it's lowest position before leaving the bed unattended - thereby the effect of fall-down / entrapment accidents is reduced.



Always apply the brakes when there is no need for moving the bed, to avoid accidents during entering or exiting the bed, and during handling of the patient.



Disconnect the plug from the mains before moving the bed. The cable must be kept clear off the floor and the castors during transport.
Make sure that the mains cord cannot be jammed or in any other way damaged, when the bed is used.



Always leave the bed in the lowest position. Otherwise there is a risk of entrapment due to accidental lowering of the mattress support. A person under the bed can be seriously injured during height adjustments.


- The *Invacare® Alegio™* is CE- marked in accordance with directive 93/42/EEC concerning Medical Apparatus.
- The motors and the control box of the *Invacare® Alegio™* are approved according to EN 60601:1996-03.
- The *Invacare® Alegio™* has undergone a risk analysis according to EN 14971.
-

The control unit, motors and hand control are protected according to IP 66. A lock cam **must** be used on the control unit, otherwise *Invacare®* France Operation SAS cannot guarantee the protection.

Max. load: 170 kg – max. patient weight: 135 kg. Max. load is stated on the identification label.

NOTE! The max. load of the bed must not be exceeded.

Operating the hand control



← **Sitting position (green button) - up/down**
(Only 4-sectioned beds)
Use the button with the symbol shown on the left.

← **Adjustment of backrest - up/down**
Use the button with the symbol shown on the right.

← **Adjustment of the thigh section - up/down**
(Only 4-sectioned beds)
Use the button with the symbol shown on the right.

← **Height adjustment of mattress section - up/down**
NB! Not applicable to beds with manual leg adjustment.


2. Operating the *Invacare® Alegio™*

Adjustment of leg section (only on 3-sectioned beds)




Operate the leg section by lifting the mattress handle.

Up: Lift the mattress handle on the leg section.

Down: Lift the mattress handle on the leg section all the way up and then lower the leg section.



Operating the braking castors

<p>Operating castors without central braking system</p> <p>When the bed is positioned correctly, at least one castor at the head end and one castor at the foot end must be locked.</p> <p>1) Braking: Step on the pedal.</p> <p>2) Releasing the brake: Step on the release pedal.</p>	
<p>Operating castors with central braking system</p> <p>When the bed is positioned correctly, the castors must be locked.</p> <p>1) Braking: Step on the bow.</p> <p>2) Releasing the brake: Lift the bow, until the brake is released.</p>	
<p>Operating the steerable castor</p> <p>The AlegioTM with central braking system may be equipped with a steerable castor.</p> <p>The steerable castor is operated by means of the central braking pedal.</p> <p>1) Activating the steering: If the brake is in neutral, step on the green pedal.</p> <p>2) Deactivating the steering: Step on the red pedal, until the brake is in neutral position.</p>	

Set-off from the castors might under special conditions appear at different types of absorbing floor covering - including untreated or badly treated floors. In matters of doubt, *Invacare*[®] recommends to place a suitable kind of protection between the castors and the floor.

3. Safety guidelines



Disconnect the plug from the mains before moving the bed. The cables must be kept clear off the floor and the castors during transport. Make sure that the mains cord cannot be jammed or in any other way damaged, when the bed is used.



The brakes of the bed must be activated during nursing of a patient in the bed and when the bed is adjusted.



The main plug must be accessible any time to disconnect the device from mains missing. We recommend to mount the mains cable on the hook for this purpose, see following picture.



Adjust all mattress support sections to a horizontal position before transporting an assembled bed. Hold the top of the bed end with both hands while the bed is pushed/pulled.



Caution! Please ensure that the pipe pins are fully inserted through both holes of the tubular section, and that the spring clip is fully engaged prior to operating the lifting mechanism on the bed. Furthermore, please ensure that all four plastic bushings are intact and located correctly between shear arm/base and shear arm/top frame.
NB! Do not pull the clamp part of the pipe pin during dismantling – the clamp may be deformed, thus being unable to lock properly.



Risk of entrapment of fingers during assembly and operation of the bed end and side rail.




A minimum of 2 persons are required to release a mattress support.



There is a risk of entrapment while operating the side rails.



Position the lifting pole in such a way that the handle extends inwards across the bed. If the handle has been turned away from the bed while the lifting pole is being used, the bed can tip.



We recommend to mount the mains cable on the hook for this purpose, see picture below.



Adjust all mattress support sections to a horizontal position before transporting an assembled bed. Hold the top of the bed end with both hands while the bed is pushed/pulled.

Only personnel who have been trained or instructed by *Invacare*® may perform the service work described in this manual in the chapters 12 and 13.


If the functions of the bed change, check the bed according to chapter 13.

The *Invacare*® **Alegio**™ must be stored in a room with a humidity of 10 - 80%, RH within a temperature of < 0 - 50°C.

4. Dismounting/mounting the *Invacare*® **Alegio**™

Dismounting the *Invacare*® **Alegio**™

1. Brake the bed and bring it to its lowest position and then about 5 cm up.
2. Remove the accessories: Siderails, bed ends and lifting pole.
3. Dismount the finger screws at both sides (see illustration A).
4. Dismount the cable from the backrest/leg section (see illustration B).
5. Dismount the top frame at the head end and place it in a vertical position (see illustration C and D).
6. Lift the leg section and dismount the pipe pins (see illustration E).
7. Dismount the top frame at the foot end and place it in a vertical position (see illustration F).
8. Dismount the shear arm (see illustration G, H, I, J and K).



A)



Finger screw


B)



C)



D)




E)




F)




G)



H)



I)**J)****K)****L)**

Mounting the Invacare®Alegio™

1. Mount the shear arm (see illustration K, J, I, H and G).
2. Mount the leg section and the pipe pins, ensuring that the spring clip is securely latched over the end of the pin (see illustration F and E).
3. Mount the head section (see illustration D and C).
4. Push the head section into the leg section (see illustration L).
5. Mount the cable for the backrest motor/leg section (see illustration B).
6. Mount the finger screws in both sides (see illustration A).
7. Mount the accessories: Siderails, bed ends and lifting pole.



Caution! Please ensure that the pipe pins are fully inserted through both holes of the tubular section, and that the spring clip is fully engaged prior to operating the lifting mechanism on the bed. Furthermore, please ensure that all four plastic bushings are intact and located correctly between shear arm/base and shear arm/top frame.

NB! Do not pull the clamp part of the pipe pin during dismantling – the clamp may be deformed, thus being unable to lock properly.

Mounting of the pipe pin between shear arm and base

The pipe pin must be mounted with the opening upwards and locked.




5. Removal of cables from the control unit


Raise the bed to highest position and remove the mains cord from the socket outlet.



Use a small tool (such as screwdriver) to release the pawl on the locking cam.
Hold the locking cam while releasing the remaining pawls.
With all pawls released, the locking cam can be removed from the control unit.



Motor connections are indicated on the label behind the cables.




6. Mounting/dismounting the accessories



Risk of entrapment of fingers during assembly and operation of the bed end and side rail.


Mounting the Vibeke bed end and Verso side rail

1. Lower the bed ends into the U-profiles.
2. Place the side rail with the release system in the leg end, and tighten both finger screws.




Mounting the Kirsty bed end and Nina side rail

1. The bed end is mounted as shown on illustration A - push the bed end all the way down, and secure it with the two finger screws (see illustration B).



2. Raise the mattress support to approximately 1/3 from lowest position.
3. Remove the side rail stopping bracket and finger screws (see illustration C).
4. Remove the protective tape from the locking dowels in each end of the side rail.
5. Install one end at a time. Lift the lower wooden rail with one hand while guiding the metal bow into the bed end bar (see illustration D).
6. Press the locking pin of the side rail with one finger.
7. Raise the side rail until the locking pin engages with an audible click, thus preventing the side rail from falling down.
8. Re-install the side rail stopping bracket and finger screws (see illustration C).



Mounting the lifting pole

Lower the lifting pole into the lifting pole tube. The lifting pole MUST be fastened with a finger screw.

7. Emergency lowering of the backrest and/or thigh section



A minimum of 2 persons are required to release a mattress support.

Remove the plug from the mains before emergency lowering of the mattress support.

In an emergency, the mattress sections are released by pulling out the cotter pin from the motors.

Two persons hold the mattress section in locked position.

One of them pulls out the cotter pin. Both slowly lower the mattress section until it is complete down.



8. Operating the accessories

To open (fold the siderail):

Pull the locking button and start the folding movement with one hand on the upper horizontal bar. Release the locking button during the folding movement and handling the siderail until its lower position (under the mattress).



Take care to don't squeeze body of patient during this step.
To avoid all finger pinch out during the manipulation, fingers must be positioned in the two lateral dug areas designed for that purpose, in locking button (grey zone, next page).

To close (unfold the siderail): Pull the siderail by the upper bar and rising until it engage in the locking system.



Take care to don't squeeze body of patient during this step.

Operating the Nina side rail

Up: Pull up the top wooden siderail bar, until the locking pins engages with an audible click.

Down: Lift the top wooden siderail bar while pressing the two locking pins together. Lower the siderail.

Operating Line and Britt side rail


The bed end may be equipped to lock the siderail at half height as well as in the top position. Normally the side rail can only be locked in top position.

Up: Pull up the top wooden side rail bar, until the locking pin locks with an audible click.

Down: Lift the top wooden side rail bar and press the two locking rings together. Lower the side rail.



Raise



Release



Lower



There is a risk of entrapment while operating the side rails.

Adjusting the height of the lifting pole handle

Loosen the cord as shown in ill. A. The lifting handle can now be adjusted to the desired height. Press the cord together as shown in ill. B and check that the cord is secured in the cord lock by pulling the handle.



Position the lifting pole in such a way that the handle extends inwards across the bed. If the handle has been turned away from the bed while the lifting pole is being used, the bed can tip.













A



B

Max. load of the lifting pole: 80 kg.

9. Accessories

Article	
Nina wooden side rail - 1 pair Britt wooden side rail - 1 pair Line wooden side rail - 1 pair	
Verso II steel side rail - 1 pair	
Kirsty wooden bed end with slide bars - 1 pair	
Kirsty wooden bed end - 1 pair	
Vibeke wooden bed end (high) - 1 pair	
Lifting pole	
Hand control support	
IV drip rod	
Support for IV drip rod	
IR handset	

Please only use original spare parts.

Spare parts lists and extra user manuals for the *Invacare® Alegio™* can be ordered from *Invacare®*.

10. Cleaning

The *Invacare® Alegio™* does not tolerate cleaning in automatic washing plants or using water jet based cleaning equipment.

The bed is washed down using a sponge, cloth or brush. Use ordinary disinfecting detergents.

Dry the bed after cleaning.

Never use acids, alkalines or solvents such as acetone or cellulose thinner.

The hand control, motors and the control box can be wiped with a moist cloth (the water temperature must be below 30° C).

When cleaning: Bring the backrest to its top position. Bring the bed to its top position. Disconnect the power cord by pulling the plug from the socket outlet before cleaning the bed. The risk of entrapment of fingers is minimized, as accidental elevation of parts of the bed is not possible.

11. Maintenance and check-ups

Only personnel who have received the necessary instruction or training may perform service and maintenance on the *Invacare® Alegio™*.

After 3 months of use the following must be checked:

- Tightening of the finger screws at the inserts in the middle of the bed.
- Fastening of the side rail - locking system and moving system.

With normal operation, the first service inspection is required after 2 years and thereafter every second year, according to the maintenance chart shown in chapter 12.

Please note:

The mattress support must be supported during service inspections to prevent accidental lowering.

Motors, control unit and hand control

These parts are serviced by exchanging the faulty part.

Motors, control unit and hand control must be regularly cleaned from dust and dirt and must be inspected for mechanical damage or breakage.

Inspect anchor points, cables, piston rod, casing and plugs and check the correct functioning of the motor.

12. Maintenance chart

Only personnel who have received the necessary instruction or training may perform service and maintenance on the <i>Invacare® Alegio™</i>				
S/N (located on mattress support): _____				
What to check for:	Date:			
Electrical safety check in accordance with the values in the standard EN 60601-1				
Check the side rails' mounting and locking/movement.				
Check mounting and braking of wheels.				
Check height adjustment motor - suspension, and performance.				
Check backrest and leg rest motor - suspension and performance.				
Check that cables and plugs are undamaged.				
Check rastofix fitting and its function.				
Check weldings.				
Have damaged coating repaired.				
Lubrication performed: 1. Points of rotation in mattress support and base frame, with oil. 2. All of the motors' tension rods, with oil. Lubricate with medically clean oil, e.g. KEN-WO 50, order no.: 813239. NB! The wooden side rails' gliding system must not be lubricated with oil - otherwise the wooden bars will move sluggishly.				

A service contract can be made in countries where *Invacare®* has its own sales company. Furthermore, *Invacare®* offers courses in service and maintenance of the *Invacare® Alegio™*.

13. Disposal

This product has been supplied from an environmentally aware manufacturer that complies with the Waste Electrical and Electronic Equipment (WEEE) Directive 2002/96/CE.


This product may contain substances that could be harmful to the environment if disposed of in places (landfills) that are not appropriate according to legislation.

The 'crossed out wheelee bin' symbol is placed on this product to encourage you to recycle wherever possible.

Please be environmentally responsible and recycle this product through your recycling facility at its end of life.

14. Trouble-shooting the electrical system


Only personnel who have received the necessary instruction or training may perform service and maintenance on the Invacare® **Alegio**™.




15. Technical specifications

All measurements are given in cm. All angles are stated in degrees. All measurements and angles are stated without tolerances.

Invacare® reserves the right to change the stated measurements and angles.



3-sectioned



4-sectioned

Transportation and storage conditions:


Temperature between -30° and 40° C

Humidity between 10 and 80 % RH

Pressure between 700 and 1.100 hPa

16. Electrical data

Voltage supply: 230 V.
 Max. current input: 1 A.
 Motor voltage: 24 V.
 Intermittent (periodic motor operation): 2/18 minutes.
 Protection class: IP 66.
 Insulation class: II, type B.

Alternating current: 

Direct current: 

Sound level < 45 dB.

The bed is not provided with a mains switch, so the mains plug is the only separation from the mains.


The patient is not separated from the ground and the chassis.



Double insulated.



Max. load (SWL) (Patient + mattress + side rail + lifting pole + other equipment).



The product should be reused where possible.



Equipotentiality localisation



17. Weights

Mattress support 3-section	31,5 kg
Mattress support 4-section	36,0 kg
Top frame – head section 3-4-section	19,5 kg
Top frame – leg section 3-4-section.....	12,0 kg
Top frame – leg section 4-section	16,5 kg
Base and shear arm.....	32,5 kg
Base.....	13,5 kg
Shear arm	15,0 kg
Bed ends Vibeke (per pcs.)	6,0 kg
Side rails Verso II (per pcs.).....	9,3 kg
Side rails Nina (per pcs.).....	8,0 kg
Lifting pole.....	7,5 kg


18. Electro Magnetic Compliance (EMC)

The Medical Bed is intended for use in the electromagnetic environment specified below. The user of the Medical Bed should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group I	The Medical Bed uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The Medical Bed is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

Guidance and manufacturer's declaration – electromagnetic immunity			
The Medical Bed is intended for use in the electromagnetic environment specified below. The user of the Medical Bed should assure that it is used in such an environment.			
IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV differential mode ± 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % UT (>95 % dip in UT) for 0,5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95 % dip in UT) for 5 s	<5 % UT (>95 % dip in UT) for 0,5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95 % dip in UT) for 5 s	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Medical bed Alegio requires continued operation during power mains interruptions, it is recommended that the Medical Bed Alegio be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	0,3 A/m	The power frequency magnetic field should be at a characteristic level of a typical commercial or hospital environment.
NOTE UT is the a.c. mains voltage prior to application of the test level.			

Guidance and manufacturer's declaration – electromagnetic immunity

The Medical bed is intended for use in the electromagnetic environment specified below. The user of the Medical Bed should assure that it is used in such an electromagnetic environment.

IMMUNITY test	IEC 60601 TEST LEVEL	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 Vrms	1 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the Medical Bed, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 3,5 \sqrt{P}$
	150 kHz to 80 MHz outside ISM bands ^a	1 Vrms	
Radiated RF IEC 61000-4-3	10 Vrms	1 Vrms	$d = 1,2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2,3 \sqrt{P}$ 800 MHz to 2,5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). ^b Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^c should be less than the compliance level in each frequency range. ^d Interference may occur in the vicinity of equipment marked with the following symbol: 
	150 kHz to 80 MHz in ISM bands ^a	10 V/m	
IEC 61000-4-3	80 MHz to 2,5 GHz	10 V/m	

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a The ISM (industrial, scientific and medical) bands between 150 kHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz.

^b The compliance levels in the ISM frequency bands between 150 kHz and 80 MHz and in the frequency range 80 MHz to 2,5 GHz are intended to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas. For this reason, an additional factor of 10/3 has been incorporated into the formulae used in calculating the recommended separation distance for transmitters in these frequency ranges.

^c Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Medical Bed is used exceeds the applicable RF compliance level above, the Medical Bed should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Medical Bed.

^d Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 1 V/m.

**Recommended separation distances between
portable and mobile RF communications equipment and the Model 006**

The Medical Bed is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The user of the Medical Bed can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Medical Bed as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter		
	150 kHz to 80 MHz $d = 1,2\sqrt{P}$	80 MHz to 800 MHz $d = 1,2\sqrt{P}$	800 MHz to 2,5 GHz $d = 2,3\sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.



Customer Sales and Service

Denmark

INVACARE A/S

Sdr. Ringvej 37
DK-2605 Brøndby
Phone: +45 36 90 00 00
Fax: +45 36 90 00 01
www.invacare.dk
denmark@invacare.com

Sweden&Finland

INVACARE AB

Fagerstagatan 9 / Box 66
S-163 91 Spånga
Phone: +46 8 761 70 90
Fax: +46 8 761 81 08
www.invacare.se
sweden@invacare.com

Norway&Iceland

INVACARE AS

Grensesvingen 9
Postbox 6230 Etterstad
N-0603 Oslo
Phone: +47 22 57 95 00
Fax: +47 22 57 95 01
www.invacare.no
norway@invacare.com

Spain

INVACARE S.A.

C/Areny S/N
Poligon Industrial de Celrà
E-17460 Celrà (Girona)
Phone: +34 972 49 32 00
Fax: +34 972 49 32 20
www.invacare.es
contactsp@invacare.com

Switzerland

MOBITEC REHAB AG

Benkenstrasse 260
CH-4108 Witterswil
Phone: +41 61 487 70 80
Fax: +41 61 487 70 81
switzerland@invacare.com

Invacare® France Operations SAS

Ident. no.: 1441391

Version: K

Date: 12.2011

Belgium&Luxemburg

INVACARE N.V.

Autobaan 22
B-8210 Loppem, Brügge
Phone: +32 50 83 10 10
Fax: +32 50 83 10 11
www.invacare.be
belgium@invacare.com

Netherlands

INVACARE B.V.

Celsiusstraat 46
NL-6716 BZ Ede
Phone: +31 318 695 757
Fax: +31 318 695 758
www.invacare.com
nederland@invacare.com
csede@invacare.com

Germany

INVACARE® GmbH

Alemannenstrasse 10
D-88316 Isny
Phone: +49 75 62 7 00 0
Fax: +49 75 62 7 00 66
www.invacare.de
kontakt@invacare.com

Portugal

INVACARE Lda

Rua Estrada Velha 949
P-4465-784 Leça do Balio
Phone: +351 225 1059 46/47
Fax: +351 225 1057 39
www.invacare.pt
portugal@invacare.com

Australia

INVACARE Australia Pty Ltd

1 Lenton Place, North Rocks
NSW 2151
Phone: +61 2 8839 5333
Fax: +61 2 8839 5353
www.invacare.com.au
sales@invacare.com.au

3rd party certified

according to

EN ISO 9001

EN ISO 13485

France

INVACARE Poirier S.A.S

Route de St. Roch
F-37230 Fondettes
Phone: +33 2 47 62 64 66
Fax: +33 2 47 42 12 24
www.invacare.fr
contactfr@invacare.com

Italy

INVACARE MECC SAN S.R.L.

Via dei Pini 62
I-36016 Thiene (VI)
Phone: +39 0445 38 00 59
Fax: +39 0445 38 00 34
www.invacare.it
italia@invacare.com

United Kingdom&Ireland

INVACARE LTD

Pencoed Technology Park,
Pencoed
UK-Bridgend CF35 5AQ
Phone: +44 1 656 776 200
Fax: +44 1 656 776 201
www.invacare.co.uk
uk@invacare.com
ireland@invacare.com

New Zealand

INVACARE NZ

4 Westfield Place,
Mt. Wellington
Auckland
Phone: +64 9 917 3939
Fax: +64 9 917 3957
www.invacare.co.nz
sales@invacare.co.nz

Austria

MOBITECMOBILITÄTSHILFENGmbH

Herzog Odilostrasse 101
A-5310 Mondsee
Phone: +43 6232 5535 0
Fax: +43 6232 5535 4
www.mobitec-austria.com
austria@invacare.com

Manufacturer:

INVACARE Lda

Rua Estrada Velha 949

P-4465-784 Leça do Balio

www.invacare.pt