Invacare® Medley Ergo™

Medley Ergo, Medley Ergo Low

en	Bed User Manual.	3
de	Bett Gebrauchsanweisung	21
es	Cama Manual del usuario	39
fi	Vuode Käyttöohje	57
fr	Lit Manuel d'utilisation.	75
it	Letto Manuale d'uso	93
nl	Bed Gebruiksaanwijzing.	111
pt	Cama Manual de utilização	131

This manual MUST be given to the user of the product.
BEFORE using this product, read this manual and save for future
reference.

Yes, you can.®

© 2016 Invacare® Corporation

All rights reserved. Reproduction, duplication or modification in whole or in part is prohibited without prior written permission from Invacare. Trademarks are identified by ™ and ®. All trademarks are owned by or licensed to Invacare Corporation or its subsidiaries unless otherwise noted.

This bed has been tested and approved according to IEC 60601-2-52.

This bed has undergone a risk analysis according to EN ISO 14971.

Invacare® is certified according to ISO 9001 and ISO 13485.

Invacare® reserves the right to alter product specifications without further notice.

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Contents

This manual MUST be given to the user of the product.
BEFORE using this product, read this manual and save for future reference.

I General	4
1.1 Introduction	4
1.1.1 Symbols in this manual	4
1.2 Intended Use	4
1.3 Service life.	4
1.4 Warranty	4
1.5 Limitation of liability.	4
2 Safety	5
2.1 General safety information	5
2.2 Mattresses.	5
2.3 Labels and symbols on the product.	5
2.3.1 Product label	5
2.3.2 Other labels and symbols	6
3 Setup	7
3.1 General Safety Information.	7
3.2 Main parts of the bed	7
3.3 Assembling the bed	7
3.4 Install Rastofix	7
3.5 Control Box	8
3.6 Wiring	8
3.6.1 Bed with 24 V external power supply.	8
3.7 Installing Aria and Bella side rail	8
3.8 Installing the mattress support extension.	9
3.9 Mattress retainers	9
4 Operating the bed.	10
4.1 General Safety Information.	10
4.2 Hand control.	10
4.2.1 Locking function	10
4.3 Castors and brakes	10
4.3.1 Castor brake	10
4.4 Lifting pole	11
4.4.1 Place the lifting pole	11
4.4.2 Adjusting the handle height	11
4.5 Emergency release of a mattress support section.	11
4.6 Transport brackets	11
5 Maintenance	13
5.1 General servicing information	13
5.1.1 Checklist – Servicing.	13
5.2 Cleaning and disinfection	13
5.2.1 Cleaning methods.	13
5.3 Lubrication	13
5.4 Troubleshooting	14
6 After Use	15
6.1 Waste disposal	15
6.2 Reconditioning.	15
7 Technical data	16
7.1 Dimensions bed.	16
7.2 Weights	16
7.3 Mattress dimensions	16
7.4 Electrical data	16
7.5 Environmental conditions.	17
7.6 Electromagnetic compliance (EMC).	17

I General

I.1 Introduction

This user manual contains important information about the handling of the product. In order to ensure safety when using the product, read the user manual carefully and follow the safety instructions.

To ensure correct use, the bed must be tested and adjusted by qualified personnel.

All references to left and right are based on a person lying on his back in the bed, with his head in the head end.

If a problem should arise in connection with the delivered product, please contact your *Invacare*® dealer. An address list is shown on the back side at this manual.

I.1.1 Symbols in this manual

Symbols and signal words are used in this manual and apply to hazards or unsafe practices which could result in personal injury or property damage. See the information below for definitions of the signal words.

WARNING

Indicates a hazardous situation that could result in serious injury or death if it is not avoided.

CAUTION

Indicates a hazardous situation that could result in minor or slight injury if it is not avoided.

IMPORTANT

Indicates a hazardous situation that could result in damage to property if it is not avoided.

Tips and Recommendations

Gives useful tips, recommendations and information for efficient, trouble-free use.

This product complies with Directive 93/42/EEC concerning medical devices. The launch date of this product is stated in the CE declaration of conformity.

Manufacturer of the product.

I.2 Intended Use

The bed has been developed for domestic care and long term care:

- *Application environment 3*; Long-term care in a medical area where medical supervision is required and monitoring is provided if necessary and medical electrical equipment used in medical procedures may be provided to help maintain or improve the condition of the patient.
- *Application environment 4*; Care provided in a domestic area where medical electrical equipment is used to alleviate or compensate for an injury, disability or disease.

- The bed is intended for indoor use only.
- The bed is intended for adult users, having a physical size equal to or more than 146 cm, a weight equal to or more than 40 kg and a body mass index (BMI) equal to or more than 17.
- The bed is not intended for transportation of users. It is mobile within a room with a user in it. Castors are lockable.
- The bed is not intended for psychiatric patients.
- Maximum user weight and safe working load are listed on the product label and in the technical data section in this manual.

WARNING!

Any other or incorrect use could lead to hazardous situations.

Invacare accepts no liability for any use, change or assembly of the product, other than stated in this user manual.

I.3 Service life

The expected service life of this product is five years when used daily and in accordance with the safety instructions, maintenance intervals and correct use, stated in this manual. The effective service life can vary according to frequency and intensity of use.

I.4 Warranty

The warranty covers all material and production defects for two years from the date of delivery, provided it can be demonstrated that such defects were present before delivery. All manufacturing faults or defects must be promptly reported.

Invacare® may repair the fault or replace the component. The warranty provided by *Invacare*® does not cover additional costs (transport, packaging, labour, sundry expenses, etc. are for the customer's account).

The warranty does not cover:

- Damage caused during transport that is not directly reported to the forwarder at the moment of delivery.
- Repairs performed by unauthorized centers and personnel.
- Parts subject to normal wear.
- Malicious damages or damaged caused by improper use of the bed.

I.5 Limitation of liability

Invacare accepts no liability for damage arising from:

- Non-compliance with the user manual
- Incorrect use
- Natural wear and tear
- Incorrect assembly or set-up by the purchaser or a third party
- Technical modifications
- Unauthorized modifications and/or use of unsuitable spare parts

2 Safety

2.1 General safety information

WARNING!

Risk of entrapment / suffocation

- There's a risk of entrapment / suffocation between mattress support, side rail and bed end or between moving parts and objects placed nearby the bed.
- The bed must not be used by persons under 12 years of age, or by persons with a body size equivalent to an average 12 years old or smaller.
 - The bed, in combination with side rails must not be used by persons having a physical size less than 146 cm, a weight less than 40 kg or a body mass index (BMI) less than 17.
 - Due to mattress compression, an increased risk may occur over time. Periodically monitor gaps between the bed, mattress and/or side rail. Replace mattress if the gaps may lead to entrapment.

WARNING!

Risk of slipping through the openings

- The bed fulfils all requirements regarding maximum distances. However, it is possible that persons with small body dimensions slip through the openings between the side rails or through the opening between the side rail and the mattress support.
- Pay special attention, if the bed is used for the care of persons with small body dimensions.

WARNING!

Risk due to electromagnetic interference

- Electromagnetic interference between the bed and other electrical products can occur.
- To reduce or eliminate such electromagnetic interference, increase the distance between the bed and the products or switch them off.

This medical bed can be used together with medical electrical equipment connected to the heart (intracardially) or blood vessels (intravascularly) provided that following points are respected:

- The bed should be equipped with means for potential equalization connection marked out by a symbol shown in the back of this manual.
- Medical electrical equipment should not be fixed on the bed's metallic accessories such as side rails, lifting pole, drip rod, bed ends, etc.
- The medical electrical equipment power supply cord should be kept clear of the accessories or any other moving parts of the bed.

WARNING!

Risk of injury or damage to property

- Do not roll the castors over the main power cord.
- Do not bring main power cord into moving parts.
- Disconnect the plug from the mains before moving the bed.
- Make sure that no wires (mains or from other equipment) are jammed or damaged, when the bed is used.
- Keep bed components and accessories at least 30 cm away from a heated surface and not in direct sunlight.

CAUTION!

- There is a risk of entrapment of fingers in the bed moving parts.
- Pay attention to your fingers.

CAUTION!

- For a user entering or exiting the bed, always lower the bed to an appropriate height. The backrest can be used as a support. Make sure, the thigh- and leg section is horizontal to avoid overloading the mattress support.
- Always lower the bed to the lowest position before leaving the patient in the bed unattended.
- Make sure that there is nothing under, over, or near the bed that can obstruct the height adjustment, like for example furniture, lifts or window frames.

2.2 Mattresses

WARNING!

Safety aspects regarding combination of side rails and mattresses:


- To get the highest possible safety level, when using side rails on the bed, the minimum and maximum measures for mattresses, must be respected.
- For correct mattress measures see mattress table in chapter 7 *Technical data*, page 16.

WARNING!

Risk of entrapment and/or suffocation

- The user could get trapped and/or suffocate, if the horizontal space, between the mattress side and the inside of the side rail, is too big. Follow the minimum width (and length) of mattresses in combination with a side rail, as stated in the mattress table in chapter 7 *Technical data*, page 16.
- Be aware that using very thick or soft mattresses (low density), or a combination of these, increases the risk.


WARNING!

Risk of falling

- The user can fall over the edge and get seriously injured, if the vertical distance \textcircled{A} between the top of the mattress and the edge of the side rail/bed end, is too short. See image above.
- Always keep a minimum distance \textcircled{A} of 22 cm.
 - Follow the maximum mattress height in combination with the side rail as stated in the mattress table in chapter 7 *Technical data*, page 16.




2.3 Labels and symbols on the product

2.3.1 Product label

The product label is placed on the frame of the bed and contains the main product information, including technical data.

Symbols on the product label

	Serial Number
	Reference Number
	Manufacturer Address
	Manufacturing Date
	Max. User Weight
	Max. Safe Working Load
	CLASS II equipment
	Type B Applied Part
	WEEE conform
	This product complies with Directive 93/42/EEC concerning medical devices.

Abbreviations for technical data:

- Iin = Incoming Current
- Uin = Incoming Voltage
- Int. = Intermittence
- AC = Alternating Current
- Max = maximum
- min = minutes

For more information about technical data, refer to 7 Technical data, page 16.

2.3.2 Other labels and symbols

	Refer to User Manual
---	----------------------

Label - user and mattress sizes

	<p>Definition of min. weight, min. height and min. body mass index of an adult user</p> <p>Refer to user documentation for the correct mattress measures.</p>
--	---

Label on side rails

	<p>Indicates the minimum space between side rail and the head end of the bed. (removable side rails only).</p>
--	--

3 Setup

3.1 General Safety Information

When you receive the bed, check the packaging. If the packaging shows any signs of damage upon delivery, contact the shipping company.

WARNING!

Risk of injury or damage to property

The assembly of the bed and installation of accessories must be done by authorized or trained personnel.


- Follow instructions carefully. If you have any question regarding the assembly, contact your local dealer or Invacare representative.
- If the bed shows any sign of damages, do not use the bed. Contact your local dealer.
- The electrical equipment of the bed must not be dismantled or combined with other electrical equipment.
- After each assembly, check that all fittings are properly tightened and that all parts have the correct function.
-

IMPORTANT!

- To prevent condensation, the bed should not be used until it has reached a temperature of 10-40° C.

3.2 Main parts of the bed


- (A) Bed ends, 2 pcs
 - (B) Mattress support, backrest section
 - (C) Mattress support, thigh/leg section
 - (D) Castor, 4 pcs
 - (E) Motor, mattress support backrest section
 - (F) Motor, mattress support thigh/leg section
 - (G) Mattress retainers, 4 pcs
 - (H) Motor, height adjustment
- Hand control (not shown in image)


Side rails are not delivered in the main configuration. For more information about side rails and other accessories, contact your Invacare representative.

3.3 Assembling the bed

1. Push the inserts on the upper half of the mattress support into the side tubes on the lower half of the mattress support.

2. Tighten with two thumb screws.
3. Make sure that the four mattress retainers are positioned pointing upwards.
4. Turn the locking rings on the bed end to position OPEN.

5. Latch the brackets of the mattress support either in high position (A) or low position (B) to the bed end and press them firmly into position. (It's normal that there's a small space between the upper slot and the bracket)
6. Turn the locking rings into Position LOCK.

CAUTION!


The locking rings ensure that the mattress support can't be accidentally lifted off the bed ends.

- Make sure that all locking rings are locked and in the upper slot as shown in the picture.
- Make sure that all brackets are attached in the same position (high or low).


3.4 Install Rastofix

(optional)


1. Press bottom of Rastofix at the foot end of the bed.

2. Squeeze upper end of Rastofix together.

3. Rotate Rastofix between tabs and release.

3.5 Control Box

The control box is latched to the backrest motor.

Connect both plugs from the hi/lo motors (head end and foot end) to the control box placed underneath the backrest.

2.

Guide the hi/lo motor wire from the bed foot end AND the wire from the thigh/leg section motor through the thigh/leg section motor pipe pin. Make sure the pipe pin is properly locked.

3. Connect the control box to the mains.
4. Run the thigh/leg section to its highest positions .
- 5.

Place the wire hi/lo motor in the foot end on the hook at the bed end.

6. Run the backrest to its highest position
- 7.

Place the wire from the hi/lo motor in the head end, on the two hooks at the bed end and attach the wire in the backrest motor pipe pin.

8.

Install the safety clip over the plugs in the control box.


9.

Attach the mains wire in the strain relief that is attached to the bed frame.

10. Before taking the bed into use, check that all parts run smoothly and that no motor cables are hanging on the floor or being squeezed while moving bed sections.

3.6.1 Bed with 24 V external power supply

1. Fasten the short cable (A) to the control box (a “click” must be heard)
2. Attach the cable (as described in step 7 in the Wiring section).
3. Connect the short cable (B) to the transformer cable (C).
4. Connect the transformer 24 V (D) to the 230 V main socket.

3.7 Installing Aria and Bella side rail

WARNING!

Risk of personal damage


There is a risk of entrapment or squeezing, while assembling or disassembling the side rail.

- Follow instructions carefully.
- After assembly, check that the side rail has the correct function.

The installation of Bella wooden side rail and Aria steel side rail is identical.

Assembling the gliding system

1.


Push the pawl (A) from the backside, through the glider (B) and make sure it locks with an audible click.

2. Do the same for all four pawls in each glider.


Attaching the gliding system and installing the side rail

1. Raise the bed to 1/3 of full height.
- 2.


Loosen the screw at the bottom of the guides on one of the bed ends and push the glider halfway up into the guide until it attaches with an audible click.

3.


Push the upper side rail onto the two upper pawls.

4.

In the other end of the side rail bar; push the two upper pawls into the end of the bar.


5.

Push the attachment into the guide until it attaches with an audible click (=properly engaged in the lower position).

6. Mount the lower side rail bar in both ends as in step 4.
7. Push the side rail in both ends, all the way up until the upper bar locks in the upper position.

8.


Retighten the screw under the lower side rail bar, in each corner of the bed ends.


Before taking the bed into use; Make sure that the side rail is properly locked with the screw and that the bars run smoothly.

3.8 Installing the mattress support extension

(optional)

1. Remove the foot bed end.
2. Insert both sides of the mattress support extension into the frame and secure with nuts and bolts.
3. Reinstall the foot bed end.
4. Attach mattress deck extender to the mattress deck.

3.9 Mattress retainers

WARNING! **Risk of injury**

If the mattress retainers are removed the mattress will not follow the movement, when profiling the sections. The mattress can slip sideways and cause the user to fall out of the bed or get trapped in the bed.

- Always use mattress retainers and be careful to put them back after any adjustments.
- Always use the original type of retainers.
- After assembly of mattress support all four mattress retainers must point upwards.

4 Operating the bed

4.1 General Safety Information

WARNING!

Risk of personal injury and damage to property.

- The bed must be placed so that the height adjustment is not obstructed by, for example, lifts or furniture.
- Take care that no body parts are being squeezed between fixed parts (such as side rails, bed ends etc) and moving parts.
- The hand control must not be used by children.

WARNING!

Risk of injury

If the bed is used by confused or restless users or users with spasms:

- either, when the bed is equipped with a lockable hand control, lock the hand control functions,
- or make sure that the hand control is out of reach for the user.

4.2 Hand control


The hand control can be equipped with either two, three or four buttons to operate the electrical functions of the bed. The hand control can also be equipped with a locking function.

Backrest section

1. Up: press left side of the button (▲).
2. Down: press right side of the button (▼).

Height adjustment

1. Up: press left side of the button (▲).
2. Down: press right side of the button (▼).

Thigh section

1. Up: press left side of the button (▲).
2. Down: press right side of the button (▼).

Tilt function


CAUTION!

Risk of fatal injury


Tilting with head end down can have fatal effect on users who are sensitive for increased blood pressure in the upper part of the body.

- Only medically trained persons must operate the tilt function.
- Always perform a medical evaluation, before tilting the bed with the user's lower extremities positioned higher than the heart.
- The tilt function is NOT a Trendelenburger function and must not be used for medical treatment.

1. Head end up: press left side of the button (▲).
2. Foot end up: press right side of the button (▼).


If the bed is equipped with anti-tilt function only, tilting with the foot end up is not possible:

1. Head end up: press left side of the button (▲).
2. Horizontal position: bring the bed either all the way up or all the way down by using the height adjustment button.

4.2.1 Locking function

The locking function prevents usage of certain function buttons. Regardless of how many buttons your hand control has, they can all be locked individually.

1. Insert key Ⓐ in the key hole, below the wanted function.
2. To lock, turn the key clockwise.
3. To unlock, turn the key anti clockwise.

4.3 Castors and brakes

CAUTION!

Risk of trapping/squeezing

- All brakes are foot-operated.
- Do not release the brake with the fingers.

CAUTION!

Risk of injury


User can fall, while entering or exiting the bed, if brakes are not locked.

- Always lock the brakes before the user is moving in or out of bed or when attending to the user.
- At least one castor at the head end and one castor in the foot end must be locked.

Castors might make marks on different types of absorbing floor surfaces - such as untreated or badly treated floors. To prevent marks, *Invacare®* recommends to place a suitable kind of protection between the castors and the floor.

4.3.1 Castor brake

Locking the brake - Step on the outer part Ⓑ of the brake pedal

Unlocking the brake - Step on the inner part Ⓐ of the brake pedal

4.4 Lifting pole

WARNING! **Risk of injury**

The bed can tip if the handle is used, while the lifting pole is turned away from the bed.

- The lifting pole always has to be positioned with the handle hanging over the bed area.
- Do not exceed the maximum load of the lifting pole; 80 kg.

4.4.1 Place the lifting pole

The lifting pole can be placed either on the left or the right side of the head end of the bed.

1. Insert the lifting pole into the lifting pole tube and fix it with the finger screw.


It's not required to tighten the finger screw, in case you want the lifting pole to swing away to the side of the bed.

4.4.2 Adjusting the handle height

The handle height should always be adjusted to the user's need.

- 1.
- 2.

1. Hold the cord out while sliding the plastic cord lock (A) upwards or downwards until the handle reaches the preferred height.
2. Lock the handle by pressing back the cord in the lock (A) and pull the handle downwards.

IMPORTANT!

After adjusting the handle height:

- Check that the two cords above the cord lock are parallel and inside the cord lock.
- Make sure the cord is properly locked by pulling the handle hard.

4.5 Emergency release of a mattress support section

In case of a power- or motor failure, an emergency release of the back, thigh or leg section could be necessary. An emergency release of the height adjustment is NOT possible.

CAUTION!

Risk of injury


- A minimum of two persons is required for an emergency release of a mattress support section.
- When releasing a mattress support section, it might lower fast. Do not reach under the mattress support while lowering it.

IMPORTANT!

- Before an emergency release of the mattress support, remove the plug from the mains socket.

1. Both persons hold the mattress section.
2. One of them locates the motor in question and pulls out the safety pin.


3. Both persons slowly lower the mattress section until it is completely down.


4.6 Transport brackets

Use the transport brackets to store and transport the bed:


- 1.


- 2.


3.


4.

5.

6.

5 Maintenance

5.1 General servicing information

- !** **IMPORTANT!**
- The mattress support must be supported during service inspections, to prevent accidental lowering.
 - Only personnel who have received the necessary instructions or training may perform service and maintenance of the bed.
 - After reconditioning the bed, or if bed functions change, service must be carried out according to the check list.

A service contract can be made in the countries, where Invacare® has its own sales company. In certain countries Invacare® offers courses in service and maintenance of the bed. Spare parts lists and additional user manuals are available from Invacare®.

Before Use

- Ensure that all manual and electrical parts functions correctly and are in a secure state.
- Check, by raising and lowering the bed, that the shear arms run smoothly in the guideways.

After three months

- Ensure that all manual and electrical parts are functioning, and tighten bolts, screws, nuts, etc.

Every year

- We recommend a safety test comprising the motors' performance and mechanical state.

Every second year

- We strongly recommend a service according to the below checklist after two years of normal use and then every second year.

 Motors, hand control and control units are serviced by exchanging the faulty part.

5.1.1 Checklist – Servicing

Check the following checkpoints:

- All parts of the bed are visually intact (no plastic deformation and/or wear and tear of welded joints).
- All screws are tightened.
- Circlips, cotter pins and plastic fixing ring are properly locked and intact.
- All motors are running without failures (regular speed; low noise).
- Mains cable and plug are intact and not squeezed.
- All wires are correctly wired and undamaged.
- All plugs are undamaged.
- Cable insulation and actuator housing have no cracks or damages.
- Side rails are fixed properly and the locking/opening is functioning.
- Caster fittings are tightened.
- Caster brakes are locking properly.

5.2 Cleaning and disinfection

- !** **IMPORTANT!**
- Wrong fluids or methods can harm or damage your product.
- Follow instructions carefully for either Non-washable or Washable components.
 - Never use corrosive fluids (alkalines, acid, cellulose thinner, acetone etc)
 - Never use a solvent that changes the structure of the plastic or dissolves the attached labels.
 - Always make sure that the bed is carefully dried before taking it into use again.

5.2.1 Cleaning methods

Electrical components

- !** **IMPORTANT!**
- Non-washable electronics can not withstand high temperatures.
- Do not wash or dry in higher temperatures than 40 °C.

Method: Wipe off with a wet cloth or soft brush.

Max. temperature: 40 °C

Solvent/chemicals: Water

Metal components

Method: Wipe off with a wet cloth or soft brush. Water may be pressurized, but not high pressure or steam.

Max. temperature: 40 °C

Solvent/chemicals: Household detergent or soap and water, 6-8 pH

Wood (including textile straps on side rails, if existing)

Method: Wipe off with a wet cloth or soft brush.

Max. temperature: 40 °C

Solvent/chemicals: Household detergent or soap and water, 6-8 pH


Textiles (including upholstery and mattresses)

See attached label on each product.

5.3 Lubrication

We recommend lubricating the bed according to the following table:

Part of bed	Lubrication method
Points of rotation in mattress support and base frame	Oil (medically clean)
Motor attachment points on mattress support	Oil (medically clean)

 For ordering of correct oil, contact your Invacare dealer.

5.4 Troubleshooting

Symptom	Possible cause	Remedy
Bed section does not move	End of stroke reached	Operate opposite button
Full bed: No electric function does work	Bed not plugged in	Plug in mains cable
	Hand control not connected	Ensure hand control correctly connected to control unit.
	Hand control defective	Call dealer / technician to replace hand control
	Control unit defective	Call dealer / technician to replace control unit
	Cables entangled or pinched	Call dealer / technician to replace cables
Side rail does not open or close	Locking mechanism defective	Call dealer / technician to repair side rail

6 After Use

6.1 Waste disposal

 Waste disposal/recycle must comply with the laws and regulations for waste handling in each country.

Invacare® is continuously working towards ensuring that the company's impact on environment, locally and globally, is reduced to a minimum.

We comply with the current environment legislation (e.g. WEEE and RoHS directive).

We only use REACH compliant materials and components.

- All wooden parts must be dismantled and sent for incineration.

- All electric parts must be dismantled and be disposed of as electric components.
- Plastic parts must be sent for incineration or recycling.
- Steel parts and castors must be disposed of as waste metals.

! **IMPORTANT!**
 • Accumulator back-up
 – Old accumulators must be returned to *Invacare®* or recycled as car batteries

6.2 Reconditioning

This product is suitable for reuse. To recondition the product for a new user, carry out the following actions:


- Inspection according to service plan
- Cleaning and disinfection

For detailed information about inspection, cleaning and disinfection, refer to 5 Maintenance, page 13.

7 Technical data

7.1 Dimensions bed

All measurements are stated in cm. All angles are stated in degrees.
All measurements and angles are stated without tolerances.

Height	Medley Ergo	
	Low	High
A	33 – 73 cm	40 – 80 cm
B	160 – 201 cm	166 – 207 cm
C	15 – 55 cm	22 – 62 cm

Height	Medley Ergo Low	
	Low	High
A	21 – 61 cm	28 – 68 cm
B	146 – 187 cm	153 – 194 cm
C	2 – 42 cm	9 – 49 cm

7.2 Weights

	Medley Ergo	Medley Ergo Low
Max. patient weight (provided that the weight of the mattress and the accessories do not exceed 35 kg)	145 kg	
Max. safe working load (patient + accessories)	180 kg	
Total weight*	64,0 kg	
Upper part of mattress support with wood slats (= heaviest part)	17,0 kg	
Upper part of mattress support, with steel slats	16,0 kg	
Lower part of mattress support with wood slats	16,0 kg	
Lower part of mattress support, with steel slats	15,0 kg	
Bed end (1 piece)	14,7 kg	13,7 kg
Lifting pole	4,2 kg	
Shipping brackets	2,0 kg	

*Complete bed, incl. mattress support (upper and lower), steel slats, main cord, mattress retainers and hand control

7.3 Mattress dimensions

Allowed mattress sizes, depending on the side rail

Side rail	Mattress size		
	Height [cm]	Width [cm]	Length [cm]
Bella /Aria high ¹⁾	6 – 12	85,5 – 90	195 – 208
Bella /Aria low ²⁾	12 – 18	85,5 – 90	195 – 208
Verso II	8 – 16	85,5 – 95	195 – 208
Scala Basic 2	10 – 16	85,5 – 90	195 – 208
Scala Basic Plus	10 – 16	85,5 – 90	195 – 208
Scala Medium	10 – 20	85,5 – 90	195 – 208
Scala Decubi 2	23 – 29	85,5 – 90	195 – 208

¹⁾ Mattress support mounted in high position.

²⁾ Mattress support mounted in low position.

7.4 Electrical data

Voltage supply	Uin 230 Voltage, AC, 50–60 Hz (AC = Alternating current)
Maximum current input	Iin max. 2 Ampere
Intermittent (periodic motor operation)	10 % 2 min / 18 min
Insulation class	CLASS II equipment

	Type B Applied Part Applied Part complying with the specified requirements for protection against electrical shock according to IEC 60601-1.
Sound level	45–50 dB
Degree of protection	The control unit, external power supply, motors and hand control are protected according to IPx4. (IPx4: The system is protected against water splashed from any direction.)

 The bed has no isolator (main switch). Disconnect the bed by unplugging it from the mains socket.

7.5 Environmental conditions

	Storage and transportation	Operation
Temperature	-10°C to +50°C	+5°C to +40°C
Relative humidity	20% to 75%	
Atmospheric pressure	800 hPa to 1060 hPa	

 Be aware that when a bed has been stored under low temperatures, it must be adjusted to operating conditions before use.

7.6 Electromagnetic compliance (EMC)

Guidance and manufacturer’s declaration – electromagnetic emission

The medical bed is intended for use in the electromagnetic environment specified below. The customer or the user of the bed should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11 (partly)	Group 1	The medical bed uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11 (partly)	Class A	The medical bed is suitable for use in all establishments including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class B	
Voltage fluctuations /flicker emissions IEC 61000-3-3	Complies	

Guidance and manufacturer’s declaration – electromagnetic immunity

The medical bed is intended for use in the electromagnetic environment specified below. The customer or the user of the bed should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment — guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrostatic transient / burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	± 1 kV differential mode ± 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	< 5% U_T (>95% dip in U_T) for 0,5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles < 5% U_T (>95% dip in U_T) for 5 sec	< 5% U_T (>95% dip in U_T) for 0,5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles < 5% U_T (>95% dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the medical bed requires continued operation during power mains interruptions, it is recommended that the medical bed be powered from an un-interruptible power supply or a battery. U_T is the a. c. mains voltage prior to application of the test level.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment — guidance
Conducted RF IEC 61000-4-6	3 V	3 V	<p>Portable and mobile RF communications equipment should be used no closer to any part of the medical bed including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance:</p> $d = \left[\frac{3,5}{V_1} \right] \sqrt{P}$ $d = \left[\frac{3,5}{E_1} \right] \sqrt{P}$ $d = \left[\frac{7}{E_1} \right] \sqrt{P}$ <p>80 MHz to 800 MHz</p> <p>800 MHz to 2,5 GHz</p> <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).^b</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
Radiated RF IEC 61000-4-3	3 V/m	3 V/m	

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the medical bed is used exceeds the applicable RF compliance level above, the medical bed should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the medical bed.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than [V1] V/m.

At 80 MHz and 800 MHz, the higher frequency range applies.

Recommended separation distances between portable and mobile RF communications equipment and the medical bed

The medical bed is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the medical bed can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the medical bed as recommended below, according to the maximum output power of the communications equipment

Rated maximum output of transmitter [W]	Separation distance according to frequency of transmitter [m]		
	150 kHz to 80 MHz $d = \left[\frac{3,5}{V_1} \right] \sqrt{P}$	80 MHz to 800 MHz $d = \left[\frac{3,5}{E_1} \right] \sqrt{P}$	800 MHz to 2,5 GHz $d = \left[\frac{7}{E_1} \right] \sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,37	0,37	0,74
1	1,17	1,17	2,33
10	3,69	3,69	7,38
100	11,67	11,67	23,33

For transmitters rated at a maximum output power not listed above the recommended separation, distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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Invacare Poirier SAS