Invacare® Matrix® Libra®

| | | |
|----|---|----|
| en | Cushion User Manual | 1 |
| de | Kissen Gebrauchsanweisung | 4 |
| fr | Coussin Manuel d'utilisation | 8 |
| it | Cuscino Manuale d'uso | 11 |
| es | Cojín Manual del usuario | 14 |
| pt | Almofada Manual de utilização | 17 |
| nl | Kussen Gebruiksaanwijzing | 20 |
| sv | Dyna Bruksanvisning | 23 |
| no | Pute Bruksanvisning | 26 |
| da | Pude Brugsanvisning | 29 |
| fi | Pehmuste Käyttöohje | 32 |

I General

I.1 Images


Fig. 1


Fig. 2


Fig. 3


Fig. 4


Fig. 5

This manual MUST be given to the user of the product.
BEFORE using this product, read this manual and save for future reference.

Yes, you can.®

1.2 General information

This user manual contains important information about the handling of the product. In order to ensure safety when using the product, read the user manual carefully and follow the safety instructions. For further information contact Invacare® in your country (addresses see back page of this manual).

1.3 Symbols in this user manual

In this User Manual warnings are indicated by symbols. The warning symbols are accompanied by a heading that indicates the severity of the danger.

WARNING

Indicates a potentially hazardous situation which if not avoided could result in death or serious injury.

CAUTION

Indicates a potentially hazardous situation which if not avoided could result in product damage, minor injury or both.

IMPORTANT

Indicates a hazardous situation which if not avoided could result in damage to the product.

Gives useful tips, recommendations and information for efficient, trouble-free use.

This product complies with the directive 93/42/EEC for medical products. The launch date for this product is specified in the CE declaration of conformity.

Manufacturer

1.4 Warranty

We provide a manufacturer's warranty in accordance with our General Terms and Conditions of Business in the respective countries. Guarantee claims can only be made through the provider from whom the appliance was obtained.

1.5 Standards and regulations

Quality is fundamental to the company's operation, working within the disciplines of ISO 13485 and ISO 9001.

The Invacare® Matrx® Libra® cushion features the CE mark, in compliance with the Medical Device Directive 93/42/EEC Class I. Invacare® is continuously working towards ensuring that the company's impact on the environment, locally and globally, is reduced to a minimum. We use only REACH compliant materials and components.

The foam and cover used to manufacture the Invacare® Matrx® Libra® cushion are fire safety tested and certified in accordance with EN 1021-1 & -2.

For further information please contact Invacare® in your country (addresses see back page of this manual).

1.6 Product components (Fig. 1)

The cushion **A** consists of following components:

| | |
|----------|--|
| B | HR foam base |
| C | Dual-layer clinical fluid sac (silicone) |
| D | Inner Platilon liner |
| E | Outer cover |
| | Velcro® hook adapter (not shown) |

1.7 Intended use

The Matrx® Libra cushion provides highest level skin protection for users with very high risk of skin breakdown. This product is suitable for users with decubitus ulcers provided it is used in alignment with the prescribing therapist's/physician's clinical judgment and they are following an appropriate decubitus ulcer therapy program. The product is also suitable for long-term wheelchair users (for whom the

risk of decubitus development is very high). There are no known clinical contraindications for use of this product. Normal use is defined by Invacare Ltd as when the product is placed without additional covers or padding between the user and the support surface. Replace the cover in case of damage resulting in reduced tightness or elasticity of the protective barrier.

1.8 Product description

The best way to avoid problems related to pressure sores is to understand their causes and your role in a skin management program. Your therapist and physician should be consulted if you have questions regarding individual limitations and needs. All cushions should be selected carefully. Working with your therapist and physician is the best way to ensure that a cushion choice matches your individual needs. As the needs of the individual become more complex, cushion evaluation becomes more important.

The Invacare® Matrx® Libra® **A** uses a combination of a moulded foam base **B** and dual layer clinical fluid sac **C** to provide the highest level of skin protection for users with very high risk of skin breakdown.

Clinical Fluid Sac

The clinical fluid sac utilizes silicone in a unique dual-layer design placed in the ischial/sacral areas to help prevent bottoming out.

Inner Liner and Outer Cover

The inner liner **D** is moisture resistant to protect the foam. The standard outer cover **E** is water resistant and breathable. Regular cleaning and inspection of the covers is recommended, see section 5 Maintenance, page 3.

2 Safety

2.1 Safety information

WARNING!

– Do not use this product or any available optional equipment without first completely reading and understanding the user manual supplied. Invacare® product manuals are available on your local Invacare web site or at your local dealer. If you are unable to understand the warnings, cautions or instructions please contact a health care professional, dealer or technical personnel before attempting to use this equipment – otherwise, injury or damage may occur.

WARNING!

Invacare® products are specifically designed and manufactured for use in conjunction with Invacare® accessories. Accessories designed by other manufacturers have not been tested by Invacare® and are not recommended for use with Invacare® products.


IMPORTANT!

The information contained in this document is subject to change without notice.

– Check all parts for shipping damage and test before using.
– In case of damage, do not use.
– Contact Invacare® / Carrier for further instructions.

2.2 Symbols on the product

| | | | |
|--|----------------------|--|--|
| | Refer to user manual | | Do not bleach |
| | CE conform | | Do not iron |
| | Tumble dry low heat | | Machine wash temperature max. 80 °C (outer cover only) |

| | | | |
|--|----------------------|---|-----------------------|
|  | Do not dry clean |  | Orientation label |
|  | Do not pierce or cut |  | Do not put near flame |
|  | Manufacturing date |  | Maximum user weight* |

* Maximum user weight as per section 7.1 General data, page 4.

3 Setup

3.1 Installing the cushion

1. Position the cushion so the material identification tags (not shown) face the rear.
2. Align the hook and loop strips (not shown) on the bottom of the cushion with those on the wheelchair seat.

 If your wheelchair seat has loop strips, use the enclosed hook adapters to secure the cushion to the wheelchair seat.

3. Position the cushion onto the wheelchair as far back as possible.
4. Press down to ensure the hook and loop strips are firmly engaged.

4 Usage

4.1 Safety information

WARNING!

Risk of injury

- Installing a cushion on a wheelchair may affect the center of gravity of the wheelchair and may cause the wheelchair to become unstable, potentially resulting in injury.
- After any adjustments, repair or service and before use, make sure that all attaching component parts are secure.
- Ensure the seat cushion is properly secured to the wheelchair before using. Otherwise injury or damage may occur.

WARNING!

Risk of injury

- Skin condition should be checked very frequently after the installation of any new seating system.
- Your therapist and physician should be consulted if you have any questions regarding individual limitations and needs.
- Working with your therapist, physician, and equipment supplier is the best way to ensure that a seating choice matches your individual needs.
- As the needs of the individual become more complex, the seating evaluation becomes more important.
- Sitting on extremely cold or hot surfaces can cause skin damage.
- Avoid sharp objects or exposure to excessive heat or open flame.

WARNING!

- Avoid the use of additional covers or padding between the user and the cushion, as this could reduce its effectiveness.

4.2 Using the cushion

1. Remove all packaging before use.
2. Place the cushion on wheelchair or chair with grey non-slip base facing down.

4.3 Accessories (Fig. 2)

Optional Matrix® Libra® accessories are available separately and allow the cushion to be adjusted where necessary to meet the specific positioning needs of the end user.

| | |
|---|--|
| Ⓕ | Abductor |
| Ⓒ | Clinical fluid supplement sac |
| Ⓗ | Lateral side wedge |
| Ⓘ | Cut to size anterior (front) wedge (Ⓙ) |
| ⓵ | Anterior (front) wedge |
| Ⓚ | RF/LR tapered wedge |
| Ⓛ | LF/RR tapered wedge |

Installing the abductor support (Fig. 3)

1. Remove the outer cover and inner liner.
2. Align the hook strip (Ⓜ) on the abductor support (Ⓕ) with the loop (Ⓝ) strip on the top of the foam base (Ⓑ) at the desired depth.
3. Press down to secure the abductor support to the foam base.
4. Re-install the inner liner and outer cover.

Installing the clinical fluid supplement sac (Fig. 4)

1. Remove outer cover and inner liner and unfasten the hook and loop strips from the side and/or rear of the clinical fluid sac (Ⓒ) where you wish to install the clinical fluid supplement sac(s) (Ⓒ).
2. Locate the small pocket/slot (Ⓓ) on the underside of the clinical fluid sac.
3. Install the clinical fluid supplement sac(s) inside the slot/pocket and adjust the clinical fluid supplement sac(s) to the desired position inside the clinical fluid sac.
4. Refasten the clinical fluid sac hook and loop strips and re-install the inner liner and outer cover.

Installing an anterior wedge, lateral wedge and/or tapered wedge (Fig. 5)

Foam wedges can be trimmed to size as needed.

1. Remove the outer cover and inner liner and position the wedge(s) in the desired position on the underside of the foam base (Ⓑ) as shown:
 - Lateral wedge (Ⓗ)
 - Anterior wedge (Ⓘ) and cut to size anterior wedge
 - RF/LR tapered wedge (Ⓚ) (LF/RR tapered wedge is not shown)
2. Using the hook and loop strips, secure wedge(s) to the underside of the foam base.
3. Re-install the inner liner and outer cover.

5 Maintenance

5.1 Inspection

WARNING!

- Do not continue to use this product if problems are discovered. Corrective maintenance can be performed at or arranged through your equipment supplier.

1. Visually inspect components for deformation, wear or compression and replace if necessary.

5.2 Cleaning and care

WARNING!

- To avoid infectious exposure, inspect cushion thoroughly and properly, clean and disinfect any materials soiled with blood or other bodily fluids.

**CAUTION!****Risk of damage to the product**

- 1% Chlorine Solution used on a regular basis can diminish the life of the cover if not rinsed and dried properly.
- Do not use cleaning agents containing phenols, bleaches or other abrasive materials.
 - Do not use granules.
 - Keep clear of open heat sources.

**IMPORTANT!**

- All cleaning agents and disinfectants used must be effective, compatible with one another and must protect the materials they are used to clean.
- For further information on decontamination in health care environments, please refer to 'The National Institute for Clinical Excellence' guidelines on Infection Control www.nice.org.uk/CGI39 and your local infection control policy.

Cleaning (light soilage) - Outer and inner cover

1. Wipe down the outer cover with a 0.1% Chlorine Solution (1,000 ppm).
2. Rinse the outer cover with clean water using a single use nonabrasive cloth.
3. Dry the outer cover thoroughly.

The inner liner should not be removed from the foam base.

4. If necessary, remove the outer cover and repeat steps 1-3 to clean the inner cover.

Cleaning (heavy soilage) - Outer cover only

Where the outer cover is badly soiled, we recommend cleaning at 80 °C in the washing machine. Large spillages of bodily fluids should be absorbed and removed with paper towels first, before washing the cover.

1. Clean up all spillages of bodily fluids i.e. blood, urine, faeces, sputum, wound exudate and all other bodily secretions as soon as possible using a 1% Chlorine Solution (10,000 ppm).
2. Rinse the outer cover with clean water using a single use nonabrasive cloth.
3. Dry cover thoroughly before re-fitting to the cushion.

Replacing covers

1. Carefully remove the cover by unzipping the zipper found on the side of the cushion.
2. Place new cover onto the foam base.

**IMPORTANT!**

- Ensure that the corners of the foam base are positioned correctly into the corners of the cover.

6 After Use**6.1 Storage****IMPORTANT!**

- Store cushions in a dry environment.
- Store cushions within a protective cover.
- Store items on clean, dry, off-flooring free from sharp edges to avoid any possible damage.
- Never store other items on top of a cushion.
- Do not store cushions next to radiators or other heating devices.
- Protect cushions from direct sunlight.

6.2 Re-use

A cleaning record must be kept as part of cleaning the system.

The product is suitable for repeated use. The number of times it can be used depends on how often and in which way the product is used.

It is recommended that the cushion is checked (foam and cover) for strike-through (this may include fluid ingress, stains, rips or damage) after the release of each patient by a suitably qualified and competent person.

1. Before reuse, clean the product thoroughly, see section 5.2 Cleaning and care, page 3.

6.3 Disposal

The disposal and recycling of used devices and packaging must comply with the applicable local legal regulations.

1. Ensure that the cushion is cleaned prior to disposal to avoid any risk of contamination.

7 Technical data**7.1 General data**

| Product | Matrx® Libra® |
|--|--|
| Total height [mm] | 89 |
| Total width [mm] | 360 – 560 |
| Total depth [mm] | 310 – 560 |
| Maximum user weight [kg] | 136 (for total widths from 360 – 510 mm) 227 (for total widths from 530 – 560 mm) |
| Minimum user weight [kg] | No minimum user weight |
| Weight complete cushion [kg] ¹⁾ | 1.6 – 3.5 |

¹⁾ Weights can vary depending on size of cushion ordered.

Size charts and local requirements / certification available upon request from Invacare®.

7.2 Materials

| | |
|--|--|
| Foam cushion | High resilient Polyurethane foam |
| Clinical fluid sac | Polyether based thermoplastic polyurethane film |
| Inner liner | Polyurethane ether film |
| Standard outer cover | Knitted polyamide, coated with an anti-microbial polyurethane coating. Substrate treated with a fluorocarbon to give water/stain resistance. |
| Spacer cover | 100% knitted polyester |
| All cushion components are Latex-free. | |

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1. Verwijder de buitenbekleding en de binnenvoering, en maak de haak en lusbanden los van de zij- en/of achterkant van de klinische vochtzak © waar u de aanvullende klinische vochtzak(ken) wilt installeren ©.
Invacare New Zealand Ltd
Invacare UK Operations Limited