

Invacare® Spectra® XTR2

Freedom of movement



The popular powerchair from Invacare providing comfort for day long use combined with good manoeuvrability and performance. It offers different styles of seating units which easily attach to the wheelbase to accommodate changes in postural needs. The chair can dismantle into manageable parts for transportation or storage.

Key Features

- Choice of seating to suit your needs
- Cable free battery system for quick and easier maintenance
- Rear suspension for added comfort
- Narrow base making it ideal for indoor use
- Easily dismantles for transportation

Technical data

	Standard	Flex II	Firm
Seat Width	44 - 48cm	39 - 53cm	34 - 48cm
Seat Depth	41, 46, 51cm	43 - 51cm	30 - 48cm
Overall chair width	59 or 62cm*	59 or 62cm*	59 or 62cm*
Max. user weight	130kg	130kg	130kg
Batteries	50Ah	50Ah	50Ah
Speed	6 or 10km/h	6 or 10km/h	6 or 10km/h
Range	30km**	30km**	30km**

* Wheelbase 63cm, differs according to wheel type

** According to ISO7176-4 2008, and depends on driving conditions.

Invacare® Mirage®

Easier way to mobility



This basic indoor/outdoor powerchair features a folding chassis that easily dismantles for transportation or storage. With additional options of tension adjustable upholstery, cushion seating, and backrest recline, the Mirage offers optimum support and comfort. It comes with adjustable armrests and swingaway legrests.

Key Features

- Choice of seat widths
- A range of options and accessories
- Manual wheel locks for safety
- 40Ah batteries
- Choose from two speeds either 3 or 6km/h
- Colours: Transparent Red, Transparent Blue, Metallic Green, Metallic Silver, Transparent Purple & Black

Technical data

	Standard
Seat Width	35 - 50cm
Seat Depth	40, 44, 50cm
Overall chair width	53-68cm*
Max. user weight	120kg
Batteries	40Ah
Speed	3 or 6km/h
Range	25km**

* Dependent on seat adjustment

** Based on 6km/h motors, according to ISO7176

Invacare® Dragon®

Compact and versatile



This cost effective powerchair provides good manoeuvrability both indoors or out. With its compact, stylish base, it features the latest controls, adjustable seat, suspension, rear luggage carrier and 6km/h motors, plus additional options to customise to your needs. This versatile powerchair is ideal for everyday use.

Key Features

- Choice of 3 seating types
- Manual settable seat tilt for optimum posture and comfort
- Tension adjustable back for additional support where needed
- Rear luggage carrier for those essential items
- Colour: Blue

Technical data

	Standard	Contour	Comfort
Seat Width	44 - 49cm	44 - 49cm	44 - 49cm
Seat Depth	41, 46, 51cm	41, 46, 51cm	41, 46, 51cm
Overall chair width	61cm**	61cm**	61cm**
Max. user weight	140kg	140kg	140kg
Batteries	60Ah	60Ah	60Ah
Speed	6km/h	6km/h	6km/h
Range	30km**	30km**	30km**

* Dependent on seat adjustment

** According to ISO7176-4 2008

Invacare® Pronto® M61™

A fresh approach to mobility




Designed around a compact base, the Pronto M61 provides excellent manoeuvrability especially in confined spaces. It features a centre wheel drive for a stable and safe drive with an electric seat riser to improve accessibility in and around the home or when out and about. The Pronto also has Sure Step technology for a smoother drive over obstacles and low thresholds.

Technical data

Choice of 16" or 18" wide seat Speed: 6.4 km/h Batteries: 2x 30Ah
Max. user weight: 136Kg Colours: Deep Blue Pearl or Candy Red Pearl

Invacare Authorised Dealer



© 2011 Invacare® International Sarl. All rights reserved. All information quoted is believed to be correct at time of print. Invacare® reserves the right to alter product specifications without prior consultation.

Invacare Powerchair leaflet part no. 1551519
www.invacare.co.uk



Powerchair Range 2011



Yes, you can.™

Invacare® Storm® 4

Simply the best fit with easy adjustments



This stylish and revolutionary powerchair has all the comfort and power you need for day long use. Suitable for both outdoors and indoors, its key feature is adjustability allowing the chair to fit around you. With its seat riser, spring suspension and new lighting system, Storm 4 has it all. Available in a choice of colours.

Key Features

- Highly configurable and adapts to a range of needs
- Adjustable seat widths from 38 to 55cm
- Manual tilt in space
- Multi adjustable arms and legrests
- Seat riser providing 25cm height adjustment
- Spring suspension for a comfortable drive

Technical data

	Standard	Easy Adapt	Recaro	Onyx
Seat Width	41, 46, 51cm	38 - 53cm	49 - 53cm	40 - 55cm
Seat Depth	38 - 53cm	38 - 53cm	46 - 51cm	40 - 45cm
Overall chair width	63cm*	63cm*	63cm*	63cm*
Max. user weight	150kg	150kg	150kg	150kg
Batteries	70Ah	70Ah	70Ah	70Ah
Speed	6 or 10km/h	6 or 10km/h	6 or 10km/h	6 or 10km/h
Range	36km**	36km**	36km**	36km**

* Wheelbase 63cm, differs according to wheel type

** Based on Storm 4 with 6km/h motors, according to ISO7176 and depends on driving conditions

Invacare® Typhoon® II

The Power of motion



This high performance mid-wheel drive powerchair offers outstanding performance for indoor or outdoor manoeuvrability. It features Invacare's Sure Step system for improved stability especially over uneven terrain, as well as the innovative 'walking beam' which lifts the castors when in contact with an obstacle providing a safe drive.

Key Features

- Offers a range of seating types
- Manual tilt in space
- Innovative walking beam for improved driving over obstacles
- Spring suspension for a comfortable drive
- Choice of colours: Cherry Red, Pearl Night Blue, Black

Technical data

	Standard	Contour	Recaro	Flex
Seat Width	44 - 48cm	44 - 53cm	49 - 53cm	39 - 53cm
Seat Depth	41, 46, 51cm	41, 46, 51cm	46 - 51cm	38, 43, 48cm
Overall chair width	65 - 69cm*	65 - 69cm*	65 - 69cm*	65 - 69cm*
Max. user weight	150kg	150kg	150kg	150kg
Batteries	60Ah	60Ah	60Ah	60Ah
Speed	6, 10, 12km/h	6, 10, 12km/h	6, 10, 12km/h	6, 10, 12km/h
Range	52km**	52km**	52km**	52km**

* Wheelbase differs according to wheel type

** Based on 6km/h motors, according to ISO7176-4 2008, and depends on driving conditions

Invacare® TDX® SP & SPNB

Get more than mobility



Available in two versions, the TDX SP offers a unique drive experience with its sleek chassis design and highly manoeuvrable mid-wheel drive that provides good handling abilities. Both versions are big on style and performance featuring powerful motors, spring suspension and a choice of seating. The TDX SP NB model has a compact wheelbase ideal for confined spaces.

Key Features

- High performance with powerful and efficient motors
- Manual seat tilt for optimum positioning
- Stability lock keeps all 6 wheels on the ground
- Spring suspension for added comfort
- Includes Sure Step for smooth travel over obstacles

Technical data

	Standard	Contour	Flex	TDX SP NB
Seat Width	44 - 48cm	44 - 53cm	39 - 53cm	39 - 53cm
Seat Depth	41, 46, 51cm	41, 46, 51cm	38, 43, 48cm	38 - 51cm
Overall chair width	64cm*	64cm*	64cm*	61cm*
Max. user weight	180kg	180kg	180kg	150kg
Batteries	73Ah	73Ah	73Ah	50Ah
Speed	8km/h	8km/h	8km/h	6km/h
Range	39km**	39km**	39km**	26km**

* Depends on seat adjustment

** According to ISO7176-4 2008, and depends on driving conditions

Invacare® FDX®

Superior performance - indoors or out



FDX is a highly manoeuvrable front wheel drive powerchair for the active person. It features Invacare's unique G-Trac technology designed to provide directional stability outdoors, whether at high speed or across uneven ground. FDX has a new configurable seat concept, and offers low maintenance with easy access to all parts.

Key Features

- Includes a seat lift of 30cm, including powered tilt function
- Rear suspension for a safe and comfortable drive
- G-Trac fitted as standard
- A range of controls options available
- Available in 2 colours: Metallic Silver or Electric Red

Technical data

	Standard	Easy Adapt	Onyx	Recaro
Seat Width	38 - 53cm	38 - 53cm	40 - 55cm	49 - 53cm
Seat Depth	41, 46, 51cm	38 - 53cm	38 - 53cm	46 - 51cm
Overall chair width	59.5cm*	59.5cm*	59.5cm*	59.5cm*
Max. user weight	136kg	136kg	136kg	136kg
Batteries	60Ah	60Ah	60Ah	60Ah
Speed	6 or 10km/h	6 or 10km/h	6 or 10km/h	6 or 10km/h
Range	32km**	32km**	32km**	32km**

* Based on 6km/h motors with 12" wheels

** According to ISO7176-4 2008, and depends on driving conditions

Invacare® Kite™

The freedom to do... the power to go!



This highly configurable powerchair offers great performance in a compact design. It features Invacare's patented Dual Swing Technology which boosts the traction for improved driving and comfort. With a host of options such as controls, seat tilt, and backrest recline, the Kite raises the standard for innovation and technology.

Key Features

- Compact base ideal for indoor or outdoor use
- Dual Swing Technology for improved traction
- Shock Absorbers provides a smoother ride over rough terrain
- Includes 4-point tie down for transportation
- Powerful 60Ah batteries

Technical data

	Standard	Flex 2 Seat
Seat Width	44 - 48cm	39 - 53cm
Seat Depth	41, 46, 51cm	43 - 51cm
Overall chair width	59.5cm*	63cm*
Max. user weight	160kg	160kg
Batteries	60/73Ah	60/73Ah
Speed	6 or 10km/h	6 or 10km/h
Range	25-35km**	25-35km**

* Wheelbase differs according to wheel type

** According to ISO7176-4 2008, with 6km/h motors range is 27 to 35km, with 10km/h motors range is 25 to 31km, and depends on driving conditions