



Get Perfect Support - Remain Active!

The new *Invacare* **XLT** is the ideal active, rigid chair for day to day use. Designed to be perfectly adjusted to fit the body the **XLT** offers individuals wide ranging levels of support, enabling and even improving independence and mobility.



Yes, you can.®



Top Configurability made Simple!

The *Invacare XLT* has been designed with the individual in mind. Measured, constructed and fitted to suit the precise requirements of each person, the **XLT** is the ideal chair to help regain, retain and improve both mobility and independence.

Complex chair-specific solutions and spare parts are a thing of the past with the new **XLT**. Benefit from the extensive and readily available range of Invacare options which can be fitted to the **XLT** without any hassle.

Stability, Rigidity and Security for an Active Life

Made entirely from titanium the **XLT** gives stability and rigidity without compromising on weight. In combination with the wide range of options, features and configuration possibilities offered with the **XLT**, the client will benefit from a chair which offers the maximum security and support needed when going about daily activities.

Proven Quality

Invacare, together with *Rea*[®] have made no compromises when it comes to the new **XLT**. Designed and built at the Rea factory in Sweden, individuals will

benefit from the strong tradition of both *Invacare* and *Rea* quality fused to create one ultimately unique chair.

Features and Options



Armrests

Wide range of armrests for best support.



Anti tipper device

The Anti-tipper device can easily be folded away by the foot. The device is also angled outwards, to give space for the feet for the attendant.



Standard Wheel

With black Marathon Evolution Plus tyres for design.

Features and Options



Front Caster Attachments
For stability and performance.



Titanium foot tube
To improve comfort as well as lightness.

Swing Away
For easy transfer and manoeuvrability as well as comfort.



Tyres
① solid, ② pneumatic profiled or puncture proof, ③ pneumatic non marking, ④ low profile, ⑤ pneumatic profiled or puncture proof, ⑥ black marathon evolution.



Frontwheels
① 75 x 24 hard (3"), ② 100 x 27 medium hard (4"), ③ 125 x 27 medium hard (5"), ④ 140 x 37 soft (5,5"), ⑤ 150 x 27 medium hard (6"), ⑥ 150 x 30 soft (6"), ⑦ 150 x 32 pneumatic (6"), ⑧ 180 x 40 pneumatic (7"), ⑨ 180 x 40 soft (7")

Technical data



Invacare XLT Active	355 - 480 mm	400 / 450 / 500 mm	400 - 520 mm	200 - 490 mm	380 - 505 mm	0 - 14°
Invacare XLT Dynamic	355 - 500 mm	400 / 450 / 500 mm	400 - 500 mm	200 - 490 mm	355 - 485 mm	0 - 14°
Invacare XLT Swing	330 - 505 mm	400 / 450 / 500 mm	390 - 520 mm	200 - 490 mm	280 - 510 mm	0 - 14°

in increments of 15 mm on 40 mm adj.



Invacare XLT Active	± 10°	SW + 200 mm	610 - 1050 mm	830 - 1030 mm	9.8 kg	135 kg
Invacare XLT Dynamic	± 10°	SW + 200 mm	610 - 1050 mm	730 - 880 mm	9.8 kg	135 kg
Invacare XLT Swing	± 10°	SW + 200 mm	610 - 1050 mm	820 - 1230 mm	11.3 kg	135 kg



Invacare XLT Active	6 kg
Invacare XLT Dynamic	6 kg
Invacare XLT Swing	6 kg



Crash-tested
ISO 7176-19

The **Invacare XLT** complies with ISO 7176-19 specifications. Note that special conditions apply regarding accessories. For detailed information regarding vehicular transportation and for more pre-sale information for this wheelchair, please consult the following pages on our website: www.invacare-rea.com under Rea wheelchairs/ Technical downloads.

Frame colours



Blue

Red

Black



Grey

Blasted Titanium*
not for Swing

Orange*

* Available as special orders

Upholstery



Black Nylon

Black Jemima

Invacare Limited
Pencoed Technology Park - Pencoed -
Bridgend CF35 5AQ - United Kingdom
Tel: +44 1 656 776 222
Fax: +44 1656 776 220
E-mail: uk@invacare.com
Sales Order E-mail:
ordersuk@invacare.com
www.invacare.co.uk

Invacare Ireland Ltd
Unit 5 - Seatown Business Campus
Seatown Road - SWORDS - County
Dublin - Ireland
Tel: +353 1 8107084
Fax +353 1 8107085
Email: Ireland@invacare.com
www.invacare.ie

© 2011 Invacare International Sàrl
All rights reserved.

All information quoted is believed to be correct at time of print.
Invacare reserves the right to alter product specifications without prior consultation.

XLT - UK - 10/2011 - 1504606



Yes, you can.®