

küschall[®]
designed for life




SET IN
MOTION



KÜSCHALL K-SERIES*attract*

KÜSCHALL K-SERIES *attract*

The küschall K-Series attract features all the hallmarks you can expect from a precision engineered Küschall product. Based on the iconic K-Series, this sleek and stylish chair features high levels of adjustment and exceptional driving performance with its rigid energy efficient design. The K-Series attract is lightweight with a minimalist design and a durable frame boasting a five year warranty* that offers a cost effective mobility solution.



A tapered frame, a variety of front and rear seat heights, centre of gravity adjustment, backrest angle and footrest height settings are just some of the product features available. It's folding backrest and quick release wheels makes the chair lightweight and easy to transport. This active chair offers the perfect combination of quality and reliability from over 35 years of Küschall innovation. Available in 4 durable, matt colours.



Technical data

	Frame angle: 75° / 90° tapered		Seat width: 360 – 480 mm in increments of 20 mm		Seat depth: 375 – 475 mm in increments of 25 mm		Seat height: front: 480 – 510 mm rear: 400 – 480 mm in increments of 10 mm
	Knee-to-heel length: 290 - 480 mm in increments of 10 mm		Backrest height: 315– 480 mm in increments of 15 mm		Backrest angle: 74°/78°/82°/86°/90°		Max. safe slope: 7°
	Reversing width: 920 - 1230 mm		Total width: SW + 170 mm Total width with 3° camber: SW + 220 mm		Total length: 75°: approx. 830 mm 90°: approx. 760 mm		Total weight unloaded: approx. 10.3 kg
	Transport weight w/out rear wheels: approx. 7 kg		Max. user weight: 120 kg		Crash Test approved** Fulfills the requirements of ISO 7176-19		

Invacare Ltd. | Pencoed Technology Park, Pencoed, Bridgend CF35 5AQ
Phone +44 (0) 1656 776222 | Fax +44 (0) 1656 776220 | www.invacare.co.uk

* side frame only
** Crash test – Back height over 405mm

küschall is a registered trademark. Every effort has been made to ensure that the contents of this publication are fully up-to-date at the time of going to print. As part of its ongoing improvement of products, Küschall AG reserves the right to modify existing models at any time. Any use of this publication or of parts thereof as well as any reproduction of images must have the written consent of Küschall AG. Patents pending K-Series attract_UK_2014 - 1582252 Copyright © 12/2014, Küschall AG, Switzerland - All rights reserved