



*Invacare*<sup>®</sup>

# ***Essential Care***

A cost effective solution for comfort and pressure redistribution



The **Essential Care** mattress is an entry level pressure redistributing support surface for the community, residential healthcare environment. Its profile cut surface provides a cost effective solution for comfort and pressure redistribution to meet the needs of individuals 'at risk' or 'high risk' of developing pressure ulcers.

- **Profile cut**
- **Multi-stretch polyurethane fabric cover**
- **Zipped around 3 sides for ease of inspection and laundering**



## Features and Options



### Cover



Multi-stretch polyurethane fabric cover.

### Profile cut foam

Improves comfort and ensures effective pressure redistribution.

## Technical data

For more information about this product, including the product's user manual, please visit [www.invacare.co.uk](http://www.invacare.co.uk).

	Size range	Max. user weight
<b>Essential Care</b>	 2000 x 880 x 150 mm 2050 x 880 x 150 mm	 140 kg (22 stone)

1. Essential nursing care is pivotal in pressure ulcer prevention. This mattress will positively contribute to the outcome of a pressure ulcer prevention care plan. Education, clinical judgement and action based planning based on vulnerability are fundamental factors in the prevention of pressure ulcers. A range of assessment scales can be used as a formal method of assessing risk from pressure ulcer development, and should be used in conjunction with an informal assessment (informed nursing judgement). Informal assessment is considered to be of greater importance and clinical value.

2. The risk category in this instance is based on the Waterlow Scale.

Mattresses is manufactured to:  
 BS EN 597-1 Cig, BS EN 597-2 Match and BS 7177: Crib Source 5

### Colours

Blue



### Warranty

1 year

Invacare Limited  
 Pencoed Technology Park - Pencoed  
 Bridgend CF35 5AQ - United Kingdom  
 Tel: +44 1656 776 222  
 Fax: +44 1656 776 220  
 E-mail: [uk@invacare.com](mailto:uk@invacare.com)  
 Sales Order E-mail:  
[ordersuk@invacare.com](mailto:ordersuk@invacare.com)  
[www.invacare.co.uk](http://www.invacare.co.uk)

Invacare Ireland Ltd  
 Unit 5 - Seatown Business Campus  
 Seatown Road - SWORDS -  
 County Dublin - Ireland  
 Tel. +353 1 8107084  
 Fax +353 1 8107085  
 Email: [Ireland@invacare.com](mailto:Ireland@invacare.com)  
[www.invacare.ie](http://www.invacare.ie)

Trusted Solutions,  
 Passionate People

