

# Invacare® | *Heelpad*



## **Pressure Redistributing Heelpad**

The *Invacare Heelpad* is designed to reduce pressure on the vulnerable heel area. The controlled volume of gel in the sacs allows pressure displacement and even weight distribution. Pressure on a patient's heels are kept to a minimum when the heels are in contact with the gel sacs. The *Invacare Heelpad* was designed to be used in conjunction with the *Invacare Odstock Wedge*, but can also be used independently as an effective pressure reducing device.



**Yes, you can.®**



Features and Options

**Designed to reduce pressure on the heels**

The *Invacare Heelpad* enables vulnerable heels to be immersed into independent gel sacs, offering significant pressure reduction and comfort to clients.

**Enhanced combined performance**

The *Heelpad* offers enhanced performance, when used in conjunction with the *Odstock Wedge*, by providing improved support for the lower limbs, and pressure redistribution for the heels.



**Single and Double Heelpads available**

Single and Double formats available to meet specific patient requirements.



**Heelpads supplied with three platilon covers**






Platilon protective covers are disposable, and help to avoid contamination of the heelpad.



**Designed to be used in conjunction with the Invacare Odstock Wedge**

Used in conjunction, these two products combine to significantly reduce pressure on the vulnerable heel area.

**Technical data**

					
<b>Heelpad</b>	Width: 62 cm	Length: 48 cm	Height: 5 cm, 3.2 cm	Weight: 3.58 kg	Cleaning: 80°

Manufactured to European Standards:  
 Interior Foam: BS 5852: 1982 Crib 5, BSEN 597-1,2  
 Cover: BS 7175: 1989 Crib 5

Invacare Limited  
 Pencoed Technology Park - Pencoed -  
 Bridgend CF35 5AQ - United Kingdom  
 Tel: +44 1 656 776 222  
 Fax: +44 1656 776 220  
 E-mail: [uk@invacare.com](mailto:uk@invacare.com)  
 Sales Order E-mail:  
[ordersuk@invacare.com](mailto:ordersuk@invacare.com)  
[www.invacare.co.uk](http://www.invacare.co.uk)

Invacare Ireland Ltd  
 Unit 5 - Seatown Business Campus  
 Seatown Road - SWORDS - County  
 Dublin - Ireland  
 Tel. +353 1 8107084  
 Fax +353 1 8107085  
 Email: [Ireland@invacare.com](mailto:Ireland@invacare.com)  
[www.invacare.ie](http://www.invacare.ie)

© 2014 Invacare International Sàrl  
 All rights reserved.

All information quoted is believed to be correct at time of print.  
 Invacare reserves the right to alter product specifications without prior consultation.

**Colour**



Fuchsia

**Warranty 1 Year**

Softform Heelpad - UK - 09/2014  
 1505431



**Yes, you can.®**