



Yes, you can.®

ACTION® 5 – FOR AN ACTIVE LIFESTYLE

COLOUR OPTIONS

- Dark Red
- Blue pop
- Anthracite Black
- Pearl White
- Electric Green
- Hypnotic Purple
- Pearl Grey
- Matt Black
- Intense Yellow
- Vibrant Orange

TOP FEATURES

- » A unique, H-Shaped folding frame allows compact folding
- » Extreme, low rolling resistance for easy propulsion
- » A patented, stepless, angle adjustable backrest system with forward folding function



ACTION® 4NG/HD – TALLER OR HEAVIER REQUIREMENTS

COLOUR OPTIONS

- Anthracite Black
- Matt Black
- Slate Grey
- Hypnotic Purple
- Blue Pop
- Electric Green
- Matt Red
- Pearl White

TOP FEATURES

- » A stepless, adjustable folding backrest allowing for easier storage and transportation
- » A market leading, lowest total weight of 14kg and transportation weight of 8kg



NEW

ACTION® 3NG LIGHT – LIGHTWEIGHT FOR AN EXCEPTIONAL DRIVE

COLOUR OPTIONS

- Anthracite Black
- Matt Black
- Slate Grey
- Hypnotic Purple
- Blue Pop
- Electric Green
- Matt Red
- Pearl White

TOP FEATURES

- » From as light as 13.0kg
- » Lightweight rear wheels
- » Quality performance with effortless propulsion



ACTION® 3NG – HIGHLY CONFIGURABLE FOR CHANGING NEEDS

COLOUR OPTIONS

- Anthracite Black
- Matt Black
- Slate Grey
- Hypnotic Purple
- Blue Pop
- Electric Green
- Matt Red
- Pearl White

TOP FEATURES

- » Folding forward backrest for easier transportation and storage
- » Reclining backrest, which allows repositioning in one continuous movement
- » A market leading, lowest total weight of 13kg and transportation weight of 7.1kg



ACTION® 2NG – BASIC / SHORT TERM USE

COLOUR OPTIONS

- Anthracite Black

TOP FEATURES

- » Adjustable centre of gravity
- » Aluminium, lightweight frame, for easy transportation



ACTION® 3JR – PAEDIATRIC REQUIREMENTS

COLOUR OPTIONS

- Pearl White
- Matt Black
- Matt Red
- Intense Yellow
- Polar Blue
- Candy Pink

TOP FEATURES

- » A telescopic cross brace which optimises the chair to the changing shape of the child as they grow
- » A range of customisation options designed for kids
- » Adjustable seat position and centre of gravity settings to optimise propelling



THE INVACARE ACTION FAMILY FOR EVERY GENERATION BY DAY BY NIGHT

