



Aquatec[®] | ***Pure Bidet***

Promoting independence and dignity



The *Aquatec Pure Bidet* is ideal for those who experience difficulty with **personal hygiene**, or for those simply seeking **added luxury** at home. The modern looking bidet **fits easily** and **discreetly** to existing toilets and provides a gentle yet effective **cleansing solution** using a warm water spray. Key features include:

- **Easy to operate with intuitive control**
- **Compact, self contained technology**
- **Installed to existing toilets**

The *Aquatec Pure Bidet* has full Water Regulations Advisory Scheme (WRAS) approval and also conforms to the directive 93/42/EEC for medical devices.

High functionality

- **Intuitive** hand control with clear buttons
- **Adjustable** water temperature up to 35°, with **aerated** spray and variable jet settings
- **Personal** spray head positioning for both posterior and feminine settings
- Adjustable **warm air dryer**
- Temperature controlled **heated seat up to 38°**
- **LED** light to enable 24 hour use
- **Soft close** seat and lid



Inbuilt technology

- **Continuous** water heating for **instant** warm water at all times
- Automatic spray head cleansing before and after use to **maintain cleanliness** of the system
- Equipped with a **thermal protector** to prevent overheating
- Seat sensor activates **automatic shut down** when bidet is not in use
- **Self-diagnosis** feature to automatically detect malfunction
- Residual-current device (RCD) for **added safety** and to provide additional peace of mind



Fully compliant to fit existing toilets

- Uses **standard pan fixing points** for easy fitting ¹
- Straightforward to install ²
- WRAS approved to ensure safe connection to public water pipes
- Conforms to DIN EN 1717 and Directive 93/42/EEC for medical devices

1. An adapter plate is available for non standard pan sizes
2. Must be installed by a competent professional



Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please visit www.invacare.co.uk.

Total width	389 mm	Type of protection, entire device	IP X4
Total length	515 mm	Battery type remote control	2 x AA (LR6) / 1.5 V
Total depth	156 mm	Water supply pressure	79 – 735 kPa
Product weight	3.9 kg	Water supply temperature	3 – 35 °C
Max. user weight	125 kg	Water flow rate CLEANSING	0.7 l/min
Length of power cord	1.4 m	Water flow rate BIDET	0.85 l/min
Power supply	220 – 240 V AC, 50/60 Hz	Water pressure adjustment	3 levels
Max. power consumption	1725 – 2053 W	Warranty	2 years

Invacare Limited
 Pencoed Technology Park - Pencoed
 Bridgend CF35 5AQ - United Kingdom
 Tel: +44 1 656 776 222
 Fax: +44 1656 776 220
 E-mail: uk@invacare.com
 Sales Order E-mail: ordersuk@invacare.com
www.invacare.co.uk

Invacare Ireland Ltd
 Unit 5 - Seatown Business Campus
 Seatown Road - SWORDS - County Dublin - Ireland
 Tel. +353 1 8107084
 Fax +353 1 8107085
 Email: Ireland@invacare.com
www.invacare.ie

