

# Birdie®/Birdie® Compact



A straightforward solution for safe and comfortable transfers



With maximum space for the client, both the *Invacare Birdie* and *Birdie Compact* offer a truly comfortable lift and transfer to or from beds, chairs or even the floor. Ideal for use in domestic and residential care environments, the hoists are designed to fold and unfold easily without the need for tools.



### Comfort and ease of use combined

The high boom maximises space for the client even when the hoist is in its highest position. This gives an ergonomic lift and transfer as the client can be easily rotated 360°, which makes positioning much easier for the carer. The optimisation of space allows increased space in front of the actuator, reducing the risk of injury to the client's knees and feet.

### Ease of handling: simple to fold, unfold and dismantle without tools

Handling the hoist has been made as smooth and easy as possible as the hoist can easily be folded and unfolded without the need for tools. Hoists take up minimal space when stored and are easy to push and transport. If required, the mast and boom can also be easily detached from the base for transportation (no tools required).



### TÜV-certification for maximum security

For total security, the folding mechanism has a safety system. The *Birdie* and *Birdie Compact* are CE-certified in line with the Medical Device Directive 93/42/ECC. To ensure the hoists comply with EN ISO 10535 (the European standard for hoists), the hoists have been rigorously tested both by Invacare and the German TÜV institute, an accredited external body.

### Features and options



#### Long boom - high user comfort

The client can easily be rotated and positioned. Reduced risk of injury to the knees.



#### Easy to fold and dismantle

Ease of transport and storage. No tools required when folding / unfolding or when dismantling the hoist.



#### Snap lock

Optimal security.



#### Emergency stop, lowering and integrated charger


Manual or electric emergency lowering (compact version only electric lowering). Control box is supplied with a detachable or fixed battery.

## Features and options



### Clip on spreader bar

Easy to change the spreader bar – just twist the safety catch, press to open - and then detach/attach the spreader bar.



### Mounted spreader bar

Allows correct and safe positioning of the spreader bar. A locking carabiner offers exceptional security.



### Manual leg spread

Simply step on the pedal to adjust the leg spread.



### Handle for leg spread

Handle can easily be attached to operate the manual leg spread. (Optional).



### Electrical leg spread

Available for *Invacare Birdie* models only. (See pricelist).



### 2 point quick release spreader bars

Easy to attach and detach. Available widths: 45cm and 55cm.



### Push bar and hand control

Superior ergonomics for carers.



### Castors

100 mm castors (Rear braking castor is shown).



### Scale

If needed, a Homecare scale is available. Approved, according to medical standard EN 45501.



### 4 point quick release spreader bars

Easy to attach and detach. Available widths: 45 cm and 55 cm.

## Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please visit [www.invacare.co.uk](http://www.invacare.co.uk).

	Overall length of base	Outside width: legs closed	Total product weight	Lifting range	Height of legs	Inside width: legs opened	Max. user weight
<b>Birdie</b>	1240 mm	660mm	42 kg	460 - 1710 mm	100 mm	1025mm	180 kg
<b>Birdie Compact</b>	1090 mm	540mm	36 kg	540 - 1630 mm	100 mm	860mm	150 kg

	Number of cycles	Battery capacity	Turning diameter
<b>Birdie</b>	40 cycles	2,9 Ah	1400 mm
<b>Birdie Compact</b>	40 cycles	2,9 Ah	1070 mm

All measurements are taken with 100m castors and 2 point spreader bar.

## Frame colour

Light grey  
Dark grey



Please note, colors may vary slightly from those displayed above.

## Warranty

5 year warranty:  
Frame  
Actuator  
Control box

2 year warranty:  
Hand control

1 year warranty:  
Battery

The *Invacare Birdie* and *Birdie Compact* are part of the Safe Patient Handling range.



Safe  
Patient  
Handling

Invacare Limited  
Pencoed Technology Park - Pencoed  
Bridgend CF35 5AQ - United Kingdom  
Tel: +44 1 656 776 222  
Fax: +44 1656 776 220  
E-mail: [uk@invacare.com](mailto:uk@invacare.com)  
Sales Order E-mail: [ordersuk@invacare.com](mailto:ordersuk@invacare.com)  
[www.invacare.co.uk](http://www.invacare.co.uk)

Invacare Ireland Ltd  
Unit 5 - Seatown Business Campus  
Seatown Road - SWORDS  
County Dublin - Ireland  
Tel. +353 1 8107084  
Fax +353 1 8107085  
Email: [Ireland@invacare.com](mailto:Ireland@invacare.com)  
[www.invacare.ie](http://www.invacare.ie)



© 2017 Invacare International Sàrl. All rights reserved. All information quoted is believed to be correct at time of print.

Invacare reserves the right to alter product specifications without prior consultation. Birdie - UK - 03/2018 - 1504431