



Invacare®

Comet® PRO SPORT

Power and performance in
a stylish package



Whether you're popping to the local shops, taking longer trips into town, or looking to go further afield, the all new **Comet PRO SPORT** will take you the distance in comfort, safety and style. Its powerful outdoor capabilities, combined with advanced safety features and options, ensures a speedy, powerful and reliable drive. With its exceptional driving features and modern aesthetics, the **Comet PRO SPORT** offers power and performance in a stylish package.

Features and options



Advanced suspension

All new suspension system for a smooth and comfortable drive over all terrains.



Sport wheels

13" Sport wheel with pneumatic black tyres giving a sporty look.



Rear view mirrors

Newly designed sleek, sports rear view mirrors.



Technical data

For more comprehensive information about this product, including the product's user manual, please visit www.invacare.co.uk

Backrest width	Seat depth	Seat to floor height	Backrest height	Driving unit width	Total length	Empty weight	Max. driving distance	Max. driving speed
510 mm	460 mm	Min. 440 mm Max. 510 mm	630 mm	665 - 685 mm	1450 mm	136 kg	55 km (34 miles)	12.8 km/h (8MPH)
The heaviest part	Max. user weight	Battery capacity	Turning diameter	Max. climbable obstacle height	Ground clearance	Max. safe slope	Motor capacity	
								1. Peak performance
approx. 63 kg	160 kg (25st)	2 x 12 V/75 Ah (C20) AGM	2750 mm	100 mm	100 mm	10 °	550 W / 1500 W ¹	

Colour of shrouds

Classic

Marble White Sapphire Blue Jade Green Diamond Silver Jasper Red




Carmine Red Iconic Chrome Electric Blue Luxor Beige Crystal Copper



Hydrographics

Premium

Polygon shades of blue Polygon shades of red Hexagon shades of white



Discover all your options at www.invacarescooter.com



Invacare Limited
Pencoed Technology Park - Pencoed
Bridgend CF35 5AQ - United Kingdom
Tel: +44 1656 776 222
Fax: +44 1656 776 220
E-mail: uk@invacare.com
Sales Order E-mail:
ordersuk@invacare.com
www.invacare.co.uk

Invacare Ireland Ltd
Unit 5 - Seatown Business Campus
Seatown Road - SWORDS -
County Dublin - Ireland
Tel. +353 1 8107084
Fax +353 1 8107085
Email: Ireland@invacare.com
www.invacare.ie



© 2017 Invacare International Sàrl. All rights reserved. All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation.
Comet Pro Sport - UK - 11/17 - 1629268