

Invacare® | **Flexipad**



A range of versatile support pads

Invacare Flexipads are highly versatile pads filled with silicone gel, and can be used to position and support limbs, joints, the torso or the head.



Yes, you can.®

Improves positioning

Invacare **Flexipads** can be used to position limbs for treatment, whilst recovering from an injury, or to assist in maintaining a more comfortable position to alleviate pain.

Improves support

Easy to use and adjust to position, **Flexipads**, can provide supplementary short-term support for patients experiencing difficulty in maintaining posture.



Features and Options



Single Flexipad

Single gel sac one side with foam the other - Fitted with a strap for fastening to chairs, the Invacare **Flexipad** increases the surface area under elbows, arm or back.



Double Flexipad

Two gel sacs wrapped in a visco elastic (memory) foam – highly versatile in a number of applications, including supporting large bony prominences.



Flexipad Ultra

Two gel sacs wrapped in a stand high density foam - A lighter, more flexible format, available with or without straps.

Technical data

					
Flexipad	Width: 22 cm	Length: 33 cm	Height: 5 cm, 3.2 cm	Weight: 0.8 kg	Cleaning: 80°

Manufactured to European Standards:
Cover: BS 7175: 1989 (Crib 5)

Invacare Limited
Pencoed Technology Park - Pencoed -
Bridgend CF35 5AQ - United Kingdom
Tel: +44 1 656 776 222
Fax: +44 1656 776 220
E-mail: uk@invacare.com
Sales Order E-mail:
ordersuk@invacare.com
www.invacare.co.uk

Invacare Ireland Ltd
Unit 5 - Seatown Business Campus
Seatown Road - SWORDS - County
Dublin - Ireland
Tel. +353 1 8107084
Fax +353 1 8107085
Email: Ireland@invacare.com
www.invacare.ie

© 2014 Invacare International Sàrl
All rights reserved.

All information quoted is believed to be correct at time of print.
Invacare reserves the right to alter product specifications without prior consultation.

Colours



Black Smooth

Warranty 1 Year

Flexipad - UK - 09/2014 - 1505438



Yes, you can.®