



# Invacare Robin

## Product Data Sheet



### Product Range Available

- Robin – Unique strap system without the use of a spreader bar
- Robin Mover – With an integrated traverse motor

### Standard Features

- Two strap design gives unsurpassed stability during the transfers
- Smooth soft start function makes positioning gentle
- Two speed motor
- Absence of a spreader provides the client with a sense of freedom reducing the fear of being close to obtrusive parts of the lifting equipment
- Lowering, raising and stop functions as standard
- Strap inlet self-adjusts automatically to size of client
- Hook self-locks when weighted
- Cable disengages to avoid damage if unit is pulled along via hand control

### Technical Specification

#### Robin & Robin Mover

- Product weight 7kg
- Maximum user weight 200kg (32 stone)
- Warranty 3 years

**Robin Mover** - Product weight 9.5kg

### Supporting Documentation

- Invacare's Business Development Managers can provide support with client assessments and provide training for OT's and carers

### Manufacturers Contact Details

Invacare Ltd, Unit 4 Pencoed Technology Park, Bridgend, CF35 5AQ

Tel: 01656 776222 E-mail: [UK@invacare.com](mailto:UK@invacare.com)

Web: [www.invacare.co.uk](http://www.invacare.co.uk) / [www.thinkindependence.co.uk](http://www.thinkindependence.co.uk)



safepatienthandling