



Aquatec® | **Pure Bidet**

Promoting independence and dignity



The *Aquatec Pure Bidet* is ideal for those who experience difficulty with **personal hygiene**, or for those simply seeking **added luxury** at home. The modern looking bidet **fits easily** and **discreetly** to existing toilets and provides a gentle yet effective **cleansing solution** using a warm water spray. Key features include:

- **Easy to operate with intuitive control**
- **Compact, self contained technology**
- **Installed to existing toilets**

The *Aquatec Pure Bidet* has full Water Regulations Advisory Scheme (WRAS) approval and also conforms to the directive 93/42/EEC for medical devices.

### High functionality

- **Intuitive** hand control with clear buttons
- **Adjustable** water temperature up to 35°, with **aerated** spray and variable jet settings
- **Personal** spray head positioning for both posterior and feminine settings
- Adjustable **warm air dryer**
- Temperature controlled **heated seat up to 38°**
- **LED** light to enable 24 hour use
- **Soft close** seat and lid



### Inbuilt technology

- **Continuous** water heating for **instant** warm water at all times
- Automatic spray head cleansing before and after use to **maintain cleanliness** of the system
- Equipped with a **thermal protector** to prevent overheating
- Seat sensor activates **automatic shut down** when bidet is not in use
- **Self-diagnosis** feature to automatically detect malfunction
- Residual-current device (RCD) for **added safety** and to provide additional peace of mind



### Fully compliant to fit existing toilets

- Uses **standard pan fixing points** for easy fitting <sup>1</sup>
- Straightforward to install <sup>2</sup>
- WRAS approved to ensure safe connection to public water pipes
- Conforms to DIN EN 1717 and Directive 93/42/EEC for medical devices

1. An adapter plate is available for non standard pan sizes
2. Must be installed by a competent professional



## Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please visit [www.invacare.co.uk](http://www.invacare.co.uk).

Total width	389 mm	Type of protection, entire device	IP X4
Total length	515 mm	Battery type remote control	2 x AA (LR6) / 1.5 V
Total depth	156 mm	Water supply pressure	79 – 735 kPa
Product weight	3.9 kg	Water supply temperature	3 – 35 °C
Max. user weight	125 kg	Water flow rate CLEANSING	0.7 l/min
Length of power cord	1.4 m	Water flow rate BIDET	0.85 l/min
Power supply	220 – 240 V AC, 50/60 Hz	Water pressure adjustment	3 levels
Max. power consumption	1725 – 2053 W	Warranty	2 years

Invacare Limited  
 Pencoed Technology Park - Pencoed  
 Bridgend CF35 5AQ - United Kingdom  
 Tel: +44 1 656 776 222  
 Fax: +44 1656 776 220  
 E-mail: [uk@invacare.com](mailto:uk@invacare.com)  
 Sales Order E-mail: [ordersuk@invacare.com](mailto:ordersuk@invacare.com)  
[www.invacare.co.uk](http://www.invacare.co.uk)

Invacare Ireland Ltd  
 Unit 5 - Seatown Business Campus  
 Seatown Road - SWORDS - County Dublin - Ireland  
 Tel. +353 1 8107084  
 Fax +353 1 8107085  
 Email: [Ireland@invacare.com](mailto:Ireland@invacare.com)  
[www.invacare.ie](http://www.invacare.ie)

