



*Superlight &
Simply Extreme*



With the name The KSL, küschall went back to its roots, combining the K, a symbol for küschall quality and design, with the letters "SL", for Super Light.

With its sleek and dynamic design and high quality materials, The KSL ensures the best driving performance, rigidity and stability combined with incredible lightness. The total weight is below 7kg.

The customization possibilities of The KSL allows a person the ability to express their own individual style and flair. The most active of lifestyles can be accommodated with The KSL.

TECHNICAL DATA



Seat width:
340 – 440 mm
in increments of
20 mm



Seat depth:
385 – 485 mm
in increments of
25 mm



Seat height:
front: 450 – 530 mm
rear: 400 – 480 mm
in increments of
10 mm



Knee-to-heel length:
350 – 500 mm
in increments of
10 mm



Backrest height:
270 – 420 mm
in increments of
15 mm



Centre of gravity:
80 mm/ 95 mm/ 110 mm/
125 mm/ 140 mm



Backrest angle:
74°/78°/82°/86°/90°



Max. safe slope:
7°



Reversing width:
920 - 1230 mm



Total width:
SW + 170 mm



Total length:
75°: approx. 800 mm
90°: approx. 730 mm



Max. user weight:
100 kg



Total weight
unloaded:
approx. 6.8 kg



Transport weight
w/out rear wheels:
approx. 4.3 kg

LEGEND

- CONFIGURATION:** This wheelchair has an aluminium tapered frame in the new frame colour 'atomic violet' combined with a CPS chrome, standard carbon camber bar and clothes-guards, aluminium foot-tube, mini push-handles, active brakes, 4" Starec front wheels and High Performance rear wheels with Schwalbe tyres and polished stainless steel handrims.
- BACKREST:** The KSL is available either with a fixed backrest welded to the angle of your choice or with a foldable backrest. The latter option offers easier transportation.
- FRAME:** An optional tapered frame is available on The KSL. With a reduced front width, this option allows better leg support and improves the overall responsiveness and manoeuvrability of the wheelchair.
- DESIGN:** The minimal frame of The KSL combined with its sleek design elements not only make it one of the lightest wheelchairs on the market, but also give it an elegant designer aesthetic.
- BACKREST UPHOLSTERY:** with a choice of four different backrest upholsteries (standard, light, tension Velcro® adjustable and Aero), The KSL offers a configuration to suit a range of postural requirements.
- CUSTOMIZE YOUR WHEELCHAIR:** For an extra personal touch, a name or words of choice can be engraved on the castor housing cap. A truly unique feature.

United Kingdom: Invacare Ltd. | Pencoed Technology Park, Pencoed, Bridgend CF35 5AQ | Phone +44 (0) 1656 776222 | Fax +44 (0) 1656 776220 | www.invacare.co.uk
Ireland: Invacare Ltd. | Unit 5, Seatown Business Campus | Seatown Road | Swords, Co. Dublin | Phone +353 (0)1 810 7084 | Fax +353 (0)1 810 7085 | www.invacare.ie

küschall is a registered trademark. Every effort has been made to ensure that the contents of this publication are fully up-to-date at the time of going to print. As part of its ongoing improvement of products, Küschall AG reserves the right to modify existing models at any time. Any use of this publication or of parts thereof as well as any reproduction of images must have the written consent of Küschall AG. Patents pending

The KSL_UK_2014_1526588

Copyright © 10/2014, Küschall AG, Switzerland - All rights reserved



2



3



4



5



6



THE KSL