

Effective  
infection  
control features  
now added



Invacare®

# Softform® Trolley Mattress

Pressure Redistributing Trolley Mattresses



The *Invacare Softform Trolley mattress range* has been designed to specifically meet the demands of the modern hospital environment. Suitable for those at 'High Risk/Very High Risk'<sup>1</sup> of developing a pressure ulcer, the high specification pressure redistributing foam mattress offers a comfortable support surface and interim pressure area care for patients waiting for treatment. The mattress is available in a wide range of sizes and options for use on a variety of trolley frames and is also suitable for X-ray and theatre use.

### Effective Pressure Ulcer Prevention

The foam core has independent geometrically-cut surface cells which helps to reduce shear and friction, increases the overall support surface area for the vulnerable patient and thereby reducing peak pressures.

### Multi stretch PU cover

The cover fabric provides excellent recovery and durability and helps to reduce shear and friction forces. The base is a tough coated polyurethane, which protects the mattress against potential trolley damage from trolley frames.

### Improved performance

To meet infection control policies, the cover is fully welded with a specially designed zip, which is welded in place and covered by a deep flap on all 4 sides.

## Features and Options



### Infection Control

A safety flap for concealment of the zip spans the perimeter of the mattress ensuring the interior of the trolley mattress is protected.

### Buckle and Straps

This free of charge optional extra prevents the mattress from slipping whilst on the trolley, ensuring patient safety.

### Waterproof zip

The trolley mattress now has a **new waterproof zip** allowing for easy inspection, to ensure the foam core is free of damage and potential contamination. The zip teeth interlock, making it impenetrable to fluid ingress and fully protects the internal foam of the trolley mattress.

### Customisation Options

The TM60, trolley mattress is available with, octagonal corners and/or velcro fastenings, making it suitable for a range of platforms and requirements.

**Also available with a Visco Elastic (Memory) foam insert if required.**

## Technical data

For more information about this product, including the product's user manual, or see [www.invacare.co.uk](http://www.invacare.co.uk)

	Width	Length	Height	Total product weight	Max. user weight	Max. cleaning temperature	Colour
<b>TM60 Trolley Mattress</b>	560 - 700 mm	1830 - 1930 mm	120 mm	7 kg	22 stone/140 kg	80°	Blue

Manufactured to European Standards: BS 5852: Crib 5 and BS 7175: Crib 5.

## Warranty

1 year

<sup>1</sup>/ Essential nursing care is pivotal in pressure ulcer prevention. This mattress will positively contribute to the outcome of a pressure ulcer prevention care plan. Education, clinical judgement and action based planning based on vulnerability are fundamental factors in the prevention of pressure ulcers. A range of assessment scales can be used as a formal method of assessing risk from pressure ulcer development, and should be used in conjunction with an informal assessment (informed nursing judgement). Informal assessment is considered to be of greater importance and clinical value. Risk assessment used in this case is the Waterlow scale

Invacare Limited  
Pencoed Technology Park - Pencoed  
Bridgend CF35 5AQ - United Kingdom  
Tel: +44 1656 776 222  
Fax: +44 1656 776 220  
E-mail: [uk@invacare.com](mailto:uk@invacare.com)  
Sales Order E-mail:  
[ordersuk@invacare.com](mailto:ordersuk@invacare.com)  
[www.invacare.co.uk](http://www.invacare.co.uk)

Invacare Ireland Ltd  
Unit 5 - Seatown Business Campus  
Seatown Road - SWORDS -  
County Dublin - Ireland  
Tel. +353 1 8107084  
Fax +353 1 8107085  
Email: [Ireland@invacare.com](mailto:Ireland@invacare.com)  
[www.invacare.ie](http://www.invacare.ie)

