

Invacare® | **Leg Trough**



Shaped leg support - provides comfort and promotes recovery

The *Invacare Leg Trough* has been designed in consultation with clinicians to provide maximum support for the whole leg. The **Leg Trough** elevates the heels and allows them to be suspended over the edge of the trough. This provides maximum pressure relief for the vulnerable heel area.





Features and Options

Provides comfort and support to promote recovery

The leg trough offers optimum support and protection for the whole leg, reducing pressure and promoting recovery.

Relieves pressure on vulnerable heels

Elevates the heels and allows them to be suspended over the edge of the trough, providing pressure relief for the vulnerable heel area.



Heavy duty water-resistant cover

A two-way stretch, vapour permeable and anti-static polyurethane (PU) water-resistant cover.



Cover includes a zip and is removable

The cover can be unzipped and removed from the foam core to allow laundering of the cover.



Available with or without straps

The adjustable side squeeze clip buckles provide additional positioning for the leg.

Part numbers for Leg Trough range with straps:

- Small - MLT/S
- Medium - MLT/M
- Large - MLT/L



Part number for Leg Trough range without straps:

- Small - 1575812
- Medium - 1575813
- Large - 1575829

Technical data

Leg Trough



Length:
Small: 510mm
Medium: 610mm
Large: 660mm



Width:
110 mm



Height:
115 mm



Weight:
0.8 kg



Cleaning: 80°

Manufactured to European Standards:
Interior Foam: BS 5852: 1982 Crib 5, BSEN597-1,2
Cover: BS7175: 1989 Crib 5

Invacare Limited
Pencoed Technology Park - Pencoed -
Bridgend CF35 5AQ - United Kingdom
Tel: +44 1 656 776 222
Fax: +44 1656 776 220
E-mail: uk@invacare.com
Sales Order E-mail:
ordersuk@invacare.com
www.invacare.co.uk

Invacare Ireland Ltd
Unit 5 - Seatown Business Campus
Seatown Road - SWORDS - County
Dublin - Ireland
Tel. +353 1 8107084
Fax +353 1 8107085
Email: Ireland@invacare.com
www.invacare.ie

© 2013 Invacare International Sàrl
All rights reserved.

All information quoted is believed to be correct at time of print.
Invacare reserves the right to alter product specifications without prior consultation.

Leg Trough - UK - 09/2014 -

Colour



Blue

Warranty 1 Year



Yes, you can.®