

Invacare<sup>®</sup>

# **Essential Basic**

Cost effective comfort mattress



The **Essential Basic** mattress offers the client a comfortable and supportive surface. This is an ideal entry level mattress for individuals who do not require therapy for pressure ulcer prevention.

- **Multi-stretch vapor-permeable, water resistant polyurethane fabric cover.**
- **Zippered around 3 sides for ease of inspection and cleaning.**
- **Meets fire safety standards.**

## Features and Options



### Entry Level



Ideal for clients who do not require therapy for pressure ulcer prevention.

### Cover

Multi-stretch polyurethane fabric cover.

## Technical data

For more information about this product, including the product's user manual, please visit [www.invacare.co.uk](http://www.invacare.co.uk).

	Size range	Max. user weight	Warranty
<b>Essential Basic</b>	 2000 x 880 x 150 mm 2050 x 880 x 150 mm	 140 kg (22 stone)	1 year

## Colour

Blue



Mattresses is manufactured to fire safety standard:  
BS EN 597-1 Cig, BS EN 597-2 Match and BS 7177: Crib Source 5

Invacare Limited  
Pencoed Technology Park - Pencoed  
Bridgend CF35 5AQ - United Kingdom  
Tel: +44 1656 776 222  
Fax: +44 1656 776 220  
E-mail: [uk@invacare.com](mailto:uk@invacare.com)  
Sales Order E-mail:  
[ordersuk@invacare.com](mailto:ordersuk@invacare.com)  
[www.invacare.co.uk](http://www.invacare.co.uk)

Invacare Ireland Ltd  
Unit 5 - Seatown Business Campus  
Seatown Road - SWORDS -  
County Dublin - Ireland  
Tel. +353 1 8107084  
Fax +353 1 8107085  
Email: [Ireland@invacare.com](mailto:Ireland@invacare.com)  
[www.invacare.ie](http://www.invacare.ie)

Trusted Solutions,  
Passionate People



© 2017 Invacare International Sàrl. All rights reserved. All information quoted is believed to be correct at time of print.  
Invacare reserves the right to alter product specifications without prior consultation. Essential Basic - UK - 10/2017 - 1629269