



Invacare TDX[®]2 Sprint

Are you ready to play ball?



Yes, you can.[®]

Beat the competition and make a tactical choice with Invacare's football powerchair **TDX[®] 2 Sprint**



TDX[®]2 Sprint

Take control of your game

Comfort

- ▶ Focus on your position with Invacare's highly adaptable **Modulite** seating system, featuring **Matrix**[®] positioning backs and cushions.

Control

- ▶ Take control of your game with the intuitive **DX2** control system.
- ▶ Strike harder with the **TDX2 Sprint** mid wheel drive football chair.

Confidence

- ▶ Interchangeable shrouds now featuring 7 modern colours, also allows for customisation to team colours.
- ▶ Quick and easy to adapt to suit any player's needs.
- ▶ Improve your game with better weight distribution and posture with Modulite seating.

Safety

- ▶ Optional H-Belt for enhanced safety when attacking the opposition.
- ▶ Rear bumper as standard.
- ▶ Optional front bumper.

Manoeuvrability

- ▶ Out manoeuvre the competition with the compact and stable wheel base of the **TDX2 Sprint** featuring stability lock and **SureStep**[®] technology.



Invacare® **Modulite™**

An evolution in seating designed around you



Invacare's latest product for powerchair football draws all its features from the **TDX® SP2**, such as Invacare's patented **Surestep®** suspension technology and stability lock. This ensures all 6 wheels stay firmly on the ground. Plus it features the highly adaptable **Modulite™** seating that can accommodate all shapes and sizes in one system to ensure you can give your best performance.



1. With Inrigged armrest.
2. 50mm armrest adjustment.
3. Ground to seat plate without cushion.
4. Measured from seat plate.
5. Depending on configuration.
6. Depending on seat position.
7. Width of wheel base.

Invacare Limited

Pencoed technology Park, Bridgend, CF35 5AQ
Tel: +44 1656 776 222 Fax: +44 1656 776 220
Email: uk@invacare.com

www.invacare.co.uk

Part No. 1586751 © 2015 Invacare International Sàrl. All rights reserved.



Follow us on Twitter:
[@Think_Mobility](https://twitter.com/Think_Mobility)



Follow us on YouTube:
Invacare Ltd - UK



Like us on Facebook:
Invacare UK



Yes, you can.®