

küschall[®]
designed for life



KÜSCHALL CHAMPION SK

KÜSCHALL CHAMPION SK

THE MOST FOLDABLE - SMALL BUT SMART

Sleek, smart styling meets compactness and practicality. The active rigid wheelchair that folds is now even smaller.

Having the benefit of an everyday wheelchair with such small folded dimensions makes light work of all aspects of travel. It's the perfect companion in your busy day to day commuting and essential for any holiday destination.

Evolved from our highly popular Champion, the SK features a folding front frame, the result is a wheelchair that is lightweight, easy to handle and convenient to transport.

Exceptional driving performance, rigidity and quality are also key features, ensuring as well as being practical you can be sure it delivers the strength and durability needed for an active lifestyle.



FEATURES AND BENEFITS

1. CONFIGURABLE: Smart, practical and with a wide range of options, accessories and colour choices the Champion SK can be configured to suit your exact needs, style and personality.

2. SMALL DIMENSIONS: Making it a practical choice.

Dimensions folded:

Height at rear 280mm

Height at front 330mm

Length of chair 670mm

Width of chair 220mm

Dimensions are based on a 400mm x 400mm seat width and depth

3. EASY FOLDING: Simply fold the main frame, then the backrest and finally tuck in the front frames for compact and easy transportation or travel.

4. ADJUSTABLE: As well as being small, practical and configurable, this wheelchair is also adjustable to ensure the perfect fit around you. The centre of gravity, backrest height and backrest angle are just some of the adjustments possible.

TECHNICAL DATA



Frame angle

75°/85°

Aluminium folding front frame



Seat width

360 – 460 mm
in increments of 20 mm



Seat depth

340 – 480 mm
in increments of 20 mm



Seat height

front: 450 – 540 mm
rear: 390 – 490 mm
stepless adjustable



Knee-to-heel length

380 – 500 mm
stepless adjustable



Backrest height

300 – 465 mm
in increments of 15 mm



Backrest angle

76°/80.5°/85°/89.5°/94°



Max. safe slope

7°



Reversing width

920 – 1230 mm



Total width

SW + 160 mm



Total folded dimensions

Rear 280mm
Front 330mm
Length 670mm
Width 220mm



Total length

75°: approx. 930 mm
85°: approx. 850 mm



Total weight unloaded
approx. 9.2 kg



Transport weight w/out rear wheels
approx. 6.7 kg



Max. user weight
120 kg

Küschall AG | Benkenstrasse 260 | CH-4108 Witterswil | Switzerland | Phone +41 (0) 61 487 70 70 | Fax +41 (0) 61 487 70 79 | www.kuschall.com

UNITED KINGDOM: Invacare Ltd. | Pencoed Technology Park, Pencoed, Bridgend CF35 5AQ | Phone +44 (0) 1656 776222 | Fax +44 (0) 1656 776220 | www.invacare.co.uk

IRELAND: Invacare Ltd. | Unit 5, Seatown Business Campus | Seatown Road | Swords, Co. Dublin | Phone +353 (0)1 810 7084 | Fax +353 (0)1 810 7085 | www.invacare.ie

küschall is a registered trademark. Every effort has been made to ensure that the contents of this publication are fully up-to-date at the time of going to print. As part of its ongoing improvement of products, Küschall AG reserves the right to modify existing models at any time. Any use of this publication or of parts thereof as well as any reproduction of images must have the written consent of Küschall AG. Patents pending

Champion SK_UK_2017

Copyright © 01/2017, Küschall AG, Switzerland - All rights reserved

1627814