

Invacare® | *Odstock Wedge*



Pressure redistributing wedge

The *Invacare Odstock Wedge* was developed in partnership with clinicians to assist in the reduction of interface pressures on the sacrum and heels without significantly increasing pressure on the thighs or calves. It is designed to improve and maintain patient posture in the recumbent position, and to reduce the tendency for the client to slip down the bed, thereby reducing shearing forces on the sacrum. The *Odstock Wedge* is ideal for improving posture for those with back discomfort.



Yes, you can.®



Features and Options

Reduces pressure on the heels

Clinical research has effectively demonstrated that the degree of knee flexion necessary to achieve substantial levels of pressure relief is relatively small. The *Invacare Odstock Wedge* was designed by clinicians to reduce pressure on the heels through re-positioning the hip and elevating the knee.

Reduces lumbar lordosis

By tilting the pelvis, the **Odstock Wedge** also reduces lumbar lordosis – desirable following back surgery and ideal for improving the sleeping posture of patients with back problems.



High frequency welded seams

Ensures the core is protected from fluid ingress, increasing product integrity.



High quality construction

The multi stretch cover, is waterproof and vapour-permeable. The foam core is constructed from, high density foam.



Designed to be used in conjunction with the *Invacare Heelpad*

Used in conjunction, these two products combine to significantly reduce pressure on the vulnerable heel area.

Technical data



Odstock Wedge Width: 57 cm



Length: 86.5 cm



Height: 11.5 cm



Weight: 1.5 kg



Cleaning: 80°

Refer to instructions for use and relevant infection control guidelines.

Manufactured to European Standards:
Interior Foam: BS 5852: Crib 5,
Cover: BS7175: Crib 5

Invacare Limited
Pencoed Technology Park - Pencoed -
Bridgend CF35 5AQ - United Kingdom
Tel: +44 1 656 776 222
Fax: +44 1656 776 220
E-mail: uk@invacare.com
Sales Order E-mail:
ordersuk@invacare.com
www.invacare.co.uk

Invacare Ireland Ltd
Unit 5 - Seatown Business Campus
Seatown Road - SWORDS - County
Dublin - Ireland
Tel. +353 1 8107084
Fax +353 1 8107085
Email: Ireland@invacare.com
www.invacare.ie

© 2014 Invacare International Sàrl
All rights reserved.

All information quoted is believed to be correct at time of print.
Invacare reserves the right to alter product specifications without prior consultation.

Odstock Wedge - UK - 09/2014 - 1505395

Colours



Fuchsia (Standard)



Blue (Option)

Warranty: 1 year



Yes, you can.®