

One carer  
operation\*



Invacare® | **Roze**

Versatile and safe - an ideal solution for  
Residential and Nursing care



The *Invacare Roze* stand assist offers an alternative to a passive lift when a degree of weight bearing ability remains. Ideal for rehabilitation, degenerative conditions and clients whose condition varies from day to day, the *Roze* can provide dignity and a greater degree of independence. With the latest control box technology, the *Roze* offers a 'Soft Start' and 'Soft Stop' function that provides a more comfortable transfer. Combining a range of adjustable features with a lightweight aluminum construction, the *Roze* is equally at home in Residential and Acute Care environments.

### A versatile stand assist

The *Invacare Roze* is a versatile stand assist lifter that supports the transfer of individuals who can only bear a minimal amount of their own weight. Ideal for transfers from a seated to seated or seated to standing position, the small footprint of the *Roze* and wide leg opening offers a flexible and comfortable transfer.

## Features and Options



### Height-adjustable leg support

The padded leg support can be adjusted for different needs.



### Removable foot plate

The foot plate is horizontal and can be removed for standing or walking a patient during rehabilitation.



### Control box with diagnostic system

This intelligent monitoring system includes a service reminder, emergency operation, audio battery alarm and battery indicator.

### Safety first, with an intelligent control box

The *Roze* is fitted with a control box that has a diagnostic system: an intelligent service lifetime monitoring system that highlights the battery status, indicates number of lifts, overloading information and service condition of the lifter to ensure safe patient handling.



## Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please visit [www.invacare.co.uk](http://www.invacare.co.uk).

	Lifting range	Inside width: legs opened	Outside width: legs closed	Overall length of base	Turning diameter	Height of legs	Number of cycles	Total product weight
<i>Invacare Roze</i>	940 -1610 mm	1050 mm	660 mm	1110 mm	1380 mm	120 mm	100 - 200 lifts per full charge	53 kg

	Max. user weight
<i>Invacare Roze</i>	200 kg

Frame colours	
Moderate white	Blue dusk

Warranty
5 years warranty: Frame, control box & actuator
2 years warranty: Hand control
1 year warranty: Battery

The *Invacare Roze* is part of the Safe Patient Handling range.  
\* A full risk assessment must be carried out before using equipment with one carer.



Invacare Limited  
Pencoed Technology Park - Pencoed  
Bridgend CF35 5AQ - United Kingdom  
Tel: +44 1 656 776 222  
Fax: +44 1656 776 220  
E-mail: [uk@invacare.com](mailto:uk@invacare.com)  
Sales Order E-mail: [ordersuk@invacare.com](mailto:ordersuk@invacare.com)  
[www.invacare.co.uk](http://www.invacare.co.uk)

Invacare Ireland Ltd  
Unit 5 - Seatown Business Campus  
Seatown Road - SWORDS - County Dublin - Ireland  
Tel. +353 1 8107084  
Fax +353 1 8107085  
Email: [Ireland@invacare.com](mailto:Ireland@invacare.com)  
[www.invacare.ie](http://www.invacare.ie)

