

One carer
operation*



Invacare® | **Roze**

Versatile and safe - an ideal solution for
Residential and Nursing care



The *Invacare Roze* stand assist offers an alternative to a passive lift when a degree of weight bearing ability remains. Ideal for rehabilitation, degenerative conditions and clients whose condition varies from day to day, the *Roze* can provide dignity and a greater degree of independence. With the latest control box technology, the *Roze* offers a 'Soft Start' and 'Soft Stop' function that provides a more comfortable transfer. Combining a range of adjustable features with a lightweight aluminum construction, the *Roze* is equally at home in Residential and Acute Care environments.

A versatile stand assist

The *Invacare Roze* is a versatile stand assist lifter that supports the transfer of individuals who can only bear a minimal amount of their own weight. Ideal for transfers from a seated to seated or seated to standing position, the small footprint of the *Roze* and wide leg opening offers a flexible and comfortable transfer.

Features and Options



Height-adjustable leg support

The padded leg support can be adjusted for different needs.



Removable foot plate

The foot plate is horizontal and can be removed for standing or walking a patient during rehabilitation.



Control box with diagnostic system

This intelligent monitoring system includes a service reminder, emergency operation, audio battery alarm and battery indicator.

Safety first, with an intelligent control box

The *Roze* is fitted with a control box that has a diagnostic system: an intelligent service lifetime monitoring system that highlights the battery status, indicates number of lifts, overloading information and service condition of the lifter to ensure safe patient handling.



Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please visit www.invacare.co.uk.

	Lifting range	Inside width: legs opened	Outside width: legs closed	Overall length of base	Turning diameter	Height of legs	Number of cycles	Total product weight
<i>Invacare Roze</i>	940 -1610 mm	1050 mm	660 mm	1110 mm	1380 mm	120 mm	100 - 200 lifts per full charge	53 kg

	Max. user weight
<i>Invacare Roze</i>	200 kg

Frame colours	
Moderate white	Blue dusk

Warranty
5 years warranty: Frame, control box & actuator
2 years warranty: Hand control
1 year warranty: Battery

The *Invacare Roze* is part of the Safe Patient Handling range.
* A full risk assessment must be carried out before using equipment with one carer.



Invacare Limited
Pencoed Technology Park - Pencoed
Bridgend CF35 5AQ - United Kingdom
Tel: +44 1 656 776 222
Fax: +44 1656 776 220
E-mail: uk@invacare.com
Sales Order E-mail: ordersuk@invacare.com
www.invacare.co.uk

Invacare Ireland Ltd
Unit 5 - Seatown Business Campus
Seatown Road - SWORDS - County Dublin - Ireland
Tel. +353 1 8107084
Fax +353 1 8107085
Email: Ireland@invacare.com
www.invacare.ie

