



Yes, you can:

Powerchairs Retail Prescription Form

Tel: 01656 776222 Fax: 01656 776220
 Email: ordersuk@invacare.com
 Online Spares available at: www.invacare.co.uk

JANUARY 2021
 Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

M100 AVIVA RX20 Modulite

Standard Features

- Max user weight 136kg (21 Stone)
- Crash tested ISO 7176 -Part 19
- Modulite highly adjustable seating system
- Telescopic adjustable seat widths from 330mm - 530mm (13" - 21") with an additional 50mm, (2") by arms
- Seat depths 410mm - 510mm stepless adjustment (16" - 20")
- Flex3 adjustable Backrest Plates with Laguna Backrest cushion
- Standard T-Armrests, detachable with short armpad
- Neckrest depth adjustable
- Seat lifter with powered tilt (300 mm, 30°) with CoG shift
- Choice of seat height's
- Manual Settable Recliner (in 3.75° steps) 90° - 120°
- Hanger ADE pair with angle adjustable footplate
- Remote LiNX DLX-REM211 with standard bracket
- Choice of 5 Shroud colours, 5 Rim insert colours
- 8" (200 mm) x 2.50-4 (65 mm) Castor, puncture proof solid tyre
- 14" (355 mm) x 3" (75 mm), Solid inlay black drive wheels
- CTC Suspension Control - Traction - Comfort
- 6 km/h 4 pole Durawatt motors
- 52Ah batteries
- 8A (200 - 250V) charger
- Anti tippers
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)

Standard Base price: £8,290



Please make your selection with a cross in the box ○ Free of charge ◆ standard feature POA special request

Customer weight = KG

Code **Description**

BASE MODEL*

JBC0000	<input type="checkbox"/>	AVIVA RX 20 Base Model Max user weight up to 136 kg	£8,290
---------	--------------------------	---	--------

CHASSIS WIDTH*

JBL2810	<input type="checkbox"/>	Chassis Narrow (560 mm/22")	◆
---------	--------------------------	-----------------------------	---

SEATING TYPE*			
JBC2100	<input type="checkbox"/>	Invacare® Modulite Seat Unit	◆
SEAT FRAME*			
JBC2111	<input type="checkbox"/>	Telescopic Seat Frame and Plate with Side Rails	◆
SEAT WIDTH*			
JBC2121	<input type="checkbox"/>	380 mm Seat Plate	○
JBC2122	<input type="checkbox"/>	430 mm Seat Plate	○
JBC2123	<input type="checkbox"/>	480 mm Seat Plate	○
JBC2124	<input type="checkbox"/>	530 mm Seat Plate (NC 2106)	○
Body Fit Kit			
JBC2106	<input type="checkbox"/>	Body fit Kit	○
SEAT DEPTH*			
JBC2132	<input type="checkbox"/>	Seat Depth 410 mm	○
JBC2134	<input type="checkbox"/>	Seat Depth 460 mm	○
JBC2136	<input type="checkbox"/>	Seat Depth 510 mm	○
SEAT CUSHION*			
JBC2142	<input type="checkbox"/>	Matrx Flo-tech Cushion Flo-tech Lite	◆
BACKREST SUPPORT*			
JBC2163	<input type="checkbox"/>	Backposts short (480 mm) (NC 2173)	◆
JBC2162	<input type="checkbox"/>	Backposts tall (540 mm)	○
PUSH HANDLES FOR BACKPOSTS			
JBC2108	<input type="checkbox"/>	Backposts without push handles	◆
BACKREST UPHOLSTERY / CUSHION*			
JBC2173	<input type="checkbox"/>	Flex3 adjustable Backrest Plates with Laguna Backrest cushion	◆
UPHOLSTERY TYPE*			
JBC2151	<input type="checkbox"/>	Black Fabric Upholstery	○
JBC2152	<input type="checkbox"/>	Black Black multi stretch PU coated fabric Upholstery	○
JBC2155	<input type="checkbox"/>	Black Spacer Fabric Upholstery	○
ARMREST (Pair)*			
JBC2191	<input type="checkbox"/>	T-Arm detachable, Standard	◆
ARMPAD (Pair)*			
JBC2195	<input type="checkbox"/>	Short Armpad (325 mm)	◆
HEADREST / NECKREST - Adequate headrest option must be selected*			
JBC1350	<input type="checkbox"/>	Neckrest depth adjustable	◆
SEAT LIFTER / SEAT TILT*			
JBC1094	<input type="checkbox"/>	Seat lifter with powered tilt (300 mm, 30°) with CoG shift	◆
SEAT TO GROUND HEIGHT*			
JBC2857	<input type="checkbox"/>	Seat-to-Floor Height 435 mm	○
JBC2858	<input type="checkbox"/>	Seat-to-Floor Height 460 mm	○
JBC2859	<input type="checkbox"/>	Seat-to-Floor Height 485 mm	○
BACKREST RECLINER*			
JBC0977	<input type="checkbox"/>	Powered Recliner 30°	◆
FOOTREST HANGER TYPE*			
JBC0531	<input type="checkbox"/>	Hanger ADE pair Lower leg length: 290 - 460 mm Powered Angle adjustment: 80 - 0° with calf pad	◆

FOOTPLATE TYPE*			
JBC0504	<input type="checkbox"/>	Angle Adjustable Footplate pair	◆
CONTROLS*			
JBC2203	<input type="checkbox"/>	Remote LiNX DLX-REM211 = Drive & Actuators 3 Drive functions and power module incl. 2 Actuators	◆
CONTROL BRACKETRY*			
JBC1207	<input type="checkbox"/>	Standard Bracket RH	○
JBC1208	<input type="checkbox"/>	Standard Bracket LH	○
SHROUD COLOUR*			
JBC2871	<input type="checkbox"/>	Mid range - Black/Anthracite Shroud	○
JBC2872	<input type="checkbox"/>	Mid range - Silver Shroud	○
JBC2873	<input type="checkbox"/>	Mid range - Bengal Red Shroud	○
JBC2874	<input type="checkbox"/>	Mid range - Mercury Blue Shroud	○
JBC2875	<input type="checkbox"/>	Mid range - Midnight Purple Shroud	○
JBX0299	<input type="checkbox"/>	Premium - Individual Shroud Colour (includes coloured spreader bars on Modulite back) RAL number required Extended lead lime: + 15 days	£195
RIM COLOUR*			
JBC2882	<input type="checkbox"/>	Mid range - Black Rim Inserts	○
JBC2883	<input type="checkbox"/>	Mid range - Silver Rim Inserts	○
JBC2884	<input type="checkbox"/>	Mid range - Bengal Red Rim Inserts	○
JBC2885	<input type="checkbox"/>	Mid range - Mercury Blue Rim Inserts	○
JBC2886	<input type="checkbox"/>	Mid range - Midnight Purple Rim Inserts	○
FORKS*			
JBC2632	<input type="checkbox"/>	Double sided fork Without castor rim insert (enlarges wide base to 656 mm)	◆
FRONT WHEEL *			
JBC0828	<input type="checkbox"/>	8" (200 mm) x 2.50-4 (65 mm) Castor, puncture proof solid tyre	◆
DRIVE WHEEL*			
JBC0705	<input type="checkbox"/>	DW 14" (355 mm) x 3" (75 mm), Solid inlay Dual density foam insert	◆
C.T.C SUSPENSION*			
JBC2870	<input type="checkbox"/>	CTC Suspension Control - Traction - Comfort	◆
SUSPENSION SETTING*			
JBC1910	<input type="checkbox"/>	Soft (suggested for up to 70 kg User Weight)	(NC 1094) ◆
JBC1912	<input type="checkbox"/>	Medium (suggested for up to 90 kg User Weight)	(NC 1094) ○
JBC1911	<input type="checkbox"/>	Firm (suggested for up to 110 kg User Weight)	○
JBC1915	<input type="checkbox"/>	Very firm (suggested for up to 136 kg User Weight)	○
MOTOR TYPE*			
JBC1282	<input type="checkbox"/>	6 km/h	◆
BATTERY TYPE*			
JBC1259	<input type="checkbox"/>	60 Ah Gel	◆
CHARGER*			
JBC1343	<input type="checkbox"/>	Opti Charge Charger 8A (IP44)	◆
ANTI-TIPPERS*			
JBC1205	<input type="checkbox"/>	Anti-tipper	◆

TRANSPORTATION KIT*

JBC1401

6-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19)) required for variants with combined Lift & Tilt module Adequate headrest is compulsory

**LAP BELT***

JBC2113

Lap belt with Auto Lock Flexible mount in side rails

**TOTAL PRICE**

Non Conforming Options:

JBC2106	Body fit Kit
JBC2173	Flex3 adjustable Backrest Plates with Laguna Backrest cushion
JBC1094	Seat lifter with powered tilt (300 mm, 30°) with CoG shift



For a more complex customer requirement please refer to our Invacare Unique special order form. email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: _____ Email: _____

Phone: _____ Fax: _____

Description of special request (please include relevant option code or part number):