



Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

JANUARY 2021
Customer Ref:

Customer Ref* Account No*
Contact No* Order No*
Delivery Address* Order Qty*
Order date*

M101 AVIVA RX40 Std seating

Standard Features

- Max user weight 136kg
Crash tested ISO 7176 -Part 19
One piece seat plate width 380mm (15") or 430mm (17") with standard foam cushion (50mm, 2")
seat depths 410mm - 510mm (16" - 20")
Tension adjustable backrest
Standard T-Armrests, detachable with short armpad
Manually settable Seat Tilt 0 - 15°
Manual settable backrest Recline (0° - 30°, in 3.75° steps)
80° hangers with fixed footplate
LiNX Control REM216 with standard Bracket
Choice of 5 Shroud colours, 5 Rim insert colours
9" (225 mm) x 2.80 (70 mm) Castor, puncture proof solid tyre
DW 14" (355 mm) x 3" (75 mm), Solid inlay Dual density foam insert
CTC Suspension Control - Traction - Comfort
12 km/h Ultra High Speed Package
LED Light system with Indicators, 30 Lux headlights
73,5 Ah Gel Batteries
8A (200 - 250V) charger
Anti tippers
4 Point tie down system
Lapbelt with autolock (adjustable both sides)

Standard Base price: £6,756



Please make your selection with a cross in the box ○ Free of charge ◆ standard feature POA special request

Customer weight = KG

Code Description

BASE MODEL*

Table with 2 columns: Code, Description, Price. Row: JBF0001, AVIVA RX 40 Base Model Max user weight up to 136 kg includes Performance Package: 120A Power Module + minimum 60 Ah Batteries, £6,756

CHASSIS WIDTH*

Table with 2 columns: Code, Description, Price. Row: JBL2810, Chassis Narrow (560 mm/22"), ◆

SEATING TYPE*

Table with 2 columns: Code, Description, Price. Row: JBF2105, Modulite Standard Seat Standard Seat Plate and Backrest Frame in fixed Width, ◆

SEAT FRAME*			
JBF2110	<input type="checkbox"/>	One Piece Seat Plate	◆
SEAT WIDTH*			
JBF2121	<input type="checkbox"/>	380 mm Seat Plate Can be reduced to 330 mm in combination with inrigged armrests	○
JBF2122	<input type="checkbox"/>	430 mm Seat Plate Can be reduced to 380 mm in combination with inrigged armrests	◆
SEAT DEPTH*			
JBF2132	<input type="checkbox"/>	Seat Depth 410 mm	○
JBF2134	<input type="checkbox"/>	Seat Depth 460 mm	○
JBF2136	<input type="checkbox"/>	Seat Depth 510 mm	○
SEAT CUSHION*			
JBF2141	<input type="checkbox"/>	Standard Foam Cushion 50 mm height	◆
JBF2142	<input type="checkbox"/>	Matrx Flo-tech Cushion Flo-tech Lite	◆
BACKREST SUPPORT*			
JBF2161	<input type="checkbox"/>	Backrest frame, detachable, height 480 mm For Tension Adjustable Upholstery/Matrx Backs	◆
BACKREST UPHOLSTERY / CUSHION*			
JBF2171	<input type="checkbox"/>	Tension Adjustable Upholstery nylon overlay short = 480 mm / 7 straps measured from seat plate to top of upholstery 480 mm	◆
UPHOLSTERY TYPE*			
JBF2151	<input type="checkbox"/>	Black Fabric Upholstery (Nylon w 2141; 2-way-stretch material w 2142, 2143)	◆
ARMREST (Pair)*			
JBF2191	<input type="checkbox"/>	T-Arm detachable, Standard Adjustment range from seat plate: 245 - 315 mm	◆
ARMPAD (Pair)*			
JBF2195	<input type="checkbox"/>	Short Armpad (325 mm)	◆
HEADREST / NECKREST - Adequate headrest option must be selected*			
JBF1495	<input type="checkbox"/>	Non Invacare Headrest - no mounting hardware With Backbrace/Pushbar dia. 22mm in conjunction with telescopic seat frame/backposts	◆
SEAT LIFTER / SEAT TILT*			
JBF1091	<input type="checkbox"/>	Manually settable Seat Tilt 0 - 15°	◆
SEAT T GROUND HEIGHT*			
JBF2857	<input type="checkbox"/>	Seat-to-Floor Height 435 mm	○
JBF2858	<input type="checkbox"/>	Seat-to-Floor Height 460 mm	○
JBF2859	<input type="checkbox"/>	Seat-to-Floor Height 485 mm	○
BACKREST RECLINER*			
JBF0975	<input type="checkbox"/>	Manual Settable Recliner (in 3.75° steps) 90° - 120°	◆
FOOTREST HANGER TYPE*			
JBF0526	<input type="checkbox"/>	Standard hanger 80° pair Lower leg length: 350 - 450 mm	◆
FOOTPLATE TYPE*			
JBF0512	<input type="checkbox"/>	Fixed Footplate pair	◆
CONTROLS*			
JBF2204	<input type="checkbox"/>	Remote LiNX DLX-REM216 = Drive, Actuators & Lights 3 Drive functions and power module incl. 2 Actuators & Lights	◆
CONTROL BRACKETRY*			
JBF1207	<input type="checkbox"/>	Standard Bracket RH	○
JBF1208	<input type="checkbox"/>	Standard Bracket LH	○

SHROUD COLOUR*			
JBF2871		Mid range - Black/Anthracite Shroud	○
JBF2872		Mid range - Silver Shroud	○
JBF2873		Mid range - Bengal Red Shroud	○
JBF2874		Mid range - Mercury Blue Shroud	○
JBF2875		Mid range - Midnight Purple Shroud	○
JBX0299		Premium - Individual Shroud Colour RAL number required Extended lead lime: + 15 days	(NC 2883, 2884, 2885, 2886, 2888, 2889, 2890, 2891) £195
JBF2865		Premium - Racing Green Shroud	(NC 2883, 2884, 2885, 2886, 2887, 2889, 2890, 2891) £95
JBF2866		Premium - Deep Ocean Blue Satin Shroud	(NC 2883, 2884, 2885, 2886, 2887, 2888, 2890, 2891) £95
JBF2867		Premium - Frosty Blue Shroud	(NC 2883, 2884, 2885, 2886, 2887, 2888, 2889, 2891) £95
JBF2868		Premium - Gold Shroud	(NC 2883, 2884, 2885, 2886, 2887, 2888, 2889, 2890) £95
RIM COLOUR*			
JBF2882		Mid range - Black Rim Inserts	○
JBF2883		Mid range - Silver Rim Inserts	○
JBF2884		Mid range - Bengal Red Rim Inserts	○
JBF2885		Mid range - Mercury Blue Rim Inserts	○
JBF2886		Mid range - Midnight Purple Rim Inserts	○
JBX2887		Premium - Individual Rim Insert Colour RAL number required Extended lead lime: + 15 days	(NC 2872, 2873, 2874, 2875, 2865, 2866, 2867, 2868) ○
JBF2888		Premium - Racing Green Rim Inserts Price Included with Premium Shroud Kit	(NC 2872, 2873, 2874, 2875, 0299, 2866, 2867, 2868) ○
JBF2889		Premium - Deep Ocean Blue Satin Rim Inserts Price Included with Premium Shroud Kit	(NC 2872, 2873, 2874, 2875, 0299, 2865, 2867, 2868) ○
JBF2890		Premium - Frosty Blue Rim Inserts Price Included with Premium Shroud Kit	(NC 2872, 2873, 2874, 2875, 0299, 2865, 2866, 2868) ○
JBF2891		Premium - Gold Rim Inserts RAL1035 Price Included with Premium Shroud Kit	(NC 2872, 2873, 2874, 2875, 0299, 2865, 2866, 2867) ○
FORKS*			
JBF2632		Double sided fork Without castor rim insert (enlarges wide base to 656 mm)	◆
FRONT WHEEL*			
JBF0836		9" (225 mm) x 2.80 (70 mm) Castor, puncture proof solid tyre	◆
DRIVE WHEEL*			
JBF0705		DW 14" (355 mm) x 3" (75 mm), Solid inlay Dual density foam insert	◆
C.T.C SUSPENSION*			
JBF2870		CTC Suspension Control - Traction - Comfort	◆
SUSPENSION SETTING*			
JBF1910		Soft (suggested for up to 70 kg User Weight)	◆
JBF1912		Medium (suggested for up to 90 kg User Weight)	○
JBF1911		Firm (suggested for up to 110 kg User Weight)	○
JBF1915		Very firm (suggested for up to 136 kg User Weight)	○
MOTOR TYPE*			
JBF1909		12 km/h Ultra High Speed Package Always combined with LED Light system with Indicators, 30 Lux headlights	◆
BATTERY TYPE*			
JBF1260		73,5 Ah Gel	◆
CHARGER*			
JBF1343		Opti Charge Charger 8A (IP44)	◆

ANTI-TIPPERS*

JBF1205	<input type="checkbox"/>	Anti-tipper	◆
---------	--------------------------	-------------	---

TRANSPORTATION KIT*

JBF1494	<input type="checkbox"/>	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19) Adequate headrest is compulsory	◆
---------	--------------------------	--	---

LAP BELT*

JBF1329	<input type="checkbox"/>	Lap belt with Auto Lock Standard mount on backrest	◆
---------	--------------------------	--	---

REAR VIEW MIRROR*

JBF1415	<input type="checkbox"/>	Rear view mirror RH	○
---------	--------------------------	---------------------	---

JBF1339	<input type="checkbox"/>	Rear view mirror LH	○
---------	--------------------------	---------------------	---

TOTAL PRICE	
--------------------	--

Non Conforming Options:

JBF2872	Mid range - Silver Shroud
JBF2873	Mid range - Bengal Red Shroud
JBF2874	Mid range - Mercury Blue Shroud
JBF2875	Mid range - Midnight Purple Shroud
JBX0299	Premium - Individual Shroud Colour RAL number required Extended lead lime: + 15 days
JBF2865	Premium - Racing Green Shroud
JBF2866	Premium - Deep Ocean Blue Satin Shroud
JBF2867	Premium - Frosty Blue Shroud
JBF2868	Premium - Gold Shroud
JBF2883	Mid range - Silver Rim Inserts
JBF2884	Mid range - Bengal Red Rim Inserts
JBF2885	Mid range - Mercury Blue Rim Inserts
JBF2886	Mid range - Midnight Purple Rim Inserts
JBX2887	Premium - Individual Rim Insert Colour RAL number required Extended lead lime: + 15 days
JBF2888	Premium - Racing Green Rim Inserts Price Included with Premium Shroud Kit
JBF2889	Premium - Deep Ocean Blue Satin Rim Inserts Price Included with Premium Shroud Kit
JBF2890	Premium - Frosty Blue Rim Inserts Price Included with Premium Shroud Kit
JBF2891	Premium - Gold Rim Inserts RAL1035 Price Included with Premium Shroud Kit



For a more complex customer requirement please refer to our Invacare Unique special order form. email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: _____ Email: _____

Phone: _____ Fax: _____

Description of special request (please include relevant option code or part number):