



Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

JANUARY 2021
Customer Ref:

| | |
|--------------------------|--------------------|
| Customer Ref* | Account No* |
| Contact No* | Order No* |
| Delivery Address* | Order Qty* |
| | Order date* |

M113 AVIVA RX40 HD Modulite

Standard Features

- Max user weight 160kg (25 Stone)
- Crash tested ISO 7176 -Part 19
- Modulite adjustable seating system
- Telescopic adjustable seat widths from 330mm - 530mm (13" - 21") with an additional 50mm, (2") by arms
- Seat depths 410mm - 510mm stepless adjustment (16" - 20")
- Flex3 Tension Adjustable Upholstery
- Neckrest depth adjustable
- Standard T-Armrests, detachable with short armpad
- Powered "fixed pivot" Tilt 20°
- Manual settable backrest Recline (0° - 30° in 3.75° steps)
- 80° hangers with fixed footplate
- LiNX control REM211 with standard bracket
- Choice of 5 Shroud colours, 5 Rim insert colours
- 9" (225 mm) x 2.80 (70 mm) Castor, puncture proof solid tyre
- 14" (355 mm) x 3" (75 mm), Solid inlay black drive wheels
- CTC Suspension Control - Traction - Comfort
- 6 km/h (4mph) 4 pole Durawatt motors and 60Ah batteries
- 8A (200 - 250V) charger
- Anti tippers
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)

Standard Base price: £6,481



Please make your selection with a cross in the box Free of charge standard feature POA special request

Customer weight = KG

| Code | Description | |
|-----------------------|--|-------------------------------------|
| BASE MODEL* | | |
| JBH0002 | AVIVA RX 40 HD Base Model Max user weight up to 160 kg HD Version, includes RX 40 Performance Package Wide Base only | £6,481 |
| CHASSIS WIDTH* | | |
| JBH2830 | Chassis Wide (635 mm/25") | <input checked="" type="checkbox"/> |

| SEATING TYPE* | | | |
|--|--------------------------|---|---|
| JBH2100 | <input type="checkbox"/> | Invacare® Modulite Seat Unit | ◆ |
| SEAT FRAME* | | | |
| JBH2111 | <input type="checkbox"/> | Telescopic Seat Frame and Plate with Side Rails width 380 to 530 mm Additional 50 mm by arms | ○ |
| JBH2112 | <input type="checkbox"/> | Telescopic Seat Frame + Sling Seat with Side Rails width 380 to 530 mm; + 30 mm seat height Additional 50 mm by arms | ○ |
| SEAT WIDTH* | | | |
| JBH2121 | <input type="checkbox"/> | 380 mm Seat Plate Can be reduced to 330 mm in combination with inrigged armrests | ○ |
| JBH2122 | <input type="checkbox"/> | 430 mm Seat Plate Can be reduced to 380 mm in combination with inrigged armrests | ○ |
| JBH2123 | <input type="checkbox"/> | 480 mm Seat Plate | ○ |
| JBH2124 | <input type="checkbox"/> | 530 mm Seat Plate (NC 2106) | ○ |
| Body Fit Kit | | | |
| JBH2106 | <input type="checkbox"/> | Body fit Kit Outrigging Backrest Brackets and Pivot Point (outrig: 1 inch/2,54 cm per side) | ○ |
| SEAT DEPTH* | | | |
| JBH2132 | <input type="checkbox"/> | Seat Depth 410 mm | ○ |
| JBH2134 | <input type="checkbox"/> | Seat Depth 460 mm | ○ |
| JBH2136 | <input type="checkbox"/> | Seat Depth 510 mm | ○ |
| SEAT CUSHION* | | | |
| JBH2142 | <input type="checkbox"/> | Matrx Flo-tech Cushion Flo-tech Lite | ◆ |
| BACKREST SUPPORT* | | | |
| JBH2162 | <input type="checkbox"/> | Backposts tall (540 mm) For Flex3 backrests/Tension Adjustable Upholstery/Matrx Backs, measured from seat plate to top of upholstery Width adjustable cast bracket | ◆ |
| BACKREST UPHOLSTERY / CUSHION* | | | |
| JBH2179 | <input type="checkbox"/> | Flex3 Tension Adjustable Upholstery height 540 mm (seat plate to top of upholstery) | ◆ |
| UPHOLSTERY TYPE* | | | |
| JBH2152 | <input type="checkbox"/> | Black multi stretch PU coated fabric Upholstery | ○ |
| ARMREST (Pair)* | | | |
| JBH2191 | <input type="checkbox"/> | T-Arm detachable, Standard Adjustment range from seat plate: 245 - 315 mm | ◆ |
| ARMPAD (Pair)* | | | |
| JBH2195 | <input type="checkbox"/> | Short Armpad (325 mm) | ◆ |
| HEADREST / NECKREST - Adequate headrest option must be selected* | | | |
| JBH1350 | <input type="checkbox"/> | Neckrest depth adjustable | ◆ |
| SEAT LIFTER / SEAT TILT* | | | |
| JBH1092 | <input type="checkbox"/> | Powered "fixed pivot" Tilt 20° | ◆ |
| SEAT TO GROUND HEIGHT* | | | |
| JBH2857 | <input type="checkbox"/> | Seat-to-Floor Height 435 mm | ○ |
| JBH2858 | <input type="checkbox"/> | Seat-to-Floor Height 460 mm | ○ |
| JBH2859 | <input type="checkbox"/> | Seat-to-Floor Height 485 mm | ○ |
| BACKREST RECLINER* | | | |
| JBH0975 | <input type="checkbox"/> | Manual Settable Recliner (in 3.75° steps) 90° - 120° | ◆ |
| FOOTREST HANGER TYPE* | | | |
| JBH0526 | <input type="checkbox"/> | Standard hanger 80° pair Lower leg length: 350 - 450 mm | ◆ |

| FOOTPLATE TYPE* | | | |
|---------------------|--------------------------|--|------|
| JBH0512 | <input type="checkbox"/> | Fixed Footplate pair | ◆ |
| CONTROLS* | | | |
| JBH2203 | <input type="checkbox"/> | Remote LiNX DLX-REM211 = Drive & Actuators 3 Drive functions and power module incl. 2 Actuators | ◆ |
| CONTROL BRACKETRY* | | | |
| JBH1207 | <input type="checkbox"/> | Standard Bracket RH | ○ |
| JBH1208 | <input type="checkbox"/> | Standard Bracket LH | ○ |
| SHROUD COLOUR* | | | |
| JBH2871 | <input type="checkbox"/> | Mid range - Black/Anthracite Shroud | ○ |
| JBH2872 | <input type="checkbox"/> | Mid range - Silver Shroud | ○ |
| JBH2873 | <input type="checkbox"/> | Mid range - Bengal Red Shroud | ○ |
| JBH2874 | <input type="checkbox"/> | Mid range - Mercury Blue Shroud | ○ |
| JBH2875 | <input type="checkbox"/> | Mid range - Midnight Purple Shroud | ○ |
| JBX0299 | <input type="checkbox"/> | Premium - Individual Shroud Colour (includes coloured spreader bars on Modulite back) RAL number required Extended lead lime: + 15 days | £195 |
| JBH2865 | <input type="checkbox"/> | Premium - Racing Green Shroud | £95 |
| JBH2866 | <input type="checkbox"/> | Premium - Deep Ocean Blue Satin Shroud | £95 |
| JBH2867 | <input type="checkbox"/> | Premium - Frosty Blue Shroud | £95 |
| JBH2868 | <input type="checkbox"/> | Premium - Ultimate Red Shroud | £95 |
| JBH2869 | <input type="checkbox"/> | Premium - Falcon White Shroud | £95 |
| RIM COLOUR* | | | |
| JBH2882 | <input type="checkbox"/> | Mid range - Black Rim Inserts | ○ |
| JBH2883 | <input type="checkbox"/> | Mid range - Silver Rim Inserts | ○ |
| JBH2884 | <input type="checkbox"/> | Mid range - Bengal Red Rim Inserts | ○ |
| JBH2885 | <input type="checkbox"/> | Mid range - Mercury Blue Rim Inserts | ○ |
| JBH2886 | <input type="checkbox"/> | Mid range - Midnight Purple Rim Inserts | ○ |
| JBX2887 | <input type="checkbox"/> | Premium - Individual Rim Insert Colour RAL number required Extended lead lime: + 15 days | ○ |
| JBH2888 | <input type="checkbox"/> | Premium - Racing Green Rim Inserts | ○ |
| JBH2889 | <input type="checkbox"/> | Premium - Deep Ocean Blue Satin Rim Inserts | ○ |
| JBH2890 | <input type="checkbox"/> | Premium - Frosty Blue Rim Inserts | ○ |
| JBH2891 | <input type="checkbox"/> | Premium - Ultimate Red Rim Inserts | ○ |
| JBH2892 | <input type="checkbox"/> | Premium - Falcon White Rim Inserts | ○ |
| FORKS* | | | |
| JBH2632 | <input type="checkbox"/> | Double sided fork Without castor rim insert (enlarges wide base to 656 mm) | ◆ |
| FRONT WHEEL * | | | |
| JBH0836 | <input type="checkbox"/> | 9" (225 mm) x 2.80 (70 mm) Castor, puncture proof solid tyre | ◆ |
| DRIVE WHEEL * | | | |
| JBH0705 | <input type="checkbox"/> | DW 14" (355 mm) x 3" (75 mm), Solid inlay Dual density foam insert | ◆ |
| C.T.C SUSPENSION* | | | |
| JBH2870 | <input type="checkbox"/> | CTC Suspension Control - Traction - Comfort | ◆ |
| SUSPENSION SETTING* | | | |
| JBH1916 | <input type="checkbox"/> | Suspension up to 160 kg User Weight For RX 40 HD only | ◆ |

| | | |
|----------------------------|--------------------------|---|
| MOTOR TYPE* | | |
| JBH1282 | <input type="checkbox"/> | 6 km/h |
| BATTERY TYPE* | | |
| JBH1259 | <input type="checkbox"/> | 60 Ah Gel |
| CHARGER* | | |
| JBH1343 | <input type="checkbox"/> | Opti Charge Charger 8A (IP44) |
| ANTI-TIPPERS* | | |
| JBH1205 | <input type="checkbox"/> | Anti-tipper |
| LIGHTS | | |
| JBH1254 | <input type="checkbox"/> | Without Light system |
| TRANSPORTATION KIT* | | |
| JBH1494 | <input type="checkbox"/> | 4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19) Adequate headrest is compulsory, Not compatible with Seat lifter |
| LAP BELT* | | |
| JBH2113 | <input type="checkbox"/> | Lap belt with Auto Lock Flexible mount in side rails |

| | |
|--------------------|----------------------|
| TOTAL PRICE | <input type="text"/> |
|--------------------|----------------------|

Non Conforming Options:

JBH2106 Body fit Kit Outrigging Backrest Brackets and Pivot Point (outrig: 1 inch/2,54 cm per side)



For a more complex customer requirement please refer to our Invacare Unique special order form. email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: _____ Email: _____

Phone: _____ Fax: _____

Description of special request (please include relevant option code or part number):