



**Yes, you can:**

# Powerchairs Retail Prescription Form

Tel: 01656 776222 Fax: 01656 776220  
 Email: ordersuk@invacare.com  
 Online Spares available at: www.invacare.co.uk

JANUARY 2021  
 Customer Ref:

|                          |                    |
|--------------------------|--------------------|
| <b>Customer Ref*</b>     | <b>Account No*</b> |
| <b>Contact No*</b>       | <b>Order No*</b>   |
| <b>Delivery Address*</b> | <b>Order Qty*</b>  |
|                          | <b>Order date*</b> |

## M116 AVIVA RX40 HD Modulite

### Standard Features

- Max user weight 160kg (25 Stone)
- Crash tested ISO 7176 -Part 19
- Modulite adjustable seating system
- Telescopic adjustable seat widths from 330mm - 530mm (13" - 21") with an additional 50mm, (2") by arms
- Seat depths 410mm - 510mm stepless adjustment (16" - 20")
- Flex3 adjustable Backrest Plates with Laguna Backrest cushion
- Neckrest depth adjustable
- Standard T-Armrests, detachable with short armpad
- Seat lifter with powered tilt (300 mm, 30°) with CoG shift
- Manual settable backrest Recline (0° - 30° in 3.75° steps)
- 80° hangers with fixed footplate
- LiNX control REM216 with standard bracket
- Choice of 5 Shroud colours, 5 Rim insert colours
- 9" (225 mm) x 2.80 (70 mm) Castor, puncture proof solid tyre
- 14" (355 mm) x 3" (75 mm), Solid inlay black drive wheels
- CTC Suspension Control - Traction - Comfort
- 10 km/h Package combined with LED Light system with Indicators, 30 Lux headlights
- 73Ah batteries
- 8A (200 - 250V) charger
- Anti tippers
- 6 Point tie down system
- Lapbelt with autolock (adjustable both sides)

**Standard Base price: £8,825**



Please make your selection with a cross in the box     Free of charge     standard feature    POA special request

Customer weight =                      KG

**Code**                      **Description**

**BASE MODEL\***

|         |                          |  |        |
|---------|--------------------------|--|--------|
| JBH0002 | <input type="checkbox"/> | AVIVA RX 40 HD Base Model Max user weight up to 160 kg<br>HD Version, includes RX 40 Performance Package<br>Wide Base only | £8,825 |
|---------|--------------------------|--|--------|

|   |                          |  |   |
|---|--------------------------|--|---|
| <b>CHASSIS WIDTH*</b>   |                          |  |   |
| JBH2830   | <input type="checkbox"/> | Chassis Wide (635 mm/25")  | ◆ |
| <b>SEATING TYPE*</b>  |                          |  |   |
| JBH2100   | <input type="checkbox"/> | Invacare® Modulite Seat Unit   | ◆ |
| <b>SEAT FRAME*</b>  |                          |  |   |
| JBH2111   | <input type="checkbox"/> | Telescopic Seat Frame and Plate with Side Rails width 380 to 530 mm Additional 50 mm by arms   | ○ |
| JBH2112   | <input type="checkbox"/> | Telescopic Seat Frame + Sling Seat with Side Rails width 380 to 530 mm; + 30 mm seat height Additional 50 mm by arms   | ○ |
| <b>SEAT WIDTH*</b>  |                          |  |   |
| JBH2121   | <input type="checkbox"/> | 380 mm Seat Plate Can be reduced to 330 mm in combination with inrigged armrests   | ○ |
| JBH2122   | <input type="checkbox"/> | 430 mm Seat Plate Can be reduced to 380 mm in combination with inrigged armrests   | ○ |
| JBH2123   | <input type="checkbox"/> | 480 mm Seat Plate  | ○ |
| JBH2124   | <input type="checkbox"/> | 530 mm Seat Plate (NC 2106)  | ○ |
| <b>Body Fit Kit</b>   |                          |  |   |
| JBH2106   | <input type="checkbox"/> | Body fit Kit Outrigging Backrest Brackets and Pivot Point (outrig: 1 inch/2,54 cm per side)  | ○ |
| <b>SEAT DEPTH*</b>  |                          |  |   |
| JBH2132   | <input type="checkbox"/> | Seat Depth 410 mm  | ○ |
| JBH2134   | <input type="checkbox"/> | Seat Depth 460 mm  | ○ |
| JBH2136   | <input type="checkbox"/> | Seat Depth 510 mm  | ○ |
| <b>SEAT CUSHION*</b>  |                          |  |   |
| JBH2142   | <input type="checkbox"/> | Matrx Flo-tech Cushion Flo-tech Lite   | ◆ |
| <b>BACKREST SUPPORT*</b>  |                          |  |   |
| JBH2162   | <input type="checkbox"/> | Backposts tall (540 mm) For Flex3 backrests/Tension Adjustable Upholstery/Matrx Backs, measured from seat plate to top of upholstery Width adjustable cast bracket | ◆ |
| <b>BACKREST UPHOLSTERY / CUSHION*</b>                                   |                          |  |   |
| JBH2173   | <input type="checkbox"/> | Flex3 adjustable Backrest Plates with Laguna Backrest cushion Flex3 Backrest plate width and height adjustable 560 mm - 790 mm                                     | ◆ |
| <b>UPHOLSTERY TYPE*</b>   |                          |  |   |
| JBH2151   | <input type="checkbox"/> | Black Fabric Upholstery (Nylon with option 2141; 2-way-stretch material with options 2142, 2143)   | ◆ |
| JBH2152   | <input type="checkbox"/> | Black multi stretch PU coated fabric Upholstery  | ○ |
| JBH2155   | <input type="checkbox"/> | Black Spacer Fabric Upholstery   | ○ |
| <b>ARMREST (Pair)*</b>  |                          |  |   |
| JBH2191   | <input type="checkbox"/> | T-Arm detachable, Standard Adjustment range from seat plate: 245 - 315 mm  | ◆ |
| <b>ARMPAD (Pair)*</b>   |                          |  |   |
| JBH2195   | <input type="checkbox"/> | Short Armpad (325 mm)  | ◆ |
| <b>HEADREST / NECKREST - Adequate headrest option must be selected*</b> |                          |  |   |
| JBH1350   | <input type="checkbox"/> | Neckrest depth adjustable  | ◆ |
| <b>SEAT LIFTER / SEAT TILT*</b>   |                          |  |   |
| JBH1094   | <input type="checkbox"/> | Seat lifter with powered tilt (300 mm, 30°) with CoG shift   | ◆ |
| <b>SEAT TO GROUND HEIGHT*</b>   |                          |  |   |
| JBH2857   | <input type="checkbox"/> | Seat-to-Floor Height 435 mm  | ○ |
| JBH2858   | <input type="checkbox"/> | Seat-to-Floor Height 460 mm  | ○ |
| JBH2859   | <input type="checkbox"/> | Seat-to-Floor Height 485 mm  | ○ |

| <b>BACKREST RECLINER*</b>    |                          |  |      |
|------------------------------|--------------------------|--|------|
| JBH0975                      | <input type="checkbox"/> | Manual Settable Recliner (in 3.75° steps) 90° - 120°   | ◆    |
| <b>FOOTREST HANGER TYPE*</b> |                          |  |      |
| JBH0526                      | <input type="checkbox"/> | Standard hanger 80° pair Lower leg length: 350 - 450 mm  | ◆    |
| <b>FOOTPLATE TYPE*</b>       |                          |  |      |
| JBH0512                      | <input type="checkbox"/> | Fixed Footplate pair   | ◆    |
| <b>CONTROLS*</b>             |                          |  |      |
| JBH2204                      | <input type="checkbox"/> | Remote LiNX<br>DLX-REM216 = Drive, Actuators & Lights 3 Drive functions and power module incl. 2 Actuators & Lights                              | ◆    |
| <b>CONTROL BRACKETRY*</b>    |                          |  |      |
| JBH1207                      | <input type="checkbox"/> | Standard Bracket RH  | ○    |
| JBH1208                      | <input type="checkbox"/> | Standard Bracket LH  | ○    |
| <b>SHROUD COLOUR*</b>        |                          |  |      |
| JBH2871                      | <input type="checkbox"/> | Mid range - Black/Anthracite Shroud  | ○    |
| JBH2872                      | <input type="checkbox"/> | Mid range - Silver Shroud  | ○    |
| JBH2873                      | <input type="checkbox"/> | Mid range - Bengal Red Shroud  | ○    |
| JBH2874                      | <input type="checkbox"/> | Mid range - Mercury Blue Shroud  | ○    |
| JBH2875                      | <input type="checkbox"/> | Mid range - Midnight Purple Shroud   | ○    |
| JBX0299                      | <input type="checkbox"/> | Premium - Individual Shroud Colour<br>(includes coloured spreader bars on Modulite back)<br>RAL number required<br>Extended lead lime: + 15 days | £195 |
| JBH2865                      | <input type="checkbox"/> | Premium - Racing Green Shroud  | £95  |
| JBH2866                      | <input type="checkbox"/> | Premium - Deep Ocean Blue Satin Shroud   | £95  |
| JBH2867                      | <input type="checkbox"/> | Premium - Frosty Blue Shroud   | £95  |
| JBH2868                      | <input type="checkbox"/> | Premium - Ultimate Red Shroud  | £95  |
| JBH2869                      | <input type="checkbox"/> | Premium - Falcon White Shroud  | £95  |
| <b>RIM COLOUR*</b>           |                          |  |      |
| JBH2882                      | <input type="checkbox"/> | Mid range - Black Rim Inserts  | ○    |
| JBH2883                      | <input type="checkbox"/> | Mid range - Silver Rim Inserts   | ○    |
| JBH2884                      | <input type="checkbox"/> | Mid range - Bengal Red Rim Inserts   | ○    |
| JBH2885                      | <input type="checkbox"/> | Mid range - Mercury Blue Rim Inserts   | ○    |
| JBH2886                      | <input type="checkbox"/> | Mid range - Midnight Purple Rim Inserts  | ○    |
| JBX2887                      | <input type="checkbox"/> | Premium - Individual Rim Insert Colour RAL number required<br>Extended lead lime: + 15 days  | ○    |
| JBH2888                      | <input type="checkbox"/> | Premium - Racing Green Rim Inserts   | ○    |
| JBH2889                      | <input type="checkbox"/> | Premium - Deep Ocean Blue Satin Rim Inserts  | ○    |
| JBH2890                      | <input type="checkbox"/> | Premium - Frosty Blue Rim Inserts  | ○    |
| JBH2891                      | <input type="checkbox"/> | Premium - Ultimate Red Rim Inserts   | ○    |
| JBH2892                      | <input type="checkbox"/> | Premium - Falcon White Rim Inserts   | ○    |
| <b>FORKS*</b>                |                          |  |      |
| JBH2632                      | <input type="checkbox"/> | Double sided fork Without castor rim insert<br>(enlarges wide base to 656 mm)  | ◆    |
| <b>FRONT WHEEL*</b>          |                          |  |      |
| JBH0836                      | <input type="checkbox"/> | 9" (225 mm) x 2.80 (70 mm) Castor, puncture proof solid tyre   | ◆    |
| <b>DRIVE WHEEL*</b>          |                          |  |      |
| JBH0705                      | <input type="checkbox"/> | DW 14" (355 mm) x 3" (75 mm), Solid inlay Dual density foam insert   | ◆    |

| <b>C.T.C SUSPENSION*</b>   |                          |   |   |
|----------------------------|--------------------------|---|---|
| JBH2870                    | <input type="checkbox"/> | CTC Suspension Control - Traction - Comfort   | ◆ |
| <b>SUSPENSION SETTING*</b> |                          |   |   |
| JBH1916                    | <input type="checkbox"/> | Suspension up to 160 kg User Weight For RX 40 HD only   | ◆ |
| <b>MOTOR TYPE*</b>         |                          |   |   |
| JBH1908                    | <input type="checkbox"/> | 10 km/h Package<br>Always combined with LED Light system with Indicators, 30 Lux headlights   | ◆ |
| <b>BATTERY TYPE*</b>       |                          |   |   |
| JBH1260                    | <input type="checkbox"/> | 73,5 Ah Gel   | ◆ |
| <b>CHARGER*</b>            |                          |   |   |
| JBH1343                    | <input type="checkbox"/> | Opti Charge Charger 8A (IP44)   | ◆ |
| <b>ANTI-TIPPERS*</b>       |                          |   |   |
| JBH1205                    | <input type="checkbox"/> | Anti-tipper   | ◆ |
| <b>TRANSPORTATION KIT*</b> |                          |   |   |
| JBH1401                    | <input type="checkbox"/> | 6-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19)) required for variants with combined Lift & Tilt module Adequate headrest is compulsory | ◆ |
| <b>LAP BELT*</b>           |                          |   |   |
| JBH2113                    | <input type="checkbox"/> | Lap belt with Auto Lock Flexible mount in side rails  | ◆ |
| <b>REAR VIEW MIRROR*</b>   |                          |   |   |
| JBH1415                    | <input type="checkbox"/> | Rear view mirror RH   | ○ |
| JBH1339                    | <input type="checkbox"/> | Rear view mirror LH   | ○ |

|                    |  |
|--------------------|--|
| <b>TOTAL PRICE</b> |  |
|--------------------|--|

**Non Conforming Options:**

JBH2106      Body fit Kit Outrigging Backrest Brackets and Pivot Point (outrig: 1 inch/2,54 cm per side)

---



For a more complex customer requirement please refer to our Invacare Unique special order form. email [specialsuk@invacare.com](mailto:specialsuk@invacare.com) or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: \_\_\_\_\_ Email: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Description of special request (please include relevant option code or part number):