



Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

JANUARY 2021
Customer Ref:

Customer Ref*, Account No*, Contact No*, Order No*, Delivery Address*, Order Qty*, Order date*

M118 AVIVA RX40 HD Modulite

Standard Features

- Max user weight 160kg (25 Stone)
Crash tested ISO 7176 -Part 19
Modulite adjustable seating system
Telescopic adjustable seat widths from 330mm - 530mm (13" - 21") with an additional 50mm, (2") by arms
Seat depths 410mm - 510mm stepless adjustment (16" - 20")
Flex3 adjustable Backrest Plates with Laguna Backrest cushion
Neckrest depth adjustable
Standard T-Armrests, detachable with short armpad
Powered "fixed pivot" Tilt 20°
Manual settable backrest Recline (0° - 30° in 3.75° steps)
80° hangers with fixed footplate
LiNX control REM211 with standard bracket
Choice of 5 Shroud colours, 5 Rim insert colours
9" (225 mm) x 2.80 (70 mm) Castor, puncture proof solid tyre
14" (355 mm) x 3" (75 mm), Solid inlay black drive wheels
CTC Suspension Control - Traction - Comfort
12 km/h Ultra High Speed Package combined with LED Light system with Indicators, 30 Lux headlights
73Ah batteries
8A (200 - 250V) charger
Anti tippers
4 Point tie down system
Lapbelt with autolock (adjustable both sides)

Standard Base price: £7,568



Please make your selection with a cross in the box ○ Free of charge ◆ standard feature POA special request

Customer weight = KG

Code Description

BASE MODEL*

Table with 3 columns: Code, Description, Price. Row: JBH0002, AVIVA RX 40 HD Base Model Max user weight up to 160 kg HD Version, includes RX 40 Performance Package Wide Base only, £7,568

CHASSIS WIDTH*

Table with 3 columns: Code, Description, Price. Row: JBH2830, Chassis Wide (635 mm/25"), ◆

SEATING TYPE*			
JBH2100	<input type="checkbox"/>	Invacare® Modulite Seat Unit	◆
SEAT FRAME*			
JBH2111	<input type="checkbox"/>	Telescopic Seat Frame and Plate with Side Rails width 380 to 530 mm Additional 50 mm by arms	○
JBH2112	<input type="checkbox"/>	Telescopic Seat Frame + Sling Seat with Side Rails width 380 to 530 mm; + 30 mm seat height Additional 50 mm by arms	○
SEAT WIDTH*			
JBH2121	<input type="checkbox"/>	380 mm Seat Plate Can be reduced to 330 mm in combination with inrigged armrests	○
JBH2122	<input type="checkbox"/>	430 mm Seat Plate Can be reduced to 380 mm in combination with inrigged armrests	○
JBH2123	<input type="checkbox"/>	480 mm Seat Plate	○
JBH2124	<input type="checkbox"/>	530 mm Seat Plate (NC 2106)	○
Body Fit Kit			
JBH2106	<input type="checkbox"/>	Body fit Kit Outrigging Backrest Brackets and Pivot Point (outrig: 1 inch/2,54 cm per side)	○
SEAT DEPTH*			
JBH2132	<input type="checkbox"/>	Seat Depth 410 mm	○
JBH2134	<input type="checkbox"/>	Seat Depth 460 mm	○
JBH2136	<input type="checkbox"/>	Seat Depth 510 mm	○
SEAT CUSHION*			
JBH2142	<input type="checkbox"/>	Matrx Flo-tech Cushion Flo-tech Lite	◆
BACKREST SUPPORT*			
JBH2162	<input type="checkbox"/>	Backposts tall (540 mm) For Flex3 backrests/Tension Adjustable Upholstery/Matrx Backs, measured from seat plate to top of upholstery Width adjustable cast bracket	◆
BACKREST UPHOLSTERY / CUSHION*			
JBH2173	<input type="checkbox"/>	Flex3 adjustable Backrest Plates with Laguna Backrest cushion Flex3 Backrest plate width and height adjustable 560 mm - 790 mm	◆
UPHOLSTERY TYPE*			
JBH2151	<input type="checkbox"/>	Black Fabric Upholstery (Nylon with option 2141; 2- way-stretch material with options 2142, 2143)	○
JBH2152	<input type="checkbox"/>	Black multi stretch PU coated fabric Upholstery	○
JBH2155	<input type="checkbox"/>	Black Spacer Fabric Upholstery	○
ARMREST (Pair)*			
JBH2191	<input type="checkbox"/>	T-Arm detachable, Standard Adjustment range from seat plate: 245 - 315 mm	◆
ARMPAD (Pair)*			
JBH2195	<input type="checkbox"/>	Short Armpad (325 mm)	◆
HEADREST / NECKREST - Adequate headrest option must be selected*			
JBH1350	<input type="checkbox"/>	Neckrest depth adjustable	◆
SEAT LIFTER / SEAT TILT*			
JBH1092	<input type="checkbox"/>	Powered "fixed pivot" Tilt 20°	◆
SEAT TO GROUND HEIGHT*			
JBH2857	<input type="checkbox"/>	Seat-to-Floor Height 435 mm	○
JBH2858	<input type="checkbox"/>	Seat-to-Floor Height 460 mm	○
JBH2859	<input type="checkbox"/>	Seat-to-Floor Height 485 mm	○
BACKREST RECLINER*			
JBH0975	<input type="checkbox"/>	Manual Settable Recliner (in 3.75° steps) 90° - 120°	◆

FOOTREST HANGER TYPE*			
JBH0526	<input type="checkbox"/>	Standard hanger 80° pair Lower leg length: 350 - 450 mm	◆
FOOTPLATE TYPE*			
JBH0512	<input type="checkbox"/>	Fixed Footplate pair	◆
CONTROLS*			
JBH2204	<input type="checkbox"/>	Remote LiNX DLX-REM216 = Drive, Actuators & Lights 3 Drive functions and power module incl. 2 Actuators & Lights	◆
CONTROL BRACKETRY*			
JBH1207	<input type="checkbox"/>	Standard Bracket RH	○
JBH1208	<input type="checkbox"/>	Standard Bracket LH	○
SHROUD COLOUR*			
JBH2871	<input type="checkbox"/>	Mid range - Black/Anthracite Shroud	○
JBH2872	<input type="checkbox"/>	Mid range - Silver Shroud	○
JBH2873	<input type="checkbox"/>	Mid range - Bengal Red Shroud	○
JBH2874	<input type="checkbox"/>	Mid range - Mercury Blue Shroud	○
JBH2875	<input type="checkbox"/>	Mid range - Midnight Purple Shroud	○
JBX0299	<input type="checkbox"/>	Premium - Individual Shroud Colour (includes coloured spreader bars on Modulite back) RAL number required Extended lead lime: + 15 days	£195
JBH2865	<input type="checkbox"/>	Premium - Racing Green Shroud	£95
JBH2866	<input type="checkbox"/>	Premium - Deep Ocean Blue Satin Shroud	£95
JBH2867	<input type="checkbox"/>	Premium - Frosty Blue Shroud	£95
JBH2868	<input type="checkbox"/>	Premium - Ultimate Red Shroud	£95
JBH2869	<input type="checkbox"/>	Premium - Falcon White Shroud	£95
RIM COLOUR*			
JBH2882	<input type="checkbox"/>	Mid range - Black Rim Inserts	○
JBH2883	<input type="checkbox"/>	Mid range - Silver Rim Inserts	○
JBH2884	<input type="checkbox"/>	Mid range - Bengal Red Rim Inserts	○
JBH2885	<input type="checkbox"/>	Mid range - Mercury Blue Rim Inserts	○
JBH2886	<input type="checkbox"/>	Mid range - Midnight Purple Rim Inserts	○
JBX2887	<input type="checkbox"/>	Premium - Individual Rim Insert Colour RAL number required Extended lead lime: + 15 days	○
JBH2888	<input type="checkbox"/>	Premium - Racing Green Rim Inserts	○
JBH2889	<input type="checkbox"/>	Premium - Deep Ocean Blue Satin Rim Inserts	○
JBH2890	<input type="checkbox"/>	Premium - Frosty Blue Rim Inserts	○
JBH2891	<input type="checkbox"/>	Premium - Ultimate Red Rim Inserts	○
JBH2892	<input type="checkbox"/>	Premium - Falcon White Rim Inserts	○
FORKS*			
JBH2632	<input type="checkbox"/>	Double sided fork Without castor rim insert (enlarges wide base to 656 mm)	◆
FRONT WHEEL*			
JBH0836	<input type="checkbox"/>	9" (225 mm) x 2.80 (70 mm) Castor, puncture proof solid tyre	◆
DRIVE WHEEL*			
JBH0705	<input type="checkbox"/>	DW 14" (355 mm) x 3" (75 mm), Solid inlay Dual density foam insert	◆
C.T.C SUSPENSION*			
JBH2870	<input type="checkbox"/>	CTC Suspension Control - Traction - Comfort	◆

SUSPENSION SETTING*			
JBH1916	<input type="checkbox"/>	Suspension up to 160 kg User Weight For RX 40 HD only	◆
MOTOR TYPE*			
JBH1909	<input type="checkbox"/>	12 km/h Ultra High Speed Package Always combined with LED Light system with Indicators, 30 Lux headlights	◆
BATTERY TYPE*			
JBH1260	<input type="checkbox"/>	73,5 Ah Gel	◆
CHARGER*			
JBH1343	<input type="checkbox"/>	Opti Charge Charger 8A (IP44)	◆
ANTI-TIPPERS*			
JBH1205	<input type="checkbox"/>	Anti-tipper	◆
TRANSPORTATION KIT*			
JBH1494	<input type="checkbox"/>	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19) Adequate headrest is compulsory, Not compatible with Seat lifter	◆
LAP BELT*			
JBH2113	<input type="checkbox"/>	Lap belt with Auto Lock Flexible mount in side rails	◆
REAR VIEW MIRROR*			
JBH1415	<input type="checkbox"/>	Rear view mirror RH	○
JBH1339	<input type="checkbox"/>	Rear view mirror LH	○

TOTAL PRICE	
--------------------	--

Non Conforming Options:

JBH2106 Body fit Kit Outrigging Backrest Brackets and Pivot Point (outrig: 1 inch/2,54 cm per side)



For a more complex customer requirement please refer to our Invacare Unique special order form. email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: _____ Email: _____

Phone: _____ Fax: _____

Description of special request (please include relevant option code or part number):