



Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

JANUARY 2021
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

M16 TDX SP2 NB

Standard Features

- Max User Weight 150kg (23.5 Stone)
- Crash test ISO7176 Part-19
- One Piece seat plate width 430mm (adjustable by arms 440mm - 480mm)
- Seat depth 410mm - 560mm stepless adjustment (16" - 22")
- Standard foam cushion 50mm thick
- Backrest frame, detachable, height 480mm
- Tension adjustable backrest
- Standard T-armrests, detachable with short armpad
- Push bar dia 22mm in conjunction with telescopic seat frame / backposts, No headrest - No mounting hardware
- Manual settable seat Tilt (0° - 20°)
- Manual settable backrest Recline (0° - 30° in 3.75° steps
- 80° hangers with fixed footplate
- LiNX control REM110 with standard bracket
- Selection of 10 stunning shroud and Rim colours
- Centre (driving wheels) 355mm x 76mm Black tyres
- Selection of 10 stunning shroud and Rim colours
- 6km/h (4mph) 4 pole motors and 50Ah (Gel) battery
- 8A (200 - 250V) charger
- 4 Point tie-down system (Q'Straint QLK-150 compatible)
- Lapbelt with autolock (adjustable both sides)

Standard Base price: £4,990



Please make your selection with a cross in the box Free of charge standard feature POA special request

Customer weight = KG

Code	Description	
BASE PRODUCT*		
ARM2100	<input type="checkbox"/> TDX SP2 NARROW BASE	£4,990
USER WEIGHT*		
ARM1920	<input type="checkbox"/> User weight up to 125kg	<input type="radio"/>
ARM1921	<input type="checkbox"/> User weight 126kg to 136kg	<input type="radio"/>
ARM1922	<input type="checkbox"/> User weight 137kg to 150kg	<input type="radio"/>
BASE MODEL*		
ARM0000	<input type="checkbox"/> TDX SP2 Narrow Base	<input checked="" type="checkbox"/>

SEATING TYPE*			
ARM2100	<input type="checkbox"/>	Modulite Seating Unit	◆
SEAT FRAME*			
ARM2110	<input type="checkbox"/>	One Piece Seat Plate	◆
SEAT WIDTH*			
ARM2122	<input type="checkbox"/>	430mm (17")	◆
SEAT DEPTH*			
ARM2132	<input type="checkbox"/>	410mm (16")	○
ARM2134	<input type="checkbox"/>	460mm (18")	○
ARM2136	<input type="checkbox"/>	510mm (20")	○
SEAT CUSHION*			
ARM2141	<input type="checkbox"/>	Standard Foam	◆
BACKREST SUPPORT*			
ARM2161	<input type="checkbox"/>	Detachable 480mm	◆
BACKREST UPHOLSTERY / CUSHION*			
ARM2171	<input type="checkbox"/>	Tension Adjustable Backrest	◆
UPHOLSTERY TYPE*			
ARM2151	<input type="checkbox"/>	Black fabric	◆
ARMREST*			
ARM2191	<input type="checkbox"/>	T-Armrest standard	◆
ARMPAD*			
ARM2195	<input type="checkbox"/>	Short 325mm Armpad	◆
HEADREST / NECKREST - Adequate headrest option must be selected*			
ARM1495	<input type="checkbox"/>	Push bar 22mm	◆
SEAT LIFTER / SEAT TILT*			
ARM1091	<input type="checkbox"/>	Manually settable	◆
SEAT HEIGHT SETTING FOR 20° SEAT TILT*			
ARM0354	<input type="checkbox"/>	Seat height 450mm (17.75")	○
BACKREST RECLINER*			
ARM0975	<input type="checkbox"/>	Manual settable Recliner	◆
FOOTREST HANGER TYPE*			
ARM0526	<input type="checkbox"/>	Side mounted : Standard 80° hanger	◆
FOOTPLATE TYPE*			
ARM0512	<input type="checkbox"/>	Fixed Footplate	◆
CONTROLS*			
ARM2201	<input type="checkbox"/>	REM110	◆
CONTROL BRACKETRY*			
ARM1207	<input type="checkbox"/>	Standard Bracket RH	○
ARM1208	<input type="checkbox"/>	Standard Bracket LH	○
ARM1707	<input type="checkbox"/>	Swing away RH Bracket	£108
ARM1708	<input type="checkbox"/>	Swing away LH Bracket	£108

SHROUD COLOUR*		
ARM0228	<input type="checkbox"/>	Wet Black ○
ARM2601	<input type="checkbox"/>	SpaceStation Silver ○
ARM2602	<input type="checkbox"/>	Island Blue ○
ARM2603	<input type="checkbox"/>	Grape JellyBean (Purple) ○
ARM2604	<input type="checkbox"/>	Monster Green ○
ARM2605	<input type="checkbox"/>	Tangy Orange ○
ARM2606	<input type="checkbox"/>	WhiteOut ○
ARM2607	<input type="checkbox"/>	Rockstar Red ○
ARM2608	<input type="checkbox"/>	Lights Out Black - Matt ○
ARM2609	<input type="checkbox"/>	Invacare Blue ○
RIM COLOUR*		
ARM2611	<input type="checkbox"/>	SpaceStation Silver ○
ARM2612	<input type="checkbox"/>	Island Blue ○
ARM2613	<input type="checkbox"/>	Grape JellyBean (Purple) ○
ARM2614	<input type="checkbox"/>	Monster Green ○
ARM2615	<input type="checkbox"/>	Tangy Orange ○
ARM2616	<input type="checkbox"/>	WhiteOut ○
ARM2617	<input type="checkbox"/>	Rockstar Red ○
ARM2618	<input type="checkbox"/>	Lights Out Black - Matt ○
ARM2619	<input type="checkbox"/>	Invacare Blue ○
ARM2620	<input type="checkbox"/>	Wet Black ○
ARM2621	<input type="checkbox"/>	Without rim insert ○
FORKS*		
ARM2632	<input type="checkbox"/>	Double sided ◆
CENTRE WHEELS (DRIVING WHEELS)*		
ARM0733	<input type="checkbox"/>	Black tyres with gel foam inserts (solid) ◆
MOTOR TYPE*		
ARM1282	<input type="checkbox"/>	6 km/h (4mph) ◆
BATTERY*		
ARM1257	<input type="checkbox"/>	50 Ah Gel ◆
CHARGER*		
ARM1343	<input type="checkbox"/>	OptiCharge 8A (200 - 250V) ◆
TRANSPORTATION KIT*		
ARM1450	<input type="checkbox"/>	Preparation for DAHL docking station £478
ARM1494	<input type="checkbox"/>	4-Point tie-down ◆
FIXATION BELT*		
ARM1329	<input type="checkbox"/>	Lap Belt fitted to base of backrest ◆
TOTAL PRICE		