



Yes, you can:

Powerchairs Retail Prescription Form

Tel: 01656 776222 Fax: 01656 776220
 Email: ordersuk@invacare.com
 Online Spares available at: www.invacare.co.uk

JANUARY 2021
 Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

M22 TDX SP2 HD

Standard Features

- Max user weight 150kg (28 stone)
- Crash tested ISO 7176 -Part 19
- One Piece seat plate 430mm (adjustable by arms 440mm - 480mm)
- Standard foam cushion 50mm thick
- Backrest frame, detachable, height 480mm
- Tension adjustable backrest, Black fabric upholstery
- Standard T-armrests, detachable with short armpad
- Headrest, height, depth and sideways adjustable
- Manual settable seat Tilt 0° - 20°
- Manual settable backrest Recline 0° - 30° in 3.75° steps
- 80° hangers with fixed footplate
- LiNX control REM116 standard bracket
- Selection of 10 stunning shroud and Rim colours
- Centre (driving wheels) 14" x 3" Black tyres
- 10km/h (6mph) 4 pole motors and 73.5 Ah batteries
- 8A (200 - 250V) charger
- 4-Point tie-down system (Q'Straint QLK-150 compatible)
- Light system with indicators, LED headlights and Mirror
- Lapbelt with autolock (adjustable both sides)

Standard Base price: £5,731



Q'STRAIN
QLK-150
 Compatible

Please make your selection with a cross in the box ○ Free of charge ◆ standard feature POA special request

Customer weight = KG

Code	Description	
CHAIR TYPE*		
ARJ0000	<input type="checkbox"/> TDX SP2 Heavy Duty	◆
User Weight*		
ARJ1918	<input type="checkbox"/> User weight less 100 kg	○
ARJ1919	<input type="checkbox"/> User weight 101 to 125 kg	○
ARJ1921	<input type="checkbox"/> User weight 126 to 136 kg	○
ARJ1922	<input type="checkbox"/> User weight 137 to 150 kg	○
BASE TYPE*		
ARJ2100	<input type="checkbox"/> Modulite Seating unit	£5,731

SEATING TYPE*		
ARJ2110	<input type="checkbox"/>	One piece seat plate width 430mm (adjustable by arms 440mm - 480mm) ◆
SEAT WIDTH*		
ARJ2122	<input type="checkbox"/>	Set at 430mm (17") (adjustable from 430mm to 480mm (17" - 19) with stepless adjustment on armrests, legrests and backrest) (NC 2138 + 2142) ◆
SEAT DEPTH SETTING*		
ARJ2132	<input type="checkbox"/>	410mm (16") (410mm - 510mm simple stepless adjustment from the side) ○
ARJ2134	<input type="checkbox"/>	460mm (18") (410mm - 510mm simple stepless adjustment from the side) ○
ARJ2136	<input type="checkbox"/>	510mm (20") (410mm - 510mm simple stepless adjustment from the side) ○
SEAT CUSHION*		
ARJ2141	<input type="checkbox"/>	Standard foam cushion 50mm thick ◆
BACKREST SUPPORT*		
ARJ2161	<input type="checkbox"/>	Backrest frame, detachable, height 480mm, 19" (one size, manual settings only) ◆
BACKREST UPHOLSTERY / CUSHION*		
ARJ2171	<input type="checkbox"/>	Tension Adjustable Backrest (540mm or 480mm upholstery height) Straps easily accessible and adjustable from the rear (Black fabric only) ◆
UPHOLSTERY COLOUR*		
ARJ2151	<input type="checkbox"/>	Black fabric ◆
ARMREST*		
ARJ2191	<input type="checkbox"/>	T-Armrest Standard, Range of adjustment from seat plate 245mm to 315mm ◆
ARMPAD*		
ARJ2195	<input type="checkbox"/>	Short Armpad (325mm) ◆
HEADREST / NECKREST - Adequate headrest option must be selected*		
ARJ1432	<input type="checkbox"/>	Headrest, height, depth and sideways adjustable ◆
SEAT LIFTER / SEAT TILT*		
ARJ1091	<input type="checkbox"/>	Manually settable Seat Tilt (0° - 20°) Seat Height 450mm ◆
SEAT HEIGHT SETTING FOR 20 DEGREES SEAT TILT*		
ARJ0354	<input type="checkbox"/>	Seat height set at 450mm (17.75") (from floor to front of seat plate at 0° tilt) ◆
BACKREST RECLINER*		
ARJ0975	<input type="checkbox"/>	Manual Settable Recliner (0° - 30°, in 3.75° steps) Foldable ◆
FOOTREST HANGER TYPE*		
ARJ0526	<input type="checkbox"/>	Side Mounted : Standard 80° hangers ◆
FOOTPLATE TYPE*		
ARJ0512	<input type="checkbox"/>	Fixed Footplate ◆
CONTROLS*		
ARJ2204	<input type="checkbox"/>	REM216 Drive, Actuators & Lights ◆
CONTROL BRACKETRY*		
ARJ1207	<input type="checkbox"/>	Standard Bracket RH ○
ARJ1208	<input type="checkbox"/>	Standard Bracket LH ○
ARJ1707	<input type="checkbox"/>	Swing Away RH Bracket £108
ARJ1708	<input type="checkbox"/>	Swing Away LH Bracket £108

SHROUD COLOUR*			
ARJ0228	<input type="checkbox"/>	Wet Black	<input type="radio"/>
ARJ2601	<input type="checkbox"/>	SpaceStation Silver	<input type="radio"/>
ARJ2602	<input type="checkbox"/>	Island Blue	<input type="radio"/>
ARJ2603	<input type="checkbox"/>	Grape JellyBean (Purple)	<input type="radio"/>
ARJ2604	<input type="checkbox"/>	Monster Green	<input type="radio"/>
ARJ2605	<input type="checkbox"/>	Tangy Orange	<input type="radio"/>
ARJ2606	<input type="checkbox"/>	WhiteOut	<input type="radio"/>
ARJ2607	<input type="checkbox"/>	Rockstar Red	<input type="radio"/>
ARJ2608	<input type="checkbox"/>	Lights Out Black - Matt	<input type="radio"/>
ARJ2609	<input type="checkbox"/>	Invacare Blue	<input type="radio"/>
RIM COLOUR*			
ARJ2611	<input type="checkbox"/>	SpaceStation Silver	<input type="radio"/>
ARJ2612	<input type="checkbox"/>	Island Blue	<input type="radio"/>
ARJ2613	<input type="checkbox"/>	Grape JellyBean (Purple)	<input type="radio"/>
ARJ2614	<input type="checkbox"/>	Monster Green	<input type="radio"/>
ARJ2615	<input type="checkbox"/>	Tangy Orange	<input type="radio"/>
ARJ2616	<input type="checkbox"/>	WhiteOut	<input type="radio"/>
ARJ2617	<input type="checkbox"/>	Rockstar Red	<input type="radio"/>
ARJ2618	<input type="checkbox"/>	Lights Out Black - Matt	<input type="radio"/>
ARJ2619	<input type="checkbox"/>	Invacare Blue	<input type="radio"/>
ARJ2620	<input type="checkbox"/>	Wet Black	<input type="radio"/>
ARJ2621	<input type="checkbox"/>	Without Rim insert	<input type="radio"/>
FORKS*			
ARJ2632	<input type="checkbox"/>	Double sided fork	<input type="radio"/>
CENTRE WHEELS (DRIVING WHEELS)*			
ARJ0733	<input type="checkbox"/>	355mm x 76mm (14" x 3") Black Tyres with Gel Foam Inserts	<input type="radio"/>
MOTOR TYPE*			
ARJ1908	<input type="checkbox"/>	10km/h (6mph) (4 pole) (includes LED light system and indicators)	<input checked="" type="radio"/>
BATTERY*			
ARJ1260	<input type="checkbox"/>	73.5 Ah (Gel)	<input checked="" type="radio"/>
CHARGER*			
ARJ1343	<input type="checkbox"/>	OptiCharge Charger 8A (200 - 250V)	<input checked="" type="radio"/>
TRANSPORTATION KIT*			
ARJ1450	<input type="checkbox"/>	Preparation for DAHL Docking station: Docking station for car implementation, please order docking station separately (Adequate headrest is compulsory, see options below)	£478
ARJ1494	<input type="checkbox"/>	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19:2001) (Adequate headrest is compulsory, see options below)	<input checked="" type="radio"/>
FIXATION BELT*			
ARJ2113	<input type="checkbox"/>	Lapbelt (adjustable both sides) adjustable mounting on seat rail extrusion	<input type="radio"/>
REAR VIEW MIRROR*			
ARJ1415	<input type="checkbox"/>	Rear view Mirror RH - Must be selected with 10 km/h Motor	<input type="radio"/>
ARJ1339	<input type="checkbox"/>	Rear view Mirror LH - Must be selected with 10 km/h Motor	<input type="radio"/>
TOTAL PRICE			<input type="text"/>