



Yes, you can:

Powerchairs Retail Prescription Form

Tel: 01656 776222 Fax: 01656 776220
 Email: ordersuk@invacare.com
 Online Spares available at: www.invacare.co.uk

JANUARY 2021
 Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

M24 TDX SP2 NB

Standard Features

- Max User Weight 150kg (23.5 Stone)
- Crash test ISO7176 Part-19
- Compact, manoeuvrable, CWD TDX SP2 Narrow Base
- Telescopic adjustable seat widths from 380mm - 530mm (15" - 21") with an additional 50mm, (2") by arms
- Seat depths 410mm - 560mm stepless adjustment (16" - 22")
- Matrix Flo-tech cushion
- Flex3 adjustable backrest plates with Laguna backrest cushion
- Standard T-armrests, detachable, short armpad (325mm)
- Push bar dia 22mm in conjunction with telescopic seat frame / Backposts, No headrest - No mounting hardware
- Manual settable seat Tilt (0° - 20°)
- Manual settable Recliner (0° - 30° in 3.75° steps)
- 80° hangers with fixed footplate
- LiNX control REM211 with standard bracket
- Selection of 10 stunning shroud and Rim colours
- Centre (driving wheels) 355mm x 76mm Black tyres
- 6km/h (4mph) motors and 50Ah Gel batteries
- 8A (200 - 250V) charger
- 4 Point tie-down brackets (Q'Straint QLK-150 compatible)
- Lap belt with auto lock (flexible mount in side rails)



Standard Base price: £5,443

Please make your selection with a cross in the box ○ Free of charge ◆ standard feature POA special request

Customer weight = KG

Code	Description	
BASE PRODUCT*		
ARM2100	<input type="checkbox"/> TDX SP2 NARROW BASE	£5,443
USER WEIGHT*		
ARM1920	<input type="checkbox"/> User weight up to 125kg	○
ARM1921	<input type="checkbox"/> User weight 126kg to 136kg	○
ARM1922	<input type="checkbox"/> User weight 137kg to 150kg	○

BASE MODEL*			
ARM0000	<input type="checkbox"/>	TDX SP2 Narrow Base	◆
SEATING TYPE*			
ARM2100	<input type="checkbox"/>	Modulite Seating Unit	◆
SEAT FRAME*			
ARM2111	<input type="checkbox"/>	Telescopic Seat Frame	◆
SEAT WIDTH*			
ARM2121	<input type="checkbox"/>	380mm (15")	○
ARM2122	<input type="checkbox"/>	430mm (17")	○
ARM2123	<input type="checkbox"/>	480mm (19")	○
ARM2124	<input type="checkbox"/>	530mm (21")	○
SEAT DEPTH RAIL LENGTH*			
ARM2129	<input type="checkbox"/>	Standard 410mm - 510mm	(NC 2138) ○
ARM2190	<input type="checkbox"/>	Extended 510mm - 560mm	(NC 2132, 2134) ○
SEAT DEPTH*			
ARM2132	<input type="checkbox"/>	410mm (16")	○
ARM2134	<input type="checkbox"/>	460mm (18")	○
ARM2136	<input type="checkbox"/>	510mm (20")	○
ARM2138	<input type="checkbox"/>	560mm (22")	○
SEAT CUSHION*			
ARM2142	<input type="checkbox"/>	Matrx Flo-Tech Lite	◆
BACKREST SUPPORT*			
ARM2162	<input type="checkbox"/>	Tall 540mm	◆
BACKREST UPHOLSTERY / CUSHION*			
ARM2173	<input type="checkbox"/>	Flex3 Backrest plates with Laguna cushion	◆
UPHOLSTERY TYPE*			
ARM2151	<input type="checkbox"/>	Black fabric	○
ARM2152	<input type="checkbox"/>	Black Black multi stretch PU coated fabric Upholstery	○
ARM2155	<input type="checkbox"/>	Black Spacer fabric	○
ARMREST*			
ARM2191	<input type="checkbox"/>	T-Armrest standard	◆
ARMPAD*			
ARM2195	<input type="checkbox"/>	Short 325mm Armpad	◆
HEADREST / NECKREST - Adequate headrest option must be selected*			
ARM1495	<input type="checkbox"/>	Push bar 22mm	◆
SEAT LIFTER / SEAT TILT*			
ARM1091	<input type="checkbox"/>	Manually settable	◆
SEAT HEIGHT SETTING FOR 20° SEAT TILT*			
ARM0354	<input type="checkbox"/>	Seat height 450mm (17.75")	○
BACKREST RECLINER*			
ARM0975	<input type="checkbox"/>	Manual settable Recliner	◆
FOOTREST HANGER TYPE*			
ARM0526	<input type="checkbox"/>	Side mounted : Standard 80° hanger	◆
FOOTPLATE TYPE*			
ARM0512	<input type="checkbox"/>	Fixed Footplate	◆
CONTROLS*			
ARM2203	<input type="checkbox"/>	REM211	◆

CONTROL BRACKETRY*			
ARM1207	<input type="checkbox"/>	Standard Bracket RH	○
ARM1208	<input type="checkbox"/>	Standard Bracket LH	○
ARM1707	<input type="checkbox"/>	Swing away RH Bracket	£108
ARM1708	<input type="checkbox"/>	Swing away LH Bracket	£108
SHROUD COLOUR*			
ARM0228	<input type="checkbox"/>	Wet Black	○
ARM2601	<input type="checkbox"/>	SpaceStation Silver	○
ARM2602	<input type="checkbox"/>	Island Blue	○
ARM2603	<input type="checkbox"/>	Grape JellyBean (Purple)	○
ARM2604	<input type="checkbox"/>	Monster Green	○
ARM2605	<input type="checkbox"/>	Tangy Orange	○
ARM2606	<input type="checkbox"/>	WhiteOut	○
ARM2607	<input type="checkbox"/>	Rockstar Red	○
ARM2608	<input type="checkbox"/>	Lights Out Black - Matt	○
ARM2609	<input type="checkbox"/>	Invacare Blue	○
RIM COLOUR*			
ARM2611	<input type="checkbox"/>	SpaceStation Silver	○
ARM2612	<input type="checkbox"/>	Island Blue	○
ARM2613	<input type="checkbox"/>	Grape JellyBean (Purple)	○
ARM2614	<input type="checkbox"/>	Monster Green	○
ARM2615	<input type="checkbox"/>	Tangy Orange	○
ARM2616	<input type="checkbox"/>	WhiteOut	○
ARM2617	<input type="checkbox"/>	Rockstar Red	○
ARM2618	<input type="checkbox"/>	Lights Out Black - Matt	○
ARM2619	<input type="checkbox"/>	Invacare Blue	○
ARM2620	<input type="checkbox"/>	Wet Black	○
ARM2621	<input type="checkbox"/>	Without rim insert	○
FORKS*			
ARM2632	<input type="checkbox"/>	Double sided	◆
CENTRE WHEELS (DRIVING WHEELS)*			
ARM0733	<input type="checkbox"/>	Black tyres with gel foam inserts (solid)	◆
MOTOR TYPE*			
ARM1282	<input type="checkbox"/>	6 km/h (4mph)	◆
BATTERY*			
ARM1257	<input type="checkbox"/>	50 Ah Gel	◆
CHARGER*			
ARM1343	<input type="checkbox"/>	OptiCharge 8A (200 - 250V)	◆
TRANSPORTATION KIT*			
ARM1450	<input type="checkbox"/>	Preparation for DAHL docking station	£478
ARM1494	<input type="checkbox"/>	4-Point tie-down	◆
FIXATION BELT*			
ARM1329	<input type="checkbox"/>	Lap Belt fitted to base of backrest	○
ARM2113	<input type="checkbox"/>	Lap Belt adjustable mounting on seat rail extrusion	◆

TOTAL PRICE	
--------------------	--

Customer Ref:

M24 TDX SP2 NB
Page 3 of 4Subject to Invacares standard terms and conditions.
Invacare reserve the right to alter product specification without prior notice

Order Number:

Non Conforming Options:

ARM2132 410mm (16")

ARM2134 460mm (18")

ARM2138 560mm (22")