



Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

JANUARY 2021
Customer Ref:

Customer Ref*, Contact No*, Delivery Address*, Account No*, Order No*, Order Qty*, Order date*

M48 Esprit Action 4NG Adult

Standard Features

- Max User Weight : 125kg (19.7 Stone)
Alber E-FIX 12" x 2.1/4" propulsion, 12Ah battery
Seat widths 380mm to 605mm (15" - 24 ")
Seat depths 400mm, 450mm or 500mm (16", 18" or 20")
Seat heights 460mm or 485mm (18" or 19")
Nylon seat upholstery
Half folding 7° 430mm backrest
Contour adjustable Nylon backrest canvas
80° swing in / swing out legrest hangers
Flip up footplates includes heel strap
Height adjustable, removable full armpad armrests
Choice of 6 colours
Anti-tippers with chassis jack and buckle belt

Standard Base price: £3,443



Please make your selection with a cross in the box ○ Free of charge ◆ standard feature POA special request

Customer weight = KG

Table with columns: Code, Description, Price. Rows include BASE PRODUCT* (SFD0210), MODEL* (SFD0040), CHAIR TYPE* (SFD0100), and CHASSIS* (SFD0210).

PROPULSION TYPE*			
SFD0300	<input type="checkbox"/>	Alber E-FIX 12" x 2,1/4" (E26, 125 kg) Pneumatic tyres + solid insert + bumper	◆
USER REMOTE*			
SFD0400	<input type="checkbox"/>	User RH remote	◆
SFD0410	<input type="checkbox"/>	User LH remote	○
BATTERY*			
SFD0500	<input type="checkbox"/>	Battery 12Ah	◆
POWER OPTIONS			
SFD0600	<input type="checkbox"/>	Additional battery 12Ah	£247
SEAT WIDTH*			
SFD1040	<input type="checkbox"/>	380mm (15") (NC 1240, 1250)	○
SFD1050	<input type="checkbox"/>	405mm (16") (NC 1240, 1250)	○
SFD1060	<input type="checkbox"/>	430mm (17")	○
SFD1070	<input type="checkbox"/>	455mm (18")	○
SFD1080	<input type="checkbox"/>	480mm (19")	○
SFD1090	<input type="checkbox"/>	505mm (20")	○
SFD1100	<input type="checkbox"/>	555mm (22") (NC 1210, 2105, 2125, 2126, 2130, 2140)	£89
SFD1110	<input type="checkbox"/>	605mm (24") (NC 2130, 2140)	£94
SEAT DEPTH*			
SFD1210	<input type="checkbox"/>	400mm (16") (NC 5072, 5073)	○
SFD1220	<input type="checkbox"/>	400mm to 450mm (16" - 18") (NC 5073)	○
SFD1230	<input type="checkbox"/>	450mm (18") (NC 5071, 5073)	○
SFD1240	<input type="checkbox"/>	450mm to 500mm (18" - 20") (NC 5071)	£28
SFD1250	<input type="checkbox"/>	500mm (20") (NC 5071)	○
SEAT HEIGHT*			
SFD1300	<input type="checkbox"/>	460mm (18")	○
SFD1310	<input type="checkbox"/>	485mm (19")	○
SEAT UPHOLSTERY TYPE*			
SFD1410	<input type="checkbox"/>	Nylon with hook and loop strips	◆
BACKREST TYPE*			
SFD2000	<input type="checkbox"/>	Folding Backrest 430mm height	◆
BACKREST UPHOLSTERY TYPE*			
SFD2105	<input type="checkbox"/>	Lumbar support (+ Strap)	○
SFD2110	<input type="checkbox"/>	Contour adjustable (4 + 1 Straps)	○
SFD2120	<input type="checkbox"/>	Contour adjustable (10 + 1 Straps)	○
SFD2125	<input type="checkbox"/>	Lumbar support (+ Strap)	○
SFD2126	<input type="checkbox"/>	Contour adjustable (3 + 1 Straps)	£48
SFD2130	<input type="checkbox"/>	Contour adjustable (4 + 1 Straps)	◆
SFD2140	<input type="checkbox"/>	Contour adjustable (10 + 1 Straps)	○
HANGERS RH*			
SFD3040	<input type="checkbox"/>	80° Straight swing in / swing out	◆
HANGERS LH*			
SFD3140	<input type="checkbox"/>	80° Straight, swing in / swing out	◆
FOOTPLATES RH*			
SFD3200	<input type="checkbox"/>	Flip up RH footplate	◆
FOOTPLATES LH*			
SFD3300	<input type="checkbox"/>	Flip up LH footplate	◆

ARMRESTS RH*		
SFD4010	<input type="checkbox"/>	Removable RH full length armpad ◆
ARMRESTS LH*		
SFD4110	<input type="checkbox"/>	Removable LH full length armpad ◆
SEAT CUSHION*		
SFD5022	<input type="checkbox"/>	Foam cushion 50mm ◆
SEAT CUSHION DEPTH		
SFD5071	<input type="checkbox"/>	Seat cushion depth 400mm ○
SFD5072	<input type="checkbox"/>	Seat cushion depth 450mm ○
SFD5073	<input type="checkbox"/>	Seat cushion depth 500mm ○
COLOURS*		
SFD6010	<input type="checkbox"/>	Slate Grey ○
SFD6030	<input type="checkbox"/>	Pearl White ○
SFD6031	<input type="checkbox"/>	Blue Pop ○
SFD6032	<input type="checkbox"/>	Hypnotic Purple ○
SFD6033	<input type="checkbox"/>	Red Matt ○
SFD6034	<input type="checkbox"/>	Matt Black ○
FRONT WHEELS*		
SFD7030	<input type="checkbox"/>	200mm x 50mm (8" x 2") Solid ◆
SAFETY OPTIONS*		
SFD8010	<input type="checkbox"/>	Rear anti-tippers with chassis jacks ◆
SFD8020	<input type="checkbox"/>	Buckle belt ◆
SFD8025	<input type="checkbox"/>	Transportation Kit ◆

TOTAL PRICE	<input style="width: 100px; height: 20px;" type="text"/>
--------------------	--

Non Conforming Options:

SFD1210	400mm (16")
SFD1240	450mm to 500mm (18" - 20")
SFD1250	500mm (20")
SFD2105	Lumbar support (+ Strap)
SFD2125	Lumbar support (+ Strap)
SFD2126	Contour adjustable (3 + 1 Straps)
SFD2130	Contour adjustable (4 + 1 Straps)
SFD2140	Contour adjustable (10 + 1 Straps)
SFD5071	Seat cushion depth 400mm
SFD5072	Seat cushion depth 450mm
SFD5073	Seat cushion depth 500mm



For a more complex customer requirement please refer to our Invacare Unique special order form. email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: _____ Email: _____

Phone: _____ Fax: _____

Description of special request (please include relevant option code or part number):