



Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

JANUARY 2021
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

M51 TDX SP2 NB

Standard Features

- Max User Weight 150kg (23.5 Stone)
- Crash test ISO7176 Part-19
- Compact, manoeuvrable, CWD TDX SP2 Narrow Base
- Telescopic adjustable seat widths from 380mm - 530mm (15" - 21") with an additional 50mm, (2") by arms
- Seat depths 410mm - 560mm stepless adjustment (16" - 22")
- Matrix Flo-tech cushion
- Flex3 adjustable backrest plates with Laguna backrest cushion
- Standard T-armrests, detachable with short armpad
- Elan Headrest Standard Pad
- Electric Tilt (0° - 30°) with CoG shift
- Manual settable backrest Recliner (0° - 30° in 3.75° steps)
- 80° hangers with angle and depth adjustable footplate (adapt)
- LiNX control REM216 with standard bracket
- Selection of 10 stunning shroud and Rim colours
- Centre (driving wheels) 355mm x 76mm Black tyres
- 10km/h (6mph) 4 pole motors (including light system with indicators, LED headlights and mirror)
- 8A (200 - 250V) charger
- 50Ah batteries
- Light system with indicators & LED headlights
- 4 Point tie-down system (Q'Straint QLK-150 compatible)
- Lapbelt with autolock (adjustable both sides)

Standard Base price: £7,113



Please make your selection with a cross in the box ○ Free of charge ◆ standard feature POA special request

Customer weight = KG

Code **Description**

BASE PRODUCT*

ARM2100	<input type="checkbox"/>	TDX SP2 NARROW BASE	£7,113
---------	--------------------------	---------------------	--------

USER WEIGHT*			
ARM1920	<input type="checkbox"/>	User weight up to 125kg	<input type="radio"/>
ARM1921	<input type="checkbox"/>	User weight 126kg to 136kg	<input type="radio"/>
ARM1922	<input type="checkbox"/>	User weight 137kg to 150kg	<input type="radio"/>
BASE MODEL*			
ARM0000	<input type="checkbox"/>	TDX SP2 Narrow Base	<input checked="" type="radio"/>
SEATING TYPE*			
ARM2100	<input type="checkbox"/>	Modulite Seating Unit	<input checked="" type="radio"/>
SEAT FRAME*			
ARM2111	<input type="checkbox"/>	Telescopic Seat Frame	<input checked="" type="radio"/>
SEAT WIDTH*			
ARM2121	<input type="checkbox"/>	380mm (15")	<input type="radio"/>
ARM2122	<input type="checkbox"/>	430mm (17")	<input type="radio"/>
ARM2123	<input type="checkbox"/>	480mm (19")	<input type="radio"/>
ARM2124	<input type="checkbox"/>	530mm (21")	<input type="radio"/>
SEAT DEPTH RAIL LENGTH*			
ARM2129	<input type="checkbox"/>	Standard 410mm - 510mm	<input type="radio"/>
			(NC 2138)
ARM2190	<input type="checkbox"/>	Extended 510mm - 560mm	<input type="radio"/>
			(NC 2132, 2134)
SEAT DEPTH*			
ARM2132	<input type="checkbox"/>	410mm (16")	<input type="radio"/>
ARM2134	<input type="checkbox"/>	460mm (18")	<input type="radio"/>
ARM2136	<input type="checkbox"/>	510mm (20")	<input type="radio"/>
ARM2138	<input type="checkbox"/>	560mm (22")	<input type="radio"/>
SEAT CUSHION*			
ARM2142	<input type="checkbox"/>	Matrx Flo-Tech Lite	<input checked="" type="radio"/>
BACKREST SUPPORT*			
ARM2162	<input type="checkbox"/>	Tall 540mm	<input checked="" type="radio"/>
BACKREST UPHOLSTERY / CUSHION*			
ARM2173	<input type="checkbox"/>	Flex3 Backrest plates with Laguna cushion	<input checked="" type="radio"/>
UPHOLSTERY TYPE*			
ARM2151	<input type="checkbox"/>	Black fabric	<input type="radio"/>
ARM2152	<input type="checkbox"/>	Black Black multi stretch PU coated fabric Upholstery	<input type="radio"/>
ARM2155	<input type="checkbox"/>	Black Spacer fabric	<input type="radio"/>
ARMREST*			
ARM2191	<input type="checkbox"/>	T-Armrest standard	<input checked="" type="radio"/>
ARMPAD*			
ARM2195	<input type="checkbox"/>	Short 325mm Armpad	<input checked="" type="radio"/>
HEADREST / NECKREST - Adequate headrest option must be selected*			
ARM2324	<input type="checkbox"/>	Elan standard pad 255mm (10")	<input checked="" type="radio"/>
SEAT LIFTER / SEAT TILT*			
ARM1093	<input type="checkbox"/>	Power Tilt with centre of gravity shift	<input checked="" type="radio"/>
BACKREST RECLINER*			
ARM0975	<input type="checkbox"/>	Manual settable Recliner	<input checked="" type="radio"/>
FOOTREST HANGER TYPE*			
ARM0526	<input type="checkbox"/>	Side mounted : Standard 80° hanger	<input checked="" type="radio"/>
FOOTPLATE TYPE*			
ARM0513	<input type="checkbox"/>	Angle and Depth adjustable	<input checked="" type="radio"/>

CONTROLS*			
ARM2204	<input type="checkbox"/>	REM216	◆
CONTROL BRACKETRY*			
ARM1207	<input type="checkbox"/>	Standard Bracket RH	○
ARM1208	<input type="checkbox"/>	Standard Bracket LH	○
ARM1707	<input type="checkbox"/>	Swing away RH Bracket	£108
ARM1708	<input type="checkbox"/>	Swing away LH Bracket	£108
SHROUD COLOUR*			
ARM0228	<input type="checkbox"/>	Wet Black	○
ARM2601	<input type="checkbox"/>	SpaceStation Silver	○
ARM2602	<input type="checkbox"/>	Island Blue	○
ARM2603	<input type="checkbox"/>	Grape JellyBean (Purple)	○
ARM2604	<input type="checkbox"/>	Monster Green	○
ARM2605	<input type="checkbox"/>	Tangy Orange	○
ARM2606	<input type="checkbox"/>	WhiteOut	○
ARM2607	<input type="checkbox"/>	Rockstar Red	○
ARM2608	<input type="checkbox"/>	Lights Out Black - Matt	○
ARM2609	<input type="checkbox"/>	Invacare Blue	○
RIM COLOUR*			
ARM2611	<input type="checkbox"/>	SpaceStation Silver	○
ARM2612	<input type="checkbox"/>	Island Blue	○
ARM2613	<input type="checkbox"/>	Grape JellyBean (Purple)	○
ARM2614	<input type="checkbox"/>	Monster Green	○
ARM2615	<input type="checkbox"/>	Tangy Orange	○
ARM2616	<input type="checkbox"/>	WhiteOut	○
ARM2617	<input type="checkbox"/>	Rockstar Red	○
ARM2618	<input type="checkbox"/>	Lights Out Black - Matt	○
ARM2619	<input type="checkbox"/>	Invacare Blue	○
ARM2620	<input type="checkbox"/>	Wet Black	○
ARM2621	<input type="checkbox"/>	Without rim insert	○
FORKS*			
ARM2632	<input type="checkbox"/>	Double sided	◆
CENTRE WHEELS (DRIVING WHEELS)*			
ARM0733	<input type="checkbox"/>	Black tyres with gel foam inserts (solid)	◆
MOTOR TYPE*			
ARM1908	<input type="checkbox"/>	10 km/h (6mph)	◆
BATTERY*			
ARM1257	<input type="checkbox"/>	50 Ah Gel	◆
CHARGER*			
ARM1343	<input type="checkbox"/>	OptiCharge 8A (200 - 250V)	◆
TRANSPORTATION KIT*			
ARM1450	<input type="checkbox"/>	Preparation for DAHL docking station	£478
ARM1494	<input type="checkbox"/>	4-Point tie-down	◆
FIXATION BELT*			
ARM1329	<input type="checkbox"/>	Lap Belt fitted to base of backrest	○
ARM2113	<input type="checkbox"/>	Lap Belt adjustable mounting on seat rail extrusion	◆

REAR VIEW MIRROR*

ARM1415		RH Mirror	○
ARM1339		LH Mirror	○

TOTAL PRICE			
--------------------	--	--	--

Non Conforming Options:

ARM2132 410mm (16")

ARM2134 460mm (18")

ARM2138 560mm (22")



For a more complex customer requirement please refer to our Invacare Unique special order form. email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: _____ Email: _____

Phone: _____ Fax: _____

Description of special request (please include relevant option code or part number):