



**Yes, you can:**

# Powerchairs Retail Prescription Form

Tel: 01656 776222 Fax: 01656 776220  
 Email: ordersuk@invacare.com  
 Online Spares available at: www.invacare.co.uk

JANUARY 2021  
 Customer Ref:

<b>Customer Ref*</b>	<b>Account No*</b>
<b>Contact No*</b>	<b>Order No*</b>
<b>Delivery Address*</b>	<b>Order Qty*</b>
	<b>Order date*</b>

## M61 TDX SP2 NB

### Standard Features

- Max User Weight 150kg (23.5 Stone)
- Crash test ISO7176 Part-19
- Compact, manoeuvrable, CWD TDX SP2 Narrow Base
- Telescopic adjustable seat widths from 380mm - 530mm (15" - 21") with an additional 50mm, (2") by arms
- Seat depths 410mm - 560mm stepless adjustment (16" - 22")
- Matrix Libra cushion
- Flex3 adjustable backrest plates with Laguna backrest cushion
- Flip-up armrests, (height, and depth adjustable) with standard armpad
- Elan Headrest Standard Pad
- Powered Tilt 0° - 20°
- Powered Recliner (0° - 30° in 3.75° steps)
- Manual centre mounted hanger, with angle adjustable footplates
- LiNX control REM211 with swing away bracket
- Centre (driving wheels) 355mm x 76mm Black tyres
- 6km/h (4mph) 4 pole motors and 50Ah batteries
- 8A (200 - 250V) charger
- 4 Point tie-down system (Q'Straint QLK-150 compatible)
- Lapbelt with autolock (adjustable both sides)

**Standard Base price: £7,185**



Please make your selection with a cross in the box  Free of charge  standard feature POA special request

Customer weight =                      KG

Code	Description	
<b>BASE PRODUCT*</b>		
ARM2100	<input type="checkbox"/> TDX SP2 NARROW BASE	£7,185
<b>USER WEIGHT*</b>		
ARM1920	<input type="checkbox"/> User weight up to 125kg	<input type="radio"/>
ARM1921	<input type="checkbox"/> User weight 126kg to 136kg	<input type="radio"/>
ARM1922	<input type="checkbox"/> User weight 137kg to 150kg	<input type="radio"/>
<b>BASE MODEL*</b>		
ARM0000	<input type="checkbox"/> TDX SP2 Narrow Base	<input checked="" type="checkbox"/>

<b>SEATING TYPE*</b>			
ARM2100	<input type="checkbox"/>	Modulite Seating Unit	◆
<b>SEAT FRAME*</b>			
ARM2111	<input type="checkbox"/>	Telescopic Seat Frame	◆
<b>SEAT WIDTH*</b>			
ARM2121	<input type="checkbox"/>	380mm (15")	○
ARM2122	<input type="checkbox"/>	430mm (17")	○
ARM2123	<input type="checkbox"/>	480mm (19")	○
ARM2124	<input type="checkbox"/>	530mm (21")	○
<b>SEAT DEPTH RAIL LENGTH*</b>			
ARM2129	<input type="checkbox"/>	Standard 410mm - 510mm	(NC 2138) ○
ARM2190	<input type="checkbox"/>	Extended 510mm - 560mm	(NC 2132, 2134) ○
<b>SEAT DEPTH*</b>			
ARM2132	<input type="checkbox"/>	410mm (16")	○
ARM2134	<input type="checkbox"/>	460mm (18")	○
ARM2136	<input type="checkbox"/>	510mm (20")	○
ARM2138	<input type="checkbox"/>	560mm (22")	○
<b>SEAT CUSHION*</b>			
ARM2147	<input type="checkbox"/>	Matrix Libra	◆
<b>BACKREST SUPPORT*</b>			
ARM2162	<input type="checkbox"/>	Tall 540mm	◆
<b>BACKREST UPHOLSTERY / CUSHION*</b>			
ARM2173	<input type="checkbox"/>	Flex3 Backrest plates with Laguna cushion	◆
<b>UPHOLSTERY TYPE*</b>			
ARM2151	<input type="checkbox"/>	Black fabric	○
ARM2152	<input type="checkbox"/>	Black Black multi stretch PU coated fabric Upholstery	○
ARM2155	<input type="checkbox"/>	Black Spacer fabric	○
<b>ARMREST*</b>			
ARM2193	<input type="checkbox"/>	Flip-up armrests	◆
<b>ARMPAD*</b>			
ARM2195	<input type="checkbox"/>	Short 325mm Armpad	◆
<b>HEADREST / NECKREST - Adequate headrest option must be selected*</b>			
ARM2324	<input type="checkbox"/>	Elan standard pad 255mm (10")	◆
<b>SEAT LIFTER / SEAT TILT*</b>			
ARM1092	<input type="checkbox"/>	Power Tilt (max user weight 136kg)	◆
<b>SEAT HEIGHT SETTING FOR 20° SEAT TILT*</b>			
ARM0354	<input type="checkbox"/>	Seat height 450mm (17.75")	○
ARM0364	<input type="checkbox"/>	Seat height 475mm (18.75")	○
<b>BACKREST RECLINER*</b>			
ARM0977	<input type="checkbox"/>	Powered Recliner	◆
<b>FOOTREST HANGER TYPE*</b>			
ARM0556	<input type="checkbox"/>	Centre mounted : Manual	◆
<b>FOOTPLATE TYPE*</b>			
ARM0559	<input type="checkbox"/>	Angle Adjustable Footplate Medium 140mm x 240mm (5.5" x 9.5") for Centre Mounted Hanger only	◆
<b>CONTROLS*</b>			
ARM2203	<input type="checkbox"/>	REM211	◆

<b>CONTROL BRACKETRY*</b>			
ARM1707	<input type="checkbox"/>	Swing away RH Bracket	<input type="radio"/>
ARM1708	<input type="checkbox"/>	Swing away LH Bracket	<input type="radio"/>
<b>SHROUD COLOUR*</b>			
ARM0228	<input type="checkbox"/>	Wet Black	<input type="radio"/>
ARM2601	<input type="checkbox"/>	SpaceStation Silver	<input type="radio"/>
ARM2602	<input type="checkbox"/>	Island Blue	<input type="radio"/>
ARM2603	<input type="checkbox"/>	Grape JellyBean (Purple)	<input type="radio"/>
ARM2604	<input type="checkbox"/>	Monster Green	<input type="radio"/>
ARM2605	<input type="checkbox"/>	Tangy Orange	<input type="radio"/>
ARM2606	<input type="checkbox"/>	WhiteOut	<input type="radio"/>
ARM2607	<input type="checkbox"/>	Rockstar Red	<input type="radio"/>
ARM2608	<input type="checkbox"/>	Lights Out Black - Matt	<input type="radio"/>
ARM2609	<input type="checkbox"/>	Invacare Blue	<input type="radio"/>
<b>RIM COLOUR*</b>			
ARM2611	<input type="checkbox"/>	SpaceStation Silver	<input type="radio"/>
ARM2612	<input type="checkbox"/>	Island Blue	<input type="radio"/>
ARM2613	<input type="checkbox"/>	Grape JellyBean (Purple)	<input type="radio"/>
ARM2614	<input type="checkbox"/>	Monster Green	<input type="radio"/>
ARM2615	<input type="checkbox"/>	Tangy Orange	<input type="radio"/>
ARM2616	<input type="checkbox"/>	WhiteOut	<input type="radio"/>
ARM2617	<input type="checkbox"/>	Rockstar Red	<input type="radio"/>
ARM2618	<input type="checkbox"/>	Lights Out Black - Matt	<input type="radio"/>
ARM2619	<input type="checkbox"/>	Invacare Blue	<input type="radio"/>
ARM2620	<input type="checkbox"/>	Wet Black	<input type="radio"/>
ARM2621	<input type="checkbox"/>	Without rim insert	<input type="radio"/>
<b>FORKS*</b>			
ARM2632	<input type="checkbox"/>	Double sided	<input type="checkbox"/>
<b>CENTRE WHEELS (DRIVING WHEELS)*</b>			
ARM0733	<input type="checkbox"/>	Black tyres with gel foam inserts (solid)	<input type="checkbox"/>
<b>MOTOR TYPE*</b>			
ARM1282	<input type="checkbox"/>	6 km/h (4mph)	<input type="checkbox"/>
<b>BATTERY*</b>			
ARM1257	<input type="checkbox"/>	50 Ah Gel	<input type="checkbox"/>
<b>CHARGER*</b>			
ARM1343	<input type="checkbox"/>	OptiCharge 8A (200 - 250V)	<input type="checkbox"/>
<b>TRANSPORTATION KIT*</b>			
ARM1450	<input type="checkbox"/>	Preparation for DAHL docking station	£478
ARM1494	<input type="checkbox"/>	4-Point tie-down	<input type="checkbox"/>
<b>FIXATION BELT*</b>			
ARM1329	<input type="checkbox"/>	Lap Belt fitted to base of backrest	<input type="radio"/>
ARM2113	<input type="checkbox"/>	Lap Belt adjustable mounting on seat rail extrusion	<input type="checkbox"/>

<b>TOTAL PRICE</b>	
--------------------	--

**Non Conforming Options:**

ARM2132 410mm (16")

ARM2134 460mm (18")

ARM2138 560mm (22")



For a more complex customer requirement please refer to our Invacare Unique special order form. email [specialsuk@invacare.com](mailto:specialsuk@invacare.com) or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: \_\_\_\_\_ Email: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Description of special request (please include relevant option code or part number):