



**Yes, you can:**

# Powerchairs Retail Prescription Form

Tel: 01656 776222 Fax: 01656 776220  
 Email: ordersuk@invacare.com  
 Online Spares available at: www.invacare.co.uk

JANUARY 2021  
 Customer Ref:

<b>Customer Ref*</b>	<b>Account No*</b>
<b>Contact No*</b>	<b>Order No*</b>
<b>Delivery Address*</b>	<b>Order Qty*</b>
	<b>Order date*</b>

## M66 Esprit Action 4NG Junior

### Standard Features

- Max User Weight : 80kg (12.6 Stone)
- Alber E-FIX 12" x 2.1/4" propulsion, 12Ah battery
- Seat width 305mm to 380mm (12" - 15")
- Seat depth 325mm to 375mm (13" - 15") adjustable (set at 350mm)
- Seat height 460mm (18")
- Nylon seat upholstery
- 50mm foam cushion
- Half folding 7° 430mm backrest
- Standard Nylon sling type backrest upholstery
- 80° swing in / swing out legrest hangers
- Flip up angle and depth adjustable footplates includes heel strap
- Height adjustable, removable armrests
- Choice of 7 colours
- Anti-tippers with chassis jack and buckle belt

**Standard Base price: £3,443**



Please make your selection with a cross in the box    ○ Free of charge    ◆ standard feature    POA special request

Customer weight =                      KG

Code	Description	
<b>BASE PRODUCT*</b>		
SFE0200	<input type="checkbox"/> ESPRIT ACTION 4NG JUNIOR	£3,443
<b>MODEL*</b>		
SFE0050	<input type="checkbox"/> Esprit Action 4NG - Junior Version	◆
<b>CHAIR TYPE*</b>		
SFE0100	<input type="checkbox"/> Esprit Action 4NG (NC 0200)	◆
<b>CHASSIS*</b>		
SFE0200	<input type="checkbox"/> Single Crossbar	◆

<b>PROPULSION TYPE*</b>			
SFE0300	<input type="checkbox"/>	Alber E-FIX 12" x 2,1/4" (E26, 125 kg) Pneumatic tyres + solid insert + bumper	◆
<b>USER REMOTE*</b>			
SFE0400	<input type="checkbox"/>	User remote control RH	◆
SFE0410	<input type="checkbox"/>	User remote control LH	○
<b>BATTERY*</b>			
SFE0500	<input type="checkbox"/>	Battery 12Ah	◆
<b>POWER OPTIONS</b>			
SFE0600	<input type="checkbox"/>	Additional Battery 12Ah	£247
<b>SEAT WIDTH*</b>			
SFE1000	<input type="checkbox"/>	305mm (12")	○
SFE1010	<input type="checkbox"/>	330mm (13")	○
SFE1020	<input type="checkbox"/>	355mm (14")	○
SFE1030	<input type="checkbox"/>	380mm (15")	○
<b>SEAT DEPTH*</b>			
SFE1200	<input type="checkbox"/>	325mm to 375mm (13" - 15")	◆
<b>SEAT HEIGHT*</b>			
SFE1300	<input type="checkbox"/>	460mm (18")	◆
<b>SEAT UPHOLSTERY TYPE*</b>			
SFE1400	<input type="checkbox"/>	Nylon with hook and loop strips, width & depth adjustable	◆
<b>BACKREST TYPE*</b>			
SFE2000	<input type="checkbox"/>	Folding backrest 430mm Height	◆
<b>BACKREST UPHOLSTERY TYPE*</b>			
SFE2100	<input type="checkbox"/>	Width & height adjustable nylon sling type	◆
SFE2120	<input type="checkbox"/>	Contour adjustable (10 + 1 Straps) padded shaped cover	○
<b>HANGERS RH*</b>			
SFE3000	<input type="checkbox"/>	80° short RH swing in / swing out (NC 3100)	◆
<b>HANGERS LH*</b>			
SFE3100	<input type="checkbox"/>	80° short LH swing in / swing out	◆
<b>FOOTPLATES RH*</b>			
SFE3210	<input type="checkbox"/>	Flip up angle and depth adjustable RH footplate	◆
<b>FOOTPLATES LH*</b>			
SFE3310	<input type="checkbox"/>	Flip up angle and depth adjustable LH footplate	◆
<b>ARMRESTS RH*</b>			
SFE4010	<input type="checkbox"/>	Removable RH full length armpad	◆
<b>ARMRESTS LH*</b>			
SFE4110	<input type="checkbox"/>	Removable LH full length armpad	◆
<b>SEATING OPTIONS*</b>			
SFE5021	<input type="checkbox"/>	Essential cushion depth 375mm	◆

COLOURS*			
SFE6030	<input type="checkbox"/>	Pearl White	<input type="radio"/>
SFE6031	<input type="checkbox"/>	Blue Pop	<input type="radio"/>
SFE6032	<input type="checkbox"/>	Hypnotic Purple	<input type="radio"/>
SFE6033	<input type="checkbox"/>	Red Matt	<input type="radio"/>
SFE6034	<input type="checkbox"/>	Matt Black	<input type="radio"/>
SFE6035	<input type="checkbox"/>	Polar Blue	<input type="radio"/>
SFE6041	<input type="checkbox"/>	Candy Pink	<input type="radio"/>
FRONT WHEELS*			
SFE7010	<input type="checkbox"/>	175mm x 50mm (7" x 2") Solid	<input type="checkbox"/>
SAFETY OPTIONS*			
SFE8010	<input type="checkbox"/>	Rear anti-tippers with chassis jacks	<input type="checkbox"/>
SFE8020	<input type="checkbox"/>	Buckle belt	<input type="checkbox"/>
SFE8025	<input type="checkbox"/>	Transportation Kit	<input type="checkbox"/>

	<b>TOTAL PRICE</b>	
--	--------------------	--

**Non Conforming Options:**

SFE0200          ESPRIT ACTION 4NG JUNIOR

SFE0200          Single Crossbar

SFE3100          80° short LH swing in / swing out



For a more complex customer requirement please refer to our Invacare Unique special order form. email [specialsuk@invacare.com](mailto:specialsuk@invacare.com) or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: \_\_\_\_\_ Email: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Description of special request (please include relevant option code or part number):