



Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

JANUARY 2021
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

M76 Storm 4 Modulite

Standard Features

- Max user weight 150Kg (23.5 stone)
- Crash tested ISO 7176 -Part 19
- Modulite highly adjustable seating system
- Telescopic adjustable Seat Widths from 330mm - 530mm (13" - 21") with an additional 50mm, (2") by arms
- Seat Depths 410mm - 510mm Stepless Adjustment (16" - 20")
- Matrix Flo-tech Lite Cushion
- Tension Adjustable Backrest
- Standard T-Armrests, detachable with Short armpad
- Manual settable seat Tilt 0° - 20°
- Manual settable Recliner (0° - 30° in 3.75° steps)
- VARI F hangers with Fixed Footplate
- LiNX Control REM216 with standard bracket
- Choice of 6 Colours
- Puncture Proof Castors
- Puncture Proof Drive Wheels
- 73.5 Ah Batteries
- 12km/h (8mph) Motors
- 8A (200 - 250V) Charger
- LED light system with indicators and Mirror
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)
- Anti-Tippers
- Rear wheel Suspension

Standard Base price: £6,677



Please make your selection with a cross in the box ○ Free of charge ◆ standard feature POA special request

Customer weight = KG

Code	Description	
CHAIR TYPE*		
ATB0000	<input type="checkbox"/> Storm 4 Modulite	£6,677
SEATING TYPE*		
ATB2100	<input type="checkbox"/> Modulite Seating Unit (Flex3 seating unit)	◆
SEATING OPTIONS*		
ATB2111	<input type="checkbox"/> Telescopic Seat Frame with metal seat plate (with stepless adjustment on armrests, legrests and backrest)	◆

SEAT WIDTH*			
ATBXXXX	<input type="checkbox"/>	Seat width 330mm (13") seat width achieved on a 380mm seat frame (2121) with inrigged armrests (2183)	<input type="radio"/>
ATB2121	<input type="checkbox"/>	Set at 380mm (adjustable from 380mm to 430mm (15" - 17") with stepless adjustment on armrests, legrests and backrest)	<input type="radio"/>
ATB2122	<input type="checkbox"/>	Set at 430mm (adjustable from 430mm to 480mm (17" - 19) with stepless adjustment on armrests, legrests and backrest)	<input type="radio"/>
ATB2123	<input type="checkbox"/>	Set at 480mm (adjustable from 480mm to 530mm (19" - 21") with stepless adjustment on armrests, legrests and backrest)	<input type="radio"/>
ATB2124	<input type="checkbox"/>	Set at 530mm (adjustable from 530mm to 580mm (21" - 23") with stepless adjustment on armrests, legrests and backrest)	<input type="radio"/>
SEAT DEPTH*			
ATB2132	<input type="checkbox"/>	410mm (16") (adjustable 410mm - 510mm stepless adjustment)	<input type="radio"/>
ATB2134	<input type="checkbox"/>	460mm (18") (adjustable 410mm - 510mm stepless adjustment)	<input type="radio"/>
ATB2136	<input type="checkbox"/>	510mm (20") (adjustable 410mm - 510mm stepless adjustment)	<input type="radio"/>
SEAT CUSHION*			
ATB2142	<input type="checkbox"/>	Matrx Flo-tech Lite	<input checked="" type="checkbox"/>
BACKREST SUPPORT*			
ATB2162	<input type="checkbox"/>	Backposts 540mm (Tall) must be selected for Flex3 Laguna backrest	<input checked="" type="checkbox"/>
PUSH HANDLES FOR BACKPOSTS*			
ATB2107	<input type="checkbox"/>	Backposts WITH push handles	<input checked="" type="checkbox"/>
BACKREST UPHOLSTERY / CUSHION*			
ATB2171	<input type="checkbox"/>	Basic Tension Adjustable Backrest (540mm or 480mm upholstery height) Straps easily accessible and adjustable from the rear, Black fabric only	<input checked="" type="checkbox"/>
UPHOLSTREY TYPE*			
ATB2151	<input type="checkbox"/>	Black fabric	<input checked="" type="checkbox"/>
ARMREST (PAIR)*			
ATB2191	<input type="checkbox"/>	T-Armrest Standard, Range of adjustment from seat plate 245mm to 315mm	<input checked="" type="checkbox"/>
ARMPAD (PAIR)*			
ATB2195	<input type="checkbox"/>	Short Armpad (325mm)	<input checked="" type="checkbox"/>
HEADREST / NECKREST - Adequate headrest option must be selected*			
ATB1432	<input type="checkbox"/>	Headrest, height, depth and sideways adjustable	<input checked="" type="checkbox"/>
SEAT HEIGHT AND SEAT RISER*			
ATB1023	<input type="checkbox"/>	Fixed Seat Height 450mm (17.75")	<input type="radio"/>
ATB1025	<input type="checkbox"/>	Fixed Seat Height 480mm (19")	<input type="radio"/>
SEAT TILT*			
ATB1002	<input type="checkbox"/>	Manual Seat Tilt by spindle	<input checked="" type="checkbox"/>
BACKREST RECLINER*			
ATB0975	<input type="checkbox"/>	Manual Settable Recliner (0° - 30°, in 3.75° steps) Foldable	<input checked="" type="checkbox"/>
FOOTREST HANGER TYPE (PAIR)*			
ATB0503	<input type="checkbox"/>	Vari F : Angle adjustable with tools 0° - 70°, length 290mm - 460mm	<input checked="" type="checkbox"/>
FOOTPLATE TYPE*			
ATB0512	<input type="checkbox"/>	Fixed Footplate	<input checked="" type="checkbox"/>

CONTROLS*			
ATB2204	<input type="checkbox"/>	REM216 Drive, Actuators & Lights	◆
CONTROL BRACKETRY*			
ATB1207	<input type="checkbox"/>	Standard Bracket RH	○
ATB1208	<input type="checkbox"/>	Standard Bracket LH	○
ATB1707	<input type="checkbox"/>	Swing Away RH Bracket	£108
ATB1708	<input type="checkbox"/>	Swing Away LH Bracket	£108
SHROUD COLOUR*			
ATB0202	<input type="checkbox"/>	Blue	○
ATB0228	<input type="checkbox"/>	Black	○
ATB0220	<input type="checkbox"/>	Champagne	○
ATB0204	<input type="checkbox"/>	Transparent Cherry Red	○
ATB0222	<input type="checkbox"/>	White	○
ATB0216	<input type="checkbox"/>	Olive green	○
FRONT WHEEL*			
ATB0853	<input type="checkbox"/>	3.00 - 6 Grey Puncture Proof (sealed)	◆
REAR WHEEL*			
ATB0773	<input type="checkbox"/>	3.00 - 8 Grey Puncture Proof (sealed)	◆
MOTOR TYPE*			
ATB1909	<input type="checkbox"/>	12 km/h (8mph) (includes LED light system with indicators)	◆
BATTERY*			
ATB1364	<input type="checkbox"/>	70 Ah Gel	◆
SUSPENSION*			
ATB1240	<input type="checkbox"/>	Suspension low (up to 80kg user weight)	○
ATB1241	<input type="checkbox"/>	Suspension middle (from 70kg up to 100kg user weight)	○
ATB1242	<input type="checkbox"/>	Suspension hard (from 90kg up to 150kg user weight)	○
G TRAC*			
ATB1215	<input type="checkbox"/>	G-Trac (ensures the powerchair takes the intended direction as indicated by the user control, without the need for continual correction by the client on uneven or sloping terrain)	£673
CHARGER*			
ATB1343	<input type="checkbox"/>	OptiCharge Charger 8A (200 - 250V)	◆
ANTI TIPPER*			
ATB1205	<input type="checkbox"/>	Anti-tipper (standard version)	◆
TRANSPORTATION KIT*			
ATB1494	<input type="checkbox"/>	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19:2001) (Adequate headrest is compulsory, see options below)	◆
FIXATION BELT*			
ATB1329	<input type="checkbox"/>	Lap Belt with Auto Lock	◆
REAR VIEW MIRROR*			
ATB1415	<input type="checkbox"/>	Rear view Mirror RH - Must be selected with 10 km/h Motor	○
ATB1339	<input type="checkbox"/>	Rear view Mirror LH - Must be selected with 10 km/h Motor	○
TOTAL PRICE			

Non Conforming Options:



For a more complex customer requirement please refer to our Invacare Unique special order form. email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: _____ Email: _____

Phone: _____ Fax: _____

Description of special request (please include relevant option code or part number):