



Yes, you can:

Powerchairs Retail Prescription Form

Tel: 01656 776222 Fax: 01656 776220
 Email: ordersuk@invacare.com
 Online Spares available at: www.invacare.co.uk

JANUARY 2021
 Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

M73 - M74 TDX SP2 NB Express

Standard Features

- Max user weight 150kg (23.5 Stone)
- Crash tested ISO 7176 -Part 19
- Telescopic adjustable seat width 430mm (17") (adjustable 440mm - 480mm)
- Seat depth set at 460mm (adjustable 410mm - 510mm)
- Flex3 adjustable backrest, Matrix Flo-tech Lite cushion
- Standard T-armrests, Detachable with long armpad
- Powered Tilt (0°-20°)
- Manual settable backrest Recline (0°- 30° in 3.75° steps)
- Standard 80° hangers with fixed Footplates
- LiNX REM110 control unit, Standard RH Bracket (can be swapped to LH)
- 6" Black castors, 14" solid Black drive wheels
- Sure Step and Stability Lock
- Matt Black body shrouds, Black rims
- 6 km/h (4mph) motors and 50Ah Batteries, 8Ah charger
- 4 Point tie down system, lapbelt
-
-
- Express shipment! - orders will be shipped same day if ordered before 2pm, next day if ordered after 2pm!

Standard Base price: £6,331



Please make your selection with a cross in the box Free of charge standard feature POA special request

Customer weight = KG

Code	Description	
TDX SP2 NB EXPRESS STOCK		
1650310	TDX SP2 NB Express Spec 3, see standard features, M73	£5,785
1650312	TDX SP2 NB Express Spec 4, see standard features with additions of Rem216, 10km/h (6mph) lights, indicators and mirror. M74	£6,331

TOTAL PRICE	
--------------------	--

Non Conforming Options:



For a more complex customer requirement please refer to our Invacare Unique special order form. email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: _____ Email: _____

Phone: _____ Fax: _____

Description of special request (please include relevant option code or part number):