



Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

JANUARY 2021
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

M75 TDX SP2 Low rider

Standard Features

- Max user weight 160kg (25.2 Stone)
- Crash tested ISO7176-Part 19
- Compact, manoeuvrable, CWD TDX SP2
- Telescopic adjustable seat widths from 380mm - 530mm (15" - 21") with an additional 50mm, (2") by arms
- Seat depths 410mm - 560mm stepless adjustment (16" - 22")
- Matrix Flo-tech Lite cushion
- Without push handles
- Flex 3 Adjustable backrest plates with laguna backrest cushion
- Elan Headrest
- Standard T-armrests, detachable with short armpad
- Seat Lifter with powered Tilt (300mm, 30°) CoG Shift
- Choice of seat height either 420mm or 403mm
- Manual settable Recliner 0° - 30° in 3.75° steps
- Standard 80° Hnagers with Angle and Depth Adjustable footplates(adapt)
- LiNX control REM211 with swing away bracket
- Selection of 10 stunning shroud and Rim colours
- Centre (driving wheels) 14" X 3" Black tyres
- Suspension and sure step technology
- 6km/h (4mph) 4 Pole motors
- 8A (200 - 250V) charger
- 4-Point tie-down system
- Lapbelt with autolock (adjustable both sides)

Standard Base price: £7,323



Q'STRAIT.
QLK-150
Compatible

Please make your selection with a cross in the box ○ Free of charge ◆ standard feature POA special request

Customer weight = **KG**

Code	Description	
BASE PRODUCT*		
ARQ0071	TDX SP2 Low Rider 73ah Batteries 420mm seat to ground height	○
ARQ0072	TDX SP2 Low Rider 60ah Batteries 403mm seat to ground height	○
User Weight*		
ARQ1918	User weight less 100 kg	○
ARQ1919	User weight 101 to 125 kg	○
ARQ1921	User weight 126 to 136 kg	○
ARQ1922	User weight 137 to 150 kg	○
ARQ1923	User weight 151 to 160 kg	○
SEATING TYPE*		
ARQ2100	Modulite Seating Unit	£7,323
SEAT FRAME*		
ARQ2111	Telescopic Seat Frame with metal seat plate	◆
SEAT WIDTH*		
ARQ2122	Set at 430mm (17") (adjustable from 430mm to 480mm (17" - 19))	○
ARQ2123	Set at 480mm(19") (adjustable from 480mm to 530mm (19" - 21"))	○
ARQ2124	Set at 530mm (21") (adjustable from 530mm to 580mm (21" - 23"))	○
SEAT DEPTH RAIL LENGTH*		
ARQ2129	Standard Depth Package (16" - 20", 410mm - 510mm)	○
ARQ2130	Extended Depth Package (20" - 22", 510mm - 560mm) (NC 2134)	○
SEAT DEPTH*		
ARQ2132	410mm (16") (410mm - 510mm)	○
ARQ2134	460mm (18") (410mm - 510mm)	○
ARQ2136	510mm (20") (410mm - 510mm)	○
ARQ2138	560mm (22") (510mm - 560mm)	○
SEAT CUSHION*		
ARQ2142	Matrx Flo-tech Lite	◆
BACKREST SUPPORT*		
ARQ2162	Backposts 540mm (Tall) must be selected for Flex3 adjustable backrest plates	◆
PUSH HANDLES FOR BACKPOSTS*		
ARQ2108	Backposts WITHOUT push handles	◆
BACKREST UPHOLSTERY / CUSHION*		
ARQ2173	Flex3 Adjustable Backrest Plates with Laguna Backrest cushion	◆
UPHOLSTERY TYPE*		
ARQ2151	Black fabric	○
ARQ2152	Black multi stretch PU coated fabric	○
ARQ2155	Black Spacer fabric	○
ARMREST*		
ARQ2191	T-Armrest Standard, Range of adjustment	◆

ARMPAD*		
ARQ2195	<input type="checkbox"/>	Short Armpad (325mm) ◆
HEADREST / NECKREST*		
ARQ2324	<input type="checkbox"/>	Elan Standard Pad 355mm (10") wide ◆
SEAT LIFTER / SEAT TILT*		
ARQ1094	<input type="checkbox"/>	Seat lifter with Electrical Tilt ◆
BACKREST RECLINER*		
ARQ0975	<input type="checkbox"/>	Manual Settable Recliner ◆
FOOTREST HANGER TYPE*		
ARQ0526	<input type="checkbox"/>	Side Mounted : Standard 80° hangers ◆
FOOTPLATE TYPE*		
ARQ0513	<input type="checkbox"/>	Angle and depth adjustable Footplate (Adapt) ◆
CONTROLS*		
ARQ2203	<input type="checkbox"/>	REM211 Drive & Actuators ◆
CONTROL BRACKETRY*		
ARQ1707	<input type="checkbox"/>	Swing Away RH Bracket ○
ARQ1708	<input type="checkbox"/>	Swing Away LH Bracket ○
SHROUD COLOUR*		
ARQ0228	<input type="checkbox"/>	Wet Black ○
ARQ2601	<input type="checkbox"/>	SpaceStation Silver ○
ARQ2602	<input type="checkbox"/>	Island Blue ○
ARQ2603	<input type="checkbox"/>	Grape JellyBean (Purple) ○
ARQ2604	<input type="checkbox"/>	Monster Green ○
ARQ2605	<input type="checkbox"/>	Tangy Orange ○
ARQ2606	<input type="checkbox"/>	WhiteOut ○
ARQ2607	<input type="checkbox"/>	Rockstar Red ○
ARQ2608	<input type="checkbox"/>	Lights Out Black - Matt ○
ARQ2609	<input type="checkbox"/>	Invacare Blue ○
RIM COLOUR*		
ARQ2611	<input type="checkbox"/>	SpaceStation Silver ○
ARQ2612	<input type="checkbox"/>	Island Blue ○
ARQ2613	<input type="checkbox"/>	Grape JellyBean (Purple) ○
ARQ2614	<input type="checkbox"/>	Monster Green ○
ARQ2615	<input type="checkbox"/>	Tangy Orange ○
ARQ2616	<input type="checkbox"/>	WhiteOut ○
ARQ2617	<input type="checkbox"/>	Rockstar Red ○
ARQ2618	<input type="checkbox"/>	Lights Out Black - Matt ○
ARQ2619	<input type="checkbox"/>	Invacare Blue ○
ARQ2620	<input type="checkbox"/>	Wet Black ○
ARQ2621	<input type="checkbox"/>	Without Rim insert ○
FORKS*		
ARQ2632	<input type="checkbox"/>	Double sided fork ◆
CENTRE WHEELS (DRIVING WHEELS)*		
ARQ0733	<input type="checkbox"/>	355mm x 76mm (14" x 3") Black Tyres with Gel Foam Inserts ◆
MOTOR TYPE*		
ARQ1282	<input type="checkbox"/>	6 km/h (4mph) (4 pole) ◆

CHARGER*		
ARQ1343	<input type="checkbox"/>	OptiCharge Charger 8A (200 - 250V) ◆
LIGHTS*		
ARQ1254	<input type="checkbox"/>	Without Light system ◆
TRANSPORTATION KIT*		
ARQ1494	<input type="checkbox"/>	4-Point tie-down ◆
FIXATION BELT*		
ARQ1329	<input type="checkbox"/>	Lapbelt (adjustable both sides) ○
ARQ2113	<input type="checkbox"/>	Lapbelt (adjustable both sides) ◆

TOTAL PRICE	
--------------------	--

Non Conforming Options:

ARQ2134 460mm (18") (410mm - 510mm)



For a more complex customer requirement please refer to our Invacare Unique special order form. email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: _____ Email: _____

Phone: _____ Fax: _____

Description of special request (please include relevant option code or part number):