



Yes, you can:

Powerchairs Retail Prescription Form

Tel: 01656 776222 Fax: 01656 776220
 Email: ordersuk@invacare.com
 Online Spares available at: www.invacare.co.uk

JANUARY 2021
 Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

M77 TDX SP2 HD

Standard Features

- Max user weight 180kg (28 stone)
- Crash tested ISO 7176 -Part 19
- Modulite Seating System
- Matrix Flo-tech Cushion
- Backrest frame, detachable, height 480mm
- Flex 3 adjustable Backrest plates with Laguna backrest cushion
- Standard T-armrests, detachable with short armpad
- Headrest, height, depth and sideways adjustable
- Manual settable seat Tilt 0° - 20°
- Manual settable backrest Recline 0° - 30° in 3.75° steps
- 80° hangers with fixed footplate
- LiNX control REM116 standard bracket
- Selection of 10 stunning shroud and Rim colours
- Centre (driving wheels) 14" x 3" Black tyres
- 8km/h (5mph) 4 pole motors and 73.5 Ah batteries
- 8A (200 - 250V) charger
- 4-Point tie-down system (Q'Straint QLK-150 compatible)
- Light system with indicators, LED headlights and Mirror
- Lapbelt with autolock (adjustable both sides)

Standard Base price: £7,767



**Q'STRAIN
QLK-150
Compatible**

Please make your selection with a cross in the box ○ Free of charge ◆ standard feature POA special request

Customer weight = KG

Code	Description	
CHAIR TYPE*		
ARJ0000	TDX SP2 Heavy Duty	◆
User Weight*		
ARJ1918	User weight less 100 kg	○
ARJ1919	User weight 101 to 125 kg	○
ARJ1921	User weight 126 to 136 kg	○
ARJ1922	User weight 137 to 150 kg	○
ARJ1923	User weight 151 to 160 kg	○
ARJ1924	User weight 161 to 180 kg	○

BASE TYPE*		
ARJ2100	<input type="checkbox"/>	Modulite Seating unit £7,767
SEAT WIDTH*		
ARJ2121	<input type="checkbox"/>	Set at 380mm (15") (adjustable from 380mm to 430mm (15" - 17") with stepless adjustment on armrests, legrests and backrest) (NC 2138 + 2142) ○
ARJ2122	<input type="checkbox"/>	Set at 430mm (17") (adjustable from 430mm to 480mm (17" - 19") with stepless adjustment on armrests, legrests and backrest) (NC 2138 + 2142) ○
ARJ2123	<input type="checkbox"/>	Set at 480mm (19") (adjustable from 480mm to 530mm (19" - 21") with stepless adjustment on armrests, legrests and backrest) (NC 2138 + 2142) ○
ARJ2124	<input type="checkbox"/>	Set at 530mm (21") (adjustable from 530mm to 580mm (21" - 23") with stepless adjustment on armrests, legrests and backrest) ○
SEAT DEPTH RAIL LENGTH*		
ARJ2129	<input type="checkbox"/>	Standard Depth Package (16" - 20", 410mm - 510mm simple stepless adjustment from the side) (NC 2138) ○
ARJ2130	<input type="checkbox"/>	Extended Depth Package (20" - 22", 510mm - 560mm simple stepless adjustment from the side) (NC 2132, 2134) ○
SEAT DEPTH SETTING*		
ARJ2132	<input type="checkbox"/>	410mm (16") (410mm - 510mm simple stepless adjustment from the side) ○
ARJ2134	<input type="checkbox"/>	460mm (18") (410mm - 510mm simple stepless adjustment from the side) ○
ARJ2136	<input type="checkbox"/>	510mm (20") (410mm - 510mm simple stepless adjustment from the side) ○
ARJ2138	<input type="checkbox"/>	560mm (22") (510mm - 560mm simple stepless adjustment from the side) (Cushion (option 2142) is only available with Seat Width (option 2124) 530mm) ○
SEAT CUSHION*		
ARJ2142	<input type="checkbox"/>	Matrx Flo-tech Lite ◆
BACKREST SUPPORT*		
ARJ2162	<input type="checkbox"/>	Backposts 540mm (Tall) must be selected for Flex3 adjustable backrest plates ◆
PUSH HANDLES FOR BACKPOSTS*		
ARJ2107	<input type="checkbox"/>	Backposts WITH push handles ◆
BACKREST UPHOLSTERY / CUSHION*		
ARJ2173	<input type="checkbox"/>	Flex3 adjustable Backrest Plates with Laguna Backrest cushion ◆
UPHOLSTERY COLOUR*		
ARJ2151	<input type="checkbox"/>	Black fabric ○
ARJ2152	<input type="checkbox"/>	Black Black multi stretch PU coated fabric Upholstery ○
ARJ2155	<input type="checkbox"/>	Black Spacer fabric ○
ARMREST*		
ARJ2191	<input type="checkbox"/>	T-Armrest Standard, Range of adjustment from seat plate 245mm to 315mm ◆
ARMPAD*		
ARJ2195	<input type="checkbox"/>	Short Armpad (325mm) ◆
HEADREST / NECKREST - Adequate headrest option must be selected*		
ARJ2324	<input type="checkbox"/>	Elan Standard Pad 355mm (10") wide with standard cover complete with Elan Headrest Hardware ◆
SEAT LIFTER / SEAT TILT*		
ARJ1095	<input type="checkbox"/>	Powered Basic with increased Seat Height 0° to 20° Seat height 465mm (Increase in seat height to facilitate rearward adjustment of the seating) ◆

SEAT HEIGHT SETTING FOR 20 DEGREES SEAT TILT		
ARJ0354	<input type="checkbox"/>	Seat height set at 450mm (17.75") (from floor to front of seat plate at 0° tilt) ○
ARJ0364	<input type="checkbox"/>	Seat height set at 475mm (18.75") (from floor to front of seat plate at 0° tilt) ○
BACKREST RECLINER*		
ARJ0975	<input type="checkbox"/>	Manual Settable Recliner (0° - 30°, in 3.75° steps) Foldable ◆
FOOTREST HANGER TYPE*		
ARJ0503	<input type="checkbox"/>	Side Mounted : Vari F : Angle adjustable with tools 0° - 70°, length 290mm - 460mm ◆
FOOTPLATE TYPE*		
ARJ0512	<input type="checkbox"/>	Fixed Footplate ◆
CONTROLS*		
ARJ2204	<input type="checkbox"/>	REM216 Drive, Actuators & Lights ◆
CONTROL BRACKETRY*		
ARJ1207	<input type="checkbox"/>	Standard Bracket RH ○
ARJ1208	<input type="checkbox"/>	Standard Bracket LH ○
ARJ1707	<input type="checkbox"/>	Swing Away RH Bracket £108
ARJ1708	<input type="checkbox"/>	Swing Away LH Bracket £108
SHROUD COLOUR*		
ARJ0228	<input type="checkbox"/>	Wet Black ○
ARJ2601	<input type="checkbox"/>	SpaceStation Silver ○
ARJ2602	<input type="checkbox"/>	Island Blue ○
ARJ2603	<input type="checkbox"/>	Grape JellyBean (Purple) ○
ARJ2604	<input type="checkbox"/>	Monster Green ○
ARJ2605	<input type="checkbox"/>	Tangy Orange ○
ARJ2606	<input type="checkbox"/>	WhiteOut ○
ARJ2607	<input type="checkbox"/>	Rockstar Red ○
ARJ2608	<input type="checkbox"/>	Lights Out Black - Matt ○
ARJ2609	<input type="checkbox"/>	Invacare Blue ○
RIM COLOUR*		
ARJ2611	<input type="checkbox"/>	SpaceStation Silver ○
ARJ2612	<input type="checkbox"/>	Island Blue ○
ARJ2613	<input type="checkbox"/>	Grape JellyBean (Purple) ○
ARJ2614	<input type="checkbox"/>	Monster Green ○
ARJ2615	<input type="checkbox"/>	Tangy Orange ○
ARJ2616	<input type="checkbox"/>	WhiteOut ○
ARJ2617	<input type="checkbox"/>	Rockstar Red ○
ARJ2618	<input type="checkbox"/>	Lights Out Black - Matt ○
ARJ2619	<input type="checkbox"/>	Invacare Blue ○
ARJ2620	<input type="checkbox"/>	Wet Black ○
ARJ2621	<input type="checkbox"/>	Without Rim insert ○
FORKS*		
ARJ2632	<input type="checkbox"/>	Double sided fork ○
CENTRE WHEELS (DRIVING WHEELS)*		
ARJ0733	<input type="checkbox"/>	355mm x 76mm (14" x 3") Black Tyres with Gel Foam Inserts ○

MOTOR TYPE*		
ARJ1284	<input type="checkbox"/>	8 km/h (5mph) (4 pole) ◆
BATTERY*		
ARJ1260	<input type="checkbox"/>	73.5 Ah (Gel) ◆
CHARGER*		
ARJ1343	<input type="checkbox"/>	OptiCharge Charger 8A (200 - 250V) ◆
LIGHTS*		
ARJ1437	<input type="checkbox"/>	Light system with Indicators & LED Headlights ◆
TRANSPORTATION KIT*		
ARJ1450	<input type="checkbox"/>	Preparation for DAHL Docking station: Docking station for car implementation, please order docking station separately (Adequate headrest is compulsory, see options below) £478
ARJ1494	<input type="checkbox"/>	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19:2001) (Adequate headrest is compulsory, see options below) ◆
FIXATION BELT*		
ARJ2113	<input type="checkbox"/>	Lapbelt (adjustable both sides) adjustable mounting on seat rail extrusion ○
REAR VIEW MIRROR*		
ARJ1415	<input type="checkbox"/>	Rear view Mirror RH - Must be selected with 10 km/h Motor £34
ARJ1339	<input type="checkbox"/>	Rear view Mirror LH - Must be selected with 10 km/h Motor £34
		TOTAL PRICE

Non Conforming Options:

ARJ2132	410mm (16") (410mm - 510mm simple stepless adjustment from the side)
ARJ2134	460mm (18") (410mm - 510mm simple stepless adjustment from the side)
ARJ2138	560mm (22") (510mm - 560mm simple stepless adjustment from the side) (Cushion (option 2142) is only available with Seat Width (option 2124) 530mm)



For a more complex customer requirement please refer to our Invacare Unique special order form. email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: _____ Email: _____

Phone: _____ Fax: _____

Description of special request (please include relevant option code or part number):