



**Yes, you can:**

# Powerchairs Retail Prescription Form

Tel: 01656 776222 Fax: 01656 776220  
 Email: ordersuk@invacare.com  
 Online Spares available at: www.invacare.co.uk

JANUARY 2021  
 Customer Ref:

<b>Customer Ref*</b>	<b>Account No*</b>
<b>Contact No*</b>	<b>Order No*</b>
<b>Delivery Address*</b>	<b>Order Qty*</b>
	<b>Order date*</b>

## M80 TDX SP2 HD

### Standard Features

- Max user weight 150kg (28 stone)
- Crash tested ISO 7176 -Part 19
- Modulite Seating System
- Matrix Flo-tech Cushion
- Backrest frame, detachable, height 480mm
- Flex 3 adjustable Backrest plates with Laguna backrest cushion
- Standard T-armrests, detachable with short armpad
- Headrest, height, depth and sideways adjustable
- Manual settable seat Tilt 0° - 20°
- Manual settable backrest Recline 0° - 30° in 3.75° steps
- 80° hangers with fixed footplate
- LiNX control REM116 standard bracket
- Selection of 10 stunning shroud and Rim colours
- Centre (driving wheels) 14" x 3" Black tyres
- 10km/h (6mph) 4 pole motors and 73.5 Ah batteries
- 8A (200 - 250V) charger
- 4-Point tie-down system (Q'Straint QLK-150 compatible)
- Light system with indicators, LED headlights and Mirror
- Lapbelt with autolock (adjustable both sides)

**Standard Base price: £7,767**



**Q'STRAIN**  
**QLK-150**  
 Compatible

Please make your selection with a cross in the box    ○ Free of charge    ◆ standard feature    POA special request

Customer weight =                      KG

Code	Description	
<b>CHAIR TYPE*</b>		
ARJ0000	<input type="checkbox"/> TDX SP2 Heavy Duty	◆
<b>User Weight*</b>		
ARJ1918	<input type="checkbox"/> User weight less 100 kg	○
ARJ1919	<input type="checkbox"/> User weight 101 to 125 kg	○
ARJ1921	<input type="checkbox"/> User weight 126 to 136 kg	○
ARJ1922	<input type="checkbox"/> User weight 137 to 150 kg	○
<b>BASE TYPE*</b>		
ARJ2100	<input type="checkbox"/> Modulite Seating unit	£7,767

<b>SEAT WIDTH*</b>			
ARJ2121	<input type="checkbox"/>	Set at 380mm (15") (adjustable from 380mm to 430mm (15" - 17") with stepless adjustment on armrests, legrests and backrest) (NC 2138 + 2142)	<input type="radio"/>
ARJ2122	<input type="checkbox"/>	Set at 430mm (17") (adjustable from 430mm to 480mm (17" - 19) with stepless adjustment on armrests, legrests and backrest) (NC 2138 + 2142)	<input type="radio"/>
ARJ2123	<input type="checkbox"/>	Set at 480mm (19") (adjustable from 480mm to 530mm (19" - 21") with stepless adjustment on armrests, legrests and backrest) (NC 2138 + 2142)	<input type="radio"/>
ARJ2124	<input type="checkbox"/>	Set at 530mm (21") (adjustable from 530mm to 580mm (21" - 23") with stepless adjustment on armrests, legrests and backrest)	<input type="radio"/>
<b>SEAT DEPTH RAIL LENGTH*</b>			
ARJ2129	<input type="checkbox"/>	Standard Depth Package (16" - 20", 410mm - 510mm simple stepless adjustment from the side) (NC 2138)	<input type="radio"/>
ARJ2130	<input type="checkbox"/>	Extended Depth Package (20" - 22", 510mm - 560mm simple stepless adjustment from the side) (NC 2132, 2134)	<input type="radio"/>
<b>SEAT DEPTH SETTING*</b>			
ARJ2132	<input type="checkbox"/>	410mm (16") (410mm - 510mm simple stepless adjustment from the side)	<input type="radio"/>
ARJ2134	<input type="checkbox"/>	460mm (18") (410mm - 510mm simple stepless adjustment from the side)	<input type="radio"/>
ARJ2136	<input type="checkbox"/>	510mm (20") (410mm - 510mm simple stepless adjustment from the side)	<input type="radio"/>
ARJ2138	<input type="checkbox"/>	560mm (22") (510mm - 560mm simple stepless adjustment from the side) (Cushion (option 2142) is only available with Seat Width (option 2124) 530mm)	<input type="radio"/>
<b>SEAT CUSHION*</b>			
ARJ2142	<input type="checkbox"/>	Matrx Flo-tech Lite	<input type="checkbox"/>
<b>BACKREST SUPPORT*</b>			
ARJ2162	<input type="checkbox"/>	Backposts 540mm (Tall) must be selected for Flex3 adjustable backrest plates	<input type="checkbox"/>
<b>PUSH HANDLES FOR BACKPOSTS*</b>			
ARJ2107	<input type="checkbox"/>	Backposts WITH push handles	<input type="checkbox"/>
<b>BACKREST UPHOLSTERY / CUSHION*</b>			
ARJ2173	<input type="checkbox"/>	Flex3 adjustable Backrest Plates with Laguna Backrest cushion	<input type="checkbox"/>
<b>UPHOLSTERY COLOUR*</b>			
ARJ2151	<input type="checkbox"/>	Black fabric	<input type="radio"/>
ARJ2152	<input type="checkbox"/>	Black Black multi stretch PU coated fabric Upholstery	<input type="radio"/>
ARJ2155	<input type="checkbox"/>	Black Spacer fabric	<input type="radio"/>
<b>ARMREST*</b>			
ARJ2191	<input type="checkbox"/>	T-Armrest Standard, Range of adjustment from seat plate 245mm to 315mm	<input type="checkbox"/>
<b>ARMPAD*</b>			
ARJ2195	<input type="checkbox"/>	Short Armpad (325mm)	<input type="checkbox"/>
<b>HEADREST / NECKREST - Adequate headrest option must be selected*</b>			
ARJ2324	<input type="checkbox"/>	Elan Standard Pad 355mm (10") wide with standard cover complete with Elan Headrest Hardware	<input type="checkbox"/>
<b>SEAT LIFTER / SEAT TILT*</b>			
ARJ1094	<input type="checkbox"/>	Seat lifter with Electrical Tilt (300mm, 30°) with CoG shift Seat Height 450mm - 480mm (Also choose option 1265)	<input type="checkbox"/>
<b>SEAT HEIGHT SETTING FOR 20 DEGREES SEAT TILT</b>			
ARZ1003	<input type="checkbox"/>	Seat height set at 460mm (18") (from floor to front of seat plate at 0° tilt)	<input type="checkbox"/>

<b>BACKREST RECLINER*</b>			
ARJ0975	<input type="checkbox"/>	Manual Settable Recliner (0° - 30°, in 3.75° steps) Foldable	◆
<b>FOOTREST HANGER TYPE*</b>			
ARJ0556	<input type="checkbox"/>	Centre Mounted: Manual : Hanger, centre mounted, angle settable (0° - 90°) : knee to heel length 390 - 510mm, (including calf pad and angle adjustable footplate)	◆
<b>FOOTPLATE TYPE*</b>			
ARJ0559	<input type="checkbox"/>	Angle Adjustable Footplate for Centre Mounted Hanger	◆
<b>CONTROLS*</b>			
ARJ2204	<input type="checkbox"/>	REM216 Drive, Actuators & Lights	◆
<b>CONTROL BRACKETRY*</b>			
ARJ1207	<input type="checkbox"/>	Standard Bracket RH	○
ARJ1208	<input type="checkbox"/>	Standard Bracket LH	○
ARJ1707	<input type="checkbox"/>	Swing Away RH Bracket	£108
ARJ1708	<input type="checkbox"/>	Swing Away LH Bracket	£108
<b>SHROUD COLOUR*</b>			
ARJ0228	<input type="checkbox"/>	Wet Black	○
ARJ2601	<input type="checkbox"/>	SpaceStation Silver	○
ARJ2602	<input type="checkbox"/>	Island Blue	○
ARJ2603	<input type="checkbox"/>	Grape JellyBean (Purple)	○
ARJ2604	<input type="checkbox"/>	Monster Green	○
ARJ2605	<input type="checkbox"/>	Tangy Orange	○
ARJ2606	<input type="checkbox"/>	WhiteOut	○
ARJ2607	<input type="checkbox"/>	Rockstar Red	○
ARJ2608	<input type="checkbox"/>	Lights Out Black - Matt	○
ARJ2609	<input type="checkbox"/>	Invacare Blue	○
<b>RIM COLOUR*</b>			
ARJ2611	<input type="checkbox"/>	SpaceStation Silver	○
ARJ2612	<input type="checkbox"/>	Island Blue	○
ARJ2613	<input type="checkbox"/>	Grape JellyBean (Purple)	○
ARJ2614	<input type="checkbox"/>	Monster Green	○
ARJ2615	<input type="checkbox"/>	Tangy Orange	○
ARJ2616	<input type="checkbox"/>	WhiteOut	○
ARJ2617	<input type="checkbox"/>	Rockstar Red	○
ARJ2618	<input type="checkbox"/>	Lights Out Black - Matt	○
ARJ2619	<input type="checkbox"/>	Invacare Blue	○
ARJ2620	<input type="checkbox"/>	Wet Black	○
ARJ2621	<input type="checkbox"/>	Without Rim insert	○
<b>FORKS*</b>			
ARJ2632	<input type="checkbox"/>	Double sided fork	○
<b>CENTRE WHEELS (DRIVING WHEELS)*</b>			
ARJ0733	<input type="checkbox"/>	355mm x 76mm (14" x 3") Black Tyres with Gel Foam Inserts	○
<b>MOTOR TYPE*</b>			
ARJ1224	<input type="checkbox"/>	12 km/h (8mph) (includes light system and indicators) G-Trac recommended	◆

<b>BATTERY*</b>		
ARJ1260	<input type="checkbox"/>	73.5 Ah (Gel) ◆
<b>CHARGER*</b>		
ARJ1343	<input type="checkbox"/>	OptiCharge Charger 8A (200 - 250V) ◆
<b>TRANSPORTATION KIT*</b>		
ARJ1450	<input type="checkbox"/>	Preparation for DAHL Docking station: Docking station for car implementation, please order docking station separately (Adequate headrest is compulsory, see options below) £478
ARJ1494	<input type="checkbox"/>	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19:2001) (Adequate headrest is compulsory, see options below) ◆
<b>FIXATION BELT*</b>		
ARJ2113	<input type="checkbox"/>	Lapbelt (adjustable both sides) adjustable mounting on seat rail extrusion ○
<b>REAR VIEW MIRROR*</b>		
ARJ1415	<input type="checkbox"/>	Rear view Mirror RH - Must be selected with 10 km/h Motor ○
ARJ1339	<input type="checkbox"/>	Rear view Mirror LH - Must be selected with 10 km/h Motor ○

<b>TOTAL PRICE</b>	<input type="text"/>
--------------------	----------------------

**Non Conforming Options:**

ARJ2132	410mm (16") (410mm - 510mm simple stepless adjustment from the side)
ARJ2134	460mm (18") (410mm - 510mm simple stepless adjustment from the side)
ARJ2138	560mm (22") (510mm - 560mm simple stepless adjustment from the side) (Cushion (option 2142) is only available with Seat Width (option 2124) 530mm)



For a more complex customer requirement please refer to our Invacare Unique special order form. email [specialsuk@invacare.com](mailto:specialsuk@invacare.com) or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: \_\_\_\_\_ Email: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Description of special request (please include relevant option code or part number):