



Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

JANUARY 2021
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

M81 Storm 4 Xplore

Standard Features

- Max user weight 150Kg (23.5 stone)
- Crash tested ISO7176-Part 19
- Modulite highly adjustable seating system
- Telescopic adjustable Seat Widths from 330mm to 530mm (13" - 21") with an additional 50mm, (2") by arms
- Seat depths from 410mm, 510mm (16" - 20")
- Flex3 adjustable backrest
- Matrix Flo-tech Lite cushion
- Backposts 540mm (Tall)
- Standard T-armrests, detachable with short armpads
- Manual settable seat Tilt
- Manual settable backrest Recline
- Angle adjustable legrest hangers (VARI F) with Fixed footplates
- LINX REM216 control with standard bracket
- Puncture proof (sealed) drive wheels and front castors
- 12 km/h (8mph) includes LED light system with indicators) Mirror
- 73.5 Ah Gel batteries
- 8A (200 - 250V) Charger
- LED Lights and Indicators
- Lapbelt (adjustable both sides)

Standard Base price: £7,160



Please make your selection with a cross in the box Free of charge standard feature POA special request

Customer weight = KG

Code	Description	
CHAIR TYPE*		
ATE0017	Storm 4 X-plore	£7,160
ATZ0017	Special Storm 4 X-plore	POA
SEATING TYPE*		
ATE2100	Invacare Modulite 'Flex3'	◆
SEAT FRAME*		
ATE2111	Telescopic Seat Frame with metal seat plate	◆

SEAT WIDTH*			
ATEXXXX	<input type="checkbox"/>	Seat width 330mm (13")	<input type="radio"/>
ATE2121	<input type="checkbox"/>	Set at 380mm	<input type="radio"/>
ATE2122	<input type="checkbox"/>	Set at 430mm	<input type="radio"/>
ATE2123	<input type="checkbox"/>	Set at 480mm	<input type="radio"/>
ATE2124	<input type="checkbox"/>	Set at 530mm	<input type="radio"/>
SEAT DEPTH*			
ATE2132	<input type="checkbox"/>	410mm (16")	<input type="radio"/>
ATE2134	<input type="checkbox"/>	460mm (18")	<input type="radio"/>
ATE2136	<input type="checkbox"/>	510mm (20")	<input type="radio"/>
SEAT CUSHION*			
ATE2142	<input type="checkbox"/>	Matrx Flo-tech Lite	<input checked="" type="radio"/>
BACKREST SUPPORT*			
ATE2162	<input type="checkbox"/>	Backposts 540mm (Tall)	<input checked="" type="radio"/>
PUSH HANDLES FOR BACKPOSTS*			
ATE2107	<input type="checkbox"/>	Backposts WITH push handles	<input checked="" type="radio"/>
BACKREST UPHOLSTREY / CUSHION*			
ATE2173	<input type="checkbox"/>	Flex3 adjustable Backrest Plates with Laguna Backrest cushion	<input checked="" type="radio"/>
UPHOLSTERY COLOUR*			
ATE2151	<input type="checkbox"/>	Black fabric	<input type="radio"/>
ATE2152	<input type="checkbox"/>	Black multi stretch PU coated fabric	<input type="radio"/>
ATE2155	<input type="checkbox"/>	Black Spacer fabric	<input type="radio"/>
ARMREST*			
ATE2191	<input type="checkbox"/>	T-Armrest Standard	<input checked="" type="radio"/>
ARMPAD*			
ATE2195	<input type="checkbox"/>	Short Armpad (325mm)	<input checked="" type="radio"/>
HEADREST / NECKREST - Adequate headrest option must be selected*			
ATE1495	<input type="checkbox"/>	Push bar	<input checked="" type="radio"/>
SEAT HEIGHT SETTING & SEAT RISER*			
ATE1023	<input type="checkbox"/>	Fixed Seat Height 450mm	<input checked="" type="radio"/>
ATE1025	<input type="checkbox"/>	Fixed Seat Height 480mm	<input type="radio"/>
SEAT TILT*			
ATE1002	<input type="checkbox"/>	Manual Seat Tilt by spindle	<input checked="" type="radio"/>
BACKREST RECLINER*			
ATE0975	<input type="checkbox"/>	Manual Settable Recliner (0° - 30°, in 3.75° steps) Foldable	<input checked="" type="radio"/>
FOOTREST HANGER TYPE (PAIR)*			
ATE0503	<input type="checkbox"/>	Vari F : Angle adjustable with tools 0° - 70°, length 290mm - 460mm	<input checked="" type="radio"/>
FOOTPLATE TYPE*			
ATE0512	<input type="checkbox"/>	Fixed Footplate	<input checked="" type="radio"/>
CONTROLS*			
ATE2204	<input type="checkbox"/>	REM216 Drive, Actuators & Lights	<input checked="" type="radio"/>

CONTROL BRACKETRY*			
ATE1207	<input type="checkbox"/>	Standard Bracket RH	○
ATE1208	<input type="checkbox"/>	Standard Bracket LH	○
ATE1707	<input type="checkbox"/>	Swing Away RH Bracket	£108
ATE1708	<input type="checkbox"/>	Swing Away LH Bracket	£108
SHROUD COLOUR*			
ATE0202	<input type="checkbox"/>	Blue	○
ATE0228	<input type="checkbox"/>	Black	○
ATE0220	<input type="checkbox"/>	Champagne	○
ATE0204	<input type="checkbox"/>	Transparent Cherry Red	○
ATE0222	<input type="checkbox"/>	White	○
ATE0216	<input type="checkbox"/>	Olive green	○
CASTOR WHEEL*			
ATE0853	<input type="checkbox"/>	3.00 - 6 Grey Puncture Proof (sealed)	◆
DRIVE WHEEL*			
ATE0773	<input type="checkbox"/>	3.00 - 8 Grey Puncture Proof (sealed)	◆
MOTOR TYPE*			
ATE1909	<input type="checkbox"/>	12 km/h (8mph) includes LED light system with indicators	◆
BATTERY*			
ATE1364	<input type="checkbox"/>	73.5 Ah GEL	◆
REAR SUSPENSION*			
ATE1233	<input type="checkbox"/>	Hitec spring X-plore	○
FRONT SUSPENSION*			
ATE1238	<input type="checkbox"/>	Front suspension (suggested for up to 100kg)	○
ATE1239	<input type="checkbox"/>	Front suspension (suggested for above 100kg)	○
CHARGER*			
ATE1343	<input type="checkbox"/>	OptiCharge Charger 8A (200 - 250V)	◆
ANTI-TIPPERS*			
ATE1205	<input type="checkbox"/>	Anti-tipper (standard version)	◆
TRANSPORTATION*			
ATE1494	<input type="checkbox"/>	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19:2001) (Adequate headrest is compulsory).	◆
FIXATION BELT*			
ATE1329	<input type="checkbox"/>	Lapbelt (adjustable both sides) fitted to base of backrest	◆
REAR VIEW MIRROR*			
ATE1415	<input type="checkbox"/>	Rear view Mirror RH - Must be selected with 10 km/h Motor	○
ATE1339	<input type="checkbox"/>	Rear view Mirror LH - Must be selected with 10 km/h Motor	○
TOTAL PRICE			

Non Conforming Options:



For a more complex customer requirement please refer to our Invacare Unique special order form. email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: _____ Email: _____

Phone: _____ Fax: _____

Description of special request (please include relevant option code or part number):