



**Yes, you can:**

# Powerchairs Retail Prescription Form

Tel: 01656 776222 Fax: 01656 776220  
 Email: ordersuk@invacare.com  
 Online Spares available at: www.invacare.co.uk

JANUARY 2021  
 Customer Ref:

<b>Customer Ref*</b>	<b>Account No*</b>
<b>Contact No*</b>	<b>Order No*</b>
<b>Delivery Address*</b>	<b>Order Qty*</b>
	<b>Order date*</b>

## M95 AVIVA RX20 Modulite

### Standard Features

- Max user weight 136kg (21 Stone)
- Crash tested ISO 7176 -Part 19
- Modulite highly adjustable seating system
- Telescopic adjustable seat widths from 330mm - 530mm (13" - 21") with an additional 50mm, (2") by arms
- Seat depths 410mm - 510mm stepless adjustment (16" - 20")
- Tension Adjustable Upholstery, nylon overlay short = 480 mm / 7 straps
- Standard T-Armrests, detachable with short armpad
- Neckrest depth adjustable
- Powered "fixed pivot" Tilt 20°
- Choice of seat height's
- Powered Recliner 30°
- 80° hangers with fixed footplate
- Remote LiNX DLX-REM211 with standard bracket
- Choice of 5 Shroud colours, 5 Rim insert colours
- 8" (200 mm) x 2.50-4 (65 mm) Castor, puncture proof solid tyre
- 14" (355 mm) x 3" (75 mm), Solid inlay black drive wheels
- CTC Suspension Control - Traction - Comfort
- 6 km/h 4 pole Durawatt motors
- 52Ah batteries
- 8A (200 - 250V) charger
- Anti tippers
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)

**Standard Base price: £5,953**



Please make your selection with a cross in the box    ○ Free of charge    ◆ standard feature    POA special request

Customer weight =                      KG

**Code**                      **Description**

**BASE MODEL\***

JBC0000	<input type="checkbox"/>	AVIVA RX 20 Base Model Max user weight up to 136 kg	£5,953
---------	--------------------------	---	--------

**CHASSIS WIDTH\***

JBL2810	<input type="checkbox"/>	Chassis Narrow (560 mm/22")	◆
---------	--------------------------	-----------------------------	---

**SEATING TYPE\***

JBC2100	<input type="checkbox"/>	Invacare® Modulite Seat Unit	◆
---------	--------------------------	------------------------------	---

<b>SEAT FRAME*</b>			
JBC2111	<input type="checkbox"/>	Telescopic Seat Frame and Plate with Side Rails	◆
<b>SEAT WIDTH*</b>			
JBC2121	<input type="checkbox"/>	380 mm Seat Plate	○
JBC2122	<input type="checkbox"/>	430 mm Seat Plate	○
JBC2123	<input type="checkbox"/>	480 mm Seat Plate	○
JBC2124	<input type="checkbox"/>	530 mm Seat Plate (NC 2106)	○
<b>Body Fit Kit</b>			
JBC2106	<input type="checkbox"/>	Body fit Kit (NC 2171)	○
<b>SEAT DEPTH*</b>			
JBC2132	<input type="checkbox"/>	Seat Depth 410 mm	○
JBC2134	<input type="checkbox"/>	Seat Depth 460 mm	○
JBC2136	<input type="checkbox"/>	Seat Depth 510 mm	○
<b>SEAT CUSHION*</b>			
JBC2142	<input type="checkbox"/>	Matrx Flo-tech Cushion Flo-tech Lite	◆
<b>BACKREST SUPPORT*</b>			
JBC2163	<input type="checkbox"/>	Backposts short (480 mm)	◆
JBC2162	<input type="checkbox"/>	Backposts tall (540 mm)	○
<b>PUSH HANDLES FOR BACKPOSTS</b>			
JBC2108	<input type="checkbox"/>	Backposts without push handles	◆
<b>BACKREST UPHOLSTERY / CUSHION*</b>			
JBC2171	<input type="checkbox"/>	Tension Adjustable Upholstery (NC 2106)	◆
<b>UPHOLSTERY TYPE*</b>			
JBC2151	<input type="checkbox"/>	Black Fabric Upholstery	◆
<b>ARMREST (Pair)*</b>			
JBC2191	<input type="checkbox"/>	T-Arm detachable, Standard	◆
<b>ARMPAD (Pair)*</b>			
JBC2195	<input type="checkbox"/>	Short Armpad (325 mm)	◆
<b>HEADREST / NECKREST - Adequate headrest option must be selected*</b>			
JBC1350	<input type="checkbox"/>	Neckrest depth adjustable	◆
<b>SEAT LIFTER / SEAT TILT*</b>			
JBC1092	<input type="checkbox"/>	Powered "fixed pivot" Tilt 20°	◆
<b>SEAT TO GROUND HEIGHT*</b>			
JBC2857	<input type="checkbox"/>	Seat-to-Floor Height 435 mm	○
JBC2858	<input type="checkbox"/>	Seat-to-Floor Height 460 mm	○
JBC2859	<input type="checkbox"/>	Seat-to-Floor Height 485 mm	○
<b>BACKREST RECLINER*</b>			
JBC0977	<input type="checkbox"/>	Powered Recliner 30°	◆
<b>FOOTREST HANGER TYPE*</b>			
JBC0526	<input type="checkbox"/>	Standard hanger 80° pair Lower leg length: 350 - 450 mm	◆
<b>FOOTPLATE TYPE*</b>			
JBC0512	<input type="checkbox"/>	Fixed Footplate pair	◆
<b>CONTROLS*</b>			
JBC2203	<input type="checkbox"/>	Remote LiNX DLX-REM211 = Drive & Actuators 3 Drive functions and power module incl. 2 Actuators	◆

<b>CONTROL BRACKETRY*</b>			
JBC1207	<input type="checkbox"/>	Standard Bracket RH	<input type="radio"/>
JBC1208	<input type="checkbox"/>	Standard Bracket LH	<input type="radio"/>
<b>SHROUD COLOUR*</b>			
JBC2871	<input type="checkbox"/>	Mid range - Black/Anthracite Shroud	<input type="radio"/>
JBC2872	<input type="checkbox"/>	Mid range - Silver Shroud	<input type="radio"/>
JBC2873	<input type="checkbox"/>	Mid range - Bengal Red Shroud	<input type="radio"/>
JBC2874	<input type="checkbox"/>	Mid range - Mercury Blue Shroud	<input type="radio"/>
JBC2875	<input type="checkbox"/>	Mid range - Midnight Purple Shroud	<input type="radio"/>
JBX0299	<input type="checkbox"/>	Premium - Individual Shroud Colour (includes coloured spreader bars on Modulite back) RAL number required Extended lead lime: + 15 days	£195
<b>RIM COLOUR*</b>			
JBC2882	<input type="checkbox"/>	Mid range - Black Rim Inserts	<input type="radio"/>
JBC2883	<input type="checkbox"/>	Mid range - Silver Rim Inserts	<input type="radio"/>
JBC2884	<input type="checkbox"/>	Mid range - Bengal Red Rim Inserts	<input type="radio"/>
JBC2885	<input type="checkbox"/>	Mid range - Mercury Blue Rim Inserts	<input type="radio"/>
JBC2886	<input type="checkbox"/>	Mid range - Midnight Purple Rim Inserts	<input type="radio"/>
<b>FORKS*</b>			
JBC2632	<input type="checkbox"/>	Double sided fork Without castor rim insert (enlarges wide base to 656 mm)	<input checked="" type="radio"/>
<b>FRONT WHEEL*</b>			
JBC0828	<input type="checkbox"/>	8" (200 mm) x 2.50-4 (65 mm) Castor, puncture proof solid tyre	<input checked="" type="radio"/>
<b>DRIVE WHEEL*</b>			
JBC0705	<input type="checkbox"/>	DW 14" (355 mm) x 3" (75 mm), Solid inlay Dual density foam insert	<input checked="" type="radio"/>
<b>C.T.C SUSPENSION*</b>			
JBC2870	<input type="checkbox"/>	CTC Suspension Control - Traction - Comfort	<input checked="" type="radio"/>
<b>SUSPENSION SETTING*</b>			
JBC1910	<input type="checkbox"/>	Soft (suggested for up to 70 kg User Weight)	<input checked="" type="radio"/>
JBC1912	<input type="checkbox"/>	Medium (suggested for up to 90 kg User Weight)	<input type="radio"/>
JBC1911	<input type="checkbox"/>	Firm (suggested for up to 110 kg User Weight)	<input type="radio"/>
JBC1915	<input type="checkbox"/>	Very firm (suggested for up to 136 kg User Weight)	<input type="radio"/>
<b>MOTOR TYPE*</b>			
JBC1282	<input type="checkbox"/>	6 km/h	<input checked="" type="radio"/>
<b>BATTERY TYPE*</b>			
JBC2860	<input type="checkbox"/>	52 Ah AGM (Top-Gel)	<input checked="" type="radio"/>
<b>CHARGER*</b>			
JBC1343	<input type="checkbox"/>	Opti Charge Charger 8A (IP44)	<input checked="" type="radio"/>
<b>ANTI-TIPPERS*</b>			
JBC1205	<input type="checkbox"/>	Anti-tipper	<input checked="" type="radio"/>
<b>TRANSPORTATION KIT*</b>			
JBC1494	<input type="checkbox"/>	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19) Adequate headrest is compulsory	<input checked="" type="radio"/>
<b>LAP BELT*</b>			
JBC2113	<input type="checkbox"/>	Lap belt with Auto Lock Flexible mount in side rails	<input checked="" type="radio"/>

<b>TOTAL PRICE</b>	
--------------------	--

**Non Conforming Options:**

JBC2106      Body fit Kit

JBC2171      Tension Adjustable Upholstery



For a more complex customer requirement please refer to our Invacare Unique special order form. email [specialsuk@invacare.com](mailto:specialsuk@invacare.com) or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: \_\_\_\_\_ Email: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Description of special request (please include relevant option code or part number):