



**Yes, you can:**

Tel: 01656 776222 Fax: 01656 776220  
Email: ordersuk@invacare.com  
Online Spares available at: www.invacare.co.uk

JANUARY 2021  
Customer Ref:

|                          |                    |
|--------------------------|--------------------|
| <b>Customer Ref*</b>     | <b>Account No*</b> |
| <b>Contact No*</b>       | <b>Order No*</b>   |
| <b>Delivery Address*</b> | <b>Order Qty*</b>  |
|                          | <b>Order date*</b> |

**M97 AVIVA RX20 Modulite**

**Standard Features**

- Max user weight 136kg (21 Stone)
- Crash tested ISO 7176 -Part 19
- Modulite highly adjustable seating system
- Telescopic adjustable seat widths from 330mm - 530mm (13" - 21") with an additional 50mm, (2") by arms
- Seat depths 410mm - 510mm stepless adjustment (16" - 20")
- Flex3 adjustable Backrest Plates with Laguna Backrest cushion
- Standard T-Armrests, detachable with short armpad
- Neckrest depth adjustable
- Powered "fixed pivot" Tilt 20°
- Choice of seat height's
- Powered Recliner 30°
- 80° hangers with fixed footplate
- Remote LiNX DLX-REM211 with standard bracket
- Choice of 5 Shroud colours, 5 Rim insert colours
- 8" (200 mm) x 2.50-4 (65 mm) Castor, puncture proof solid tyre
- 14" (355 mm) x 3" (75 mm), Solid inlay black drive wheels
- CTC Suspension Control - Traction - Comfort
- 6 km/h 4 pole Durawatt motors
- 52Ah batteries
- 8A (200 - 250V) charger
- Anti tippers
- Kerb Rider With park position feature
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)

**Standard Base price: £6,516**



Please make your selection with a cross in the box    ○ Free of charge    ◆ standard feature    POA special request

Customer weight =                      KG

**Code**                      **Description**

**BASE MODEL\***

|         |                          |   |        |
|---------|--------------------------|---|--------|
| JBC0000 | <input type="checkbox"/> | AVIVA RX 20 Base Model Max user weight up to 136 kg | £6,516 |
|---------|--------------------------|---|--------|

**CHASSIS WIDTH\***

|         |                          |                             |   |
|---------|--------------------------|-----------------------------|---|
| JBL2810 | <input type="checkbox"/> | Chassis Narrow (560 mm/22") | ◆ |
|---------|--------------------------|-----------------------------|---|

|   |                          |   |   |
|---|--------------------------|---|---|
| <b>SEATING TYPE*</b>  |                          |   |   |
| JBC2100   | <input type="checkbox"/> | Invacare® Modulite Seat Unit                                  | ◆ |
| <b>SEAT FRAME*</b>  |                          |   |   |
| JBC2111   | <input type="checkbox"/> | Telescopic Seat Frame and Plate with Side Rails               | ◆ |
| <b>SEAT WIDTH*</b>  |                          |   |   |
| JBC2121   | <input type="checkbox"/> | 380 mm Seat Plate   | ○ |
| JBC2122   | <input type="checkbox"/> | 430 mm Seat Plate   | ○ |
| JBC2123   | <input type="checkbox"/> | 480 mm Seat Plate   | ○ |
| JBC2124   | <input type="checkbox"/> | 530 mm Seat Plate (NC 2106)                                   | ○ |
| <b>Body Fit Kit</b>   |                          |   |   |
| JBC2106   | <input type="checkbox"/> | Body fit Kit  | ○ |
| <b>SEAT DEPTH*</b>  |                          |   |   |
| JBC2132   | <input type="checkbox"/> | Seat Depth 410 mm   | ○ |
| JBC2134   | <input type="checkbox"/> | Seat Depth 460 mm   | ○ |
| JBC2136   | <input type="checkbox"/> | Seat Depth 510 mm   | ○ |
| <b>SEAT CUSHION*</b>  |                          |   |   |
| JBC2142   | <input type="checkbox"/> | Matrx Flo-tech Cushion Flo-tech Lite                          | ◆ |
| <b>BACKREST SUPPORT*</b>  |                          |   |   |
| JBC2163   | <input type="checkbox"/> | Backposts short (480 mm) (NC 2173)                            | ◆ |
| JBC2162   | <input type="checkbox"/> | Backposts tall (540 mm)                                       | ○ |
| <b>PUSH HANDLES FOR BACKPOSTS</b>                                       |                          |   |   |
| JBC2108   | <input type="checkbox"/> | Backposts without push handles                                | ◆ |
| <b>BACKREST UPHOLSTERY / CUSHION*</b>                                   |                          |   |   |
| JBC2173   | <input type="checkbox"/> | Flex3 adjustable Backrest Plates with Laguna Backrest cushion | ◆ |
| <b>UPHOLSTERY TYPE*</b>   |                          |   |   |
| JBC2151   | <input type="checkbox"/> | Black Fabric Upholstery                                       | ○ |
| JBC2152   | <input type="checkbox"/> | Black Black multi stretch PU coated fabric Upholstery         | ○ |
| JBC2155   | <input type="checkbox"/> | Black Spacer Fabric Upholstery                                | ○ |
| <b>ARMREST (Pair)*</b>  |                          |   |   |
| JBC2191   | <input type="checkbox"/> | T-Arm detachable, Standard                                    | ◆ |
| <b>ARMPAD (Pair)*</b>   |                          |   |   |
| JBC2195   | <input type="checkbox"/> | Short Armpad (325 mm)   | ◆ |
| <b>HEADREST / NECKREST - Adequate headrest option must be selected*</b> |                          |   |   |
| JBC1350   | <input type="checkbox"/> | Neckrest depth adjustable                                     | ◆ |
| <b>SEAT LIFTER / SEAT TILT*</b>   |                          |   |   |
| JBC1092   | <input type="checkbox"/> | Powered "fixed pivot" Tilt 20°                                | ◆ |
| <b>SEAT TO GROUND HEIGHT*</b>   |                          |   |   |
| JBC2857   | <input type="checkbox"/> | Seat-to-Floor Height 435 mm                                   | ○ |
| JBC2858   | <input type="checkbox"/> | Seat-to-Floor Height 460 mm                                   | ○ |
| JBC2859   | <input type="checkbox"/> | Seat-to-Floor Height 485 mm                                   | ○ |
| <b>BACKREST RECLINER*</b>   |                          |   |   |
| JBC0977   | <input type="checkbox"/> | Powered Recliner 30°  | ◆ |
| <b>FOOTREST HANGER TYPE*</b>  |                          |   |   |
| JBC0526   | <input type="checkbox"/> | Standard hanger 80° pair Lower leg length: 350 - 450 mm       | ◆ |
| <b>FOOTPLATE TYPE*</b>  |                          |   |   |
| JBC0512   | <input type="checkbox"/> | Fixed Footplate pair  | ◆ |

| <b>CONTROLS*</b>           |                          |  |      |
|----------------------------|--------------------------|--|------|
| JBC2203                    | <input type="checkbox"/> | Remote LiNX<br>DLX-REM211 = Drive & Actuators 3 Drive functions<br>and power module incl. 2 Actuators  | ◆    |
| <b>CONTROL BRACKETRY*</b>  |                          |  |      |
| JBC1207                    | <input type="checkbox"/> | Standard Bracket RH  | ○    |
| JBC1208                    | <input type="checkbox"/> | Standard Bracket LH  | ○    |
| <b>SHROUD COLOUR*</b>      |                          |  |      |
| JBC2871                    | <input type="checkbox"/> | Mid range - Black/Anthracite Shroud  | ○    |
| JBC2872                    | <input type="checkbox"/> | Mid range - Silver Shroud  | ○    |
| JBC2873                    | <input type="checkbox"/> | Mid range - Bengal Red Shroud  | ○    |
| JBC2874                    | <input type="checkbox"/> | Mid range - Mercury Blue Shroud  | ○    |
| JBC2875                    | <input type="checkbox"/> | Mid range - Midnight Purple Shroud   | ○    |
| JBX0299                    | <input type="checkbox"/> | Premium - Individual Shroud Colour<br>(includes coloured spreader bars on Modulite back)<br>RAL number required<br>Extended lead lime: + 15 days | £195 |
| <b>RIM COLOUR*</b>         |                          |  |      |
| JBC2882                    | <input type="checkbox"/> | Mid range - Black Rim Inserts  | ○    |
| JBC2883                    | <input type="checkbox"/> | Mid range - Silver Rim Inserts   | ○    |
| JBC2884                    | <input type="checkbox"/> | Mid range - Bengal Red Rim Inserts   | ○    |
| JBC2885                    | <input type="checkbox"/> | Mid range - Mercury Blue Rim Inserts   | ○    |
| JBC2886                    | <input type="checkbox"/> | Mid range - Midnight Purple Rim Inserts  | ○    |
| <b>FORKS*</b>              |                          |  |      |
| JBC2632                    | <input type="checkbox"/> | Double sided fork Without castor rim insert<br>(enlarges wide base to 656 mm)  | ◆    |
| <b>FRONT WHEEL*</b>        |                          |  |      |
| JBC0828                    | <input type="checkbox"/> | 8" (200 mm) x 2.50-4 (65 mm) Castor, puncture proof<br>solid tyre  | ◆    |
| <b>DRIVE WHEEL*</b>        |                          |  |      |
| JBC0705                    | <input type="checkbox"/> | DW 14" (355 mm) x 3" (75 mm), Solid inlay Dual<br>density foam insert  | ◆    |
| <b>C.T.C SUSPENSION*</b>   |                          |  |      |
| JBC2870                    | <input type="checkbox"/> | CTC Suspension Control - Traction - Comfort  | ◆    |
| <b>SUSPENSION SETTING*</b> |                          |  |      |
| JBC1910                    | <input type="checkbox"/> | Soft (suggested for up to 70 kg User Weight)   | ◆    |
| JBC1912                    | <input type="checkbox"/> | Medium (suggested for up to 90 kg User Weight)   | ○    |
| JBC1911                    | <input type="checkbox"/> | Firm (suggested for up to 110 kg User Weight)  | ○    |
| JBC1915                    | <input type="checkbox"/> | Very firm (suggested for up to 136 kg User Weight)   | ○    |
| <b>MOTOR TYPE*</b>         |                          |  |      |
| JBC1282                    | <input type="checkbox"/> | 6 km/h   | ◆    |
| <b>BATTERY TYPE*</b>       |                          |  |      |
| JBC2860                    | <input type="checkbox"/> | 52 Ah AGM (Top-Gel)  | ◆    |
| <b>CHARGER*</b>            |                          |  |      |
| JBC1343                    | <input type="checkbox"/> | Opti Charge Charger 8A (IP44)  | ◆    |
| <b>ANTI-TIPPERS*</b>       |                          |  |      |
| JBC1205                    | <input type="checkbox"/> | Anti-tipper  | ◆    |
| <b>KERB RIDER</b>          |                          |  |      |
| JBC1244                    | <input type="checkbox"/> | Kerb Rider With park position feature  | ◆    |

**TRANSPORTATION KIT\***

|         |  |  |   |
|---------|--|--|---|
| JBC1494 |  | 4-Point tie-down<br>(Tested version for using wheelchair as a vehicle seat according to ISO 7176-19) Adequate headrest is compulsory | ◆ |
|---------|--|--|---|

**LAP BELT\***

|         |  |  |   |
|---------|--|--|---|
| JBC2113 |  | Lap belt with Auto Lock Flexible mount in side rails | ◆ |
|---------|--|--|---|

|                    |  |  |  |
|--------------------|--|--|--|
| <b>TOTAL PRICE</b> |  |  |  |
|--------------------|--|--|--|

**Non Conforming Options:**

JBC2106            Body fit Kit

JBC2173            Flex3 adjustable Backrest Plates with Laguna Backrest cushion



For a more complex customer requirement please refer to our Invacare Unique special order form. email [specialsuk@invacare.com](mailto:specialsuk@invacare.com) or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: \_\_\_\_\_ Email: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Description of special request (please include relevant option code or part number):