



Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
 Email: ordersuk@invacare.com
 Online Spares available at: www.invacare.co.uk

Powerchairs Retail Prescription Form

MARCH 2021
 Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

AVIVA RX40 HD

Standard Features

- Max user weight 160kg (25 Stone)
- Crash tested ISO 7176 -Part 19
- Modulite adjustable seating system
- Telescopic adjustable seat widths from 330mm - 530mm (13" - 21") with an additional 50mm, (2") by arms
- Seat depths 410mm - 510mm stepless adjustment (16" - 20")
- Tension adjustable backrest
- Standard T-Armrests, detachable with short armpad
- Manual settable seat Tilt (0° - 20°)
- Manual settable backrest Recline (0° - 30° in 3.75° steps)
- 80° hangers with fixed footplate
- LiNX control REM110 with standard bracket
- Choice of 5 Shroud colours, 5 Rim insert colours
- 9" (225 mm) x 2.80 (70 mm) Castor, puncture proof solid tyre
- 14" (355 mm) x 3" (75 mm), Solid inlay black drive wheels
- CTC Suspension Control - Traction - Comfort
- 6 km/h (4mph) 4 pole Durawatt motors and 60Ah batteries
- 8A (200 - 250V) charger
- Anti tippers
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)

Standard Base price: £6,052



Please make your selection with a cross in the box Free of charge standard feature POA special request

Customer weight = KG

Code	Description	
BASE MODEL*		
JBH0002	AVIVA RX 40 HD Base Model Max user weight up to 160 kg HD Version, includes RX 40 Performance Package Wide Base only	(NC 2193) £6,052
JBX1300	Special Extended lead time: + 10 days	POA
CHASSIS WIDTH*		
JBH2830	Chassis Wide (635 mm/25")	◆

SEATING TYPE*			
JBH2100		Invacare® Modulite Seat Unit	○
JBX2101		Invacare® Modulite HD Seat Unit On RX 40 HD only Tension Adjustable Upholstry only JBX2171	(NC 2100, 2179, 2172, 2173, 2170) ○
SEAT FRAME*			
JBH2111		Telescopic Seat Frame and Plate with Side Rails width 380 to 530 mm Additional 50 mm by arms	○
JBH2112		Telescopic Seat Frame + Sling Seat with Side Rails width 380 to 530 mm; + 30 mm seat height Additional 50 mm by arms	○
JBX2102		HD Telescopic Seat Frame and Plate with Side Rails width 530 to 550 mm	○
SEAT WIDTH*			
JBH2123		480 mm Seat Plate	(NC 2702, 2704, 2706) ○
JBH2124		530 mm Seat Plate	(NC 2174, 2175, 2176, 2177, 2178, 2702, 2704, 2706, 2106) ○
JBX2124		530 mm Seat Plate	(NC 2175, 2176, 2177, 2178, 2702, 2704, 2706, 2106, 2100) ○
JBX2104		550 mm Seat Plate	○
Body Fit Kit			
JBH2106		Body fit Kit Outrigging Backrest Brackets and Pivot Point (outrig: 1 inch/2,54 cm per side)	○
SEAT DEPTH*			
JBH2132		Seat Depth 410 mm	(NC 2196, 2197, 2198) ○
JBH2134		Seat Depth 460 mm	○
JBH2136		Seat Depth 510 mm	○
JBX2136		Seat Depth 510 mm	○
JBX2137		580 mm	○
SEAT CUSHION*			
JBH2142		Matrx Flo-tech Cushion Flo-tech Lite	£111
JBH2147		Matrx Libra Cushion, 90 mm thick Extended lead time: + 4 days	£369
JBH2140		Without seat cushion	○
JBX2148		Seat Cushion (90 mm high)	○
JBX2140		Without seat cushion	○
BACKREST SUPPORT*			
JBH2163		Backposts short (480 mm) For Tension Adjustable Upholstery/Matrx Backs, measured from seat plate to top of upholstery Width adjustable cast bracket	(NC 2172, 2173) ○
JBH2162		Backposts tall (540 mm) For Flex3 backrests/Tension Adjustable Upholstery/Matrx Backs, measured from seat plate to top of upholstery Width adjustable cast bracket	◆
JBX2162		Backposts tall (540 mm) For Flex3 backrests/Tension Adjustable Upholstery/Matrx Backs, measured from seat plate to top of upholstery Width adjustable cast bracket	◆
PUSH HANDLES FOR BACKPOSTS			
JBH2108		Backposts without push handles	(NC 2101, 2219, 2220, 2227, 2228, 2245, 2291, 2292) ○

BACKREST UPHOLSTERY / CUSHION*			
JBH2171		Tension Adjustable Upholstery nylon overlay short = 480 mm / 7 straps measured from seat plate to top of upholstery 480 mm	(NC 2155, 2230, 2247) ◆
JBH2179		Flex3 Tension Adjustable Upholstery height 540 mm (seat plate to top of upholstery)	(NC 2151, 2192, 2193, 2194, 2199, 1446, 2229, 2155, 2230, 2247) £146
JBH2172		Flex3 adjustable Backrest Plates (without Backrest Cushion) Flex3 Backrest plate width and height adjustable 560 mm - 790 mm	(NC 2163, 2128) £238
JBH2173		Flex3 adjustable Backrest Plates with Laguna Backrest cushion Flex3 Backrest plate width and height adjustable 560 mm - 790 mm	(NC 2163, 2128) £386
JBH2174		Matrx Posture Back, height 400 mm (Mesh fabric) Corresponds to SW setting chosen (up to 136 kg max user weight)	£530
JBH2175		Matrx Posture Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen (up to 136 kg max user weight)	£530
JBH2176		Matrx Elite Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen (up to 136 kg max user weight)	£530
JBH2177		Matrx Elite TR Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen	£530
JBH2178		Matrx Elite Deep Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen (up to 136 kg max user weight)	£616
JBH2170		Without Backrest Upholstery (Back Posts/Back Frame only)	(NC 2322, 2323, 2324, 2325, 2326, 2327, 2230, 2247) ○
JBX2171		Tension Adjustable Upholstery nylon overlay short = 480 mm / 7 straps measured from seat plate to top of upholstery 480 mm MUST BE SELECTED WITH HD SEATING JBX2102	(NC 2155, 2230, 2247) ◆
UPHOLSTERY TYPE*			
JBH2151		Black Fabric Upholstery (Nylon with option 2141; 2- way-stretch material with options 2142, 2143)	◆
JBH2152		Black multi stretch PU coated fabric Upholstery	(NC 2171) ○
JBH2155		Black Spacer Fabric Upholstery	(NC 2171) ○
JBX2156		Black Nylon	◆
ARMREST (Pair)*			
JBH2191		T-Arm detachable, Standard Adjustment range from seat plate: 245 - 315 mm	◆
JBH2181		T-Arm detachable, Extended height Adjustment range from seat plate: 295 - 365 mm	£38
JBH2198		Padded Armrest with integrated cloth guard	(NC 2132, 2195, 2196, 2197, 2199) £104
JBH2193		Flip-up Armrest Up to 150 kg max user weight	(NC 1445, 1446, 2101, 2179) £61
JBH2192		Following Armrest Height adjustable from 230 - 300 mm	(NC 1445, 1446, 2179) £175
JBH2194		Following Armrest Height adjustable from 300 - 360 mm	(NC 1445, 1446, 2179) £175
JBX2191		T-Arm detachable, Standard Adjustment range from seat plate: 245 - 315 mm	○
JBX2181		T-Arm detachable, Extended height Adjustment range from seat plate: 295 - 365 mm	£38

ARMPAD (Pair)*				
JBH2195		Short Armpad (325 mm)	(NC 2198)	◆
JBH2196		Long Armpad (375 mm)	(NC 2132, 2198)	○
JBH2197		Long Ergonomic Armpad (397 mm)	(NC 2132, 2198, 2230, 2247)	○
JBH2199		Soft Armpad In combination with 10-way switch, the switch 2251 has to be mounted on the opposite side of the remote	(NC 2198)	£90
JBX2195		Short Armpad (325 mm)	(NC 2198)	◆
JBX2196		Long Armpad (375 mm)	(NC 2132, 2198)	○
JBX2197		Long Ergonomic Armpad (397 mm)	(NC 2132, 2198, 2230, 2247)	○
HIP SUPPORTS LH				
JBH2184		Hip Support LH, pad size 140 x 100 mm	(NC 2198, 2701, 2702, 2703, 2704, 2705, 2706)	£240
JBH2186		Hip Support LH, pad size 230 x 100 mm	(NC 2198, 2701, 2702, 2703, 2704, 2705, 2706)	£240
JBH2188		Hip Support LH, pad size 330 x 100 mm	(NC 2198, 2701, 2702, 2703, 2704, 2705, 2706)	£240
HIP SUPPORTS RH				
JBH2185		Hip Support RH, pad size 140 x 100 mm	(NC 2198, 2701, 2702, 2703, 2704, 2705, 2706)	£240
JBH2187		Hip Support RH, pad size 230 x 100 mm	(NC 2198, 2701, 2702, 2703, 2704, 2705, 2706)	£240
JBH2189		Hip Support RH, pad size 330 x 100 mm	(NC 2198, 2701, 2702, 2703, 2704, 2705, 2706)	£240
MOUNTING HARDWARE FOR HIP SUPPORTS				
JBH2164		Mounting hardware for hip supports (6,1" / 155 mm high)		○
JBH2165		Mounting hardware for hip supports (8,1" / 206 mm high)		○

HEADREST / NECKREST - Adequate headrest option must be selected*			
JBH1350		Neckrest depth adjustable (NC 2230, 2247)	£77
JBH1432		Headrest sideways adjustable (NC 2230, 2247)	£215
JBH1427		Headrest with cheek support (NC 2230, 2247)	£215
JBH2323		Elan Standard Pad 150 mm / 6" wide w/standard cover (ESP26) Complete with Elan Headrest Hardware (NC 2230, 2247, 2170)	£343
JBH2324		Elan Standard Pad 255 mm / 10" wide w/standard cover (ESP26) Complete with Elan Headrest Hardware (NC 2230, 2247, 2170)	£343
JBH2322		Elan Standard Pad 355 mm / 14" wide w/standard cover (ESP26) Complete with Elan Headrest Hardware (NC 2230, 2247, 2170)	£343
JBH2325		Elan Occipital Pad 230x180 mm / 9"x7" w/standard cover (EOP23) Complete with Elan Headrest Hardware (NC 2230, 2247, 2170)	£343
JBH2326		Elan Occipital Pad 305x200 mm / 12"x8" w/standard cover Complete with Elan Headrest Hardware (NC 2230, 2247, 2170)	£343
JBH2327		Elan 4 Point Headrest Pad 280x255 mm / 11"x10" w/standard Complete with Elan Headrest Hardware (NC 2230, 2247, 2170)	£466
JBH1496		Non Invacare Headrest but profile mount fitting Flex3 hard back (NC 2170, 2171, 2174, 2175, 2176, 2177, 2178, 1288, 1289)	○
JBH1495		Non Invacare Headrest - no mounting hardware With Backbrace/Pushbar dia. 22mm in conjunction with telescopic seat frame/backposts	○
JBX1350		Neckrest depth adjustable (NC 2230, 2247)	£77
JBX1432		Headrest sideways adjustable (NC 2230, 2247)	£215
JBX1427		Headrest with cheek support (NC 2230, 2247)	£215
JBX2323		Elan Standard Pad 150 mm / 6" wide w/standard cover (ESP26) Complete with Elan Headrest Hardware (NC 2230, 2247, 2170)	£343
JBX2324		Elan Standard Pad 255 mm / 10" wide w/standard cover (ESP26) Complete with Elan Headrest Hardware (NC 2230, 2247, 2170)	£343
JBX2322		Elan Standard Pad 355 mm / 14" wide w/standard cover (ESP26) Complete with Elan Headrest Hardware (NC 2230, 2247, 2170)	£343
JBX2325		Elan Occipital Pad 230x180 mm / 9"x7" w/standard cover (EOP23) Complete with Elan Headrest Hardware (NC 2230, 2247, 2170)	£343
JBX2326		Elan Occipital Pad 305x200 mm / 12"x8" w/standard cover Complete with Elan Headrest Hardware (NC 2230, 2247, 2170)	£343
JBX2327		Elan 4 Point Headrest Pad 280x255 mm / 11"x10" w/standard Complete with Elan Headrest Hardware (NC 2230, 2247, 2170)	£466
JBX1495		Non Invacare Headrest - no mounting hardware With Backbrace/Pushbar dia. 22mm in conjunction with telescopic seat frame/backposts	○
LATERAL SUPPORTS			
JBH1446		Lateral Support with extended arm, rigid pad (NC 2170, 2171, 2174, 2175, 2176, 2177, 2178, 2179, 2192, 2193, 2194)	£159
JBH1445		Lateral Support with extended arm, multi-adjustable pad (NC 2170, 2171, 2174, 2175, 2176, 2177, 2178, 2192, 2193, 2194)	£159
SEAT LIFTER / SEAT TILT*			
JBH1091		Manually settable Seat Tilt 0 - 15° (NC 1401)	○
JBH1092		Powered "fixed pivot" Tilt 20° (NC 2101, 2201)	£384
JBH1093		Powered Tilt (0° - 30°) with CoG shift (NC 2101, 2201)	£1,021
JBH1094		Seat lifter with powered tilt (300 mm, 30°) with CoG shift (NC 2101, 2201, 2220)	£1,802
JBX1091		Manually settable Seat Tilt 0 - 15° (NC 1401)	○
JBX1095		Powered Tilt (0° - 20°) (NC 2201)	£384
JBX1093		Powered Tilt (0° - 30°) with CoG shift (NC 2100, 2201)	£1,021
JBX1094		Seat lifter with powered tilt (300 mm, 30°) with CoG shift, 6 point Tie downs Must be selected (NC 2100, 2201, 2220)	£1,802

SEAT TO GROUND HEIGHT*			
JBH2857		Seat-to-Floor Height 435 mm	○
JBH2858		Seat-to-Floor Height 460 mm	○
JBH2859		Seat-to-Floor Height 485 mm	○
JBX2857		Seat-to-Floor Height 435 mm	○
JBX2858		Seat-to-Floor Height 460 mm	○
JBX2859		Seat-to-Floor Height 485 mm	○
BACKREST RECLINER*			
JBH0975		Manual Settable Recliner (in 3.75° steps) 90° - 120°	◆
JBH0977		Powered Recliner 30°	£480
JBX0975		Manual Settable Recliner (in 3.75° steps) 90° - 120°	◆
JBX0977		Powered Recliner 30° (NC 2201)	£480
CENTRE MOUNTED LEGREST*			
JBH0556		Centre mount legrest, angle settable Length: 280 - 385 mm (NC 2366, 0517, 1416)	£355
JBH0589		LNx Power centre mount legrest, length compensating Length: 340 - 410 mm (NC 2201, 1416)	£1,562
JBH2855		Without Legrests	○
FOOTPLATES / FOOT PLATFORM*			
JBH2365		Angle & width adjustable footplates S (127 W x 190 mm D / 5"x7.5")	£98
JBH0559		Angle & width adjustable footplates M (140 W x 240 mm D / 5.5"x9.5")	○
JBH0516		Angle & width adjustable footplates L (150 W x 290 mm D / 6"x11.5)	£128
JBH2366		Angle adjustable foot platform S (305 W x 215 mm D / 12"x8.5")	£257
JBH0517		Angle adjustable foot platform M (305 W x 275 mm D / 12"x10.75")	£257
JBX3001		Please select If Centre mount footplate is not required. (NC 0589)	○
FOOTREST HANGER TYPE*			
JBH0503		Hanger Vari F pair Lower leg length: 290 - 460 mm Angle adjustment: 70 - 0° with tools	◆
JBH0533		Hanger Vari A pair Lower leg length: 290 - 460 mm Angle adjustment: 70 - 0° with handwheel and calf pad	£221
JBH0531		Hanger ADE pair Lower leg length: 290 - 460 mm Powered Angle adjustment: 80 - 0° with calf pad (NC 2201, 0522)	£637
JBH0501		Without Hanger pair (NC 0504, 0513, 0522)	○
FOOTPLATE TYPE*			
JBH0512		Fixed Footplate pair (NC 0533, 0531)	◆
JBH0504		Angle Adjustable Footplate pair	£30
JBH0513		Angle and depth adjustable Footplate pair	£112
JBH0522		Angle adjustable Footboard (NC 0531, 0501)	£145
JBX3003		No Side Mounted Footplate required (NC 0503, 0533, 0531)	○

CONTROLS*			
JBH2201		Remote LiNX DLX-REM110 = Drive only 1 Drive function with Drive only power module	(NC 2837, 1092, 1093, 1094, 0977, 0589, 2231, 2232, 2230, 2247, 2257, 2233, 2240, 2258, 2219, 2245, 2270, 1908) ◆
JBH2203		Remote LiNX DLX-REM211 = Drive & Actuators 3 Drive functions and power module incl. 2 Actuators	(NC 1437, 2231, 2232, 2229, 2233, 2230, 2247, 2257, 2240, 2270, 2258, 1908, 1909) ○
JBH2204		Remote LiNX DLX-REM216 = Drive, Actuators & Lights 3 Drive functions and power module incl. 2 Actuators & Lights	(NC 2232, 2229, 2233, 2230, 2247, 2257, 2240, 2270, 2258) £144
JBH2205		Rehab Remote LiNX, touch-screen user interface DLX-REM400 Configurable jack sockets	(NC 2220) £324
JBH2206		Rehab Remote LiNX without joystick for alternative access users DLX-REM500 Requires a user input/alternative joystick: 2231, 2232, 2229, 2230, 2247, 2233	(NC 2220, 2270) £738
REM400 TOGGLE OPTION			
JBH2270		Assignable Toggle Switch Module f. REM400 (DLX-SW400TGL)	(NC 2203, 2204, 2206) £300
CONTROL BRACKETRY*			
JBH1207		Standard Bracket RH	(NC 2285, 2283, 2287, 2207, 2265, 2275, 1404) ○
JBH1208		Standard Bracket LH	(NC 2286, 2284, 2288, 2208, 2266, 2276, 1405) ○
JBH1707		Swing Away Bracket RH	(NC 2285, 2283, 2287, 2207, 2265, 2275, 1404) £126
JBH1708		Swing Away Bracket LH	(NC 2286, 2284, 2288, 2208, 2266, 2276, 1405) £126
JBH2213		Centre Mount, swing-away to RH, Nucleus	(NC 2283, 2287, 2207, 2265, 2275, 1415, 2217, 1404, 1405) £474
JBH2214		Centre Mount, swing-away to LH, Nucleus	(NC 2284, 2288, 2208, 2266, 2276, 1339, 2218, 1404, 1405) £474
JBH2227		Bracket RH, mounted for Attendant	(NC 1287, 2220, 2128, 2108, 2219, 1221, 1312) £60
JBH2228		Bracket LH, mounted for Attendant	(NC 1286, 2220, 2128, 2108, 2219, 1222, 1363) £60
JBX2239		Individual mounting option per request (Special) Requires additional information/specification	POA
DUAL CONTROL			
JBH2220		Dual Control LiNX (DLX-REM050 as secondary remote)	(NC 2273, 2108, 2205, 2206, 1221, 1222, 2227, 2228) £376
JBH2219		Compact Dual Control LiNX (DLX-ACU200)	(NC 2108, 1221, 1222, 2227, 2228) £376
INVACARE LINX IDC DUAL CONTROL			
JBH2245		Invacare LiNX IDC Dual control	(NC 2108, 2259, 2260) £714
DUAL CONTROL BRACKETRY			
JBH2259		Bracket RH, Attendant	(NC 2291, 2245, 1312) ○
JBH2260		Bracket LH, Attendant	(NC 2292, 2245, 1363) ○
INVACARE LINX IDC DUAL CONTROL BRACKETRY			
JBH1221		Bracket RH, detachable	(NC 2291, 2220, 2219, 2227, 1312) ○
JBH1222		Bracket LH, detachable	(NC 2292, 2220, 2219, 2228, 1363) ○
LINX COMPACT REMOTE			
JBH2231		Proportional Compact Remote DLX-CR400	(NC 2203, 2204) £684
JBH2232		Proportional Compact Remote Low Force DLX-CR-400LF	(NC 2203, 2204) £914
LINX COMPACT REMOTE BRACKET			
JBH2207		RH Mount	(NC 1707, 2213, 2217, 1404) ○
JBH2208		LH Mount	(NC 1708, 2214, 2218, 1405) ○
JBH2265		Centre Mount, swing-away to RH (Nucleus) If several joysticks (or Joystick + REM500) should be mounted on one Nucleus, please select special mounting option 2215!	(NC 1707, 2213, 2217, 1404, 1405, 1415) £474
JBH2266		Centre Mount, swing-away to LH (Nucleus) If several joysticks (or Joystick + REM500) should be mounted on one Nucleus, please select special mounting option 2215!	(NC 1708, 2214, 2218, 1404, 1405, 1339) £474
JBX2215		Individual mounting option per request (Special) Requires additional information/specification	POA

CHIN CONTROL			
JBH2229		Proportional Chin Joystick (ASL138) SELECTED HEADREST ALWAYS WITH STD MOUNTING HARDWARE	(NC 2201, 2203, 2204, 2230, 2247, 2233, 2257, 2240, 2179, 2128) £1,650
CHIN CONTROL BRACKETRY MANUAL			
JBH1289		Swing Away Bracket RH, side manual (ASL628) Requires Elan Headrest	(NC 1350, 1432, 1427, 2322, 1496) ○
JBH1288		Swing Away Bracket LH, side manual (ASL628) Requires Elan Headrest	(NC 1350, 1432, 1427, 2322, 1496) ○
CHIN CONTROL BRACKETRY POWERED			
JBH1287		Swing Away Bracket RH, powered	(NC 2259, 1221, 2291, 1312, 1415) £2,079
JBH1286		Swing Away Bracket LH, powered	(NC 2260, 1222, 2292, 1363, 1339) £2,079
LINX HEAD CONTROL			
JBH2230		Head Control with ultrasonic proximity sensors (ATOM, ASL104) with standard wings M3	(NC 2203, 2204, 2229, 2233, 2323, 2324, 2322, 2325, 2326, 2327, 2257, 2240, 2175, 2171, 2170, 2197, 1350, 1432, 1427) £3,542
JBH2247		Head Control with ultrasonic proximity sensors (PROTON, ASL104P) with PROTON wings	(NC 2203, 2204, 2229, 2233, 2323, 2324, 2322, 2325, 2326, 2327, 2257, 2240, 2175, 2171, 2170, 2197, 1350, 1432, 1427) £4,148
LINX FOUR -SWITCH PROXIMITY ARRAY			
JBH2257		Four-Switch Proximity Array (ASL106)	(NC 2203, 2204, 2229, 2230, 2247, 2233, 2240, 2128) £5,485
JBH2275		Centre Mount, swing-away to RH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm)	(NC 1707, 2213, 2217, 1404, 1405, 1415) £474
JBH2276		Centre Mount, swing-away to LH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm)	(NC 1708, 2214, 2218, 1404, 1405, 1339) £474
JBX2269		Individual mounting option per request (Special) Requires additional information/specification	POA
LINX WIRELESS ACCESSORIES FOR HEAD CONTROL AND PROXIMITY CONTROL			
JBH2249		Wireless Triple Switch Receiver (ASL557-3)	£324
JBH2250		Wireless Mouse Emulator (ASL558)	£778
LINX MICRO EXTREMITY CONTROL BRACKET			
JBH2233		MEC Proportional Joystick (Micro Extremity Control ASL130)	(NC 2203, 2204, 2229, 2230, 2247, 2257, 2128) £3,439
LINX HEAVY DUTY REMOTE			
JBH2240		Proportional Heavy Duty Remote Joystick w mode switch (ASL133)	(NC 2203, 2204, 2229, 2230, 2247, 2257, 2128) £1,917
LINX MICRO EXTREMITY AND HEAVY DUTY CONTROL BRACKET			
JBH2285		Tray Mount RH (100x150 mm)	(NC 1707, 2213, 2217, 1404, 1405, 1415) ○
JBH2286		Tray Mount LH (100x150 mm)	(NC 1708, 2214, 2218, 1404, 1405, 1339) ○
JBH2283		Centre Mount, swing-away to RH, Nucleus	(NC 1707, 2213, 2217, 1404, 1405, 1415) £474
JBH2284		Centre Mount, swing-away to LH, Nucleus	(NC 1708, 2214, 2218, 1404, 1405, 1339) £474
JBH2287		Centre Mount, swing-away to RH, Nucleus Incl. Tray ASL603, 100x150 mm	(NC 1707, 2213, 2217, 1404, 1405, 1415) £474
JBH2288		Centre Mount, swing-away to LH, Nucleus Incl. Tray ASL603, 100x150 mm	(NC 1708, 2214, 2218, 1404, 1405, 1339) £474
JBX2289		Individual mounting option per request (Special) Requires additional information/specification	POA
LINX 10-WAY SWITCH			
JBH2251		10-Way Switch DLX-SW10-A for direct Powered Seating Access	£480
LINX 10-WAY SWITCH BRACKET			
JBH2217		Bracket RH, mounted for user	(NC 2207, 2287, 2283, 2275, 2213, 2265, 1404) ○
JBH2218		Bracket LH, mounted for user	(NC 2208, 2288, 2284, 2286, 2276, 2214, 2266, 1405) ○
JBH2291		Bracket RH, mounted for attendant	(NC 2259, 1221, 1287, 2108) £60
JBH2292		Bracket LH, mounted for attendant	(NC 2260, 1222, 1288, 2108) £60
LINX USB CHARGER PORT			
JBH2253		USB Charger Port DLX-USB02-A	(NC 2702, 2703, 2704, 2705, 2706) £156

LINX USB CHARGER PORT BRACKET			
JBH2281		RH Mount	○
JBH2282		LH Mount	○
REMOTE STOP SWITCH			
JBH2258		Remote Stop Switch (ASL504) (NC 2203, 2204)	£1,000
LINX IN200 MODULE			
JBH2273		IN200 Module (Jack port for external power button + control inputs via DB9)	£229
FURTHER LINX COMPONENTS / ACCESSORIES			
JBH2300		To be enclosed only	○
JBX2301		To be connected and programmed (Special)	POA
BUDDY BUTTONS			
JBH1685		Red Buddy Button	£83
JBH1686		Yellow Buddy Button	£83
JBH1687		Green Buddy Button	£83
JBH1690		Blue Buddy Button	£83
EGG SWITCH			
JBH2254		Egg Switch, black	£83
PIKO BUTTONS			
JBH1711		Red Piko Button	£150
JBH1712		Yellow Piko Button	£150
JBH1713		Green Piko Button	£150
JBH1714		Orange Piko Button	£150
JBH1715		Blue Piko Button	£150
JBH1716		Black Piko Button	£150
BODYPOINT JOYSTICK KNOBS			
JBH2721		3 inch U-Shaped, Joystick Handle, Bodypoint PC101A	£59
JBH2722		4 inch U-Shaped, Joystick Handle, Bodypoint PC102A	£59
JBH2723		Domed Rubber, Joystick Handle, Bodypoint PC107A	£59
JBH2724		Mushroom Grooved, Joystick Handle, Bodypoint PC110A	£59
5:1 SWITCH ADAPTOR			
JBH1681		5:1 Switch Adapter (9-pin sub D to 5 jack ports)	£114
SHROUD COLOUR*			
JBH2871		Mid range - Black/Anthracite Shroud	○
JBH2872		Mid range - Silver Shroud	○
JBH2873		Mid range - Bengal Red Shroud	○
JBH2874		Mid range - Mercury Blue Shroud	○
JBH2875		Mid range - Midnight Purple Shroud	○
JBX0299		Premium - Individual Shroud Colour (includes coloured spreader bars on Modulite back) RAL number required Extended lead lime: + 15 days	£203
JBH2865		Premium - Racing Green Shroud	£99
JBH2866		Premium - Deep Ocean Blue Satin Shroud	£99
JBH2867		Premium - Frosty Blue Shroud	£99
JBH2868		Premium - Ultimate Red Shroud	£99
JBH2869		Premium - Falcon White Shroud	£99

RIM COLOUR*			
JBH2882		Mid range - Black Rim Inserts	○
JBH2883		Mid range - Silver Rim Inserts	○
JBH2884		Mid range - Bengal Red Rim Inserts	○
JBH2885		Mid range - Mercury Blue Rim Inserts	○
JBH2886		Mid range - Midnight Purple Rim Inserts	○
JBX2887		Premium - Individual Rim Insert Colour RAL number required Extended lead lime: + 15 days	○
JBH2888		Premium - Racing Green Rim Inserts	○
JBH2889		Premium - Deep Ocean Blue Satin Rim Inserts	○
JBH2890		Premium - Frosty Blue Rim Inserts	○
JBH2891		Premium - Ultimate Red Rim Inserts	○
JBH2892		Premium - Falcon White Rim Inserts	○
FORKS*			
JBH2631		Single sided fork	£215
JBH2632		Double sided fork Without castor rim insert (enlarges wide base to 656 mm)	◆
FRONT WHEEL*			
JBH0835		9" (225 mm) x 2.80 (70 mm) Castor, pneumatic	○
JBH0836		9" (225 mm) x 2.80 (70 mm) Castor, puncture proof solid tyre	◆
CASTOR LOCKS			
JBH1490		Castor Lock	£188
DRIVE WHEEL*			
JBH0704		DW 14" (355 mm) x 3" (75 mm), Pneumatic	○
JBH0705		DW 14" (355 mm) x 3" (75 mm), Solid inlay Dual density foam insert	◆
C.T.C SUSPENSION*			
JBH2870		CTC Suspension Control - Traction - Comfort	◆
SUSPENSION SETTING*			
JBH1916		Suspension up to 160 kg User Weight For RX 40 HD only	◆
MOTOR TYPE*			
JBH1294		Indoor setting, reduced speed approx. 3 km/h	○
JBH1282		6 km/h	◆
JBH1908		10 km/h Package Always combined with LED Light system with Indicators, 30 Lux headlights	(NC 2203, 1437, 1254, 2837) £517
JBH1909		12 km/h Ultra High Speed Package Always combined with LED Light system with Indicators, 30 Lux headlights	(NC 2201, 2203, 1437, 1254, 2837) £621
BATTERY TYPE*			
JBH1259		60 Ah Gel	◆
JBH1260		73,5 Ah Gel	(NC 1450) £87
G TRAC			
JBH1215		G-Trac Module	£308
CHARGER*			
JBH1343		Opti Charge Charger 8A (IP44)	◆
JBH1313		Opti Charge Charger 10A (IP44)	£144
ANTI-TIPPERS*			
JBH1205		Anti-tipper	◆

LIGHTS			
JBH1437		LED Light system with Indicators 30 Lux headlights (NC 2203, 1908, 1909)	£456
JBH2837		LED Light system extra bright with indicators 60 Lux headlights (NC 2203)	£135
JBH1254		Without Light system (NC 1909)	○
TRANSPORTATION KIT*			
JBH1450		Preparation for DAHL-Docking station (Adequate headrest is compulsory) Includes a mounting bracket in the battery box Docking Station for car implementation has to be ordered separately with DAHL Engineering	£557
JBH1494		4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19) Adequate headrest is compulsory, Not compatible with Seat lifter	◆
JBH1401		6-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19)) required for variants with combined Lift & Tilt module Adequate headrest is compulsory (NC 1091, 1092, 1093)	○
EASY ACCESS FRONT TIE DOWNS			
JBH2894		Easily accessible front tie-down points, additionally mounted for 4- or 6- point-tie-down For transport worthy versions (1494/1401) only	£47
LAP BELT*			
JBH1329		Lap belt with Auto Lock Standard mount on backrest	○
JBH2113		Lap belt with Auto Lock Flexible mount in side rails	○
POSITIONING BELT			
JBH2701		BP Evoflex Pelvic Stabilizer, Rehab Latch Size L Padded length 620 mm (NC 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253)	£73
JBH2702		BP Evoflex Pelvic Stabilizer, Rehab Latch Size M Padded length 460 mm (NC 2123, 2124, 2184, 2185, 2186, 2187, 2188, 2189, 2253)	£73
JBH2703		BP 2-Point Hip Belt, Ctr-Pull, Push Butt, Cinch-Mount Size L Padded length 620 mm (NC 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253)	£75
JBH2704		BP 2-Point Hip Belt, Ctr-Pull, Push Butt, Cinch-Mount Size M Padded length 460 mm (NC 2123, 2124, 2184, 2185, 2186, 2187, 2188, 2189, 2253)	£75
JBH2705		BP 2-Point Hip Belt, Dual-Pull, Side-Release Butt. Cinch-Mount Size L Padded length 620 mm (NC 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253)	£75
JBH2706		BP 2-Point Hip Belt, Dual-Pull, Side-Release Butt. Cinch-Mount Size M Padded length 460 mm (NC 2123, 2124, 2184, 2185, 2186, 2187, 2188, 2189, 2253)	£75
REAR VIEW MIRROR			
JBH1415		Rear view mirror RH (NC 2217, 2285, 2213, 2265, 2275, 2283, 2287)	£39
JBH1339		Rear view mirror LH (NC 2218, 2286, 2214, 2266, 2276, 2284, 2288)	£39
OTHER OPTIONS			
JBH1438		Luggage rack (MAX carry capacity 5kg)	£313
JBH2128		Bag Hooks on Modulite Seat Max load up to 6 kg = impact on max user weight (NC 2172, 2173, 2227, 2228, 2229, 2233, 2240, 2257, 2291, 2292)	£73
JBH1312		Cane holder RH (NC 2227, 2259, 1221)	£39
JBH1363		Cane holder LH (NC 2228, 2260, 1222)	£39
JBH1416		Fender for Castors	£83
JBH1417		Fender for Drive wheels	£160
JBH1404		Tray, Swing Away RH (NC 1707, 2275, 2276, 2287, 2288, 2283, 2284, 2285, 2286, 2265, 2266, 2207, 2213, 2214, 2217)	£158
JBH1405		Tray, Swing Away LH (NC 1708, 2275, 2276, 2287, 2288, 2283, 2284, 2285, 2286, 2265, 2266, 2208, 2213, 2214, 2218)	£158

TOTAL PRICE	
--------------------	--

Non Conforming Options:

JBH2100	Invacare® Modulite Seat Unit
JBX2101	Invacare® Modulite HD Seat Unit On RX 40 HD only Tension Adjustable Upholstry only JBX2171
JBH2123	480 mm Seat Plate
JBH2124	530 mm Seat Plate
JBX2124	530 mm Seat Plate
JBH2106	Body fit Kit Outrigging Backrest Brackets and Pivot Point (outrig: 1 inch/2,54 cm per side)
JBH2132	Seat Depth 410 mm
JBH2163	Backposts short (480 mm) For Tension Adjustable Upholstery/Matrx Backs, measured from seat plate to top of upholstery Width adjustable cast bracket
JBH2108	Backposts without push handles
JBH2171	Tension Adjustable Upholstery nylon overlay short = 480 mm / 7 straps measured from seat plate to top of upholstery 480 mm
JBH2179	Flex3 Tension Adjustable Upholstery height 540 mm (seat plate to top of upholstery)
JBH2172	Flex3 adjustable Backrest Plates (without Backrest Cushion) Flex3 Backrest plate width and height adjustable 560 mm - 790 mm
JBH2173	Flex3 adjustable Backrest Plates with Laguna Backrest cushion Flex3 Backrest plate width and height adjustable 560 mm - 790 mm
JBH2174	Matrx Posture Back, height 400 mm (Mesh fabric) Corresponds to SW setting chosen (up to 136 kg max user weight)
JBH2175	Matrx Posture Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen (up to 136 kg max user weight)
JBH2176	Matrx Elite Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen (up to 136 kg max user weight)
JBH2177	Matrx Elite TR Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen
JBH2178	Matrx Elite Deep Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen (up to 136 kg max user weight)
JBH2170	Without Backrest Upholstery (Back Posts/Back Frame only)
JBX2171	Tension Adjustable Upholstery nylon overlay short = 480 mm / 7 straps measured from seat plate to top of upholstery 480 mm MUST BE SELECTED WITH HD SEATING JBX2102
JBH2151	Black Fabric Upholstery (Nylon with option 2141; 2-way-stretch material with options 2142, 2143)
JBH2155	Black Spacer Fabric Upholstery
JBH2198	Padded Armrest with integrated cloth guard
JBH2193	Flip-up Armrest Up to 150 kg max user weight
JBH2192	Following Armrest Height adjustable from 230 - 300 mm
JBH2194	Following Armrest Height adjustable from 300 - 360 mm
JBH2195	Short Armpad (325 mm)
JBH2196	Long Armpad (375 mm)
JBH2197	Long Ergonomic Armpad (397 mm)
JBH2199	Soft Armpad In combination with 10-way switch, the switch 2251 has to be mounted on the opposite side of the remote
JBX2195	Short Armpad (325 mm)
JBX2196	Long Armpad (375 mm)
JBX2197	Long Ergonomic Armpad (397 mm)

JBH2184	Hip Support LH, pad size 140 x 100 mm
JBH2186	Hip Support LH, pad size 230 x 100 mm
JBH2188	Hip Support LH, pad size 330 x 100 mm
JBH2185	Hip Support RH, pad size 140 x 100 mm
JBH2187	Hip Support RH, pad size 230 x 100 mm
JBH2189	Hip Support RH, pad size 330 x 100 mm
JBH1350	Neckrest depth adjustable
JBH1432	Headrest sideways adjustable
JBH1427	Headrest with cheek support
JBH2323	Elan Standard Pad 150 mm / 6" wide w/standard cover (ESP26) Complete with Elan Headrest Hardware
JBH2324	Elan Standard Pad 255 mm / 10" wide w/standard cover (ESP26) Complete with Elan Headrest Hardware
JBH2322	Elan Standard Pad 355 mm / 14" wide w/standard cover (ESP26) Complete with Elan Headrest Hardware
JBH2325	Elan Occipital Pad 230x180 mm / 9"x7" w/standard cover (EOP23) Complete with Elan Headrest Hardware
JBH2326	Elan Occipital Pad 305x200 mm / 12"x8" w/standard cover Complete with Elan Headrest Hardware
JBH2327	Elan 4 Point Headrest Pad 280x255 mm / 11"x10" w/standard Complete with Elan Headrest Hardware
JBH1496	Non Invacare Headrest but profile mount fitting Flex3 hard back
JBX1350	Neckrest depth adjustable
JBX1432	Headrest sideways adjustable
JBX1427	Headrest with cheek support
JBX2323	Elan Standard Pad 150 mm / 6" wide w/standard cover (ESP26) Complete with Elan Headrest Hardware
JBX2324	Elan Standard Pad 255 mm / 10" wide w/standard cover (ESP26) Complete with Elan Headrest Hardware
JBX2322	Elan Standard Pad 355 mm / 14" wide w/standard cover (ESP26) Complete with Elan Headrest Hardware
JBX2325	Elan Occipital Pad 230x180 mm / 9"x7" w/standard cover (EOP23) Complete with Elan Headrest Hardware
JBX2326	Elan Occipital Pad 305x200 mm / 12"x8" w/standard cover Complete with Elan Headrest Hardware
JBX2327	Elan 4 Point Headrest Pad 280x255 mm / 11"x10" w/standard Complete with Elan Headrest Hardware
JBH1446	Lateral Support with extended arm, rigid pad
JBH1445	Lateral Support with extended arm, multi-adjustable pad
JBH1091	Manually settable Seat Tilt 0 - 15°
JBH1092	Powered "fixed pivot" Tilt 20°
JBH1093	Powered Tilt (0° - 30°) with CoG shift
JBH1094	Seat lifter with powered tilt (300 mm, 30°) with CoG shift
JBX1091	Manually settable Seat Tilt 0 - 15°
JBX1093	Powered Tilt (0° - 30°) with CoG shift
JBX1094	Seat lifter with powered tilt (300 mm, 30°) with CoG shift, 6 point Tie downs Must be selected
JBH0977	Powered Recliner 30°
JBX0977	Powered Recliner 30°
JBH0589	LNx Power centre mount legrest, length compensating Length: 340 - 410 mm
JBH2366	Angle adjustable foot platform S (305 W x 215 mm D / 12"x8.5")
JBH0517	Angle adjustable foot platform M (305 W x 275 mm D / 12"x10.75")
JBH0503	Hanger Vari F pair Lower leg length: 290 - 460 mm Angle adjustment: 70 - 0° with tools
JBH0533	Hanger Vari A pair Lower leg length: 290 - 460 mm Angle adjustment: 70 - 0° with handwheel and calf pad
JBH0531	Hanger ADE pair Lower leg length: 290 - 460 mm Powered Angle adjustment: 80 - 0° with calf pad

JBH0501	Without Hanger pair
JBH0504	Angle Adjustable Footplate pair
JBH0513	Angle and depth adjustable Footplate pair
JBH0522	Angle adjustable Footboard
JBH2201	Remote LiNX DLX-REM110 = Drive only 1 Drive function with Drive only power module
JBH2203	Remote LiNX DLX-REM211 = Drive & Actuators 3 Drive functions and power module incl. 2 Actuators
JBH2204	Remote LiNX DLX-REM216 = Drive, Actuators & Lights 3 Drive functions and power module incl. 2 Actuators & Lights
JBH2205	Rehab Remote LiNX, touch-screen user interface DLX-REM400 Configurable jack sockets
JBH2206	Rehab Remote LiNX without joystick for alternative access users DLX-REM500 Requires a user input/alternative joystick: 2231, 2232, 2229, 2230, 2247, 2233
JBH2270	Assignable Toggle Switch Module f. REM400 (DLX-SW400TGL)
JBH1707	Swing Away Bracket RH
JBH1708	Swing Away Bracket LH
JBH2213	Centre Mount, swing-away to RH, Nucleus
JBH2214	Centre Mount, swing-away to LH, Nucleus
JBH2227	Bracket RH, mounted for Attendant
JBH2228	Bracket LH, mounted for Attendant
JBH2220	Dual Control LiNX (DLX-REM050 as secondary remote)
JBH2219	Compact Dual Control LiNX (DLX-ACU200)
JBH2245	Invacare LiNX IDC Dual control
JBH2259	Bracket RH, Attendant
JBH2260	Bracket LH, Attendant
JBH1221	Bracket RH, detachable
JBH1222	Bracket LH, detachable
JBH2231	Proportional Compact Remote DLX-CR400
JBH2232	Proportional Compact Remote Low Force DLX-CR-400LF
JBH2207	RH Mount
JBH2208	LH Mount
JBH2265	Centre Mount, swing-away to RH (Nucleus) If several joysticks (or Joystick + REM500) should be mounted on one Nucleus, please select special mounting option 2215!
JBH2266	Centre Mount, swing-away to LH (Nucleus) If several joysticks (or Joystick + REM500) should be mounted on one Nucleus, please select special mounting option 2215!
JBH2229	Proportional Chin Joystick (ASL138) SELECTED HEADREST ALWAYS WITH STD MOUNTING HARDWARE
JBH1289	Swing Away Bracket RH, side manual (ASL628) Requires Elan Headrest
JBH1288	Swing Away Bracket LH, side manual (ASL628) Requires Elan Headrest
JBH1287	Swing Away Bracket RH, powered
JBH1286	Swing Away Bracket LH, powered
JBH2230	Head Control with ultrasonic proximity sensors (ATOM, ASL104) with standard wings M3
JBH2247	Head Control with ultrasonic proximity sensors (PROTON, ASL104P) with PROTON wings
JBH2257	Four-Switch Proximity Array (ASL106)
JBH2275	Centre Mount, swing-away to RH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm)
JBH2276	Centre Mount, swing-away to LH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm)
JBH2233	MEC Proportional Joystick (Micro Extremity Control ASL130)

JBH2240	Proportional Heavy Duty Remote Joystick w mode switch (ASL133)
JBH2285	Tray Mount RH (100x150 mm)
JBH2286	Tray Mount LH (100x150 mm)
JBH2283	Centre Mount, swing-away to RH, Nucleus
JBH2284	Centre Mount, swing-away to LH, Nucleus
JBH2287	Centre Mount, swing-away to RH, Nucleus Incl. Tray ASL603, 100x150 mm
JBH2288	Centre Mount, swing-away to LH, Nucleus Incl. Tray ASL603, 100x150 mm
JBH2217	Bracket RH, mounted for user
JBH2218	Bracket LH, mounted for user
JBH2291	Bracket RH, mounted for attendant
JBH2292	Bracket LH, mounted for attendant
JBH2253	USB Charger Port DLX-USB02-A
JBH2258	Remote Stop Switch (ASL504)
JBH2273	IN200 Module (Jack port for external power button + control inputs via DB9)
JBH1908	10 km/h Package Always combined with LED Light system with Indicators, 30 Lux headlights
JBH1909	12 km/h Ultra High Speed Package Always combined with LED Light system with Indicators, 30 Lux headlights
JBH1437	LED Light system with Indicators 30 Lux headlights
JBH2837	LED Light system extra bright with indicators 60 Lux headlights
JBH1254	Without Light system
JBH1450	Preparation for DAHL-Docking station (Adequate headrest is compulsory) Includes a mounting bracket in the battery box Docking Station for car implementation has to be ordered separately with DAHL Engineering
JBH1401	6-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19)) required for variants with combined Lift & Tilt module Adequate headrest is compulsory
JBH2701	BP Evoflex Pelvic Stabilizer, Rehab Latch Size L Padded length 620 mm
JBH2702	BP Evoflex Pelvic Stabilizer, Rehab Latch Size M Padded length 460 mm
JBH2703	BP 2-Point Hip Belt, Ctr-Pull, Push Butt, Cinch-Mount Size L Padded length 620 mm
JBH2704	BP 2-Point Hip Belt, Ctr-Pull, Push Butt, Cinch-Mount Size M Padded length 460 mm
JBH2705	BP 2-Point Hip Belt, Dual-Pull, Side-Release Butt. Cinch-Mount Size L Padded length 620 mm
JBH2706	BP 2-Point Hip Belt, Dual-Pull, Side-Release Butt. Cinch-Mount Size M Padded length 460 mm
JBH1415	Rear view mirror RH
JBH1339	Rear view mirror LH
JBH2128	Bag Hooks on Modulite Seat Max load up to 6 kg = impact on max user weight
JBH1312	Cane holder RH
JBH1363	Cane holder LH
JBH1416	Fender for Castors
JBH1404	Tray, Swing Away RH
JBH1405	Tray, Swing Away LH



For a more complex customer requirement please refer to our Invacare Unique special order form. email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: _____ Email: _____

Phone: _____ Fax: _____

Description of special request (please include relevant option code or part number):