



Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

MARCH 2021
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Action 3NG Self-Propel

Standard Features

- Max user weight 125kg (19.7 stone)
- Crash tested ISO- 7176-Part 19
- Self Propulsion with 22" or 24" Quick release rear wheels
- Choice of 6 seat widths 380mm to 505mm (15" to 20")
- Choice of 2 seat depths 400mm or 450mm (16" or 18")
- Choice of 6 seat heights 385mm to 510mm (15" to 20")
- Black nylon upholstery
- Fixed 7° 430mm backrest
- Lumbar support padded backrest with shaped cover
- 80° Swing in / Swing out legrest hangers
- Flip up footplates
- Swing away and detachable desk armrests
- Choice of 8 colours (with 7 special colours)
- Push or pull to lock brakes (configuration dependant)
- Transportation kit and buckle belt (adjustable both sides)
- Adjustable heel straps
- User tool kit



Standard Base price: £592

Please make your selection with a cross in the box Free of charge standard feature POA special request

Customer weight = KG

Code	Description	
PROPULSION*		
CEE0100	<input type="checkbox"/> SELF PROPULSION	£592

SEAT WIDTH*			
CEE44444		13" / 330 MM SEAT WIDTH - SELECT SEAT WIDTH OPTION CEE0180 AND IN-RIGGED ARMREST OPTION CEE0952 (lapbelt adjustment 420mm - 1520mm)	£38
CEE0180		15" / 380 MM SEAT WIDTH (lapbelt adjustment 420mm - 1520mm)	○
CEE0190		16" / 405 MM SEAT WIDTH (lapbelt adjustment 420mm - 1520mm)	○
CEE0200		17" / 430 MM SEAT WIDTH (lapbelt adjustment 420mm - 1520mm)	◆
CEE0210		18" / 455 MM SEAT WIDTH (lapbelt adjustment 420mm - 1880mm)	○
CEE0220		19" / 480 MM SEAT WIDTH (lapbelt adjustment 420mm - 1880mm)	○
CEE0230		20" / 505 MM SEAT WIDTH (lapbelt adjustment 420mm - 1880mm)	○
SEAT DEPTH*			
CEE0240		16" / 400 MM SEAT DEPTH (NC 0879)	○
CEE0250		18" / 450 MM SEAT DEPTH (NC 0878)	◆
CEE0252		16" TO 18" / 400 TO 450 MM ADJUSTABLE SEAT DEPTH (NC 0330, 0350)	£19
SEAT HEIGHT*			
CEE0260		20" / 510 MM SEAT HEIGHT (NC 0511, 0561, 0738, 1035, 1050, 1054, 1056, 1066)	○
CEE0270		19" / 485 MM SEAT HEIGHT (NC 0511, 0561, 1054)	○
CEE0280		18" / 460 MM SEAT HEIGHT	◆
CEE0290		17" / 435 MM SEAT HEIGHT (NC 1010, 1030)	○
CEE0300		16" / 410 MM SEAT HEIGHT (NC 0935, 1010, 1030, 1056, 1370, 1440, 1441)	○
CEE0310		15" / 385 MM SEAT HEIGHT , (NC 0511, 0561, 0660, 0670, 0671, 0738, 0760, 0770, 0838, 0839, 0935, 1010, 1030, 1035, 1050, 1056, 1068, 1280, 1370, 1440, 1441)	○
SEAT UPHOLSTERY*			
CEE0330		NYLON W/O HOOK AND LOOP STRIPS	◆
CEE0340		NYLON WITH HOOK AND LOOP STRIPS	○
CEE0350		PADDED NYLON	£28

BACKREST - Anti-tippers should be selected with reclining backrests 0435, 0435, 0440, 0445 and 0450*			
CEE0380		FIXED 0° BACKREST 2 HEIGHTS 16"/ 400 MM or 17" / 430 MM	(NC 0471, 0477, 0945) ○
CEE0390		FIXED 0° BACKREST 20" / 510 MM HEIGHT	(NC 0471, 0477, 0945) ○
CEE0395		FIXED 7° BACKREST HT ADJ 15" / 380 MM TO 18" / 460 MM FIX PUSH HANDLES	(NC 0474, 0477, 0478, 0481, 0905, 0906, 0907, 0940) £55
CEE0396		FIXED 7° BACKREST HEIGHT ADJUST 15" / 380 TO 18" / 460 MM HT ADJ PUSH HANDLES	(NC 0474, 0477, 0478, 0481, 0905, 0906, 0907, 0940) £74
CEE0410		FIXED 7° BACKREST 2 HEIGHTS 16" / 400 or 17" / 430 MM FIX PUSH HANDLES	(NC 0471, 0477, 0945) ◆
CEE0420		FIXED 7° BACKREST 20" / 510 MM HEIGHT FIX PUSH HANDLES	(NC 0471, 0477, 0945) ○
CEE0425		FIXED 7° BACKREST 17" / 430 MM HEIGHT HT ADJ PUSHHANDLE	(NC 0471, 0477, 0945) £40
CEE0428		FIXED 7° BACKREST ACTIVE 14" / 350 MM HT ADJ TO 16" / 400 MM FIX PUSH HANDLES	(NC 0471, 0477, 0481, 0905, 0906, 0907, 0940, 1280) ○
CEE0432		ANGLE ADJUSTABLE AND FOLD DOWN (15° to + 15°) 17"	(NC 0465, 0471, 0474, 0478, 0481, 0882, 0952) £80
CEE0435		RECLINING (RATCHET) 0/30° BACKREST 3 HEIGHTS 15" 16" 17"	(NC 0471, 0477, 0738, 0838, 0940, 0945, 1190) £162
CEE0440		RECLINING (RATCHET) 0/30° BACKREST 3 HEIGHTS 18" 19" 20"	(NC 0471, 0477, 0738, 0838, 0940, 0945, 1190) £162
CEE0445		RECLINING (GAS SPRING) 0/30° BACKREST 3 HEIGHTS 15" 16" 17"	(NC 0471, 0477, 0738, 0838, 0940, 0945, 1190) £215
CEE0450		RECLINING (GAS SPRING) 0/30° 3 HEIGHTS 18" 19" 20"	(NC 0471, 0477, 0738, 0838, 0940, 0945, 1190) £215
CEE0455		FOLDING 7° BACKREST 17"	(NC 0471, 0477, 0945) £57
CEE0460		FOLDING 7° BACKREST 20"	(NC 0471, 0477, 0481, 0945) £57
BACKREST UPHOLSTERY*			
CEE0465		W/O BACK UPH	(NC 0882) ○
CEE0471		STANDARD NYLON SLING (JEMIMA)	(NC 0882, 0905, 0906, 0907) ○
CEE0474		LUMBAR SUPPORT BI MATERIAL PADDED COVER	○
CEE0477		CONTOUR ADJUSTABLE (3+1 STRAPS) BI MATERIAL PADDED COVER	£50
CEE0478		CONTOUR ADJUSTABLE (3+1 STRAPS) BI MATERIAL PADDED COVER INCL BACKBRACE	(NC 0940, 0945) £57
CEE0481		CONTOUR ADJUSTABLE (7+1 STRAPS) BI MATERIAL PADDED COVER INCL BACKBRACE	(NC 0940, 0945) £111
LEGREST HANGER RH*			
CEE0510		80° RH SWING IN/SWING OUT	(NC 1111, 66666) ◆
CEE0511		90° RH SWING IN/SWING OUT	(NC 0600, 0936, 1010, 1030, 1035, 1050, 1056, 1111, 66666) ○
CEE0520		ELEVATING LEGREST RH SWING IN/SWING OUT	(NC 0611, 0930, 0935, 0936, 66666) £107
CEE0545		AMPUTEE LEGREST RH	(NC 0600, 0610, 0611, 0930, 0936, 1111) £123
LEGREST HANGER LH*			
CEE0560		80° LH SWING IN/SWING OUT	(NC 2222) ◆
CEE0561		90° LH SWING IN/SWING OUT	(NC 0620, 0936, 1010, 1030, 1035, 1050, 1056, 2222) ○
CEE0570		ELEVATING LEGREST LH SWING IN/SWING OUT	(NC 0611, 0930, 0935, 0936, 55555) £107
CEE0595		AMPUTEE LEGREST LH	(NC 0611, 0620, 0630, 0930, 0936, 2222) £123

FOOTREST RH*			
CEE0600		FLIP UP RH FOOTPLATE	◆
CEE0610		FLIP UP DEPTH AND ANGLE AJUSTABLE RH FOOTPLATE NB: IF SELECTED WITH A FIXED FOOTREST FOOTPLATE WILL BE ANGLE AND DEPTH ADJ AND COMPATIBLE WITH ANKLE HUGGERS. IF SELECTED WITH AN ELEVATING LEGREST FOOTPLATE WILL BE ANGLE, DEPTH AND WIDTH ADJ NOT COMPATIBLE WITH ANKLE HUGGERS.	£40
SWO1111		FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY)	£40
CEE0611		FLIP UP ONE PIECE (NC 0620, 0630, 0935, 0936)	£39
CEE66666		NO FOOTREST	○
FOOTREST LH*			
CEE0620		FLIP UP LH FOOTPLATE	◆
CEE0630		FLIP UP DEPTH AND ANGLE AJUSTABLE LH FOOTPLATE NB: IF SELECTED WITH A FIXED FOOTREST FOOTPLATE WILL BE ANGLE AND DEPTH ADJ AND COMPATIBLE WITH ANKLE HUGGERS. IF SELECTED WITH AN ELEVATING LEGREST FOOTPLATE WILL BE ANGLE, DEPTH AND WIDTH ADJ NOT COMPATIBLE WITH ANKLE HUGGERS.	£40
SWO2222		FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY)	£40
CEE55555		NO FOOTREST	○

ARMREST RH*				
CEE0640		REMOVABLE RH (N°1) DESK ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED	(NC 0760, 0770, 0780, 0790, 0795, 0796, 0838, 0839, 1480, 1483)	○
CEE0650		REMOVABLE RH (N°1) FULL LENGTH ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED	(NC 0760, 0770, 0780, 0790, 0795, 0796, 0838, 0839)	○
CEE0655		REMOVABLE RH (N°1) "RIO" Special Only For Denmark AUTOLOCK SAFETY DEVICE INCLUDED FULL LENGTH ADJ 190290 MM HEIGHT by step of 15 MM	(NC 0740, 0750, 0760, 0770, 0780, 0790, 0795, 0796, 0955, 0838, 0839, 0952, 1480, 1483)	○
CEE0656		SHORT CONCAVE ARMPAD NOT AVAILABLE YET	(NC 0760, 0770, 0780, 0790, 0795, 0796, 0838, 0839, 1480, 1483)	○
CEE0657		ROTATIVE HEMI ARMPAD WITH PALMAR SUPPORT (360°)	(NC 0838, 0839, 1480, 1483)	○
CEE0660		SWING AWAY DETACHABLE RH (N°2) DESK 240 MM HEIGHT	(NC 0740, 0750, 0780, 0790, 0795, 0796, 0838, 0839, 0952, 0955, 1480, 1483)	◆
CEE0670		SWING AWAY DETACHABLE RH (N°2) FULL LENGTH 240 MM HEIGHT	(NC 0740, 0750, 0780, 0790, 0795, 0796, 0838, 0839, 0952, 0955)	○
CEE0671		SHORT CONCAVE ARMPAD NOT AVAILABLE YET	(NC 0740, 0750, 0780, 0790, 0795, 0796, 0838, 0839, 0952, 0955, 1480, 1483)	○
CEE0672		SHORT CONCAVE ARMPAD NOT AVAILABLE YET	(NC 0740, 0750, 0760, 0770, 0795, 0796, 0838, 0839, 0952, 0955, 1480, 1483)	○
CEE0675		CONCAVE ARMPAD	(NC 0740, 0750, 0760, 0770, 0795, 0796, 0838, 0839, 0952, 0955, 1480, 1483)	○
CEE0677		ROTATIVE HEMI ARMPAD WITH PALMAR SUPPORT (360°)	(NC 0838, 0839, 0952, 1480, 1483)	○
CEE0680		SWING AWAY DETACHABLE RH (N°3) DESK ADJ. IN HEIGHT (170 to 295 MM)	(NC 0740, 0750, 0760, 0770, 0795, 0796, 0838, 0839, 0952, 0955, 1480, 1483)	£51
CEE0690		SWING AWAY DETACHABLE RH (N°3) FULL LENGTH ADJ. IN HEIGHT (170 to 295 MM) AND DEPTH	(NC 0740, 0750, 0760, 0770, 0795, 0796, 0838, 0839, 0952, 0955, 1483)	£51
CEE0695		REMOVABLE RH (N°4) DESK ADJ. IN HEIGHT (210345 MM) + SKIRTGUARD HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED	(NC 0740, 0750, 0760, 0770, 0780, 0790, 0838, 0839, 0952, 0955, 1480, 1483)	£55
CEE0696		REMOVABLE RH (N°4) FULL LENGTH ADJ. IN HEIGHT (210345 MM) AND DEPTH + SKIRTGUARD HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED	(NC 0740, 0750, 0760, 0770, 0780, 0790, 0838, 0839, 0952, 0955, 1483)	£55
CEE0732		REMOVABLE RH (N°4) ARMPAD HEMI BLACK MULTI STRETCH PU COATED FABRIC MULTI POSITION AUTOLOCK SAFETY DEVICE INCLUDED	(NC 0838, 0839, 0952, 1480, 1483)	£54
CEE0738		REMOVABLE MUDGUARD RH (N°6) (CUSHION IS RECOMMENDED FOR SEAT HEIGHT 410 AND 435 MM) HEIGHT ADJUSTABLE (155 TO 220 MM)	(NC 0740, 0750, 0760, 0770, 0780, 0790, 0795, 0796, 0832, 0839, 0952, 0955, 1066, 1200, 1201, 1480, 1483)	£32

ARMREST LH*			
CEE0740		REMOVABLE LH (N°1) DESK ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED (NC 1480, 1484)	○
CEE0750		REMOVABLE LH (N°1) FULL LENGTH ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED	○
CEE0760		SWING AWAY DETACHABLE LH (N°2) DESK 240 MM HEIGHT (NC 0952, 0955, 1480, 1484)	◆
CEE0770		SWING AWAY DETACHABLE LH (N°2) FULL LENGTH 240 MM HEIGHT (NC 0952, 0955)	○
CEE0780		SWING AWAY DETACHABLE LH (N°3) DESK ADJ. IN HEIGHT (170 to 295 MM) (NC 0952, 0955, 1480, 1484)	£51
CEE0790		SWING AWAY DETACHABLE LH (N°3) FULL LENGTH ADJ. IN HEIGHT (170 to 295 MM) AND DEPTH (NC 0952, 0955, 1484)	£51
CEE0795		REMOVABLE LH (N°4) DESK ADJ. IN HEIGHT (210345 MM) + SKIRTGUARD HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED (NC 0952, 0955, 1480, 1484)	£55
CEE0796		REMOVABLE LH (N°4) FULL LENGTH ADJ. IN HEIGHT (210345 MM) AND DEPTH + SKIRTGUARD HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED (NC 0952, 0955, 1484)	£55
CEE0832		REMOVABLE LH N°4 ARMPAD HEMI BLACK MULTI STRETCH PU COATED FABRIC MULTI POSITION ARMPAD HEMI MULTI POSITION (NC 0952, 1480, 1484)	£54
CEE0838		REMOVABLE MUDGUARD LH (N°6) (CUSHION IS RECOMMENDED FOR SEAT HEIGHT 410 AND 435 MM) HEIGHT ADJUSTABLE (155 TO 220 MM) (NC 0952, 0955, 1066, 1200, 1201, 1484)	£32
CEE0839		TUBULAR ARMREST LH, SWING AWAY PADDED TUBE WITH REMOVABLE MUDGUARD LH (N°6) (CUSHION IS RECOMMENDED FOR SEAT HEIGHT 410 AND 435 MM) (NC 0952, 0955, 1200, 1201, 1345, 1430, 1480, 1484)	○
SEAT CUSHION			
CEE0849		BASIC 50 MM FOAM CUSHION NYLON AND ANTISLIP COVER (NC 0952)	£22
SEAT CUSHION DEPTH			
CEE0878		FOR 400 MM SEAT DEPTH	○
CEE0879		FOR 450 MM SEAT DEPTH	○
OTHER SEATING OPTIONS			
CEE0882		LATERAL POSITIONING PADS (PAIR)	£56
CEE0905		NECKREST BLACK MULTI STRETCH PU COATED FABRIC (NC 0906, 0907, 0940, 0945)	£121
CEE0906		HEADREST BLACK MULTI STRETCH PU COATED FABRIC (NC 0907, 0940, 0945)	£146
CEE0907		HEADREST WITH CHEEK SUPPORT BLACK MULTI STRETCH PU COATED FABRIC (NC 0940, 0945)	£169
SWO3333		MULT ADJUSTABLE LATERALS	£262
CEE0930		CALF STRAP	○
CEE0935		LONGER FOOTREST TUBE	£28
CEE0936		SHORTER FOOTREST SUPPORT	£39
CEE0940		BACKBRACE (NC 0945)	£32
CEE0945		FOLDABLE BACKBRACE	£36
CEE0952		INRIGGED ARMREST (NC 1480, 1483, 1484)	£39
CEE0955		W/O AUTOLOCK SAFETY DEVICE FOR ARMREST N°1 ONLY	○

COLOUR*			
CEE0961		SLATE GREY (RAL7024 FINE TEXTURE)	○
CEE0964		PEARL WHITE	○
CEE0965		ELECTRIC GREEN	○
CEE0966		MATT BLACK	○
CEE0967		BLUE POP	○
CEE0968		HYPNOTIC PURPLE	○
CEE0969		RED MATT	○
CEE0980		ANTHRACITE BLACK	○
CEE0991		SPECIAL COLOUR ON DEMAND ONLY	£156
SWO4444		DARK RED SPECIAL	£156
SWO5555		PEARL GREY SPECIAL	£156
SWO6666		INTENSE YELLOW	£156
SWO7777		VIBRANT ORANGE	£156
SWO8888		POLAR BLUE	£156
SWO9999		CANDY PINK	£156
CASTOR HOUSING*			
CEE0995		STANDARD	◆
CEE0996		LIGHT WEIGHT	£26
FRONT CASTORS*			
CEE1010		8" X 1.1/4" BLACK SOLID TYRE	◆
CEE1030		8" X 2" BLACK SOLID TYRE	£26
CEE1035		6" X 1.1/4" BLACK SOLID TYRE MEDIUM HARD	£26
CEE1050		6" X 1.1/4" BLACK SOLID TYRE	£26
CEE1054		5" X 1" BLACK SOLID TYRE	£26
CEE1056		5.1/2" X 1.3/4" BLACK SOLID TYRE MEDIUM HARD	£26
REAR WHEEL SIZE*			
CEE1066		22" (NC 1092, 1260)	○
CEE1068		24" (NC 1261)	◆
REAR WHEELS*			
CEE1071		STANDARD (BASIC) (BLACK HUB, COLORLESS RIM AND ZINC PLATTED SPOKES) (NC 1092, 1280)	◆
CEE1072		SELF PROPEL REINFORCED REAR WHEEL (NC 1092, 1280)	£21
CEE1076		LIGHTWEIGHT (NC 1170, 1201, 1230, 1260, 1261, 1280)	£54
TYRES*			
CEE1082		SOLID BLACK	◆
CEE1084		PNEUMATIC PROFILE BLACK	£19
CEE1088		PNEUMATIC LOW PROFILE BLACK	○
CEE1092		PNEUMATIC SCHWALBE MARATHON PLUS EVO	£36
HANDRIM*			
CEE1140		ALUMINIUM PREANODIZED HANDRIM	◆
CEE1150		HANDRIM COATED IN BLACK ANTISLIP PLASTIC	£31
REAR WHEEL AXLE*			
CEE1170		FIXED AXLE	○
CEE1180		Q/R AXLE	◆

REAR WHEEL POSITION*			
CEE1190	<input type="checkbox"/>	STANDARD REAR WHEELS POSITION	◆
CEE1200	<input type="checkbox"/>	PASSIVE REAR WHEELS POSITION	○
CEE1201	<input type="checkbox"/>	DUAL AMPUTEE REAR WHEEL POSITION (NC 1275, 1280)	£67
PUMP			
CEE1230	<input type="checkbox"/>	PUMP AND SUPPORT	£13
CEE1235	<input type="checkbox"/>	PUMP HIGH PRESSURE	£38
CHASSIS OPTION			
CEE1241	<input type="checkbox"/>	FRONT WHEELS EXTENSION	£20
CEE1260	<input type="checkbox"/>	SPOKE GUARDS PAIR 24" BLK (NC 1261)	£72
CEE1261	<input type="checkbox"/>	SPOKE GUARDS PAIR 560 AND 515 MM (22" AND 20") CLR	£90
BRAKES*			
CEE1270	<input type="checkbox"/>	PUSH TO LOCK	◆
CEE1275	<input type="checkbox"/>	PULL TO LOCK BRAKES	○
CEE1280	<input type="checkbox"/>	ATTENDANT HUB BRAKES	£213
TRANSPORTATION KIT*			
CEE1345	<input type="checkbox"/>	TRANSPORTATION BRACKETS	◆
SAFETY OPTIONS*			
CEE1360	<input type="checkbox"/>	POSTURE BELT (NC 1361, 1362)	◆
CEE1361	<input type="checkbox"/>	BODYPOINT 2PT POSTURE BELT M (NC 1362)	£74
CEE1362	<input type="checkbox"/>	BODYPOINT 2PT POSTURE BELT L	£74
CEE1370	<input type="checkbox"/>	REAR ANTITIPPERS PAIR (NC 1375, 1440, 1441)	£33
CEE1375	<input type="checkbox"/>	REAR ANTITIPPERS PAIR WITH TIPPING AID (NC 1440, 1441)	£67
CEE1380	<input type="checkbox"/>	REFLECTORS KIT	£7
CEE1390	<input type="checkbox"/>	BRAKE EXTENSION PAIR	£13
OTHER OPTIONS			
CEE1430	<input type="checkbox"/>	CANE HOLDER	£21
CEE1440	<input type="checkbox"/>	TIPPING AID RH	£13
CEE1441	<input type="checkbox"/>	TIPPING AID LH	£13
CEE1480	<input type="checkbox"/>	ACRYLIC TRAY (NC 1483, 1484, 1485)	£107
CEE1483	<input type="checkbox"/>	HALF TABLE TRAY RH (NC 1484)	£82
CEE1484	<input type="checkbox"/>	HALF TABLE TRAY LH	£82
CEE1485	<input type="checkbox"/>	COMFORT PADDING FOR HALF TABLE TRAY	£41
SWO10002	<input type="checkbox"/>	ABDUCTED LEGRESTS (PR) 0° TO 30° - INTERNAL REFERENCE DE259618	£99

TOTAL PRICE	
--------------------	--

Non Conforming Options:

CEE0330	NYLON W/O HOOK AND LOOP STRIPS
CEE0350	PADDED NYLON
CEE0465	W/O BACK UPH
CEE0471	STANDARD NYLON SLING (JEMIMA)
CEE0474	LUMBAR SUPPORT BI MATERIAL PADDED COVER
CEE0477	CONTOUR ADJUSTABLE (3+1 STRAPS) BI MATERIAL PADDED COVER
CEE0478	CONTOUR ADJUSTABLE (3+1 STRAPS) BI MATERIAL PADDED COVER INCL BACKBRACE
CEE0481	CONTOUR ADJUSTABLE (7+1 STRAPS) BI MATERIAL PADDED COVER INCL BACKBRACE
CEE0511	90° RH SWING IN/SWING OUT
CEE0561	90° LH SWING IN/SWING OUT
CEE0600	FLIP UP RH FOOTPLATE
CEE0610	FLIP UP DEPTH AND ANGLE AJUSTABLE RH FOOTPLATE NB: IF SELECTED WITH A FIXED FOOTREST FOOTPLATE WILL BE ANGLE AND DEPTH ADJ AND COMPATIBLE WITH ANKLE HUGGERS. IF SELECTED WITH AN ELEVATING LEGREST FOOTPLATE WILL BE ANGLE, DEPTH AND WIDTH ADJ NOT COMPATIBLE WITH ANKLE HUGGERS.
SWO1111	FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY)
CEE0611	FLIP UP ONE PIECE
CEE66666	NO FOOTREST
CEE0620	FLIP UP LH FOOTPLATE
CEE0630	FLIP UP DEPTH AND ANGLE AJUSTABLE LH FOOTPLATE NB: IF SELECTED WITH A FIXED FOOTREST FOOTPLATE WILL BE ANGLE AND DEPTH ADJ AND COMPATIBLE WITH ANKLE HUGGERS. IF SELECTED WITH AN ELEVATING LEGREST FOOTPLATE WILL BE ANGLE, DEPTH AND WIDTH ADJ NOT COMPATIBLE WITH ANKLE HUGGERS.
SWO2222	FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY)
CEE55555	NO FOOTREST
CEE0660	SWING AWAY DETACHABLE RH (N°2) DESK 240 MM HEIGHT
CEE0670	SWING AWAY DETACHABLE RH (N°2) FULL LENGTH 240 MM HEIGHT
CEE0671	SHORT CONCAVE ARMPAD NOT AVAILABLE YET
CEE0738	REMOVABLE MUDGUARD RH (N°6) (CUSHION IS RECOMMENDED FOR SEAT HEIGHT 410 AND 435 MM) HEIGHT ADJUSTABLE (155 TO 220 MM)
CEE0740	REMOVABLE LH (N°1) DESK ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED
CEE0750	REMOVABLE LH (N°1) FULL LENGTH ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED
CEE0760	SWING AWAY DETACHABLE LH (N°2) DESK 240 MM HEIGHT
CEE0770	SWING AWAY DETACHABLE LH (N°2) FULL LENGTH 240 MM HEIGHT
CEE0780	SWING AWAY DETACHABLE LH (N°3) DESK ADJ. IN HEIGHT (170 to 295 MM)
CEE0790	SWING AWAY DETACHABLE LH (N°3) FULL LENGTH ADJ. IN HEIGHT (170 to 295 MM) AND DEPTH
CEE0795	REMOVABLE LH (N°4) DESK ADJ. IN HEIGHT (210345 MM) + SKIRTGUARD HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED
CEE0796	REMOVABLE LH (N°4) FULL LENGTH ADJ. IN HEIGHT (210345 MM) AND DEPTH + SKIRTGUARD HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED
CEE0832	REMOVABLE LH N°4 ARMPAD HEMI BLACK MULTI STRETCH PU COATED FABRIC MULTI POSITION ARMPAD HEMI MULTI POSITION
CEE0838	REMOVABLE MUDGUARD LH (N°6) (CUSHION IS RECOMMENDED FOR SEAT HEIGHT 410 AND 435 MM) HEIGHT ADJUSTABLE (155 TO 220 MM)
CEE0839	TUBULAR ARMREST LH, SWING AWAY PADDED TUBE WITH REMOVABLE MUDGUARD LH (N°6) (CUSHION IS RECOMMENDED FOR SEAT HEIGHT 410 AND 435 MM)
CEE0878	FOR 400 MM SEAT DEPTH
CEE0879	FOR 450 MM SEAT DEPTH
CEE0882	LATERAL POSITIONING PADS (PAIR)
CEE0905	NECKREST BLACK MULTI STRETCH PU COATED FABRIC

CEE0906	HEADREST BLACK MULTI STRETCH PU COATED FABRIC
CEE0907	HEADREST WITH CHEEK SUPPORT BLACK MULTI STRETCH PU COATED FABRIC
CEE0930	CALF STRAP
CEE0935	LONGER FOOTREST TUBE
CEE0936	SHORTER FOOTREST SUPPORT
CEE0940	BACKBRACE
CEE0945	FOLDABLE BACKBRACE
CEE0952	INRIGGED ARMREST
CEE0955	W/O AUTOLOCK SAFETY DEVICE FOR ARMREST N°1 ONLY
CEE1010	8" X 1.1/4" BLACK SOLID TYRE
CEE1030	8" X 2" BLACK SOLID TYRE
CEE1035	6" X 1.1/4" BLACK SOLID TYRE MEDIUM HARD
CEE1050	6" X 1.1/4" BLACK SOLID TYRE
CEE1054	5" X 1" BLACK SOLID TYRE
CEE1056	5.1/2" X 1.3/4" BLACK SOLID TYRE MEDIUM HARD
CEE1066	22"
CEE1068	24"
CEE1092	PNEUMATIC SCHWALBE MARATHON PLUS EVO
CEE1170	FIXED AXLE
CEE1190	STANDARD REAR WHEELS POSITION
CEE1200	PASSIVE REAR WHEELS POSITION
CEE1201	DUAL AMPUTEE REAR WHEEL POSITION
CEE1230	PUMP AND SUPPORT
CEE1260	SPOKE GUARDS PAIR 24" BLK
CEE1261	SPOKE GUARDS PAIR 560 AND 515 MM (22" AND 20") CLR
CEE1275	PULL TO LOCK BRAKES
CEE1280	ATTENDANT HUB BRAKES
CEE1345	TRANSPORTATION BRACKETS
CEE1361	BODYPOINT 2PT POSTURE BELT M
CEE1362	BODYPOINT 2PT POSTURE BELT L
CEE1370	REAR ANTITIPPERS PAIR
CEE1375	REAR ANTITIPPERS PAIR WITH TIPPING AID
CEE1430	CANE HOLDER
CEE1440	TIPPING AID RH
CEE1441	TIPPING AID LH
CEE1480	ACRYLIC TRAY
CEE1483	HALF TABLE TRAY RH
CEE1484	HALF TABLE TRAY LH
CEE1485	COMFORT PADDING FOR HALF TABLE TRAY



For a more complex customer requirement please refer to our Invacare Unique special order form. email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: _____ Email: _____

Phone: _____ Fax: _____

Description of special request (please include relevant option code or part number):