



Manual Retail Prescription Form

Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

MARCH 2021
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Action 3NG Transit

Standard Features

- Max user weight 125kg (19.7 stone)
- Crash tested ISO 7176-Part 19
- Attendant Propulsion
- Choice of 6 seat widths 380mm to 505mm (15" to 20")
- Choice of 2 seat depths 400mm or 450mm (16" or 18)
- Choice of 3 seat heights 460mm to 510mm (18" to 20")
- Black nylon upholstery
- Fixed 7° 430mm backrest
- Lumbar support padded backrest with shaped cover
- 80° Swing in / swing out legrest hangers
- Flip up footplates
- Swing away and detachable desk armrests
- Choice of 8 colours (with 7 special colours)
- Push to lock brakes
- Transportation kit and buckle belt (adjustable both sides)
- Adjustable heel straps
- User tool kit



Standard Base price: £592

Please make your selection with a cross in the box Free of charge standard feature POA special request

Customer weight = KG

Code	Description	
BASE PRODUCT*		
CEF0004	ACTION 3NG TRANSIT (NC 1074)	<input type="radio"/>
CEFSPEC	SPECIAL	POA
PROPULSION*		
CEF0110	TRANSIT PROPULSION	£592

SEAT WIDTH*			
CEF44444		13" / 330 MM SEAT WIDTH - SELECT SEAT WIDTH OPTION CEE0180 AND IN-RIGGED ARMREST OPTION CEE0952 (lapbelt adjustment 420mm - 1520mm)	£38
CEF0180		15" / 380 MM SEAT WIDTH	○
CEF0190		16" / 405 MM SEAT WIDTH	○
CEF0200		17" / 430 MM SEAT WIDTH	◆
CEF0210		18" / 455 MM SEAT WIDTH	○
CEF0220		19" / 480 MM SEAT WIDTH	○
CEF0230		20" / 505 MM SEAT WIDTH	○
SEAT DEPTH*			
CEF0240		16" / 400 MM SEAT DEPTH (NC 0879)	○
CEF0250		18" / 450 MM SEAT DEPTH (NC 0878)	◆
CEF0252		16" TO 18" / 400 TO 450 MM ADJUSTABLE SEAT DEPTH (NC 0330, 0350)	£19
SEAT HEIGHT*			
CEF0260		20" / 510 MM SEAT HEIGHT	○
CEF0270		19" / 485 MM SEAT HEIGHT	○
CEF0280		18" / 460 MM SEAT HEIGHT	○
SEAT UPHOLSTERY*			
CEF0330		NYLON W/O HOOK AND LOOP STRIPS	◆
CEF0340		NYLON WITH HOOK AND LOOP STRIPS	○
CEF0350		PADDED NYLON	£27
BACKREST - Anti-tippers should be selected with reclining backrests 0435, 0440, 0445 and 0450*			
CEF0380		FIXED 0° BACKREST 2 HEIGHTS 16" / 400 MM or 17" / 430 MM (NC 0471, 0477, 0945)	○
CEF0390		FIXED 0° BACKREST 20" / 510 MM HEIGHT (NC 0471, 0477, 0945)	○
CEF0395		FIXED 7° BACKREST HT ADJ 15" / 380 MM TO 18" / 460 MM FIX PUSH HANDLES (NC 0474, 0477, 0478, 0481, 0905, 0906, 0907, 0940)	£62
CEF0396		FIXED 7° BACKREST HEIGHT ADJUST 15" / 380 TO 18" / 460 MM HT ADJ PUSH HANDLES (NC 0474, 0477, 0478, 0481, 0905, 0906, 0907, 0940)	£74
CEF0410		FIXED 7° BACKREST 2 HEIGHTS 16" / 400 or 17" / 430 MM FIX PUSH HANDLES (NC 0471, 0477, 0945)	◆
CEF0420		FIXED 7° BACKREST 20" / 510 MM HEIGHT FIX PUSH HANDLES (NC 0471, 0477, 0945)	○
CEF0425		FIXED 7° BACKREST 17" / 430 MM HEIGHT HT ADJ PUSHHANDLE (NC 0471, 0477, 0945)	£40
CEF0432		ANGLE ADJUSTABLE AND FOLD DOWN (15° to + 15°) 17" (NC 0465, 0471, 0474, 0478, 0481, 0882, 0952)	£80
CEF0450		RECLINING (GAS SPRING) 0/30° 3 HEIGHTS 18" 19" 20" (NC 0471, 0477, 0940, 0945)	£215
CEF0455		FOLDING 7° BACKREST 17" (NC 0471, 0477, 0945)	£57
CEF0460		FOLDING 7° BACKREST 20" (NC 0471, 0477, 0481, 0945)	£57
BACKREST UPHOLSTERY*			
CEF0465		W/O BACK UPH (NC 0882)	○
CEF0471		STANDARD NYLON SLING (JEMIMA) (NC 0882, 0905, 0906, 0907)	○
CEF0474		LUMBAR SUPPORT BI MATERIAL PADDED COVER	○
CEF0477		CONTOUR ADJUSTABLE (3+1 STRAPS) BI MATERIAL PADDED COVER	£50
CEF0478		CONTOUR ADJUSTABLE (3+1 STRAPS) BI MATERIAL PADDED COVER INCL BACKBRACE (NC 0940, 0945)	£57
CEF0481		CONTOUR ADJUSTABLE (7+1 STRAPS) BI MATERIAL PADDED COVER INCL BACKBRACE (NC 0940, 0945)	£111

LEGREST HANGER RH*				
CEF0510		80° RH SWING IN/SWING OUT	(NC 1111, 66666)	◆
CEF0520		ELEVATING LEGREST RH SWING IN/SWING OUT	(NC 0611, 0930, 0935, 0936, 66666)	£107
CEF0545		AMPUTE E LEGREST RH	(NC 0600, 0610, 0611, 0930, 0936, 1111)	£123
LEGREST HANGER LH*				
CEF0560		80° LH SWING IN/SWING OUT	(NC 2222)	◆
CEF0595		AMPUTE E LEGREST LH	(NC 0611, 0620, 0630, 0930, 0936, 2222)	£123
FOOTREST RH*				
CEF0600		FLIP UP RH FOOTPLATE		◆
CEF0610		FLIP UP DEPTH AND ANGLE AJUSTABLE RH FOOTPLATE - NB: IF SELECTED WITH A FIXED FOOTREST - FOOTPLATE WILL BE ANGLE AND DEPTH ADJ AND COMPATIBLE WITH ANKLE HUGGERS. IF SELECTED WITH AN ELEVATING LEGREST - FOOTPLATE WILL BE ANGLE, DEPTH AND WIDTH ADJ - NOT COMPATIBLE WITH ANKLE HUGGERS.		£40
SWO1111		FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE - NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY)		£40
CEF0611		FLIP UP ONE PIECE	(NC 0620, 0630, 0935, 0936)	£72
CEF66666		NO FOOTREST		○
FOOTREST LH*				
CEF0620		FLIP UP LH FOOTPLATE		◆
CEF0630		FLIP UP DEPTH AND ANGLE AJUSTABLE LH FOOTPLATE - NB: IF SELECTED WITH A FIXED FOOTREST - FOOTPLATE WILL BE ANGLE AND DEPTH ADJ AND COMPATIBLE WITH ANKLE HUGGERS. IF SELECTED WITH AN ELEVATING LEGREST - FOOTPLATE WILL BE ANGLE, DEPTH AND WIDTH ADJ - NOT COMPATIBLE WITH ANKLE HUGGERS.		£40
SWO2222		FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE - NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY)		£40
CEF55555		NO FOOTREST		○
ARMREST RH*				
CEF0640		REMOVABLE RH (N°1) DESK ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED	(NC 0760, 0770, 0780, 0790, 0795, 0796, 1480, 1483)	○
CEF0650		REMOVABLE RH (N°1) FULL LENGTH ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED	(NC 0760, 0770, 0780, 0790, 0795, 0796)	○
CEF0660		SWING AWAY DETACHABLE RH (N°2) DESK 240 MM HEIGHT	(NC 0740, 0750, 0780, 0790, 0795, 0796, 0952, 0955, 1480, 1483)	◆
CEF0670		SWING AWAY DETACHABLE RH (N°2) FULL LENGTH 240 MM HEIGHT	(NC 0740, 0750, 0780, 0790, 0795, 0796, 0952, 0955)	○
CEF0680		SWING AWAY DETACHABLE RH (N°3) DESK ADJ. IN HEIGHT (170-295 MM)	(NC 0740, 0750, 0760, 0770, 0795, 0796, 0952, 0955, 1480, 1483)	£51
CEF0690		SWING AWAY DETACHABLE RH (N°3) FULL LENGTH ADJ. IN HEIGHT (170-295 MM) AND DEPTH	(NC 0740, 0750, 0760, 0770, 0795, 0796, 0952, 0955, 1483)	£51
CEF0695		REMOVABLE RH (N°4) DESK ADJ. IN HEIGHT (210-345 MM) + SKIRTGUARD HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED	(NC 0740, 0750, 0760, 0770, 0780, 0790, 0952, 0955, 1480, 1483)	£55
CEF0696		REMOVABLE RH (N°4) FULL LENGTH ADJ. IN HEIGHT (210-345 MM) AND DEPTH + SKIRTGUARD HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED	(NC 0740, 0750, 0760, 0770, 0780, 0790, 0952, 0955, 1483)	£55
CEF0732		REMOVABLE RH (N°4) ARMPAD HEMI BLACK MULTI STRETCH PU COATED FABRIC MULTI POSITION AUTOLOCK SAFETY DEVICE INCLUDED	(NC 0952, 1480, 1483)	£54

ARMREST LH*			
CEF0740		REMOVABLE LH (N°1) DESK ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED (NC 1480, 1484)	○
CEF0750		REMOVABLE LH (N°1) FULL LENGTH ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED	○
CEF0760		SWING AWAY DETACHABLE LH (N°2) DESK 240 MM HEIGHT (NC 0952, 0955, 1480, 1484)	◆
CEF0770		SWING AWAY DETACHABLE LH (N°2) FULL LENGTH 240 MM HEIGHT (NC 0952, 0955)	○
CEF0780		SWING AWAY DETACHABLE LH (N°3) DESK ADJ. IN HEIGHT (170-295 MM) (NC 0952, 0955, 1480, 1484)	£51
CEF0790		SWING AWAY DETACHABLE LH (N°3) FULL LENGTH ADJ. IN HEIGHT (170-295 MM) AND DEPTH (NC 0952, 0955, 1484)	£51
CEF0795		REMOVABLE LH (N°4) DESK ADJ. IN HEIGHT (210-345 MM) + SKIRTGUARD HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED (NC 0952, 0955, 1480, 1484)	£55
CEF0796		REMOVABLE LH (N°4) FULL LENGTH ADJ. IN HEIGHT (210-345 MM) AND DEPTH + SKIRTGUARD HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED (NC 0952, 0955, 1484)	£55
CEF0832		REMOVABLE LH N°4 ARMPAD HEMI BLACK MULTI STRETCH PU COATED FABRIC MULTI POSITION ARMPAD HEMI MULTI POSITION (NC 0952, 1480, 1484)	£54
SEAT CUSHION			
CEF0849		BASIC 50 MM FOAM CUSHION - NYLON AND ANTI-SLIP COVER (NC 0952)	£22
SEAT CUSHION DEPTH			
CEF0878		FOR 400 MM SEAT DEPTH	○
CEF0879		FOR 450 MM SEAT DEPTH	○
OTHER SEATING OPTIONS			
CEF0882		LATERAL POSITIONING PADS (PAIR)	£56
CEF0905		NECKREST BLACK MULTI STRETCH PU COATED FABRIC (NC 0906, 0907, 0940, 0945)	£122
CEF0906		HEADREST BLACK MULTI STRETCH PU COATED FABRIC (NC 0907, 0940, 0945)	£146
CEF0907		HEADREST WITH CHEEK SUPPORT BLACK MULTI STRETCH PU COATED FABRIC (NC 0940, 0945)	£169
SWO3333		MULT ADJUSTABLE LATERALS	£263
CEF0930		CALF STRAP	◆
CEF0935		LONGER FOOTREST TUBE	£28
CEF0936		SHORTER FOOTREST SUPPORT	£39
CEF0940		BACKBRACE (NC 0945)	£32
CEF0945		FOLDABLE BACKBRACE	£36
CEF0952		IN-RIGGED ARMREST (NC 1480, 1483, 1484)	£39
CEF0955		W/O AUTOLOCK SAFETY DEVICE FOR ARMREST N°1 ONLY	○

COLOUR*			
CEF0961	<input type="checkbox"/>	SLATE GREY (RAL7024 FINE TEXTURE)	○
CEF0964	<input type="checkbox"/>	PEARL WHITE	○
CEF0965	<input type="checkbox"/>	ELECTRIC GREEN	○
CEF0966	<input type="checkbox"/>	MATT BLACK	○
CEF0967	<input type="checkbox"/>	BLUE POP	○
CEF0968	<input type="checkbox"/>	HYPNOTIC PURPLE	○
CEF0969	<input type="checkbox"/>	RED MATT	○
CEF0980	<input type="checkbox"/>	ANTHRACITE BLACK	○
CEF0991	<input type="checkbox"/>	SPECIAL COLOUR ON DEMAND ONLY	£156
SWO4444	<input type="checkbox"/>	DARK RED - SPECIAL	£156
SWO5555	<input type="checkbox"/>	PEARL GREY - SPECIAL	£156
SWO6666	<input type="checkbox"/>	INTENSE YELLOW	£156
SWO7777	<input type="checkbox"/>	VIBRANT ORANGE	£156
SWO8888	<input type="checkbox"/>	POLAR BLUE	£156
SWO9999	<input type="checkbox"/>	CANDY PINK	£156
CASTOR HOUSING*			
CEF0995	<input type="checkbox"/>	STANDARD	◆
CEF0996	<input type="checkbox"/>	LIGHT WEIGHT	£27
FRONT CASTORS*			
CEF1010	<input type="checkbox"/>	8" X 1.1/4" BLACK SOLID TYRE	◆
CEF1030	<input type="checkbox"/>	8" X 2" BLACK SOLID TYRE	£26
REAR WHEEL SIZE*			
CEF1062	<input type="checkbox"/>	12"	◆
REAR WHEELS*			
CEF1073	<input type="checkbox"/>	DRUM BRAKE - TRANSIT 12" (NC 1201, 1270)	○
CEF1074	<input type="checkbox"/>	TRANSIT COMPOSITE (BLACK PLASTIC HUB AND RIM) (NC 1280)	◆
TYRES*			
CEF1082	<input type="checkbox"/>	SOLID BLACK	◆
REAR WHEEL AXLE*			
CEF1170	<input type="checkbox"/>	FIXED AXLE	◆
REAR WHEEL POSITION*			
CEF1200	<input type="checkbox"/>	PASSIVE REAR WHEELS POSITION	◆
CEF1201	<input type="checkbox"/>	DUAL AMPUTEE REAR WHEEL POSITION (NC 1280)	£67
CHASSIS OPTION			
CEF1241	<input type="checkbox"/>	FRONT WHEELS EXTENSION	£20
BRAKES*			
CEF1270	<input type="checkbox"/>	PUSH TO LOCK	◆
CEF1280	<input type="checkbox"/>	ATTENDANT HUB BRAKES	£213
TRANSPORTATION KIT*			
CEF1345	<input type="checkbox"/>	TRANSPORTATION BRACKETS	◆

SAFETY OPTIONS*

CEF1360		POSTURE BELT	(NC 1361, 1362)	◆
CEF1361		BODYPOINT 2-PT POSTURE BELT M	(NC 1362)	£74
CEF1362		BODYPOINT 2-PT POSTURE BELT L		£74
CEF1370		REAR ANTI-TIPPERS PAIR	(NC 1375, 1440, 1441)	£33
CEF1375		REAR ANTI-TIPPERS PAIR WITH TIPPING AID	(NC 1440, 1441)	£67

OTHER OPTIONS

CEF1430		CANE HOLDER		£21
CEF1440		TIPPING AID RH		£13
CEF1441		TIPPING AID LH		£13
CEF1480		ACRYLIC TRAY	(NC 1483, 1484, 1485)	£107
CEF1483		HALF TABLE TRAY RH	(NC 1484)	£83
CEF1484		HALF TABLE TRAY LH		£83
CEF1485		COMFORT PADDING FOR HALF TABLE TRAY		£41
SWO10001		ABDUCTED LEGRESTS (PR) 0° TO 30° - INTERNAL REFERENCE DE259618		£99

TOTAL PRICE

Non Conforming Options:

CEF0330	NYLON W/O HOOK AND LOOP STRIPS
CEF0350	PADDED NYLON
CEF0465	W/O BACK UPH
CEF0471	STANDARD NYLON SLING (JEMIMA)
CEF0474	LUMBAR SUPPORT BI MATERIAL PADDED COVER
CEF0477	CONTOUR ADJUSTABLE (3+1 STRAPS) BI MATERIAL PADDED COVER
CEF0478	CONTOUR ADJUSTABLE (3+1 STRAPS) BI MATERIAL PADDED COVER INCL BACKBRACE
CEF0481	CONTOUR ADJUSTABLE (7+1 STRAPS) BI MATERIAL PADDED COVER INCL BACKBRACE
CEF0600	FLIP UP RH FOOTPLATE
CEF0610	FLIP UP DEPTH AND ANGLE AJUSTABLE RH FOOTPLATE - NB: IF SELECTED WITH A FIXED FOOTREST - FOOTPLATE WILL BE ANGLE AND DEPTH ADJ AND COMPATIBLE WITH ANKLE HUGGERS. IF SELECTED WITH AN ELEVATING LEGREST - FOOTPLATE WILL BE ANGLE, DEPTH AND WIDTH ADJ - NOT COMPATIBLE WITH ANKLE HUGGERS.
SWO1111	FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE - NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY)
CEF0611	FLIP UP ONE PIECE
CEF66666	NO FOOTREST
CEF0620	FLIP UP LH FOOTPLATE
CEF0630	FLIP UP DEPTH AND ANGLE AJUSTABLE LH FOOTPLATE - NB: IF SELECTED WITH A FIXED FOOTREST - FOOTPLATE WILL BE ANGLE AND DEPTH ADJ AND COMPATIBLE WITH ANKLE HUGGERS. IF SELECTED WITH AN ELEVATING LEGREST - FOOTPLATE WILL BE ANGLE, DEPTH AND WIDTH ADJ - NOT COMPATIBLE WITH ANKLE HUGGERS.
SWO2222	FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE - NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY)
CEF0740	REMOVABLE LH (N°1) DESK ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED
CEF0750	REMOVABLE LH (N°1) FULL LENGTH ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED
CEF0760	SWING AWAY DETACHABLE LH (N°2) DESK 240 MM HEIGHT
CEF0770	SWING AWAY DETACHABLE LH (N°2) FULL LENGTH 240 MM HEIGHT
CEF0780	SWING AWAY DETACHABLE LH (N°3) DESK ADJ. IN HEIGHT (170-295 MM)
CEF0790	SWING AWAY DETACHABLE LH (N°3) FULL LENGTH ADJ. IN HEIGHT (170-295 MM) AND DEPTH
CEF0795	REMOVABLE LH (N°4) DESK ADJ. IN HEIGHT (210-345 MM) + SKIRTGUARD HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED
CEF0796	REMOVABLE LH (N°4) FULL LENGTH ADJ. IN HEIGHT (210-345 MM) AND DEPTH + SKIRTGUARD HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED
CEF0878	FOR 400 MM SEAT DEPTH
CEF0879	FOR 450 MM SEAT DEPTH
CEF0882	LATERAL POSITIONING PADS (PAIR)
CEF0905	NECKREST BLACK MULTI STRETCH PU COATED FABRIC
CEF0906	HEADREST BLACK MULTI STRETCH PU COATED FABRIC
CEF0907	HEADREST WITH CHEEK SUPPORT BLACK MULTI STRETCH PU COATED FABRIC
CEF0930	CALF STRAP
CEF0935	LONGER FOOTREST TUBE
CEF0936	SHORTER FOOTREST SUPPORT
CEF0940	BACKBRACE
CEF0945	FOLDABLE BACKBRACE
CEF0952	IN-RIGGED ARMREST
CEF0955	W/O AUTOLOCK SAFETY DEVICE FOR ARMREST N°1 ONLY
CEF1074	TRANSIT COMPOSITE (BLACK PLASTIC HUB AND RIM)
CEF1201	DUAL AMPUTEE REAR WHEEL POSITION
CEF1270	PUSH TO LOCK

CEF1280	ATTENDANT HUB BRAKES
CEF1361	BODYPOINT 2-PT POSTURE BELT M
CEF1362	BODYPOINT 2-PT POSTURE BELT L
CEF1375	REAR ANTI-TIPPERS PAIR WITH TIPPING AID
CEF1440	TIPPING AID RH
CEF1441	TIPPING AID LH
CEF1480	ACRYLIC TRAY
CEF1483	HALF TABLE TRAY RH
CEF1484	HALF TABLE TRAY LH
CEF1485	COMFORT PADDING FOR HALF TABLE TRAY



For a more complex customer requirement please refer to our Invacare Unique special order form. email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: _____ Email: _____

Phone: _____ Fax: _____

Description of special request (please include relevant option code or part number):