



Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

MARCH 2021
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Kuschall The KSL 2.0

Standard Features

- Max User Weight 100kg (16 stone)
- Active, Fully Welded, Lightweight, Rigid Wheelchair with Height Adjustable Rigid Footrest
- Fully Hydroformed Aluminium Front Frame 75° or 90° Front Angle
- 24" Self Propel Q/R Wheels with Grey Aluminium Handrims & Schwalbe Black & Grey Pneumatic Tyres
- Choice of 6 Seat Widths 340mm to 440mm (13" - 17")
- Choice of 5 Seat Depths 385mm to 485mm (15" - 19")
- Range of Front Seat Heights from 450mm to 530mm (17" - 21")
- Range of Rear Seat Heights from 400mm to 480mm (15" - 19")
- Lockdown Folding, Angle & Height Adjustable Backrest
- Black Nylon Upholstery
- Without Push Handles, Carbon Clothesguards
- Carbon Fibre Axle Tube with 3° Camber
- Black Comfort Front Castor
- Standard Brake Lever Bent
- Tool Kit



Standard Base price: £2,599

Please make your selection with a cross in the box ○ Free of charge ◆ standard feature POA special request

Customer weight = KG

Code	Description	
BASE PRODUCT*		
DSA0085	KÜSCHALL THE KSL 2.0	£2,599
DSA0000	MYKÜSCHALL Special Request - Please include quotation reference in special description at the end of the form where applicable	POA

FRAME*			
DSA1000		Aluminium Hydroformed V Style Frame Angle 75° to the ground	◆
DSA1010		Aluminium Hydroformed Tapered Frame Angle 75° to the ground	○
DSA1020		Aluminium Hydroformed V Style Frame Angle 90° to the ground	(NC 5300, 5350, 5400, 5450, 5500) ○
DSA1030		Aluminium Hydroformed Tapered Frame Angle 90° to the ground	(NC 5300, 5350, 5400, 5450, 5500) ○
SEAT WIDTH*			
DSA1315		340 mm	(NC 1600, 2720, 2730) ○
DSA1320		360 mm	○
DSA1325		380 mm	○
DSA1330		400 mm	◆
DSA1335		420 mm	○
DSA1340		440 mm	○
SEAT DEPTH*			
DSA1430		385 mm	(NC 2720, 2730) ○
DSA1440		410 mm	◆
DSA1455		435 mm	○
DSA1470		460 mm	○
DSA1485		485 mm	○
BACK HEIGHT*			
DSA1600		Back Height for 3rd party backrest. Final height depends on the height of the positioning backrest, select if Matrix/Fixed backrest is required.	(NC 2200, 2220, 2230) ○
DSA1605		270 mm	(NC 2240) ○
DSA1610		285 mm	(NC 2240) ○
DSA1615		300 mm	(NC 2240) ○
DSA1620		315 mm	(NC 2240) ○
DSA1625		330 mm	(NC 2240) ○
DSA1630		345 mm	(NC 2240) ○
DSA1635		360 mm	(NC 2240) ○
DSA1640		375 mm	(NC 2240) ○
DSA1645		390 mm	(NC 2240) ○
DSA1650		405 mm	(NC 2240) ○
DSA1655		420 mm	(NC 2240) ○
FRONT SEAT HEIGHT. Lower leg must be front height -100mm To -30mm on 75° frames and -80mm To -40mm on 90° frames*			
DSA1740		450 mm	(NC 1840, 1845, 1850, 1855, 1860, 6420) ○
DSA1745		460 mm	(NC 1845, 1850, 1855, 1860) ○
DSA1750		470 mm	(NC 1850, 1855, 1860) ○
DSA1755		480 mm	(NC 1855, 1860) ○
DSA1760		490 mm	(NC 1860) ○
DSA1765		500 mm	◆
DSA1770		510 mm	(NC 1820) ○
DSA1775		520 mm	(NC 1820, 1825) ○
DSA1780		530 mm	(NC 1820, 1825, 1830, 6410) ○

REAR SEAT HEIGHT. Rear height = "Front height -100 MM" TO "Front height -20 MM"*			
DSA1820		400 mm	(NC 6420) ○
DSA1825		410 mm	○
DSA1830		420 mm	○
DSA1835		430 mm	◆
DSA1840		440 mm	○
DSA1845		450 mm	○
DSA1850		460 mm	○
DSA1855		470 mm	○
DSA1860		480 mm	○
LOWER LEG LENGTH. Cushion height must be deducted*			
DSA1930		350 mm	○
DSA1932		360 mm	○
DSA1934		370 mm	○
DSA1936		380 mm	○
DSA1938		390 mm	○
DSA1940		400 mm	○
DSA1942		410 mm	○
DSA1944		420 mm	○
DSA1946		430 mm	◆
DSA1948		440 mm	○
DSA1950		450 mm	○
DSA1952		460 mm	○
DSA1954		470 mm	○
DSA1956		480 mm	○
DSA1958		490 mm	○
DSA1960		500 mm	○
BACKREST TYPE*			
DSA2110		Fixed Aluminium Hydroformed backrest in frame colour	◆
DSA2140		Folding Aluminium angle adjustable backrest in frame colour	○
BACKREST UPHOLSTERY*			
DSA2200		Standard Light Black Nylon upholstery (No tension adjustment)	◆
DSA2220		Tension Adjustable Backrest, with Premium Black Nylon and Reflective Trim Upholstery Cover, Back Pocket/hook and loop fastener	£118
DSA2230		Tension Adjustable Backrest, with Black Aero, Breathable Spacer Upholstery Cover, Back Pocket/hook and loop fastener	£193
DSA2240		Without Backrest Upholstery Select if fitting a Matrix/Fixed backrest, no upholstery required	○
PUSH HANDLES*			
DSA2400		Without Push Handles	◆
DSA2410		Push Handles Fixed Short/Mini Hand Grip	£115
DSA2420		Push Handles Fixed Long Hand Grip	£115
DSA2440		Push Handles Height Adjustable, Integrated in backrest	£189
DSA2450		Push Handles Foldable	£171

BACKREST ANGLE*			
DSA2500		74° angle from the seat	○
DSA2520		78° angle from the seat	○
DSA2540		82° angle from the seat	◆
DSA2560		86° angle from the seat	○
DSA2580		90° angle from the seat	○
SEAT UPHOLSTERY*			
DSA2600		Standard Nylon Seat Upholstery with one 5cm central band (Loop)	◆
DSA2630		Light Seat Upholstery, mesh material to reduce weight	£60
SEAT CUSHION			
DSA2700		Standard Seat Cushion, 60 mm Double Foam with Visco Top Layer	£57
DSA2720		Matrx Libra Cushion, Black Smooth Cover - Very High Risk	£338
DSA2730		Matrx Vi Cushion, Black Smooth Cover - High Risk	£312
FOOTREST*			
DSA3200		Footrest Tube, Rigid - Aluminium Hydroformed Profile	◆
DSA3210		Footrest Tube, Rigid - Titanium (Rounded tube)	£160
DSA3220		Carbon Footplate, Angle Adjustable (NC 3290)	£195
FOOTREST COVER			
DSA3290		Footrest Cover	£32
CALF STRAP*			
DSA3300		Calf strap Mounted on the frame with hook and loop	◆
DSA3420		Calf strap Mounted on the supporter with bracket	○
CLOTHESGUARDS/MUDGUARDS*			
DSA3410		Carbon Clothesguards, fixed between the back rest and seat rails	◆
DSA3430		Carbon Mudguards, fixed between the backrest and seat rails (NC 7230)	£220
FRAME COLOUR*			
DSA4000		Space Black	○
DSA4010		Anthracite Grey	○
DSA4020		Pure White (NC 4220)	○
DSA4030		Carmine Red	○
DSA4040		Azure Blue	○
DSA4060		Flashy Yellow (Additional 10 day lead time)	£26
DSA4070		Atomic Violet (Additional 10 day lead time)	£26
DSA4100		Black Glossy (Additional 10 day lead time)	£26
DSA4130		Silver Grey (Additional 10 day lead time)	£26
DSA4140		Royal Blue (Additional 10 day lead time)	£26
DSA4150		Ruby Red (Additional 10 day lead time)	£26
DSA4170		Dynamite Orange (Additional 10 day lead time)	£26
DSA4180		Forest Green (Additional 10 day lead time)	£26
DSA4050		Personal frame colour. Provide standard RAL colour code of your choice	£268
FRAME DECAL TYPE*			
DSA4190		Bold Decal (Full coloured transfer on frame) (NC 4300, 4310)	£26
DSA4195		Tag (Small Kuschall brand transfer on front of frame) (NC 4200, 4210, 4220, 4230, 4240, 4250)	◆

BOLD DECAL COLOUR*			
DSA4200	<input type="checkbox"/>	Black	<input type="radio"/>
DSA4210	<input type="checkbox"/>	Dark Grey	<input type="radio"/>
DSA4220	<input type="checkbox"/>	Light Grey	<input type="radio"/>
DSA4230	<input type="checkbox"/>	Red	<input type="radio"/>
DSA4240	<input type="checkbox"/>	Blue	<input type="radio"/>
DSA4250	<input type="checkbox"/>	Gold	<input type="radio"/>
DSA1111	<input type="checkbox"/>	No bold colour, select only with Option DSA4195 Tag (NC 4190)	<input type="radio"/>
TAG COLOUR*			
DSA4300	<input type="checkbox"/>	Black	<input type="radio"/>
DSA4310	<input type="checkbox"/>	Red	<input type="radio"/>
DSA2222	<input type="checkbox"/>	No tag colour, select only with Option DSA4190 Bold Decal (NC 4195)	<input type="radio"/>
FRONT CASTORS*			
DSA5000	<input type="checkbox"/>	75mm x 25mm (3" x 1") Supersport	<input type="radio"/>
DSA5050	<input type="checkbox"/>	100mm x 30mm (4" x 1.2") Everyday Black hub and tyre	<input type="radio"/>
DSA5100	<input type="checkbox"/>	100mm x 30mm (4" x 1.2") Comfort Black hub and tyre	<input type="radio"/>
DSA5150	<input type="checkbox"/>	100mm x 25mm (4" x 1") Starec Black hub and tyre	£72
DSA5200	<input type="checkbox"/>	100mm x 36mm (4" x 1.4") Aluminium Softroll Black hub and tyre	£80
DSA5250	<input type="checkbox"/>	100mm x 36mm (4" x 1.4") Aluminium Softroll Silver hub, Black tyre	£80
DSA5300	<input type="checkbox"/>	125mm x 30mm (5" x 1.2") Everyday Black hub and tyre	<input type="radio"/>
DSA5350	<input type="checkbox"/>	125mm x 30mm (5" x 1.2") Comfort Black hub and tyre	<input type="radio"/>
DSA5400	<input type="checkbox"/>	125mm x 25mm (5" x 1") Starec Black hub and tyre	£72
DSA5450	<input type="checkbox"/>	125mm x 36mm (5" x 1.4") Aluminium Softroll Black hub and tyre	£80
DSA5500	<input type="checkbox"/>	125mm x 36mm (5" x 1.4") Aluminium Softroll Silver hub, Black tyre	£80
SEAT POSITION WELDED*			
DSA6100	<input type="checkbox"/>	Position 1 = 140 mm MOST ACTIVE POSITION	<input type="radio"/>
DSA6110	<input type="checkbox"/>	Position 2 = 125 mm	<input type="radio"/>
DSA6120	<input type="checkbox"/>	Position 3 = 110 mm	◆
DSA6130	<input type="checkbox"/>	Position 4 = 95 mm	<input type="radio"/>
DSA6140	<input type="checkbox"/>	Position 5 = 80 mm LEAST ACTIVE POSITION	<input type="radio"/>
CAMBER OF CARBON AXLE*			
DSA6200	<input type="checkbox"/>	0° Camber	<input type="radio"/>
DSA6210	<input type="checkbox"/>	3° Camber	◆
DSA6220	<input type="checkbox"/>	6° Camber	<input type="radio"/>
QUAD QUICK RELEASE AXLE			
DSA6300	<input type="checkbox"/>	Quad Quick Release/Tetra Clip	£74
REAR WHEEL SIZE*			
DSA6410	<input type="checkbox"/>	24"	◆
DSA6420	<input type="checkbox"/>	25" (NC 6510, 6610, 6700, 6710, 6720, 6760, 6770, 6900, 6910)	<input type="radio"/>

REAR WHEEL TYPE*		
DSA6510		24" KÜSCHALL Starec Wheel, Black Cross Spokes, Rim & Hub. ◆
DSA6520		24", 25" High Performance Wheel Black hub and rim £156
DSA6550		24", 25" Spinergy LXX, Black hub, rim and spokes in Black £527
SWO0907		24", 25" Spinergy SPOX, Chrome hub and 18 spokes in Black. Coloured spokes are available via MyKÜSCHALL £420
SWO0909		24", 25" Spinergy LX, Chrome hub and 12 spokes in Black. Coloured spokes are available via MyKÜSCHALL £635
REAR WHEEL TYRE*		
DSA6600		24", 25" Schwalbe Rightrun Black/Grey Pneumatic tyre ◆
DSA6610		24", 25" Schwalbe Rightrun Blue with K Guard Pneumatic tyre £41
DSA6620		24", 25" Schwalbe Marathon Plus Smart Guard (Puncture Protected tyre with 5mm layer of flexible rubber) £85
DSA6630		24", 25" Schwalbe One V-Guard Puncture Resistant tyre (Cut-Resistant High-Tech Fibre) £83
HANDRIM*		
DSA6700		24" Supergripp, mounted at 3 cm distance £193
DSA6710		24" Steel, mounted at 3 cm distance £54
DSA6720		24" Titanium, long tab with 2 holes for 3 cm or 4 cm mountings, mounted at 3cm distance £166
DSA6740		24" & 25" Carbolife Gekko mounted at 3 cm distance. Unique shape for low-effort propulsion, with fully integrated rubber strip for extra grip £276
DSA6750		24" & 25" Carbolife Gekko S, mounted at 3 cm distance. Unique shape for low-effort propulsion, with fully integrated rubber strip for extra grip, smaller profile for smaller hands £276
DSA6760		24" Carbolife Curve L, mounted at 3 cm distance. Unique shape for low-effort propulsion £245
DSA6770		24" Carbolife Quadro, mounted at 3 cm distance. Unique shape for low-effort propulsion, with full anti-slip coating over the entire handrim £308
DSA6780		24", 25" Light Aluminium Grey, long tab with 2 holes for 3 cm or 4cm mountings, mounted at 3cm distance ◆
DSA6790		24", 25" Light Aluminium Black, mounted at 3 cm distance £35
SPOKEGUARD		
DSA6900		24" Transparent Spokeguard £108
DSA6910		24" Black Spokeguard £108
BRAKES*		
DSA7210		Standard Push to Lock Brakes, Lever Bent ◆
DSA7230		Performance Brakes, Push to Lock. Brake barrel flips away from tyre £66
DSA7240		Active/Scissor Brakes £66
DSA7250		Light/Compact Active, Aluminium Brakes £179
SAFETY OPTIONS		
DSA7500		Anti-Tipper Right £105
DSA7510		Anti-Tipper Left £105
DSA7520		Anti-Tipper Both sides £211
DSA7800		Standard Pelvic belt £63

OTHER OPTIONS			
DSA8010		Front Frame Protectors; 2 Piece. Removable with hook and loop	£32
DSA8500		Air Pump	£11
DSA8550		Tool Kit	◆
DSA8680		24" Mountain Bike Wheels with bag. Fitted with Specific Black Tyres, Grey Aluminium Handrims; Quick Release Axles Included	£276

TOTAL PRICE	
--------------------	--

Non Conforming Options:

DSA1600	Back Height for 3rd party backrest. Final height depends on the height of the positioning backrest, select if Matrix/Fixed backrest is required.
DSA1820	400 mm
DSA1825	410 mm
DSA1830	420 mm
DSA1840	440 mm
DSA1845	450 mm
DSA1850	460 mm
DSA1855	470 mm
DSA1860	480 mm
DSA2200	Standard Light Black Nylon upholstery (No tension adjustment)
DSA2220	Tension Adjustable Backrest, with Premium Black Nylon and Reflective Trim Upholstery Cover, Back Pocket/hook and loop fastener
DSA2230	Tension Adjustable Backrest, with Black Aero, Breathable Spacer Upholstery Cover, Back Pocket/hook and loop fastener
DSA2240	Without Backrest Upholstery Select if fitting a Matrix/Fixed backrest, no upholstery required
DSA2720	Matrix Libra Cushion, Black Smooth Cover - Very High Risk
DSA2730	Matrix Vi Cushion, Black Smooth Cover - High Risk
DSA3290	Footrest Cover
DSA4190	Bold Decal (Full coloured transfer on frame)
DSA4195	Tag (Small Kuschall brand transfer on front of frame)
DSA4200	Black
DSA4210	Dark Grey
DSA4220	Light Grey
DSA4230	Red
DSA4240	Blue
DSA4250	Gold
DSA4300	Black
DSA4310	Red
DSA5300	125mm x 30mm (5" x 1.2") Everyday Black hub and tyre
DSA5350	125mm x 30mm (5" x 1.2") Comfort Black hub and tyre
DSA5400	125mm x 25mm (5" x 1") Starec Black hub and tyre
DSA5450	125mm x 36mm (5" x 1.4") Aluminium Softroll Black hub and tyre
DSA5500	125mm x 36mm (5" x 1.4") Aluminium Softroll Silver hub, Black tyre
DSA6410	24"
DSA6420	25"
DSA6510	24" Kuschall Starec Wheel, Black Cross Spokes, Rim & Hub.
DSA6610	24", 25" Schwalbe Rightrun Blue with K Guard Pneumatic tyre
DSA6700	24" Supergripp, mounted at 3 cm distance
DSA6710	24" Steel, mounted at 3 cm distance
DSA6720	24" Titanium, long tab with 2 holes for 3 cm or 4 cm mountings, mounted at 3cm distance
DSA6760	24" Carbolife Curve L, mounted at 3 cm distance. Unique shape for low-effort propulsion
DSA6770	24" Carbolife Quadro, mounted at 3 cm distance. Unique shape for low-effort propulsion, with full anti-slip coating over the entire handrim
DSA6900	24" Transparent Spokeguard
DSA6910	24" Black Spokeguard



For a more complex customer requirement please refer to our Invacare Unique special order form. email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Fax:

Description of special request (please include relevant option code or part number):