



Manual Retail Prescription Form

Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

MARCH 2021
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Azalea Minor

Standard Features

- Max User Weight 75KG/12 Stone
- Crash Tested ISO 7176-Part 19
- Transit Propulsion with Attendant Hub Brakes and Black Solid Tyres
- Seat Width 340mm (+50mm or +100mm Width Adjustment)
- Choice of 2 Seat Heights 400mm or 450mm
- Depth Adjustable Seat 380mm-450mm
- Backrest Plate with Braced Push Handles
- Assistant Operated Gas Piston Seat Tilt +1° - +25° and Backrest Recline 30° with Tilt Lock
- Angle Adjustable Legrests with Calf Support
- Height Adjustable Wide Pad Armrests
- Pearl Silver Frame Colour
- User Manoeuvred Tyre Brakes
- Pelvic Belt
- Anti Tipplers



Standard Base price: £2,099

Please make your selection with a cross in the box ○ Free of charge ◆ standard feature POA special request

Customer weight = KG

Code	Description		
BASE PRODUCT*			
BNC0014	Azalea Minor		£2,099
BNC8000	SPECIAL REQUEST - PLEASE INCLUDE QUOTATION REFERENCE IN SPECIAL DESCRIPTION AT THE END OF THE FORM WHERE APPLICABLE		POA
PROPULSION*			
BNC0020	Transit 16" wheels with hub brakes	(NC 1055, 1062, 1074, 1100, 1120)	◆
BNC0025	Self Propel 24"/22" wheels with user brakes	(NC 0075, 0605, 0635, 1047, 1076, 1305)	○
BNC0030	Self Propel 24"/22" wheels with assistant hub brakes	(NC 0075, 0605, 0635, 1047, 1074, 1305)	£129

SEAT DIMENSIONS*				
BNC0037	<input type="checkbox"/>	Minor (- 50 mm In seat depth and - 30 mm In chassis length)		◆
SEAT WIDTH*				
BNC0038	<input type="checkbox"/>	340 mm/13" (Minor)		◆
SEAT WIDTH EXTENSION*				
BNC0065	<input type="checkbox"/>	+ 50 mm seat width pre-adjustment on armrests and legrests		○
BNC0075	<input type="checkbox"/>	+ 100 mm seat width pre-adjustment on armrests and legrests	(NC 0110)	○
BNC0085	<input type="checkbox"/>	Without seat width extension, set at 0mm adjustment but can be adjusted wider if required		○
SEAT DEPTH*				
BNC0095	<input type="checkbox"/>	380-450 mm / 15" - 18" (Azalea Minor)		◆
SEAT HEIGHT*				
BNC0110	<input type="checkbox"/>	400 mm / 16"	(NC 1062)	○
BNC0120	<input type="checkbox"/>	450 mm / 18"	(NC 1055)	○
SEAT ANGLE*				
BNC0000	<input type="checkbox"/>	Standard Seat Angle +1° to +25°		◆
BNC0127	<input type="checkbox"/>	Negative Seat Angle -1 To +23°		○
SEAT TYPE*				
BNC0130	<input type="checkbox"/>	Depth Adjustable		◆
SEAT CUSHION*				
BNC0165	<input type="checkbox"/>	Cushion Flosshape seat upholstery	(NC 1111)	£52
BNC0180	<input type="checkbox"/>	Without Seat Cushion, Seat Plate only	(NC 0201, 0202)	◆
SEAT UPHOLSTERY*				
BNC0201	<input type="checkbox"/>	Black Multi-stretch Black PU coated fabric		◆
BNC0202	<input type="checkbox"/>	Black Spacer TR36		○
BNC1111	<input type="checkbox"/>	No upholstery, select with options 0170 Multifunctional and 0180 Without seat cushion		○
BACKREST*				
BNC0217	<input type="checkbox"/>	Flex 2 Backrest Tension adjustable (Azalea Minor)	(NC 0925)	£127
BNC0218	<input type="checkbox"/>	Flex 2 Backrest Fixed wooden plate (Azalea Minor)	(NC 0240, 0905)	◆
BACKREST CUSHION*				
BNC0240	<input type="checkbox"/>	Cover for tension adjustable upholstery		○
BNC0250	<input type="checkbox"/>	Cushion, Laguna (Lateral Support)	(NC 0905)	○
BNC0260	<input type="checkbox"/>	Cushion, Mistral 2 (Waist Support)	(NC 0322, 0905)	£133
BACKREST UPHOLSTERY*				
BNC0311	<input type="checkbox"/>	Black Plush TR35		○
BNC0321	<input type="checkbox"/>	Black Multi-stretch Black PU coated fabric		◆
BNC0322	<input type="checkbox"/>	Black Spacer TR36		○
PUSH HANDLES*				
BNC0340	<input type="checkbox"/>	Push Handles braced (Azalea Minor)		◆
BNC0345	<input type="checkbox"/>	Push Bar, angle adjustable		£94
SEAT TILT AND BACKREST REGULATION*				
BNC0350	<input type="checkbox"/>	Assistant operated gas piston seat tilt and backrest recline		◆
TILT LOCK*				
BNC1263	<input type="checkbox"/>	Tilt Lock (locking of pre-set angle)		◆

LEGREST HANGER ATTACHMENT*			
BNC0375		Standard position for legrest mounting	◆
SWO0189		Standard position with 0-30° Legrest Abduction function	£92
BNC0378		90° Hanger attachment to achieve a 90° legrest position on the Elevating Legrest only (NC 0390, 0410, 0440, 0460)	○
LEGREST HANGER LH*			
BNC0380		Elevating Legrest, angle adjustable 0 to 75° swing in/swing out with knee pad and calf support LH (NC 0475, 0480, 0490, 0942, 0945)	◆
BNC0390		Legrest Hanger fixed 80 deg, swing in/swing out LH	○
BNC0410		Legrest Amputee LH (NC 0475, 0476, 0477, 0480, 0490, 0510, 0520, 0930, 0941, 0942, 0945)	£126
LEGREST HANGER RH*			
BNC0430		Elevating Legrest, angle adjustable 0 to 75° swing In/swing out with knee pad and calf support RH (NC 0475, 0540, 0550, 0942, 0945)	◆
BNC0440		Legrest Hanger fixed 80 deg, swing in/swing out RH	○
BNC0460		Legrest Amputee RH (NC 0475, 0476, 0477, 0540, 0550, 0570, 0580, 0940, 0941, 0942, 0945)	£126
FOOTREST LH*			
BNC0475		Footrest one piece, depth and angle adjustable short tube (NC 0477, 0480, 0490, 0510, 0520, 0540, 0550, 0570, 0580, 0930, 0940)	£109
BNC0476		Footrest one piece, depth and angle adjustable long tube (NC 0477, 0480, 0490, 0510, 0520, 0540, 0550, 0570, 0580, 0930, 0940)	£109
BNC0477		Footbox, angle adjustable, padded with short tube (NC 0480, 0490, 0510, 0520, 0540, 0550, 0570, 0580, 0930, 0940, 0941, 0942, 0945)	£501
BNC0480		Footrest lower part, fixed with short tube and heel strap, LH (NC 0930, 0941)	○
BNC0490		Footrest lower part, depth and angle adjustable with short tube and heel strap, LH	○
BNC0510		Footrest lower part, fixed with long tube and heel strap, LH. (NC 0930, 0941)	○
BNC0520		Footrest lower part, depth and angle adjustable with long tube and heel strap, LH	○
BNC3333		Without footrest, only select with option 0410 Amputee legrest	○
FOOTREST RH*			
BNC0540		Footrest lower part, fixed with short tube and heel strap, RH (NC 0940, 0941)	○
BNC0550		Footrest lower part, depth and angle adjustable with short tube and heel strap, RH	○
BNC0570		Footrest lower part, fixed with long tube and heel strap, RH (NC 0940, 0941)	○
BNC0580		Footrest lower part, depth and angle adjustable with long tube and heel strap, RH	○
BNC4444		Without footrest, only select with options 0460 Amputee legrest, 0475 and 0476 One piece footrests or 0477 Footbox	○
ARMREST LH*			
BNC0600		Armrest Ducktail Pad LH, height adjustable with autolock (NC 1232)	○
BNC0605		Armrest Wide Pad LH, height adjustable with autolock. Compatible with 16" wheels only (NC 1232)	◆
BNC0607		Armrests Reinforced Wide Pad, LH, height adjustable with autolock (NC 1230, 1231, 1232, 1260)	£82
BNC0610		Armrest Padded LH, height adjustable with autolock (NC 1230, 1231, 1232, 1260)	£59
BNC0615		Armrest Hemiplegic LH, height adjustable with autolock (NC 0695, 1230, 1231, 1232, 1260)	£165
BNC0620		Without Armrest LH Select when choosing BNC1232 Half Tray RH, Narrow armrest is included (NC 0695, 1230, 1231, 1260)	○

ARMREST RH*			
BNC0630		Armrest Ducktail Pad RH, height adjustable with autolock	(NC 1233) ○
BNC0635		Armrest Wide Pad RH, height adjustable with autolock. Compatible with 16" wheels only	(NC 1233) ◆
BNC0637		Armrest Reinforced Wide Pad, RH, height adjustable with autolock	(NC 1230, 1231, 1233, 1260) £82
BNC0640		Armrest Padded RH, height adjustable with autolock.	(NC 1230, 1231, 1233, 1260) £59
BNC0645		Armrest Hemiplegic, RH, height adjustable with autolock	(NC 0695, 1230, 1231, 1233, 1240, 1260) £165
BNC0650		Without Armrest RH Select when choosing BNC1233 Half Tray LH, Narrow armrest is included	(NC 0695, 1230, 1231, 1240, 1260) ○
ARMREST OPTIONS. Side support cover 0680, must be selected if 0690 Foam for side support is required			
BNC0680		Side Support Cover, removable with tools (wipeable cover for the sideguard)	(NC 0695) £47
BNC0690		Foam for side support cover (Internal cushioned pad, which is only possible in combination with 0680)	(NC 0695) ○
BNC0695		Pad for sideguard, removable with hook and loop fastening	£47
SEAT OPTIONS			
BNC0710		Incontinence Cover	£37
BNC0720		Pommel Abduction Cushion	£109
BNC0725		Tilt Scale for seat tilt	£41
BACKREST OPTIONS			
BNC0730		Headrest, Sideway Adjustable	(NC 0750, 0770) ◆
BNC0750		Headrest with cheek support, sideway adjustable	(NC 0740, 0770, 0780) £86
BNC0740		Cover for headrest	£28
BNC0770		Neckrest	○
BNC0780		Cover For Neckrest	£20
BNC0790		Handles Clinic Kit for Head-/Neckrest (Levers which replace screws for easy adjustment)	£45
BNC0900		Trunk Support Attachment to mount laterals later. If trunk supports are chosen, the attachment bracket is already included	○
BNC0905		Lateral Positioning Pads. Attaches with hook and loop between the strap and the cover of the tension adjustable upholstery	£56
BNC0910		Lumbar Support Cushion	£28
BNC0915		Tilt Scale for backrest recline	£41
BNC0920		Wedges (One Pair)	£9
BNC0925		Privacy Attachment (Material connecting backrest and seat)	£31
TRUNK SUPPORT LH			
BNC0820		Trunk Support Multi Adjustable with extended height arm LH	(NC 0900) ○
BNC0880		Trunk Support Swingaway LH	(NC 0900) £97
BNC0895		Trunk Support Multi Adjustable Swingaway LH	(NC 0900) £97
TRUNK SUPPORT RH			
BNC0830		Trunk Support Multi Adjustable with extended height arm RH	(NC 0900) ○
BNC0890		Trunk Support Swingaway RH	(NC 0900) £97
BNC0896		Trunk Support Multi Adjustable Swingaway RH	(NC 0900) £97

LEGREST OPTIONS			
BNC0930	<input type="checkbox"/>	Footplate Extender LH (NC 0941)	£31
BNC0940	<input type="checkbox"/>	Footplate Extender RH (NC 0941)	£31
BNC0941	<input type="checkbox"/>	Cover Footplate, Black Plush Fabric	£49
BNC0942	<input type="checkbox"/>	Calf Strap Padded	£67
BNC0945	<input type="checkbox"/>	Calf Strap (for One Piece Footrest)	○
CHASSIS COLOUR*			
BNC0944	<input type="checkbox"/>	Silver Pearl	◆
BNC0946	<input type="checkbox"/>	Azurite (Dark Blue)	○
BNC0947	<input type="checkbox"/>	Sand (Beige)	○
BNC0948	<input type="checkbox"/>	Electric Green	○
BNC0949	<input type="checkbox"/>	Happy Red	○
CASTOR ATTACHMENT*			
BNC0960	<input type="checkbox"/>	Fixed castor attachment	◆
FRONT CASTORS*			
BNC0982	<input type="checkbox"/>	140x45 mm (5 1/2" x 1 3/4") Medium Hard Black tyre	£18
BNC0992	<input type="checkbox"/>	150x30 mm (6" x 1 1/4") Medium Hard Black tyre	◆
REAR WHEEL SIZE*			
BNC1047	<input type="checkbox"/>	16" / 405 mm Transit rear wheel with fixed attachment (NC 1074, 1082, 1084, 1100, 1120, 1160)	◆
BNC1055	<input type="checkbox"/>	22" / 560 mm Self Propel rear wheel with QR attachment	○
BNC1062	<input type="checkbox"/>	24" / 600 mm Self Propel rear wheel with QR attachment	○
REAR WHEEL TYPE*			
BNC1074	<input type="checkbox"/>	Standard self propel with colourless hub, clear anodized rim and Stainless steel spokes	○
BNC1076	<input type="checkbox"/>	Hub Brake 16" Transit or 22" or 24" Self Propel Transit with Black Plastic Hub and Rim Self Propel with colourless hub, clear anodized rim and Stainless steel spokes	○
REAR WHEEL TYRE*			
BNC1082	<input type="checkbox"/>	Pneumatic Profile Black 1 3/8" tyres	○
BNC1084	<input type="checkbox"/>	Pneumatic Puncture proof / flatfree 1 3/8" Black tyres	£38
BNC1092	<input type="checkbox"/>	Solid Black tyres	◆
HANDRIM			
BNC1100	<input type="checkbox"/>	Aluminium handrim	◆
BNC1120	<input type="checkbox"/>	Plastic coated handrim	£45
SPOKEGUARD			
BNC1160	<input type="checkbox"/>	Frost spokeguard	£153
BRAKE*			
BNC1170	<input type="checkbox"/>	User manoeuvred tyre brakes	◆
PELVIC BELT*			
BNC1200	<input type="checkbox"/>	Pelvic belt	◆
SAFETY OPTIONS			
BNC1220	<input type="checkbox"/>	Reflector Kit (wheels and chassis)	£20
TRANSPORTATION BRACKETS*			
BNC1225	<input type="checkbox"/>	Transportation Brackets, mounted to the front of the wheelchair	◆

OTHER OPTIONS

BNC1230		Fixed Table Tray without elbow pads includes mounting kit (NC 1237, 1260)	£127
BNC1231		Fixed Table Tray with elbow Pads includes mounting kit (NC 1237, 1260)	£127
BNC1232		Half table tray, clear plastic, flip up, depth adjustable, left side (armrest with straight narrow pad included) (NC 1235, 1260)	£95
BNC1233		Half table tray, clear, plastic, flip up, depth adjustable, right side (armrest with straight narrow pad included) (NC 1235, 1260)	£95
BNC1240		Swing away table without elbow pads (NC 1235, 1260)	£140
BNC1235		Table tray padded cushion	£47
BNC1237		Half table tray padded cushion	£47
BNC1260		Mounting kit for table (if ordered separately)	£7
BNC1270		Pump	£20
BNC1280		Cane Holder	£46
BNC1290		Handles Clinic Kit for Head-/Neckrest (Levers which replace screws for easy adjustment)	£50
BNC1300		Harness attachment , must be selected in combination with Head/Neckrest	£79
BNC1305		Front castor locks, includes front tiedown brackets. Allows castors to remain parallel for easy access into confined spaces such as lifts or vehicles	£169
BNC1310		Attachment for positioning belt	£29

TOTAL PRICE	
--------------------	--

Non Conforming Options:

BNC0075	+ 100 mm seat width pre-adjustment on armrests and legrests
BNC0110	400 mm / 16"
BNC0201	Black Multi-stretch Black PU coated fabric
BNC0202	Black Spacer TR36
BNC1111	No upholstery, select with options 0170 Multifunctional and 0180 Without seat cushion
BNC0240	Cover for tension adjustable upholstery
BNC0322	Black Spacer TR36
BNC0390	Legrest Hanger fixed 80 deg, swing in/swing out LH
BNC0410	Legrest Amputee LH
BNC0440	Legrest Hanger fixed 80 deg, swing in/swing out RH
BNC0460	Legrest Amputee RH
BNC0475	Footrest one piece, depth and angle adjustable short tube
BNC0476	Footrest one piece, depth and angle adjustable long tube
BNC0477	Footbox, angle adjustable, padded with short tube
BNC0480	Footrest lower part, fixed with short tube and heel strap, LH
BNC0490	Footrest lower part, depth and angle adjustable with short tube and heel strap, LH
BNC0510	Footrest lower part, fixed with long tube and heel strap, LH.
BNC0520	Footrest lower part, depth and angle adjustable with long tube and heel strap, LH
BNC0540	Footrest lower part, fixed with short tube and heel strap, RH
BNC0550	Footrest lower part, depth and angle adjustable with short tube and heel strap, RH
BNC0570	Footrest lower part, fixed with long tube and heel strap, RH
BNC0580	Footrest lower part, depth and angle adjustable with long tube and heel strap, RH
BNC0605	Armrest Wide Pad LH, height adjustable with autolock. Compatible with 16" wheels only
BNC0635	Armrest Wide Pad RH, height adjustable with autolock. Compatible with 16" wheels only
BNC0695	Pad for sideguard, removable with hook and loop fastening
BNC0750	Headrest with cheek support, sideways adjustable
BNC0740	Cover for headrest
BNC0770	Neckrest
BNC0780	Cover For Neckrest
BNC0900	Trunk Support Attachment to mount laterals later. If trunk supports are chosen, the attachment bracket is already included
BNC0905	Lateral Positioning Pads. Attaches with hook and loop between the strap and the cover of the tension adjustable upholstery
BNC0925	Privacy Attachment (Material connecting backrest and seat)
BNC0930	Footplate Extender LH
BNC0940	Footplate Extender RH
BNC0941	Cover Footplate, Black Plush Fabric
BNC0942	Calf Strap Padded
BNC0945	Calf Strap (for One Piece Footrest)
BNC1047	16" / 405 mm Transit rear wheel with fixed attachment
BNC1055	22" / 560 mm Self Propel rear wheel with QR attachment
BNC1062	24" / 600 mm Self Propel rear wheel with QR attachment
BNC1074	Standard self propel with colourless hub, clear anodized rim and Stainless steel spokes
BNC1076	Hub Brake 16" Transit or 22" or 24" Self Propel Transit with Black Plastic Hub and Rim Self Propel with colourless hub, clear anodized rim and Stainless steel spokes

BNC1082	Pneumatic Profile Black 1 3/8" tyres
BNC1084	Pneumatic Puncture proof / flatfree 1 3/8" Black tyres
BNC1100	Aluminium handrim
BNC1120	Plastic coated handrim
BNC1160	Frost spokeguard
BNC1230	Fixed Table Tray without elbow pads includes mounting kit
BNC1231	Fixed Table Tray with elbow Pads includes mounting kit
BNC1232	Half table tray, clear plastic, flip up, depth adjustable, left side (armrest with straight narrow pad included)
BNC1233	Half table tray, clear, plastic, flip up, depth adjustable, right side (armrest with straight narrow pad included)
BNC1240	Swing away table without elbow pads
BNC1235	Table tray padded cushion
BNC1237	Half table tray padded cushion
BNC1260	Mounting kit for table (if ordered separately)
BNC1305	Front castor locks, includes front tiedown brackets. Allows castors to remain parallel for easy access into confined spaces such as lifts or vehicles



For a more complex customer requirement please refer to our Invacare Unique special order form. email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: _____ Email: _____

Phone: _____ Fax: _____

Description of special request (please include relevant option code or part number):