



Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

MARCH 2021
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Azalea Self Propel

Standard Features

- Max User Weight 135KG/21 Stone
- Crash Tested ISO 7176-Part 19
- 24" or 22" Self Propulsion with User Manoeuvred Tyre Brakes and Solid Black Tyres
- Choice of Seat Widths 390mm, 440mm or 490mm (+50mm Width Adjustment)
- Choice of 2 Seat Heights 400mm or 450mm
- Depth Adjustable Seat 430mm-500mm
- Without Seat Cushion, Supplied with Seat Board
- Flex 3 Backrest Plate, Width and Height Adjustable
- Passad 2 Backrest Cushion with Multi-stretch Black PU coated fabric Cover and Push Handles Connected
- Assistant Operated Gas Piston Seat Tilt +1° - +25°
- Backrest Recline 30° with Tilt Lock
- Angle Adjustable Legrests with Calf Support, Height and Depth Adjustable Footplates
- Height Adjustable Ducktail Pad Armrests with Autolock
- Pearl Silver Frame Colour, Solid Black Front Castors
- Quick Release Self Propel Wheels, Solid Tyres and Aluminium Handrims
- Pelvic Belt, Adjustable both sides
- Anti Tippers



Standard Base price: £2,207

Please make your selection with a cross in the box Free of charge standard feature POA special request

Customer weight = KG

Code	Description	
BASE PRODUCT*		
BND9001	Rea Azalea Self Propel	£2,207
BND8000	SPECIAL REQUEST - PLEASE INCLUDE QUOTATION REFERENCE IN SPECIAL DESCRIPTION AT THE END OF THE FORM WHERE APPLICABLE	POA
PROPULSION*		
BND0025	Self Propel 24"/22" wheels with user brakes (NC 1076)	◆
BND0030	Self Propel 24"/22" wheels with assistant hub brakes (NC 1074)	£129

SEAT WIDTH*			
BND0040		390 mm/15" (Small chassis)	○
BND0050		440 mm/17" (Medium chassis)	○
BND0060		490 mm/19" (Large chassis)	○
SEAT WIDTH EXTENSION*			
BND0065		+ 50 mm seat width pre-adjustment on armrests and legrests	○
BND0085		Without seat width extension, set at 0mm adjustment but can be adjusted wider if required	○
SEAT DEPTH*			
BND0100		430-500 mm / 17" - 20"	◆
SEAT HEIGHT*			
BND0110		400 mm / 16" (NC 0370, 1015, 1025, 1062)	○
BND0120		450 mm / 18" (NC 1055)	○
SEAT ANGLE*			
BND0000		Standard Seat Angle +1° to +25°	◆
BND0127		Negative Seat Angle -1 To +23° (Std +1 To +25°) (NC 0360)	○
SEAT TYPE*			
BND0130		Depth Adjustable	◆
SEAT CUSHION*			
BND0165		Cushion Flosshape seat upholstery (NC 1111)	£52
SWO0190		Cushion Flosshape seat upholstery Shorter depth 450mm deep. Only use with seat depth in shortest position (NC 0191, 0202, 1111)	£52
BND0170		Cushion Vicair Multifunctional (Black) (NC 0191, 0201, 0202)	£559
BND0180		Without Seat Cushion, Seat Plate only (NC 0191, 0201, 0202)	○
SEAT UPHOLSTERY*			
BND0191		Black Plush TR35	○
BND0201		Black Multi-stretch Black PU coated fabric	◆
BND0202		Black Spacer TR36	○
BND1111		No upholstery, select with options 0170 Multifunctional and 0180 Without seat cushion	○
BACKREST WIDTH EXTENSION*			
BND0205		+ 50 mm Backrest width pre-adjustment, set at +50mm but can be adjusted if required (NC 0221, 0223, 0225, 0228)	○
BND0207		Without Backrest width pre-adjustment, set at 0mm but can be adjusted if required	○
BACKREST*			
BND0210		Flex 3 Backrest, Fixed Plate (Width Adjustable) (NC 0233, 0234, 0240, 0343, 0771, 0905)	◆
BND0220		Flex 3 Backrest Tension (Width Adjustable) (NC 0233, 0234, 0343, 0771, 0925)	£94
BND0221		Matrx Elite (includes hardware mounted to frame, backrest supplied in box) (NC 0240, 0250, 0260, 0270, 0280, 0290, 0311, 0321, 0322, 0335, 0800, 0810, 0820, 0830, 0880, 0890, 0895, 0896, 0900, 0905, 0910, 0920, 0925)	£437
BND0223		Matrx PB Elite Deep (includes hardware mounted to frame, backrest supplied in box) (NC 0240, 0250, 0260, 0270, 0280, 0290, 0311, 0321, 0322, 0335, 0800, 0810, 0820, 0830, 0880, 0890, 0895, 0896, 0900, 0905, 0910, 0920, 0925)	£485
BND0225		Matrx Posture Back (includes hardware mounted to frame, backrest supplied in box) (NC 0240, 0250, 0260, 0270, 0280, 0290, 0311, 0321, 0322, 0335, 0800, 0810, 0820, 0830, 0880, 0890, 0895, 0896, 0900, 0905, 0910, 0920, 0925)	£437
BND0228		Matrx Elite TR (includes hardware mounted to frame, backrest supplied in box) (NC 0240, 0250, 0260, 0270, 0280, 0290, 0311, 0321, 0322, 0335, 0800, 0820, 0880, 0895, 0896, 0900, 0905, 0910, 0920, 0925)	£361
BND0230		Backrest Tubes Without Backrest (NC 0233, 0234, 0240, 0250, 0260, 0270, 0280, 0290, 0311, 0321, 0322, 0335, 0345, 0740, 0771, 0800, 0810, 0820, 0830, 0880, 0890, 0895, 0896, 0900, 0905, 0910, 0920, 0925)	○
MATRIX BACKREST HEIGHT			
BND0233		Matrx Backrest Height 410 mm / 16"	○
BND0234		Matrx Backrest Height 510 mm / 20"	○

BACKREST CUSHION*				
BND0240		Cover for tension adjustable upholstery	(NC 2222)	○
BND0250		Cushion, Laguna (Lateral Support)	(NC 0905, 2222)	○
BND0260		Cushion, Mistral 2 (Waist Support)	(NC 0322, 0905, 2222)	£133
BND0270		Cushion, Passad 2 (Shoulder Support)	(NC 0905, 2222)	◆
BND0280		Cushion, Shoulder High (Black PU coated fabric)	(NC 0311, 0322, 0905, 2222)	£138
BND0290		Cushion, Vicair Multifunctional (Black)	(NC 0311, 0321, 0322, 0905)	£537
BND0300		Without Backrest Cover/ Cushion	(NC 0311, 0321, 0322, 0905)	○
BACKREST UPHOLSTERY*				
BND0311		Black Plush TR35		○
BND0321		Black Multi-stretch Black PU coated fabric		◆
BND0322		Black Spacer TR36		○
BND2222		No upholstery, select with Matrix, Multifunctional and Without backrest options		○
PUSH HANDLES*				
BND0335		Push Handles connected,no head/neckrest attachment		◆
BND0343		Push Handles connected,with Rea head/neckrest attachment	(NC 0771)	£66
BND0345		Push Bar, angle adjustable		£94
SEAT TILT AND BACKREST REGULATION*				
BND0350		Assistant operated gas piston seat tilt and backrest recline	(NC 5555)	◆
BND0360		Electrical operated seat tilt and backrest recline	(NC 1263)	£673
TILT LOCK*				
BND1263		Tilt Lock (locking of pre-set angle)		◆
BND5555		No Tilt Lock, Select with 0360 Electric operated seat and backrest		○
CENTRAL LEGREST				
BND0370		Central Legrest, includes calfpad and footrest. Also select options 0420, 0470, 3333, 4444	(NC 0376, 0378, 0380, 0390, 0400, 0410, 0430, 0440, 0450, 0460, 0475, 0476, 0477, 0720, 0930, 0940, 0941, 0942, 1290)	£279
LEGREST HANGER ATTACHMENT*				
BND0375		Standard position for legrest mounting		◆
SWO0189		Standard position with 0-30° Legrest Abduction function	(NC 0370)	£92
BND0376		High position 50mm higher for leg length compensation		○
BND0378		90° Hanger attachment to achieve a 90° legrest position on the Elevating Legrest	(NC 0390, 0400, 0410, 0440, 0450, 0460)	○
LEGREST HANGER LH*				
BND0380		Elevating Legrest, angle adjustable 0 to 75° swing in/swing out with knee pad and calf support LH	(NC 0475, 0480, 0490, 0942, 0945)	◆
BND0390		Legrest Hanger fixed 80 deg, swing in/swing out LH		○
BND0400		Legrest Hanger fixed 90 deg, swing in/swing out LH	(NC 0480, 0520, 0930)	○
BND0410		Legrest Amputee LH	(NC 0475, 0476, 0477, 0480, 0490, 0520, 0930, 0941, 0942, 0945)	£126
BND0420		Without Legrest, LH	(NC 0475, 0476, 0477, 0480, 0490, 0520, 0930, 0941, 0942, 0945)	○
LEGREST HANGER RH*				
BND0430		Elevating Legrest, angle adjustable 0 to 75° swing In/swing out with knee pad and calf support RH	(NC 0475, 0540, 0550, 0942, 0945)	◆
BND0440		Legrest Hanger fixed 80 deg, swing in/swing out RH		○
BND0450		Legrest Hanger fixed 90 deg, swing in/swing out RH	(NC 0540, 0580, 0940)	○
BND0460		Legrest Amputee RH	(NC 0475, 0476, 0477, 0540, 0550, 0580, 0940, 0941, 0942, 0945)	£126
BND0470		Without Legrest, RH	(NC 0475, 0476, 0477, 0540, 0550, 0580, 0940, 0941, 0942, 0945)	○

FOOTREST LH*			
BND0475		Footrest one piece, depth and angle adjustable short tube (NC 0477, 0480, 0490, 0520, 0540, 0550, 0580, 0930, 0940)	£109
BND0476		Footrest one piece, depth and angle adjustable long tube (NC 0477, 0480, 0490, 0520, 0540, 0550, 0580, 0930, 0940)	£109
BND0477		Footbox, angle adjustable, padded with short tube (NC 0480, 0490, 0520, 0540, 0550, 0580, 0930, 0940, 0941, 0942, 0945)	£501
BND0480		Footrest lower part, fixed with short tube and heel strap, LH (NC 0930, 0941)	○
BND0490		Footrest lower part, depth and angle adjustable with short tube and heel strap, LH	○
BND0520		Footrest lower part, depth and angle adjustable with long tube and heel strap, LH	○
BND3333		Without footrest, only select with option 0370 Central legrest, 0410 Amputee legrest and 0420 Without legrest	○
FOOTREST RH*			
BND0540		Footrest lower part, fixed with short tube and heel strap, RH (NC 0940, 0941)	○
BND0550		Footrest lower part, depth and angle adjustable with short tube and heel strap, RH	○
BND0580		Footrest lower part, depth and angle adjustable with long tube and heel strap, RH	○
BND4444		Without footrest, only select with options 0370 Central legrest, 0460 Amputee legrest, 0470 Without legrest, 0475 and 0476 One piece footrests or 0477 Footbox	○
ARMREST LH*			
BND0600		Armrest Ducktail Pad LH, height adjustable with autolock (NC 1232)	◆
BND0607		Armrests Reinforced Wide Pad, LH, height adjustable with autolock (NC 1230, 1231, 1232, 1250, 1260)	£82
BND0610		Armrest Padded LH, height adjustable with autolock (NC 1230, 1231, 1232, 1250, 1260)	£59
BND0615		Armrest Hemiplegic LH, height adjustable with autolock (NC 0695, 1230, 1231, 1232, 1250, 1260)	£165
BND0620		Without Armrest LH Select when choosing BND1232 Half Tray LH, Narrow Armrest is included (NC 0695, 1230, 1231, 1250, 1260)	○
ARMREST RH*			
BND0630		Armrest Ducktail Pad RH, height adjustable with autolock (NC 1233)	◆
BND0637		Armrest Reinforced Wide Pad, RH, height adjustable with autolock (NC 1230, 1231, 1233, 1250, 1260)	£82
BND0640		Armrest Padded RH, height adjustable with autolock. (NC 1230, 1231, 1233, 1250, 1260)	£59
BND0645		Armrest Hemiplegic, RH, height adjustable with autolock (NC 0695, 1230, 1231, 1233, 1240, 1245, 1250, 1260)	£165
BND0650		Without Armrest RH Select when choosing BND1233 Half Tray RH, Narrow Armrest is included (NC 0695, 1230, 1231, 1240, 1245, 1250, 1260)	○
ARMREST OPTIONS. Side support cover 0680, must be selected if 0690 Foam for side support is required			
BND0680		Side Support Cover, removable with tools (wipeable cover for the sideguard) (NC 0695)	£47
BND0690		Foam for side support cover (Internal cushioned pad, which is only possible in combination with 0680) (NC 0695)	○
BND0695		Pad for sideguard, removable with hook and loop fastening	£47
SEAT OPTIONS			
BND0710		Incontinence Cover	£37
BND0720		Pommel Abduction Cushion (NC 0722)	£109
BND0722		2 cm Increased seat height plate	£59
BND0725		Tilt Scale for seat tilt	£41

BACKREST OPTIONS				
BND0730		Headrest, Sideway Adjustable	(NC 0750, 0770)	◆
BND0750		Headrest with cheek support, sideway adjustable	(NC 0740, 0770, 0780)	£99
BND0740		Cover for headrest		£28
BND0770		Neckrest	(NC 0771)	○
BND0771		Headrest Matrx, 250 mm Soft Pad	(NC 0780, 0790)	£171
BND0780		Cover For Neckrest		£20
BND0782		Alternative headrest/neckrest mounting bracket		£21
BND0790		Handles Clinic Kit for Head-/Neckrest (Levers which replace screws for easy adjustment)		○
BND0900		Trunk Support Attachment to mount laterals later. If trunk supports are chosen, the attachment bracket is already included		○
BND0905		Lateral Positioning Pads. Attaches with hook and loop between the strap and the cover of the tension adjustable upholstery		£56
BND0910		Lumbar Support Cushion		£28
BND0915		Tilt Scale for backrest recline		£41
BND0920		Wedges (One Pair)		£9
BND0925		Privacy Attachment (Material connecting backrest and seat)		£31
TRUNK SUPPORT LH				
BND0800		Trunk Support Multi Adjustable with standard arm LH	(NC 0900)	○
BND0820		Trunk Support Multi Adjustable with extended height arm LH	(NC 0900)	○
BND0880		Trunk Support Swingaway LH	(NC 0900)	£97
BND0895		Trunk Support Multi Adjustable Swingaway LH	(NC 0900)	£104
TRUNK SUPPORT RH				
BND0810		Trunk Support Multi Adjustable with standard arm RH	(NC 0900)	○
BND0830		Trunk Support Multi Adjustable with extended height arm RH	(NC 0900)	○
BND0890		Trunk Support Swingaway RH	(NC 0900)	£97
BND0896		Trunk Support Multi Adjustable Swingaway RH	(NC 0900)	£104
LEGREST OPTIONS				
BND0930		Footplate Extender LH	(NC 0941)	£31
BND0940		Footplate Extender RH	(NC 0941)	£31
BND0941		Cover Footplate, Black Plush Fabric		£49
BND0942		Calf Strap Padded		£67
BND0945		Calf Strap (for One piece footrest)		○
CHASSIS COLOUR*				
BND0944		Silver Pearl		◆
BND0946		Azurite (Dark Blue)		○
BND0947		Sand (Beige)		○
BND0948		Electric Green		○
BND0949		Happy Red		○
CHASSIS SAFETY				
BND0950		Chassis extender, castors extend 50mm forward for added stability		○
CASTOR ATTACHMENT*				
BND0960		Fixed castor attachment		◆

FRONT CASTORS*			
BND0992		150x30 mm (6" x 1 1/4") Medium Hard Black tyre	○
BND1015		200x32 mm (8" x 1 1/4") Black Solid Tyre	£19
BND1025		200x45 mm (8" X 1 3/4") Black Solid tyre	◆
REAR WHEEL SIZE*			
BND1055		22" / 560 mm Self Propel rear wheel with QR attachment	○
BND1062		24" / 600 mm Self Propel rear wheel with QR attachment	○
REAR WHEEL TYPE*			
BND1074		Standard Self Propel with colourless hub, clear anodized rim and Stainless steel spokes	◆
BND1076		Hub Brake 22" or 24" Self Propel with colourless hub, clear anodized rim and Stainless steel spokes	○
REAR WHEEL TYRE*			
BND1082		Pneumatic Profile Black 1 3/8" tyres	○
BND1084		Pneumatic Puncture proof / flatfree 1 3/8" Black tyres	£57
BND1092		Solid Black tyres	◆
HANDRIM*			
BND1100		Aluminium handrim	◆
BND1120		Plastic coated handrim	£45
SPOKEGUARD			
BND1160		Frost spokeguard	£153
BRAKE*			
BND1170		User manoeuvred tyre brakes	◆
BND1180		One Arm Brake LH	£116
BND1185		One Arm Brake RH	£116
PELVIC BELT*			
BND1200		Pelvic belt	◆
TRANSPORTATION BRACKETS*			
BND1225		Transportation Brackets, mounted to the front of the wheelchair	◆
SAFETY OPTIONS			
BND1220		Reflector Kit (wheels and chassis)	£20

OTHER OPTIONS

BND1230		Fixed Table Tray without elbow pads includes mounting kit (NC 1237, 1260)	£127
BND1231		Fixed Table Tray with elbow Pads includes mounting kit (NC 1237, 1260)	£127
BND1232		Half table tray, clear plastic, flip up, depth adjustable, left side (armrest with straight narrow pad included) (NC 1235, 1260)	£95
BND1233		Half table tray, clear, plastic, flip up, depth adjustable, right side (armrest with straight narrow pad included) (NC 1235, 1260)	£95
BND1240		Swing away table without elbow pads (NC 1235, 1250, 1260)	£140
BND1245		Swing away table with elbow pads (NC 1235, 1260)	£140
BND1235		Table tray padded cushion	£47
BND1237		Half table tray padded cushion	£43
BND1250		Locking device for table tray	£13
BND1260		Mounting kit for table (if ordered separately)	£7
BND1270		Pump	£20
BND1280		Cane Holder	£46
BND1290		Handles Clinic Kit for Head-/Neckrest (Levers which replace screws for easy adjustment)	£50
BND1300		Harness attachment , must be selected in combination with Head/Neckrest	£79
BND1310		Attachment for positioning belt	£29

TOTAL PRICE

Non Conforming Options:

BND0191	Black Plush TR35
BND0201	Black Multi-stretch Black PU coated fabric
BND0202	Black Spacer TR36
BND1111	No upholstery, select with options 0170 Multifunctional and 0180 Without seat cushion
BND0221	Matrx Elite (includes hardware mounted to frame, backrest supplied in box)
BND0223	Matrx PB Elite Deep (includes hardware mounted to frame, backrest supplied in box)
BND0225	Matrx Posture Back (includes hardware mounted to frame, backrest supplied in box)
BND0228	Matrx Elite TR (includes hardware mounted to frame, backrest supplied in box)
BND0233	Matrx Backrest Height 410 mm / 16"
BND0234	Matrx Backrest Height 510 mm / 20"
BND0240	Cover for tension adjustable upholstery
BND0250	Cushion, Laguna (Lateral Support)
BND0260	Cushion, Mistral 2 (Waist Support)
BND0270	Cushion, Passad 2 (Shoulder Support)
BND0280	Cushion, Shoulder High (Black PU coated fabric)
BND0290	Cushion, Vicair Multifunctional (Black)
BND0311	Black Plush TR35
BND0321	Black Multi-stretch Black PU coated fabric
BND0322	Black Spacer TR36
BND2222	No upholstery, select with Matrx, Multifunctional and Without backrest options
BND0335	Push Handles connected,no head/neckrest attachment
BND0343	Push Handles connected,with Rea head/neckrest attachment
BND0345	Push Bar, angle adjustable
BND0360	Electrical operated seat tilt and backrest recline
BND1263	Tilt Lock (locking of pre-set angle)
BND5555	No Tilt Lock, Select with 0360 Electric operated seat and backrest
BND0370	Central Legrest, includes calfpad and footrest. Also select options 0420, 0470, 3333, 4444
BND0376	High position 50mm higher for leg length compensation
BND0378	90° Hanger attachment to achieve a 90° legrest position on the Elevating Legrest
BND0380	Elevating Legrest, angle adjustable 0 to 75° swing in/swing out with knee pad and calf support LH
BND0390	Legrest Hanger fixed 80 deg, swing in/swing out LH
BND0400	Legrest Hanger fixed 90 deg, swing in/swing out LH
BND0410	Legrest Amputee LH
BND0430	Elevating Legrest, angle adjustable 0 to 75° swing In/swing out with knee pad and calf support RH
BND0440	Legrest Hanger fixed 80 deg, swing in/swing out RH
BND0450	Legrest Hanger fixed 90 deg, swing in/swing out RH
BND0460	Legrest Amputee RH
BND0475	Footrest one piece, depth and angle adjustable short tube
BND0476	Footrest one piece, depth and angle adjustable long tube
BND0477	Footbox, angle adjustable, padded with short tube
BND0480	Footrest lower part, fixed with short tube and heel strap, LH
BND0490	Footrest lower part, depth and angle adjustable with short tube and heel strap, LH
BND0520	Footrest lower part, depth and angle adjustable with long tube and heel strap, LH

BND0540	Footrest lower part, fixed with short tube and heel strap, RH
BND0550	Footrest lower part, depth and angle adjustable with short tube and heel strap, RH
BND0580	Footrest lower part, depth and angle adjustable with long tube and heel strap, RH
BND0695	Pad for sideguard, removable with hook and loop fastening
BND0720	Pommel Abduction Cushion
BND0722	2 cm Increased seat height plate
BND0750	Headrest with cheek support, sideway adjustable
BND0740	Cover for headrest
BND0770	Neckrest
BND0771	Headrest Matrix, 250 mm Soft Pad
BND0780	Cover For Neckrest
BND0790	Handles Clinic Kit for Head-/Neckrest (Levers which replace screws for easy adjustment)
BND0900	Trunk Support Attachment to mount laterals later. If trunk supports are chosen, the attachment bracket is already included
BND0905	Lateral Positioning Pads. Attaches with hook and loop between the strap and the cover of the tension adjustable upholstery
BND0910	Lumbar Support Cushion
BND0920	Wedges (One Pair)
BND0925	Privacy Attachment (Material connecting backrest and seat)
BND0800	Trunk Support Multi Adjustable with standard arm LH
BND0820	Trunk Support Multi Adjustable with extended height arm LH
BND0880	Trunk Support Swingaway LH
BND0895	Trunk Support Multi Adjustable Swingaway LH
BND0810	Trunk Support Multi Adjustable with standard arm RH
BND0830	Trunk Support Multi Adjustable with extended height arm RH
BND0890	Trunk Support Swingaway RH
BND0896	Trunk Support Multi Adjustable Swingaway RH
BND0930	Footplate Extender LH
BND0940	Footplate Extender RH
BND0941	Cover Footplate, Black Plush Fabric
BND0942	Calf Strap Padded
BND0945	Calf Strap (for One piece footrest)
BND1015	200x32 mm (8" x 1 1/4") Black Solid Tyre
BND1025	200x45 mm (8" X 1 3/4") Black Solid tyre
BND1055	22" / 560 mm Self Propel rear wheel with QR attachment
BND1062	24" / 600 mm Self Propel rear wheel with QR attachment
BND1074	Standard Self Propel with colourless hub, clear anodized rim and Stainless steel spokes
BND1076	Hub Brake 22" or 24" Self Propel with colourless hub, clear anodized rim and Stainless steel spokes
BND1230	Fixed Table Tray without elbow pads includes mounting kit
BND1231	Fixed Table Tray with elbow Pads includes mounting kit
BND1232	Half table tray, clear plastic, flip up, depth adjustable, left side (armrest with straight narrow pad included)
BND1233	Half table tray, clear, plastic, flip up, depth adjustable, right side (armrest with straight narrow pad included)
BND1240	Swing away table without elbow pads
BND1245	Swing away table with elbow pads
BND1235	Table tray padded cushion
BND1237	Half table tray padded cushion

BND1250	Locking device for table tray
BND1260	Mounting kit for table (if ordered separately)
BND1290	Handles Clinic Kit for Head-/Neckrest (Levers which replace screws for easy adjustment)



For a more complex customer requirement please refer to our Invacare Unique special order form. email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: _____ Email: _____

Phone: _____ Fax: _____

Description of special request (please include relevant option code or part number):